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Aqure

Sterile Spüllösung

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Aqure

Sterile Irrigating Solution

Sterile, physiological irrigation solution for ophthalmic procedures

Composition

1 ml solution contains:

Sodium chloride	6.40 mg
Potassium chloride	0.75 mg
Calcium chloride 2 H ₂ O	0.48 mg
Magnesium chloride 6 H ₂ O	0.30 mg
Sodium acetate 3 H ₂ O	3.90 mg
Sodium citrate 2 H ₂ O	1.70 mg

Excipients:

Water for injection, hydrochloric acid (for pH adjustment)

No buffered solution (according to EN ISO 16671)

Intended purpose

For the irrigation of the anterior chamber during cataract surgery and for other intraocular or extraocular procedures, such as therapeutic or diagnostic measures, where irrigation or retention of moisture is required.

Side-effects

Where there is pre-existing pathological damage to the corneal endothelium, any intraocular intervention – including intraocular irrigation – can lead to bullous keratopathy which is characterized by corneal oedema or clouding.

After the surgical procedure, inflammatory processes, including endophthalmitis, can occur. No connection with the use of Aqure has yet been demonstrated.

Dosage recommendations

The dosage is determined by the particular circumstances dictated by the operative, therapeutic or diagnostic procedure and should be decided by the doctor on a case-by-case basis.

Method of administration

To irrigate the eye for the purpose of surgical procedures as well as diagnostic and therapeutic measures that require irrigation or corneal wetting.

An irrigation set 'P' (pressure-resistant) with built-in air vent should be used. Bottles do not have added venting, which is why an irrigation set with air vent must be used. The instruction for use from the manufacturer of the irrigation set must be always followed. The irrigation set should be fitted onto the bottle as follows (aseptic conditions):

1. If applicable, install the provided plastic hook over the glass bottle.

2. Remove the protective cap from the bottle.
3. The rubber-stopper is secured by an aluminum collar which if desired can be removed.
4. Disinfect the area of the rubber-stopper.
5. Insert the puncture-spike of the irrigation set aseptically through the target area into the rubber-stopper.
6. Squeeze the drip chamber to allow the fluid to flow.

Note: The fluid should already be flowing before use, in order to ensure that there is no air in the tubing.

The period of use is limited by the duration of the procedure. One bottle of Aqure is intended for one patient only.

References in the literature reveal that the osmolality of the aqueous humour is higher in diabetics and intraoperative lens changes can arise. Since Aqure is iso-osmotic with normal aqueous humour, it should be used with caution in diabetics undergoing a vitrectomy. In these cases, it is advisable to match the osmolality to that of the aqueous humour.

Special instructions

- Use only for eye irrigation and for eye surgery, not for injection or infusion.
- Only use solution if it is completely clear and no minute particles are visible.
- Before use, check that the container is intact. Do not use, if sterile packaging was damaged.
- Use the solution immediately after opening the container.
- Aqure contains no preservatives. Therefore, do not store opened containers for further use. The medical device is intended solely for single use and is not suitable for reprocessing.
- Do not use Aqure as solvent for other medicinal products.
- The expiry date relates to an unopened product in its closed and properly stored packaging. Do not use beyond the expiry date given on the packaging.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer/distributor and the competent authority of the country in which the user and/or patient is established.
- Keep out of the reach and sight of children.

Storage

No special requirements.

Disposal information

With household waste or recycling:

Instructions for use, traceability labels, Protective cap, aluminium collar, bottle (empty the bottle into a sink before disposing)

With contaminated hospital products:

Irrigation set

Presentation

REF 1198	Glass bottle with 500 ml,
REF 1215	PP bottle with 250 ml and
REF 1196	PP bottle with 500 ml

Not all pack sizes may be marketed

Date of revision of the text

2025-09

Manufacturer

MTP Medical Technologies
Kunstseidenstrasse 4
01796 Pirna, Germany

Explanation of the symbols

	Consult instructions for use
	Do not use if package is damaged
	Do not re-use
	Sterilized using steam
	Single sterile barrier system
	Article number
	Batch code
	Manufacturing date
	Use-by date
	Manufacturer
	Medical Device
	Unique Device Identifier

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