

IPURE

ASPHERIC  
MONOFOCAL OPTIC

# Cataract Surgery Streamlined

High contrast and high quality<sup>1</sup> vision.  
Seamless, predictable, and efficient  
preloaded IOL delivery.<sup>2</sup>

# Seamless. Predictable. Efficient.<sup>3</sup>

## Seamless preparation. Smooth implantation.

IPure is the first and only preloaded monofocal IOL available in both a 1-piece and 3-piece design.

## Enhanced predictability

You can feel confident in IPure's proven predictability.<sup>3</sup>



In a study of 600 IOL release bench tests<sup>4</sup>

**100%** normal lens releasing behavior was seen with IPure

In a human trial comparing pre-loaded and non-preloaded IOL delivery systems<sup>3</sup>

**0%** of IPure cases required additional IOL manipulations in the eye


**VS** up to **32%** of cases with manually loaded IOL delivery systems



# Enhanced workflow efficiency and productivity

## Streamline your operating procedure

Save valuable time in the OR by eliminating time-consuming steps before, during, and after surgery.<sup>3</sup>



# Trusted high quality design

## Uniquely patented Three-Zone aspheric optic design<sup>1</sup>

Neutral spherical aberration in the ~2.4 mm central optic zone (zones 1 and 2) and -0.18 microns spherical aberration in the peripheral optic zone.

- High contrast in various lighting conditions
- Maintains natural corneal depth of focus
- Less sensitive to natural off-axis conditions, lens decentration, and corneal aberration

## IPure's meticulous lens manufacturing process

IPure's individually lathe cut and pad polished lenses deliver:

- Smooth and even optic surface designed to reduce light scattering on lens surface<sup>5</sup>
- Sharp optic edge designed to reduce posterior capsule opacification (PCO)
- Blue PMMA chemically bonded haptic tips are highly visible inside the injector for controlled folding behavior
  - Reduced sticking of materials<sup>3</sup>
  - No leading haptic/trailing haptic failure was seen in 201 IOL release bench tests<sup>6</sup>

Neutral SA:  
~2.4 mm diameter  
Zone 1: (+SA)  
Zone 2: (-SA)

Zone 3:  
-0.18 microns SA

Sharp Squared-Optic Edge

Controlled Folding Behavior



# CustomEyes®

## Streamline cataract surgery from surgical prep to IOL implantation with BVI

Visitec®



Beaver®



NuVisc™ OVD\*



IPure® IOL\*



### IPure® enhances the BVI CustomEyes® offering with value in mind

You can now choose IPure for your surgical procedure, and let BVI CustomEyes be the right surgical pack for your team.

- Reliable and knowledgeable support
- Transparent pricing
- Simplified surgical prep
- Easy order process

CustomEyes allows you to personalize your surgical packs with high quality surgical consumables from the brands you trust, including Beaver® and Visitec®.



\*NuVisc and IPure shipped separately.

Ask your BVI sales representative about efficiencies in the OR.


# IPure Specifications

Trusted aspheric monofocal IOL in a fully pre-loaded, single-use IOL injector system

## Pre-loaded IOL System



## IPURE



## IPURE 3P

Contact Your  
BVI Sales  
Representative  
to Learn More

866-906-8080

Model Names	B1PC (Clear) B1PY (Yellow)	B3PC (Clear)
Optic Material	Hydrophobic Acrylic, UV Filter	Hydrophobic Acrylic, UV Filter
Optic Design	Aspheric -0.18 microns patented lens design	Aspheric -0.18 microns patented lens design
Manufacturing	Lathe-cut and pad polished	Lathe-cut and pad polished
Haptic Material	Hydrophobic acrylic with blue PMMA chemically bonded haptics tips	PMMA
Haptic Configuration	Modified-C loop, 5° angulation	Modified-C loop, 5° angulation
Dimension (Optic/OAL)	6.0 mm/12.5 mm	6.0 mm/12.5 mm
Power	+6.0 to +30.0 D (0.5 D steps)	+6.0 to +30.0 D (0.5 D steps)
A-Constant <sup>7</sup>	118.4	118.4
Injector Type	Single Use Preloaded IOL System	Single Use Fully Preloaded IOL System
Front injector tip outer diameter	<b>1.82 mm</b>	<b>1.89 mm</b>
Optimized Optical constants	SRK/T (A): 118.5 Holladay 1 (SF): 1.52 Barrett <sup>8</sup> (LF/A): 1.61/118.48 Hill RBF <sup>9</sup> (A): 118.58 Haigis: a0= -0.542, a1= 0.161 a2= 0.204 Hoffer-Q (pACD)= 5.30	SRK/T (A): 118.6 Holladay 1 (SF): 1.54 Hagis a0=-0.093 a1=-0.023 a2:0.208 Hoffer Q(pACD): 5.30

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

### References

1. US Patent NO: US8647383. 2. Data on file, BVI, 2019. 3. Chung B, et al. Preloaded and non-preloaded intraocular lens delivery system and characteristics: human and porcine eyes trial. *Int J Ophthalmol*. 2018;11:6-11. 4. Data on file, HOYA Medical Singapore Pte. Ltd, 2012. 5. Werner L. Glistening's and surface light scattering in intraocular lenses. *J Cataract Refract Surg*. 2010;36:1398-1420. 6. Data on file, HOYA Medical Singapore Pte. Ltd, 2015. 7. The A-Constant mentioned above is presented as a guideline only for lens power calculations. It is recommended that the A-Constant measurement be customized based on the surgeon's experience and measure equipment. 8. Barrett: [http://calc.apacrs.org/barrett\\_universal2105/](http://calc.apacrs.org/barrett_universal2105/). 9. Hill RBF: <https://rbfcalculator.com/lens-constants.html>

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2020



PN 1564403-02