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**Instructions for use**

The instruction for use mentioned here is applicable to the following product numbers:

**116001 TEST CHAMBER (2 PCS)**

**1. DISCLAIMER**  
This device has been designed, produced and sold for use as described in this instruction for use. Use that does not correspond to these instructions provided by OPTIKON could harm the patient and/or damage the device itself as well as any instruments or accessories connected to it. These instructions for the use of the device must be read along with the instruction for use for any instruments and/or accessories which will be used with it. In case of doubt and/or if further information is required concerning the use of this product and related warranty conditions, please consult your local distributor or local OPTIKON Customer Service before use. This product must not be modified without the manufacturer's approval.

**2. NOTES ON SAFETY**  
**⚠ WARNING !**  
Hazard which leads to harm that can cause death or serious injuries.  
**⚠ CAUTION !**  
Hazard which leads to harm that can cause injuries with medical treatment.  
**NOTE**  
Hazard which leads to harm that can cause injuries without medical treatment.

**3. AREA OF USE**  
**3.1. INTENDED USE**  
This re-usable device is an accessory for OPTIKON phaco systems. It has been designed to simulate the fluid dynamic conditions in the eye during priming process before the phacoemulsification surgery.  
**3.2. INDICATION**  
The test chamber simulates the anterior chamber of the human eye. It is used during the priming of the phacoemulsification unit, in order to guarantee continuity of the irrigation/aspiration system circuit. Furthermore it is used to verify that the selected I/A parameters will prevent a collapse of the anterior chamber.  
**3.3. CONTRAINDICATION**  
There are no known contraindications for the product when used in accordance with its intended use.

**4. MESSAGE TO MANUFACTURER AND AUTHORITIES**  
If a serious incident occurs in relation to this medical device, the user and / or patient must notify their Competent Authority and the manufacturer and of this incident.

**5. DEVICE COMPOSITION**  
This device consists of:  
• connection tube to connect the test chamber with the silicone sleeve,  
• test chamber for simulation of the eye

**6. CONNECTION TO OTHER DEVICES**  
This product can be used only in combination with OPTIKON silicone sleeves.

**7. INSTRUCTIONS FOR ASSEMBLY AND USE**  
**NOTE**  
• This device is designed to be used at medical centres where eye surgery is performed. Its use is restricted to qualified medical personnel only. Follow the instructions shown here taking care not to compromise sterility in the operating area.  
**⚠ CAUTION !**  
Risks due to the user's level of experience  
Outcomes with use of the device can be affected by the user's level of experience.  
• The user should be adequately trained in handling the device prior to its use.

**⚠ CAUTION !**  
Risk due to non-sterile device.  
• Please follow the instructions for cleaning, disinfecting and sterilization in the document "Instruction for reprocessing" prior to first use of the device.  
**NOTE**  
• Before and during use, avoid contact with metal objects, particularly if they are sharp and/or pointed.  
• Place the test chamber always on the sleeve so that a perfect seal is obtained.  
• Do not push the test chamber to the bottom of the sleeve in order to avoid difficulty in removing it.  
• Do not use the equipment if it does not pass the priming test or the Irrigation/Aspiration balance test. Carefully check the test chamber for damage and make sure it is positioned correctly.  
• Dispose the test chamber if it shows any sign of wear or damage.

**7.1. Fill the test chamber with sterile saline solution.**  
**7.2. Place the test chamber on the sleeve and push gently until it is perfectly sealed.**  
Follow the instructions of the relevant phaco system for the priming phase and for the irrigation aspiration test to check the fluidic balance of the anterior chamber.

**8. INSTRUCTIONS FOR REPROCESSING**  
**CLEANING**  
It is the task of the healthcare organization to set up procedures based on recommendations to reduce the risk of personal exposure to pathogens and to prevent damage to the device, published by professional organizations. Cleaning of reusable devices must occur immediately after they are used. To clean the product, do the following:  
(in the cleaning department)

1. Immerse the device into a lukewarm neutral enzymatic detergent. Follow the instructions for use of the detergent manufacturer.
2. Carefully rinse the device using lukewarm water until complete removal of the detergent.

To avoid damage to device, do not use sharp tools while cleaning the device. Cleaning beyond expected life can compromise the integrity and functional features of the device.  
**DISINFECTION AND DECONTAMINATION**  
**Disinfection**  
To minimize the risk of professional exposure to pathogenic pathogens via the bloodstream, the product may be treated by immersion in peracetic acid disinfectant solutions.  
**Decontamination**  
If the product is used with TSE carriers (transmissible spongiform encephalopathy), it must be decontaminated by sterilization with a pre-vacuum cycle at 134 °C for 18 minutes. Disinfection/decontamination beyond expected life can compromise the integrity and functional features of the device.

**STERILIZATION**  
The device must be steam sterilized with a pre-vacuum cycle (ISO 17665-1) at 134 °C for at least 4 minutes. Once the device, whether or not sterilized, has been cleaned as described above must be placed in a suitable sealed envelope, which is a biological barrier, so as to ensure that sterility is maintained after the sterilization process.  
After sterilization let the product cool down into the environment.  
Do not immerse the device in liquids of any kind to accelerate the cooling.  
Do not immerse the device in liquids based on sodium hypochlorite or sodium hydroxide.  
Prevent contact of the device with metallic objects during the sterilization cycle.

**NOTE**  
• Always follow the instructions for cleaning, disinfecting and sterilization in the document "Instruction for reprocessing" as a deviant process could result in a decreased biocompatibility of the device.

**9. ADDITIONAL INFORMATION FOR REPROCESSING**  
• Machine cleaning and disinfecting is not allowed.

**10. AMBIENT CONDITIONS FOR OPERATION AND TRANSPORT**

Allowed temperature in operation condition:	from	+10°C	to	+35°C
Allowed temperature in transport condition:	from	-20°C	to	+60°C
Relative humidity in transport condition:	from	10%	to	95%

The package that contains the device must be stored in a dust-free environment, far from any heat sources, and it must not be exposed to direct sunlight.

**11. DISPOSAL**  
**⚠ CAUTION !**  
• Clean, disinfect and sterilize the device completely prior to disposal.  
• For disposal, also observe the legal provisions applicable in your country and the hygiene regulation of the doctor's office or hospital.

**12. STORING**  
Handle the test chamber with care, being careful not to let them come into contact with sharp or pointed objects. After being used for the first time, the test chamber may be stored in a sterilization container.

**13. EXPECTED LIFETIME**  
The test chamber is a very delicate object, replace it immediately if it shows any sign of damage. The test chamber must be replaced after 25 reprocessing cycles even if it does not show any sign of wear or damage.

**⚠ CAUTION !**  
• Any use beyond the given reprocessing cycles or the use of damaged and/or contaminated components is the sole responsibility of the user. Failure to observe these instructions nullifies liability of the manufacturer.

**14. SYMBOLOLOGY**

Symbol	Description
	Manufacturer (name and address of manufacturer and, date of manufacture, if applicable)
	Catalog number
	Batch code
	Keep dry
	Keep away from sunlight
	Consult instructions for use
	Non sterile
	CE Mark
	Medical Device
	2D Matrix UDI (01) UDI-DI, (10) lot number, (11) production date, (241) customer part number

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