



BVI

POD Platform

Double C-loop Design

**Have it
all**



POD Platform

bvimedical.com

Surgeons often have to weigh priorities when choosing an IOL platform: rotational stability versus ease of alignment, innovation versus proven performance. In addition to these decisions is the expectation of reliability and consistency.


Developed in collaboration with Dr. Christophe Chassain (France), BVI's unique open double C-loop haptic design is **engineered to meet all of these requirements** — resulting in a platform that delivers:

1. Optimal **rotational stability** with intraoperative maneuverability
2. Consistent **refractive performance**
3. High **refractive accuracy**

With over a decade of global adoption, the POD platform now serves as the foundation for all BVI IOLs and is **available with all BVI optic technologies**.


1. Optimal rotational stability

A recent systematic review including POD platform data from over 1,400 eyes concluded that it provides good rotational stability¹. Furthermore, a selection of peer-reviewed studies published between 2016 and 2023 — involving 610 eyes — consistently demonstrated an average rotation of 2° or less across a range of optic technologies:



1.22°
average rotation for
the PODEYE Toric lens.

n = 94 eyes




CLINICAL VALIDATION

«In our study, 97.87% of eyes had a rotation of <10 degrees with a mean absolute rotation of 1.22 ± 2.21 degrees. We consider that this rotation is not clinically relevant and, based on these values, we conclude that this lens showed a high degree of rotational stability when implanted in the capsular bag.⁸»

And ease of alignment


In addition to post-operative rotational stability, the double C-loop POD platform can be rotated **clockwise or counterclockwise**⁹ to align the toric IOL with the intended axis. In contrast, classic C-loop IOLs can only be rotated clockwise and require additional steps in the event of misalignment.¹⁰ As a result, the POD platform is recognized for its **intraoperative control**,¹¹ made possible by its symmetrical design.



Have it all: innovation and performance


2. Consistent refractive performance

A selection of peer-reviewed studies has demonstrated that the POD platform maintained excellent refractive performance across a range of optical technologies:



3. High refractive accuracy

A selection of peer-reviewed studies confirmed that the POD platform delivered long-term refractive accuracy across toric and non-toric models:




Over 90% of patients achieved a postoperative refraction within ± 0.50 D of target (both cylinder and spherical equivalent).

Distinct by Design

The POD platform delivers consistent performance across key criteria, thanks to its **symmetrical, open double C-loop** haptic architecture. This unique configuration is designed to conform gently and smoothly to the contours of the capsular bag.

Its four haptics provide a **high contact angle** — greater than standard C-loop IOLs.^{18,19} This increased haptic-capsule interface helps distribute forces across a broader surface and is designed to reduce stress on the capsule.



Contact Angle

Double C-loop Solutions



Our premium POD platform supports a wide range of optic technologies.



FINEVISION HP

TRIFOCAL OPTIC

TORIC

FINEVISION HP

TRIFOCAL OPTIC



ISOPURE SERENITY

PREMIUM MONOFOCAL IOL

TORIC

ISOPURE SERENITY

PREMIUM MONOFOCAL IOL



PODEYE

MONOFOCAL OPTIC

TORIC



PODEYE

MONOFOCAL OPTIC

References:

1. Daya S, Chassain C, Pagnoulle C. Double C-Loop Haptic Lenses Are Rotationally Stable – A Systematic Review. *Clinical Ophthalmology* 2025;19:2259–2269. 2. Poyales F, Garzon N, Rozema JJ, Romero C, de Zarate BO. Stability of a Novel Intraocular Lens Design: Comparison of Two Trifocal Lenses. *J Refract Surg.* 2016;32(6):394–402. | 3. Draschl P, Hirnschall N, Luft N, Schuschnitz S, Wiesinger J, Rigal K, Findl O. Rotational stability of 2 intraocular lenses with an identical design and different materials. *J Cataract Refract Surg.* 2017;43(2):234–238. | 4. Chassain C, Chamard C. Posterior capsule opacification, glistenings and visual outcomes: 3 years after implantation of a new hydrophobic IOL. *J Fr Ophtalmol.* 2018;41(6):513–520. | 5. Ribeiro FJ, Ferreira TB, Relha C, Esteves C, Gaspar S. Predictability of different calculators in the minimization of postoperative astigmatism after implantation of a toric intraocular lens. *Clin Ophthalmol.* 2019;13:1649–1656. | 6. Ribeiro FJ, Ferreira TB. Comparison of visual and refractive outcomes of 2 trifocal intraocular lenses. *J Cataract Refract Surg.* 2020;46(5):694–699. | 7. Ang RET. Long-term trifocal toric intraocular lens outcomes in Asian eyes after cataract surgery. *J Cataract Refract Surg.* 2023;49(8):832–839. | 8. Ang RET, Tana-Rivero P, Pastor-Pascual F, Stodulka P, Tetz M, Fischinger I. Visual and Refractive Outcomes After Bilateral Implantation of a Biconvex Aspheric Toric Monofocal Intraocular with a Double C-Loop Haptic Design. *Clin Ophthalmol.* 2023;17:2765–2776. | 9. Ang R, Tana-Rivero P, Pastor-Pascual F, Stodulka P, Slovak M, Tetz M, Fischinger I, Cazal J, Gessa-Sorroche M, Ibarz-Barbera M, Holland D, Groneberg T. One-Year Clinical Outcomes Following Aspheric Toric Monofocal IOL with a Double C-Loop Haptic Design Implantation. *Clinical Ophthalmology* 2024;18 3515–3525. | 10. Torio et al. Comparison of the Rotational Stability of Different Toric Intraocular Lens Implants. *Philipp J Ophthalmol* 2014;39:67–72. | 11. Ang RET. “PODEYE Toric Clinical Outcomes.” Presentation, BVI Advisory Board meeting, Boston 2024. | 12. Poyales F, Perez R, Lopez-Brea I, Zhou Y, Rico L, Garzon N. Comparison of Visual Performance and Patient Satisfaction Outcomes with Two Trifocal IOLs with Similar Optical Design but Different Materials. *Clin Ophthalmol.* 2020;14:3237–3247. | 13. Daya S, Espinosa Lagana M. Visual and Refractive Outcomes of a New Hydrophobic Trifocal Toric Intraocular Lens. *Clinical Ophthalmology* 2024;18 997–1007 | 14. Khoramnia R, Kretz F, Gerl M, Breyer D, Auffarth G. Long-term Clinical Outcomes After Bilateral Implantation of Two Trifocal Diffractive IOLs. *J Refract Surg.* 2023;39(12):798–807. | 15. Ang RET. Long Term Clinical Outcomes of Hydrophilic and Hydrophobic Versions of a Trifocal IOL with the Same Optical Design. *Clinical Ophthalmology* 2023;17:623–632. | 16. Chassain, M.K. Hallak, M. Lesaffre Rotational stability and clinical outcomes after implantation of a new monofocal toric intraocular lens with double C-loop design. *Journal Francais d'Ophtalmologie* 2022 0181-5512. | 17. Akahoshi. Hydrophobic Trifocal Toric Intraocular Lens Outcomes in Japanese Eyes After Cataract Surgery. *Clinical Ophthalmology* 2024, 18:2021–2031. | 18. Schickhardt SK, Labuz G, Munro DJ, Liebenwirth I, Zhang L, Fang H, Auffarth GU. In-vitro assessment of a novel intraocular lens made of crosslinked polyisobutylene. *Journal of the Mechanical Behavior of Biomedical Materials* 152 (2024). | 19. REP_503_1_2022_15.2 PODIGF Mechanical specifications.

Please check the availability of products in your market with your sales representative.

Note: The intraocular lenses are not FDA approved.

Contact Information:

www.bvimedical.com/customer-support/

