

BVI

Fluid Control



PVA Eye Spears

Spear overall handle length is 7.16 cm, or 2.82 in.

UltraClean™ LASIK PVA Eye Spears

Compressed shape, designed for fast expanding. Made from highly absorbent, virtually lint free PVA material. Designed specifically for rapid recovery of fluids during LASIK and other corneal procedures.



Qty: 100/Bx
Brand: UltraClean™

Description	SKU #
UltraClean™ Eye Spears	40810
UltraClean™ Eye Spears	40815

Visitec® LASIK PVA Expanded Eye Spear

Made from highly absorbent PVA material. Pre-expanded for convenience. Designed to provide instantaneous wicking. No need to pre-soak and fan.



Qty: 40/Bx
Brand: Visitec®

Description	SKU #
LASIK Expanded Spear	581709

Visitec® LASIK PVA Eye Spears

Made from highly absorbent PVA material. Designed for optimal tissue manipulation and controlled fluid absorption to facilitate flap alignment and adhesion.



Qty: 100/Bx
Brand: Visitec®

Description	SKU #
LASIK Eye Spear - Orange	581761

Weck-Cel® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Fast wicking action delivers quick absorption of blood and fluid from the surgical area. The flat tip provides precise application at the surgical incision.



Qty: Various/Bx
Brand: Weck-Cel®

Description	SKU #
Weck-Cel® Spears, 6 Pack, 144/Bx	0008680
Weck-Cel® Spears, 10 Pack, 180/Bx	0008685

Merocel® PVA Eye Spears

Made from exclusive Merocel® polyvinyl acetal (PVA) material. The Merocel® material is designed to be highly absorbent, fast wicking and virtually lint free.



Qty: Various/Bx
Brand: Merocel®

Description	SKU #
Merocel® Eye Spear, 10 Spears/Pouch, 180/Bx	400101
Merocel® Eye Spear, 5 Spears/Pouch, 50/Bx	405101

Visitec® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Fast wicking action delivers quick absorption of fluid from the surgical area. The flat tip provides precise application at the surgical incision. Magenta handle.



Qty: 100/Bx
Brand: Visitec®

Description	SKU #
Sponge Eye, 7 cm (2.75 in) overall length.	581089

Sponge Points, Strips and Wicks

Merocel® Eye Wicks

Made from exclusive Merocel® polyvinyl acetal (PVA) material. The Merocel® material is highly absorbent, fast wicking and virtually lint free.



Qty: Various/Bx
Brand: Merocel®

Description	SKU #
Merocel® Eye Wicks, 8 in, 2/Pouch, 20/Bx	400111

Corneal Light Shields

Merocel® Corneal Light Shield

Used to protect hinge and posterior side of flap. Used to protect the retina from operating microscope light. Made from exclusive Merocel PVA material. The Merocel® material is virtually lint free.



Qty: 20/Bx
Brand: Merocel®

Description	SKU #
Merocel® Corneal Light Shield, 8 mm (0.3 in) diameter	400106

Visitec® Corneal Light Shield

PVA absorbent pad. Soft, porous structure designed to protect the retina from operating microscope light and contains no cellulose. Upon moistening, shield becomes soft and conforms to shape of cornea. Designed to keep cornea moist.




Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Visitec® Corneal Light Shield, 8 mm (0.3 in) diameter	581062

Keracel® Corneal Light Shield

Round shield cut flat at one end to fit against the hinge of a LASIK corneal flap. Made from a patented high-density formulation of virtually lint free Merocel Polyvinyl Acetal (PVA) sponge. 10 mm diameter with straight edge.



Qty: 20/Bx
Brand: Keracel®

Description	SKU #
Keracel® Corneal Light Shield with straight edge	400121


Drapes: Full Body

Visidrape® full body drapes feature wider widths and a three-sided channel bag for maximum fluid capacity. Full body drapes can be utilized in all ophthalmic procedures.

Full Body Incise Drape

Qty: 10/Bx
Brand: Visidrape®


Description	SKU #
Visidrape® Non-Woven Incise Drape with Fluid Collection Bag	581087
Overall size:	140 cm x 244 cm (55 in x 96 in)
Incise area:	9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
Fluid collection bag:	25.4 cm x 15.2 cm (10 in x 5.98 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid



Full Body Aperture Drape

Qty: 10/Bx
Brand: Visidrape®


Description	SKU #
Visidrape® Full Body Aperture Drape	581088
Overall size:	140 cm x 244 cm (55 in x 96 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid



Full Body Incise Drape with Plastic Head and Channel Bag

Qty: 10/Bx
Brand: Visidrape®

Description	SKU #
Visidrape® Full Body Incise Drape with Plastic Head and Channel Bag	581715
Overall size:	165 cm x 254 cm (65 in x 100 in)
Incise Area:	9 cm x 11 cm (3.5 in x 4 in)
Channel Bag:	62 cm x 40 cm (24 in x 15.75 in)
Capacity:	1000 ml
Material Type:	SMS (63 in.) and Opaque (37 in.)




Drapes: Full Body

Non-Woven Full Body Aperture Drape

Qty: 10/Bx
Brand: Visidrape®

Description	SKU #
Visidrape® Non-Woven Full Body Aperture Drape	581722
Overall size:	165 cm x 254 cm (65 in x 100 in)
Aperture area:	6 cm x 4 cm (2.25 in x 1.5 in)
Material Type:	Non-Woven Wetlaid




Drapes: Half Body

Embossed Plastic Incise Drape with Fluid Collection Bag

Qty: 10/Bx
Brand: Visidrape®


Description	SKU #
Visidrape® Embossed Plastic Incise Drape with Fluid Collection Bag	581071
Overall size:	101.5 cm x 127 cm (40 in x 50 in)
Incise area:	9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6 1/2 in)
Capacity:	250 ml
Material Type:	Polyethylene Opaque Film



Embossed Plastic Aperture Drape with Fluid Collection Bag

Qty: 10/Bx
Brand: Visidrape®

Description	SKU #
Visidrape® Embossed Plastic Aperture Drape with Fluid Collection Bag	581078
Overall size:	101.5 cm x 127 cm (40 in x 50 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6 1/2 in)
Capacity:	250 ml
Material Type:	Polyethylene Opaque Film




Drapes: Half Body

Non-Woven Incise Drape with Fluid Collection Bag

Qty: 10/Bx
Brand: Visidrape®


Description	SKU #
Visidrape® Non-Woven Incise Drape with Fluid Collection Bag	581072
Overall size:	101.5 cm x 127 cm (40 in x 50 in)
Incise area:	9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6 1/2 in)
Capacity:	250 ml
Material Type:	Non-Woven Wetlaid



Non-Woven Aperture Drape with Fluid Collection Bag

Qty: 10/Bx
Brand: Visidrape®


Description	SKU #
Visidrape® Non-Woven Aperture Drape with Fluid Collection Bag	581079
Overall size:	101.5 cm x 127 cm (40 in x 50 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6 1/2 in)
Capacity:	250 ml
Material Type:	Non-Woven Wetlaid



Non-Woven Incise Drape with Channel Bag

Qty: 10/Bx
Brand: Visidrape®

Description	SKU #
Visidrape® Non-Woven Incise Drape with Channel Bag	581084
Overall size:	140 cm x 127 cm (55 in x 50 in)
Incise area:	9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
Fluid collection bag with channels:	14.6 cm x 25.4 cm (5 3/4 in x 10 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid




Drapes: Half Body

Non-Woven Aperture Drape with Channel Bag

Qty: 10/Bx
Brand: Visidrape®


Description	SKU #
Visidrape® Non-Woven Aperture Drape with Channel Bag	581086
Overall size:	140 cm x 127 cm (55 in x 50 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Fluid collection bag with channels:	14.6 cm x 25.4 cm (5 3/4 in x 10 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid



Non-Woven Half Body Incise Drape* with Channel Bag

Qty: 10/Bx
Brand: Visidrape®

Description	SKU #
Visidrape® Non-Woven Half Body Incise Drape* with Channel Bag	581723
Overall size:	165 cm x 183 cm (65 in x 72 in)
Incise area:	9 cm x 11 cm (3.5 in x 4 in)
Channel bag:	62 cm x 40 cm (24 in x 15.75 in)
Capacity:	1000 ml
Material Type:	SMS



*Without wick

Drapes: Mini

Ideal for refractive or office procedures.

Mini Plastic Incise Drape

Adhesive backing.

Qty: 10/Bx
Brand: Visidrape®

Description	SKU #
Visidrape® Mini Plastic Incise Drape	581095
Overall size:	41 cm x 41 cm (16 in x 16 in)
Incise area:	5 cm x 7.5 cm (2 in x 3 in)
Material Type:	Polyethylene Opaque Film



Drapes: Mini

Mini Plastic Aperture Drape

Adhesive backing.

Qty: 10/Bx

Brand: Visidrape®

Description		SKU #
Visidrape® Mini Plastic Aperture Drape		581096
Overall size:	41 cm x 41 cm (16 in x 16 in)	
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)	
Material Type:	Polythethylene Opaque Film	




Mini Plastic Aperture Drape

Small incise provides complete coverage of lashes, lids and brows. Converts aperture to incise drape.

Qty: 10/Bx

Brand: Visidrape®

Description		SKU #
Visidrape® Incise Drape		581700
Overall size:	10.2 x 20.4 cm (4 x 8 in)	
Adhesive area:	11 cm x 16 cm (4 in x 6 1/4 in)	
Material Type:	Polythethylene Opaque Film	



Dual Refractive Drape

Qty: 10/Bx

Brand: Visidrape®

Description		SKU #
Visidrape® Dual Refractive Drape		581710
Overall size:	40 cm x 26 cm (16 in x 10 1/4 in)	
Aperture area:	2.5 cm x 3.75 cm (1 in x 1 1/2 in)	
Material Type:	Polythethylene Opaque Film	



Eye Drains

Visitec® LASIK Eye Drain (Chayet)

Designed to reduce stromal bed hydration for a consistent ablation and for nasally hinged flaps.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
LASIK Eye Drain 10 mm	581098

Fluid Collection

Merocel® Eye Drain

Provides continuous control of fluid throughout the procedure. The wicks are made from highly absorbent, fast wicking, virtually lint free Merocel® PVA material.



400108

Qty: Various/Bx
Brand: Merocel®

Description	SKU #
Merocel® Eye Drain with 80 cc Bag Bag: 21.6 cm x 5.1 cm (8.5 in x 2.0 in) Wick: Includes one, 15.2 cm (6 in)	400108

Visitec® Fluid Collection Bags

For procedures where collection and control of fluid is designed to keep surgeon, patient and floor dry. Reduces operating room cleanup and turnaround time. Bags have an adhesive backing to adhere to surgical drape.



581005



581070

Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Fluid Collection Bag: 17 cm x 16 cm (6 3/4 x 6 1/4 in), 250 ml capacity	581005
Mini Fluid Collection with Wick Bag: 9 x 16 cm (3 1/2 x 6 1/4 in), 100 ml capacity Wick: Includes one, 0.40 x 15 cm (5/32 x 5 7/8 in)	581070

Instrument Wipes

Merocel® Instrument Wipes

Merocel® wipes are designed to provide the surgeon a wipe for safe cleaning of fine microsurgical instruments. The exclusive Merocel® material is highly absorbent, fast wicking, virtually lint free, polyvinyl acetal (PVA).



Qty: Various/Bx
Brand: Merocel®

Description	SKU #
Instrument Wipe 92.1 x 76.2 mm (3.625 x 3 in), 20/Bx	400200
Wipe with Wick 92.1 x 76.2 mm (3.625 x 3 in), 20/Bx	400210

Visitec® Instrument Wipes



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Instrument Wipe 76.2 x 76.2 mm (3 x 3 in), 2 mm thick	90900

Visitec® Visiwipe Instrument Wipes



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Instrument Wipe 7.3 x 7.3 cm (2 7/8 x 2 7/8 in)	581047
Adhesive Wipe 7.6 x 4.5 cm (3 x 1 3/4 in)	581076

Other Fluid Control

Absorbent Stick

6.6 cm overall length x 5 mm diameter (2.59 in x 3/16 in)



Qty: 250/Bx
Brand: Visitec®

Description	SKU #
Absorbent Stick 6.6 cm, 10 per pouch, 25 pouches per box.	581060

Visiswab Cotton Swabs

7.6 cm (3 in) overall length. Rounded end is made for scleral depression and pointed end may be used as a second instrument to move sutures, manipulate tissue or absorb blood and fluid. Tightly wound ends create less linting than other cotton swabs. 10 per pouch, 25 pouches per box.



Qty: 250/Bx
Brand: Visitec®

Description	SKU #
Visiswab Cotton 7.6 cm (3 in) length	581069

Fluid Control Sets

Merocel® Cataract Set



Qty: 10/Bx
Brand: Merocel®

Description	SKU #
Merocel® Cataract Kit	400113
Contents	SKU #
Merocel® eye spears, Qty/5	400101
Instrument wipe, Qty/1	400200
Corneal light shield, Qty/1	400106
80 cc Collection bag with wick, Qty/1	400108

