




# 90° Directional Laser Probes



## A New Vision on Flexibility

- 90° curved extendable tip for easy access to peripheral retina
- Hardened metal shaft reduces flexibility by 25%<sup>1</sup>
- Retractable mechanism improves control of instrument tip and reduces risk of accidental removal of cannula
- Available in standard (23G, 25G and 27G) and illuminated (23G and 25G) color coded for easy gauge identification
- Compatible with all major equipment

**VitreQ**  
Vitreoretinal Surgical Products

**Indications for Use:** VitreQ Disposable Laser Probes, 90° Directional Laser Probes and 90° Directional Illuminated Laser Probes are intended for use in vitreoretinal surgery to perform endo-ocular laser photocoagulation treatments at operating wavelengths of 500nm to 900nm and illumination at light wavelengths of 425nm to 700nm.

**Precautions:** Do not use product if package integrity has been breached or compromised; do not use product after the expiration date.

**Caution:** Federal law restricts this device to sale by or own order of a licensed physician.

**Attention:** See Instructions for Use included in product packaging for complete listing of indications and warnings.

**[bvimedical.com](http://bvimedical.com)**

# Vitrex® Vitreoretinal Surgical Products

## 90° Directional Laser Probes Designed For Improved Surgical Control


- Ergonomic handle design for improved intraocular surgical control
- Sliding knob to control the motion of instrument tip

### Gauge Size with Color Coding

23G


25G

27G



Standard Directional*		
G	Description	SKU #
●	23G Compatible with Alcon® and DORC®	LP23.A04
●	23G Compatible with Iridex® and B&L	LP23.I04
●	25G Compatible with Alcon® and DORC®	LP25.A04
●	25G Compatible with Iridex® and B&L	LP25.I04
●	27G Compatible with Alcon® and DORC®	LP27.A04
●	27G Compatible with Iridex® and B&L	LP27.I04

Illuminated Directional*		
G	Description	SKU #
●	23G Compatible with Alcon® and DORC®	LP23.A05
●	23G Compatible with Iridex® and B&L	LP23.I05
●	25G Compatible with Alcon® and DORC®	LP25.A05
●	25G Compatible with Iridex® and B&L	LP25.I05

Light Source Adaptors*	
	
Description	SKU #
Constellation® Compatible Light Source Adaptor	LF00.R01
Stellaris® Compatible Light Source Adaptor	LF00.R02
DORC® Compatible Light Source Adaptor	LF00.R05

\* All product names, trademarks and registered trademarks are property of their respective owners.

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

**Vitrex**

Vitreoretinal Surgical Products

Tel: 866-906-8080

Email: customersupport@bvimedical.com



**BVI**  
bvimedical.com

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2021

1516224-03