



# Malosa® Single-Use Instruments

High Quality | Efficient | Cost-Effective

United States

**Malosa®**  
Single-Use Instruments

[bvimedical.com](http://bvimedical.com)

# Malosa<sup>®</sup> Single-Use Instruments

---

- **Eliminate Reprocessing and Sterilisation Costs**

No added burden on facilities to re-sterilise or repair damaged instruments, helping to increase operational efficiencies

- **Improve Patient Safety and Decrease Contamination**

By utilising a brand new instrument for each surgery, a patient's risk for infection or cross contamination is decreased



Scheib, MD, FEBO, C. Single-use instruments deserve consideration for eye surgery. Ocular Surgery News Europe Edition. December 2016

Southworth, P.M. Infections and exposures: reported incidents associated with unsuccessful decontamination of reusable surgical instruments. The Journal of Hospital Infection. September 17, 2014

Tarantola MD, R.; Shuler MD, PhD, M.; Kaushal, MD PhD, S.; Nagpal, MS, DO, FRCS(UK), M.; Disposable Instruments are Cost Effective and as Good as Reusable. Retina Today. Jan/Feb 2015



# Table of Contents

---

## **Forceps**

|                                    |    |
|------------------------------------|----|
| Capsulorhexis Forceps              | 4  |
| Fixation Forceps                   | 5  |
| Tying Forceps                      | 6  |
| IOL Forceps                        | 6  |
| General Forceps and Fixation Rings | 7  |
| Dressing Forceps                   | 8  |
| Needle-Holders                     | 8  |
| Choppers                           | 9  |
| Double-Ended Choppers              | 9  |
| Manipulators                       | 10 |

## **Scissors**

|                     |    |
|---------------------|----|
| Sprung Scissors     | 11 |
| Dissection Scissors | 12 |
| Dressing Scissors   | 12 |

## **Speculums**

|                                   |    |
|-----------------------------------|----|
| Adjustable Speculums              | 13 |
| Sprung Wire Speculums             | 13 |
| Paediatric and Neonatal Speculums | 13 |
| Depressors                        | 14 |
| Lacrimonal                        | 14 |
| Markers and Gauges                | 15 |
| Clamps                            | 16 |
| Curettes                          | 16 |
| Flap-Lifters and Corneal Markers  | 17 |
| SMILE®                            | 18 |
| Spatulas and Wells                | 18 |



# Capsulorhexis Forceps

The Malosa® range of single-use Capsulorhexis Forceps offers tip and blade designs for performing continuous curvilinear capsulotomy and supporting sub - 1.8 mm microincision (MICS) procedures.



MICS Capsulorhexis - Curved/Sharp



MICS Capsulorhexis - Straight/Sharp




Next Generation Capsulorhexis Forceps Grasping



| Capsulorhexis Forceps                          | Description                                                                                               | SKU #     |
|------------------------------------------------|-----------------------------------------------------------------------------------------------------------|-----------|
| MICS Capsulorhexis - Curved/Sharp              | 0.4 mm x 0.3 mm Sharp tips, 12 mm Curved direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle   | MMSU1567S |
| MICS Capsulorhexis - Straight/Sharp            | 0.4 mm x 0.3 mm Sharp tips, 12 mm Straight direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle | MMSU1568S |
| Next Generation Capsulorhexis Forceps Sharp    | Sharp curved tips, Round handle                                                                           | 581488    |
| Next Generation Capsulorhexis Forceps Grasping | Grasping curved tips, Round handle                                                                        | 581486    |



# Fixation Forceps



Toothed Forceps (Bonn)




Notched Forceps (Hoskins)



| Toothed Fixation Forceps          | Description                                             | SKU #     |
|-----------------------------------|---------------------------------------------------------|-----------|
| Fine Toothed Forceps Short (Bonn) | 0.3 mm 1 x 2 Teeth, 5 mm Tying platforms, 52 mm Handle  | MMSU1123S |
| Toothed Forceps (Bonn)            | 0.45 mm 1 x 2 Teeth, 5 mm Tying platforms, 80 mm Handle | MMSU1124S |
| Colibri Toothed Forceps Short     | 0.4 mm 1 x 2 Teeth, 5 mm Tying platforms, 52 mm Handle  | MMSU1129S |

| Notched Fixation Forceps             | Description                                             | SKU #     |
|--------------------------------------|---------------------------------------------------------|-----------|
| Fine Notched Forceps Short (Hoskins) | 0.2 mm Notched tips, 5 mm Tying platforms, 52 mm Handle | MMSU1114S |
| Notched Forceps (Hoskins)            | 0.3 mm Notched tips, 5 mm Tying platforms, 80 mm Handle | MMSU1115S |
| Colibri Notched Forceps Short        | 0.2 mm Notched tips, 5 mm Tying platforms, 52 mm Handle | MMSU1209S |



Colibri Toothed Forceps Short



## Tying Forceps



| Tying Forceps                    | Description                                                                                    | SKU #     |
|----------------------------------|------------------------------------------------------------------------------------------------|-----------|
| Straight Tying Forceps           | 0.35 mm Tips, 5 mm Tying platforms, 80 mm Handle                                               | MMSU1125S |
| Curved Tying Forceps             | 0.25 mm Tips, 5 mm Curved tying platforms, 80 mm Handle                                        | MMSU1210S |
| Angled Tying Forceps (McPherson) | 0.35 mm Tip, 5 mm Tying platforms, 8 mm Shaft angled 45 degree, 80 mm Handle                   | MMSU1127S |
| Waqar Suture Removal Forceps     | 0.3 mm x 0.65 mm Pick angled 90 degree, 5 mm Tying platforms, Curved 28 mm shaft, 80 mm Handle | MMSU1583S |

## IOL Forceps



| IOL Forceps          | Description                                              | SKU #     |
|----------------------|----------------------------------------------------------|-----------|
| Lens Loading Forceps | 0.6 mm x 8 mm Tips, 45 Degree angled shaft, 80 mm Handle | MMSU1146S |

# General Forceps and Fixation Ring



| General Forceps                 | Description                                             | SKU #     |
|---------------------------------|---------------------------------------------------------|-----------|
| Moorfields Forceps              | 0.6 mm tips, 11.5 mm Serrated jaws, 105 mm Long         | MMSU0277S |
| St Martins Forceps              | 0.9 mm x 0.26 mm 1 x 2 Teeth, 52 mm Handle              | MMSU0279S |
| Epilation Forceps               | 5 mm x 7 mm Flat plane tips, 70 mm Handle               | MMSU0085S |
| Jewellers Forceps Fine Straight | 7 mm Pointed tips, 80 mm Handle                         | MMSU0286S |
| Jewellers Forceps Fine Angled   | 7 mm, 45 Degree angled pointed tips, 80 mm Handle       | MMSU0283S |
| Contact Lens Forceps            | 0.6 mm x 0.4 mm, 9 mm Curved platforms, 70 mm Handles   | MMSU1284S |
| Adson Forceps, Toothed          | 1.35 mm x 1.8 mm 1 x 2 Teeth, 122 mm Long               | MMSU0238S |
| Adson Forceps, Non-Toothed      | 0.85 mm x 1 mm Tips, 14.5 mm Serrated jaws, 122 mm Long | MMSU0241S |

## Fixation Ring



| Fixation Ring | Description                                                                                                                                   | SKU #     |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| Fixation Ring | 0.5 mm Thick 13 mm internal diameter ring, 1.5 mm x 0.6 mm Teeth at 3.5 mm intervals, 1.3 mm x 0.4 mm Teeth at 0.9 mm intervals, 75 mm Handle | MMSU1495S |

## Dressing Forceps



| Dressing Forceps | Description                                             | SKU #     |
|------------------|---------------------------------------------------------|-----------|
| Straight Rampley | 11 mm x 23 mm Serrated tips, 180 mm Long, Nylon plastic | MMSU0068S |
| Rampley Short    | 17 mm x 24 mm Serrated tips, 125 mm Long, Nylon plastic | MMSU1578S |

## Needle-Holders



Corneal Needle-Holder



Barraquer Needle-Holder



Troutman Needle-Holder

| Needle-Holders          | Description                                                                | SKU #     |
|-------------------------|----------------------------------------------------------------------------|-----------|
| Corneal Needle-Holder   | 9 mm Fine curved jaws, 0.45 mm Tips, 106 mm Long                           | MMSU1110S |
| Barraquer Needle-Holder | 0.42 mm x 0.45 mm Tips, 4 mm Fine curved jaws with 5 mm shaft, 138 mm Long | MMSU3015S |
| Troutman Needle-Holder  | 4 mm Fine curved jaws, 0.55 mm Tips, 76 mm Long handle                     | MMSU1243S |



Troutman Needle-Holder



# Choppers

The Malosa® range of single-use choppers offers a range and combination of second instruments to suit all techniques of phacoemulsification.



**Claw Chopper**  
MMSU1105S



**Phaco Chopper**  
MMSU1108S




**Greens Chopper**  
Right Handed  
MMSU1119S

# Double-Ended Choppers



**Phaco Chopper |**  
Mushroom Manipulator  
MMSU1221S



**Horizontal Chopper |**  
Kuglen Hook  
MMSU1231S

# Manipulators

Our single-use manipulators offer a range and combination of second instruments to suit all surgical techniques.



**Govan Manipulator**  
MMSU1102S



**Drysdale Manipulator**  
MMSU1144S



**Sinskey Hook 0.2 mm**  
MMSU1103S



**Mushroom Manipulator**  
MMSU1179S



**Lebuisson Manipulator**  
MMSU1458S




**Lester Hook**  
MMSU1207S



**ICL Manipulator**  
MMSU1153S



**Micro Manipulator**  
581441



**Sinskey Hook**  
581442



# Scissors

Sprung and direct-action ophthalmic scissors for all surgical specialties.



Westcott Scissors



Vannas Scissors



## Sprung Scissors

| Sprung Scissors                      | Description                                                                   | SKU #     |
|--------------------------------------|-------------------------------------------------------------------------------|-----------|
| Westcott Scissors Curved Pointed     | 0.25 mm x 0.3 mm Pointed tips, 12 mm Curved blades, 116 mm Long               | MMSU1211S |
| Westcott Scissors Curved Dophin Nose | 0.5 mm x 0.6 mm Rounded tips, 22 mm Curved 'Dolphin Nose' blades, 120 mm Long | MMSU1213S |
| Westcott Scissors Straight Pointed   | Pointed tips, 17 mm Straight blades, 120 mm Long                              | MMSU1214S |
| Vannas Scissors Curved               | Pointed tips on 8 mm curved blades, 1.5 mm when closed, 34 mm Handle          | MMSU1140S |
| Vannas Scissors Straight             | Pointed tips on 8 mm straight blades, 1.5 mm when closed, 34 mm Handle        | MMSU1141S |

# Scissors (Continued)

## Dissection Scissors



Iris Scissors  
Curved Pointed



Stevens Tenotomy  
Straight Rounded

| Dissection Scissors               | Description                                                        | SKU #     |
|-----------------------------------|--------------------------------------------------------------------|-----------|
| Stevens Tenotomy Straight Rounded | 1.5 mm x 1.5 mm Rounded tips, Straight blades, 110 mm Long         | MMSU0214S |
| Iris Scissors Curved Pointed      | 0.1 mm x 0.7 mm Pointed tips, 25 mm Curved blades, 115 mm Long     | MMSU0244S |
| Iris Scissors Straight Pointed    | 0.1 mm x 0.5 mm Pointed tips, 32.5 mm Straight blades, 120 mm Long | MMSU0249S |

## Dressing Scissors




Dressing Scissors

| Dressing Scissors                  | Description                                                               | SKU #     |
|------------------------------------|---------------------------------------------------------------------------|-----------|
| Dressing Scissors Straight Pointed | 0.3 mm Pointed tips, 24.5 mm Straight blades, 102 mm Long, Plastic handle | MMSU0454S |
| Dressing Scissors Straight Rounded | 0.8 mm Rounded tips, 24.5 mm Straight blades, 102 mm Long, Plastic handle | MMSU0455S |




# Speculums

## Adjustable Speculums



Mechanical Speculum Solid Blade



Mechanical Speculum Open Blade

| Adjustable Speculums            | Description                                                                                                                 | SKU #     |
|---------------------------------|-----------------------------------------------------------------------------------------------------------------------------|-----------|
| Mechanical Speculum Solid Blade | 6 mm x 14.5 mm Solid curved polycarbonate blades, 0.7 mm Thick, 1 mm Wire, Polycarbonate body and screw mechanical speculum | MMSU1488S |
| Mechanical Speculum Open Blade  | 5 mm x 14 mm Open blades, 1 mm Wire, Polycarbonate body and screw mechanical speculum                                       | MMSU1290S |

## Sprung Wire Speculums



Barraquer Speculum Flat Blade



Kratz-Barraquer Speculum Curved



Barraquer Speculum Solid Curved Blade



Barraquer Speculum Solid 18 mm Blade

| Sprung Wire Speculums                   | Description                                                                               | SKU #     |
|-----------------------------------------|-------------------------------------------------------------------------------------------|-----------|
| Barraquer Speculum Flat Blade           | 6.5 mm x 14.5 mm Straight closed blades, 1 mm Wire, 30 mm Wide                            | MMSU1275S |
| Barraquer Speculum Solid Curved Blade   | 6 mm x 14.5 mm Solid curved polycarbonate blades, 0.7 mm Thick, 0.9 mm Wire, 30 mm Wide   | MMSU1474S |
| Barraquer Speculum Solid Straight Blade | 7 mm x 14.5 mm Solid straight polycarbonate blades, 0.8 mm Thick, 0.9 mm Wire, 30 mm Wide | MMSU1475S |
| Barraquer Speculum Solid 18 mm Blade    | 6 mm x 14.5 mm Solid curved polycarbonate blades, 0.8 mm Thick, 1 mm Wire                 | MMSU1479S |
| Kratz-Barraquer Speculum                | 6.5 mm x 14 mm Open blades, 0.9 mm Wire, 30 mm Wide                                       | MMSU1252S |
| Kratz-Barraquer Speculum Temporal       | 6.5 mm x 14 mm Open blades, 0.9 mm Wire, 30 mm Wide                                       | MMSU1253S |
| Kratz-Barraquer Speculum Curved         | 6.5 mm x 14.5 mm Open curved blades, 0.9 mm Wire, 30 mm Wide                              | MMSU1287S |

## Paediatric and Neonatal Speculums

Neonatal Speculum Small



Paediatric Speculum




| Paediatric and Neonatal Speculums | Description                                               | SKU #     |
|-----------------------------------|-----------------------------------------------------------|-----------|
| Neonatal Speculum Small           | 3 mm x 7 mm Closed blades, 0.8 mm Wire, 25 mm Wide        | MMSU1262S |
| Paediatric Speculum               | 5.5 mm x 13.5 mm Closed blades, 0.9 mm Wire, 31.5 mm Wide | MMSU1278S |

# Depressors




| Depressors                              | Description                                                             | SKU #     |
|-----------------------------------------|-------------------------------------------------------------------------|-----------|
| Schocket Double-Ended Scleral Depressor | 2.7 mm x 6.5 mm T shaped tip, 5.5 mm x 2.2 mm Rounded tip, 68 mm Handle | MMSU1405S |
| Ball Tip Scleral Depressor              | 3 mm Ball tip, 21.5 mm Straight shaft, 50 mm Handle                     | MMSU1408S |


# Lacrimal



Punctal Dilator



Lacrimal Probe



Golf Club Spud

| Lacrimal                       | Description                                                                 | SKU #     |
|--------------------------------|-----------------------------------------------------------------------------|-----------|
| Punctum Dilator/Seeker         | 0.45 mm x 21.5 mm Tip, 5 mm Pointed tip, 18 mm Straight shaft, 77 mm Handle | MMSU1193S |
| Punctum Dilator Fine/Medium    | 0.4 mm x 21.5 mm Tip, 0.6 mm x 18 mm Tip, 77 mm Handle                      | MMSU1181S |
| Punctum Dilator/Lacrimal Probe | 0.28 mm x 25 mm Dilator, 0.7 mm x 44.5 mm Probe, 70 mm Handle               | MMSU1419S |
| Golf Club Spud                 | 0.1 mm x 0.8 mm x 2 mm Oval profile, 18 mm Straight shaft, 75 mm Handle     | MMSU1417S |


# Markers and Gauges



Braunstein IVT Marker



Pointed IVT Marker



Double Ended  
IVT Guide Pad



Intravitreal Injection Guide



| Markers and Gauges                     | Description                                                                                                                                             | SKU #     |
|----------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| Intravitreal Injection Guide           | Polycarbonate lash guard with stainless-steel swage guide tube, Triangular base plate with 3 stabilizing studs, 39 mm, 45 Degree angled wishbone handle | MMSU1577S |
| Braunstein IVT Marker 3.5   4 mm Green | 2 mm x 1.5 mm Tips with 3.5 mm spread, 2 mm x 1.5 mm Tips with 4 mm spread, Polycarbonate, 108 mm Long                                                  | MMSU1276S |
| Pointed IVT Marker 3.5   4 mm Blue     | 2 mm x 0.55 mm Tips with 3.5 mm spread, 2 mm x 0.65 mm Tips with 4 mm spread, Polycarbonate, 108 mm Long                                                | MMSU1277S |
| Double Ended IVT Guide Pad 3.5   4 mm  | 7 mm Pad with 3.5 guide, 0.8 mm Pad with 4 mm guide, 6 mm 45 Degree angled shafts, 90 mm Long                                                           | MMSU1285S |
| Braunstein IVT Marker 3.5   4 mm       | 2 mm x 1.5 mm Tips with 3.5 mm spread, 2 mm x 1.5 mm Tips with 4 mm spread, 90 mm Long                                                                  | MMSU1229S |
| Mendez Gauge                           | 2.5 mm Polycarbonate ring with 13 mm ID, 0 through 180 Degree marks in increments of 10, 75 mm Handle                                                   | MMSU1152S |



Intravitreal Injection Guide

# Clamps and Curettes



## Clamps

| Clamps                | Description                                                    | SKU #     |
|-----------------------|----------------------------------------------------------------|-----------|
| Chalazion Clamp Round | 12 mm Round profile, 9 mm Round aperture, 90 mm Long           | MMSU1150S |
| Chalazion Clamp Oval  | 13 mm x 20 mm Oval profile, 10 mm x 17 mm Aperture, 90 mm Long | MMSU1149S |

## Curettes



| Curettes                             | Description                                                          | SKU #     |
|--------------------------------------|----------------------------------------------------------------------|-----------|
| Chalazion Curette Double 1 mm   2 mm | 1 mm Curette, 1 mm Depth, 2 mm Curette, 1.35 mm Depth, 128 mm Long   | MMSU1157S |
| Chalazion Curette Double 2 mm   3 mm | 2 mm Curette, 1.35 mm Depth, 3 mm Curette, 1.6 mm Depth, 128 mm Long | MMSU1158S |



# Flap-Lifters and Corneal Markers



LASIK Flap Elevator



LASIK Flap Separator



Y Shape  
Corneal Marker

| Flap-Lifters and Corneal Markers          | Description                                                                                                            | SKU #     |
|-------------------------------------------|------------------------------------------------------------------------------------------------------------------------|-----------|
| LASIK Flap Elevator Curved 8 mm           | 0.2 mm x 0.15 mm Tip, Curved 8 mm tapered shaft at 90 degrees, 75 mm Handle                                            | MMSU1177S |
| LASIK Flap Separator and Elevator, Binder | 0.1 mm x 0.3 mm Rounded tip 0.7 mm long, 10 mm Curved 60 degree angled shaft, 75 mm Handle                             | MMSU1206S |
| Y Shape Corneal Marker                    | 2 mm x 8 mm Marking blades, 75 mm Handle                                                                               | MMSU1134S |
| LASIK Flap Separator/Elevator (Seibel)    | 0.16 mm x 0.28 mm x 1 mm Tip, 20 mm Straight shaft, 0.5 mm x 10.5 mm pointed hook, 52 mm Handle                        | MMSU1172S |
| LASIK Flap Separator/Elevator Straight    | 0.15 mm x 0.3 mm x 8 mm Straight separator at 90 degree to shaft, 0.15 mm x 0.3 mm x 1 mm Pocketing hook, 52 mm Handle | MMSU1295S |
| LASIK Flap-Lifting Forceps                | 0.3 mm x 11.5 mm Tips, 45 Degree angled shaft, 80 mm Handle                                                            | MMSU1169S |



LASIK Flap-Lifting Forceps

## SMILE®



| SMILE® Cornea                 | Description                                                                                                     | SKU #     |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------|-----------|
| Reinstein Lenticule Separator | 0.22 mm x 2 mm Pocketing hook, 1.2 mm x 2 mm Spatula tip, Curved separator at 90 degrees to shaft, 52 mm Handle | MMSU1297S |
| IOL Removal Forceps           | 0.1 mm x 0.1 mm Ridges along 4 mm of 15 mm shaft, Angled at 30 degree, 52 mm Handle                             | MMSU1460S |

\*SMILE® is a trademark of Zeiss

## Spatulas and Wells



MMSU1187

MMSU1113

MMSU1196

| Spatulas and Wells               | Description                                                                         | SKU #     |
|----------------------------------|-------------------------------------------------------------------------------------|-----------|
| Corneal Epithelium Spatula Flat  | 5 mm Spatula with 0.9 mm flat blade, 0.8 mm x 1.7 mm Round bladed tip, 75 mm Handle | MMSU1113S |
| Corneal Epithelium Spatula Round | 2.65 mm x 5.5 mm x 0.35 mm Rounded blade, 75 mm Handle                              | MMSU1187S |
| Epithelium Well 9 mm Trephining  | 3 mm x 9 mm ID well, Angled 45 degrees, 75 mm Handle                                | MMSU1196S |
| Epithelium Well                  | 3 mm x 8.5 mm ID well, Angled 45 degrees, 75 mm Handle                              | MMSU1197S |

# Index

| SKU       | DESCRIPTION                            | PG | SKU       | DESCRIPTION                                    | PG |
|-----------|----------------------------------------|----|-----------|------------------------------------------------|----|
| MMSU0068S | Straight Rampley                       | 8  | MMSU1196S | Epithelium Well 9 mm Trephining                | 18 |
| MMSU0085S | Epilation Forceps                      | 7  | MMSU1197S | Epithelium Well 8.5 mm                         | 18 |
| MMSU0214S | Stevens Tenotomy Straight Rounded      | 12 | MMSU1206S | LASIK Flap Separator and Elevator, Binder      | 17 |
| MMSU0238S | Adson Forceps, Toothed                 | 7  | MMSU1207S | Lester Hook                                    | 10 |
| MMSU0241S | Adson Forceps, Non-Toothed             | 7  | MMSU1209S | Colibri Notched Forceps Short                  | 5  |
| MMSU0244S | Iris Scissors Curved Pointed           | 12 | MMSU1210S | Curved Tying Forceps                           | 6  |
| MMSU0249S | Iris Scissors Straight Pointed         | 12 | MMSU1211S | Westcott Scissors Curved Pointed               | 11 |
| MMSU0277S | Moorfields Forceps                     | 7  | MMSU1213S | Westcott Scissors Curved Dolphin Nose          | 11 |
| MMSU0279S | St Martins Forceps                     | 7  | MMSU1214S | Westcott Scissors Straight Pointed             | 11 |
| MMSU0283S | Jewellers Forceps Fine Angled          | 7  | MMSU1221S | Phaco Chopper   Mushroom Manipulator           | 9  |
| MMSU0284S | Jewellers Forceps Curved/Plain         | 7  | MMSU1229S | Braunstein IVT Marker 3.5   4 mm               | 15 |
| MMSU0286S | Jewellers Forceps Fine Straight        | 7  | MMSU1231S | Horizontal Chopper/Kuglen Hook                 | 9  |
| MMSU0454S | Dressing Scissors Straight Pointed     | 12 | MMSU1243S | Troutman Needle-Holder                         | 8  |
| MMSU0455S | Dressing Scissors Straight Rounded     | 12 | MMSU1252S | Kratz-Barraquer Speculum                       | 13 |
| MMSU1102S | Govan Manipulator                      | 10 | MMSU1253S | Kratz-Barraquer Speculum Temporal              | 13 |
| MMSU1103S | Sinskey Hook 0.2 mm                    | 10 | MMSU1262S | Neonatal Speculum Small                        | 13 |
| MMSU1105S | Claw Chopper                           | 9  | MMSU1275S | Barraquer Speculum Flat Blade                  | 13 |
| MMSU1108S | Phaco Chopper                          | 9  | MMSU1276S | Braunstein IVT Marker 3.5   4 mm Green         | 15 |
| MMSU1110S | Corneal Needle-Holder                  | 8  | MMSU1277S | Pointed IVT Marker 3.5   4 mm Blue             | 15 |
| MMSU1113S | Corneal Epithelium Spatula Flat        | 18 | MMSU1278S | Paediatric Speculum                            | 13 |
| MMSU1114S | Fine Notched Forceps Short (Hoskins)   | 5  | MMSU1284S | Contact Lens Forceps                           | 7  |
| MMSU1115S | Notched Forceps (Hoskins)              | 5  | MMSU1285S | Double Ended IVT Guide Pad 3.5   4 mm          | 15 |
| MMSU1119S | Greens Chopper Right Handed            | 9  | MMSU1287S | Kratz-Barraquer Speculum Curved                | 13 |
| MMSU1123S | Fine Toothed Forceps Short (Bonn)      | 5  | MMSU1290S | Mechanical Speculum Open Blade                 | 13 |
| MMSU1124S | Toothed Forceps (Bonn)                 | 5  | MMSU1295S | LASIK Flap Separator/Elevator Straight         | 17 |
| MMSU1125S | Straight Tying Forceps                 | 6  | MMSU1297S | Reinstein Lenticule Separator                  | 18 |
| MMSU1127S | Angled Tying Forceps (McPherson)       | 6  | MMSU1405S | Schocket Double-Ended Scleral Depressor        | 14 |
| MMSU1129S | Colibri Toothed Forceps Short          | 5  | MMSU1408S | Ball Tip Scleral Depressor                     | 14 |
| MMSU1134S | Y Shape Corneal Marker                 | 17 | MMSU1417S | Golf Club Spud                                 | 14 |
| MMSU1140S | Vannas Scissors Curved                 | 11 | MMSU1419S | Punctum Dilator/Lacrimal Probe                 | 14 |
| MMSU1141S | Vannas Scissors Straight               | 11 | MMSU1458S | Lebuisson Manipulator                          | 10 |
| MMSU1144S | Drysdale Manipulator                   | 10 | MMSU1474S | Barraquer Speculum Solid Curved Blade          | 13 |
| MMSU1146S | Lens Loading Forceps                   | 6  | MMSU1475S | Barraquer Speculum Solid Straight Blade        | 13 |
| MMSU1149S | Chalazion Clamp Oval                   | 16 | MMSU1479S | Barraquer Speculum Solid 18 mm Blade           | 13 |
| MMSU1150S | Chalazion Clamp Round Grooved          | 16 | MMSU1488S | Mechanical Speculum Solid Blade                | 13 |
| MMSU1152S | Mendez Gauge                           | 15 | MMSU1495S | Fixation Ring                                  | 7  |
| MMSU1153S | ICL Manipulator                        | 10 | MMSU1567S | MICS Capsulorhexis - Curved/Sharp              | 4  |
| MMSU1157S | Chalazion Curette Double 1 mm   2 mm   | 16 | MMSU1568S | MICS Capsulorhexis - Straight/Sharp            | 4  |
| MMSU1158S | Chalazion Curette Double 2 mm   3 mm   | 16 | MMSU1577S | Intravitreal Injection Guide                   | 15 |
| MMSU1169S | LASIK Flap-Lifting Forceps             | 17 | MMSU1578S | Rampley Short                                  | 8  |
| MMSU1172S | LASIK Flap Separator/Elevator (Seibel) | 17 | MMSU1583S | Waqar Suture Removal Forceps                   | 6  |
| MMSU1177S | LASIK Flap Elevator Curved 8 mm        | 17 | MMSU3015S | Barraquer Needle-Holder                        | 8  |
| MMSU1179S | Mushroom Manipulator                   | 10 | 581441    | Micro Manipulator                              | 10 |
| MMSU1181S | Punctum Dilator Fine/Medium            | 14 | 581442    | Sinskey Hook                                   | 10 |
| MMSU1187S | Corneal Epithelium Spatula Round       | 18 | 581486    | Next Generation Capsulorhexis Forceps Grasping | 4  |
| MMSU1193S | Punctum Dilator/Seeker                 | 14 | 581488    | Next Generation Capsulorhexis Forceps Sharp    | 4  |





## Malosa®

Single-Use Instruments

### Questions:

T: 866-906-8080

E: [customersupport@bvimedical.com](mailto:customersupport@bvimedical.com)



[bvimedical.com](http://bvimedical.com)