


BVI Wet-Field® Eraser®

Precise, Targeted Coagulation

BVI offers a range of specialty Eraser® Instruments, multiple options for bipolar pencils, for vitreoretinal procedures and glaucoma filtration. Designed for targeted and discrete coagulation with minimal tissue trauma and minimal peripheral cellular necrosis.




Contoured shafts are designed to conform to ocular anatomy while maximizing tip visualization, range of motion, and target approach versatility.



Wet-Field® Eraser®


BVI provides clinicians a wide portfolio of bipolar pencils ranging from 18 to 25 gauge, with either straight and angled shafts. Tip options include beveled, blunt, fine, and tapered blunt or fine to provide support for many ophthalmic surgical procedures and conditions. With different shaft and tip options, clinicians can choose what will be applicable for their patient's eye anatomy and anticipated surgical outcome.

Wet-Field® Hemostatic Eraser® Bipolar 18G




| Description | SKU # |
|--|--------|
| 18G, Bevel tip, Straight, 10/Bx | 221250 |
| 18G, Blunt tip, Angled 45 degrees, 10/Bx | 221251 |

Wet-Field® Hemostatic Eraser® Bipolar 20G



| Description | SKU # |
|---|-----------|
| 20G, Blunt tip, Straight, 10/Bx | 221260 |
| 20G, Blunt tip, Angled 135 degrees, 10/Bx | 221260135 |
| 20G, Blunt tip, Angled 150 degrees, 10/Bx | 221260150 |


Wet-Field® Hemostatic Eraser® Bipolar 20-23G



| Description | SKU # |
|--|--------|
| 20-23G, Tapered, Fine tip, Straight, 5/Bx | 221265 |
| 20-23G, Tapered, Blunt tip, Straight, 5/Bx | 221266 |


Wet-Field® Hemostatic Eraser® Bipolar 25G

Designed for fine tip hemostasis without compromising power delivery.




| Description | SKU # |
|---|-----------|
| 25G, Fine tip, Straight, 5/Bx | 221267 |
| 25G, Fine tip, Angled 135 degrees, 5/Bx | 221267135 |
| 25G, Fine tip, Angled 150 degrees, 5/Bx | 221267150 |

Wet-Field® Hemostatic Eraser® Bipolar 23G



| Description | SKU # |
|-----------------------------------|-----------|
| 23G, Fine tip, Straight | 221268 |
| 23G, Fine tip, Angled 135 degrees | 221268135 |
| 23G, Fine tip, Angled 150 degrees | 221268150 |

Single-Use Bipolar Cable



| Description | SKU # |
|--------------|----------|
| Shrouded pin | 22155010 |

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

Tel: 1.866.906.8080 Fax: 1.866.906.4304



BVI, BVI Logo and all other trademarks are property of BVI © 2019

1517867-01