

BVI

40 Years of Innovation and Quality Cannulas

Supporting All Areas of Ophthalmic Surgery

- A market leader of ophthalmic cannulas for over 40 years
- Full range of single-use cannulas developed to support all ophthalmic surgical procedures
- A history of innovation to meet your diverse surgical needs



Visitec®
Cannulas

bvimedical.com

Irrigation

Anterior Chamber Cannula

Angled 45 degrees, 6 mm from end.




581284

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	581642
27G x 7/8 in, 0.4 x 22 mm	581284
30G x 7/8 in, 0.3 x 22 mm	581285

Viscoflow Cannula

Angled 45 Degrees, 9 mm from end.



585009

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
23G x 7/8 in, 0.6 x 22 mm	585009
25G x 7/8 in, 0.5 x 22 mm	585045
27G x 7/8 in, 0.4 x 22 mm	585026
30G x 7/8 in, 0.3 x 22 mm	585046

Anterior Chamber (Rycroft)


Angled 45 degrees, 4 mm from end.



585025

Qty: 10/Bx
Brand: Visitec®


Description	SKU #
21G x 7/8 in, 0.8 x 22 mm	581281
23G x 7/8 in, 0.6 x 22 mm	585025
25G x 7/8 in, 0.5 x 22 mm	581639
27G x 7/8 in, 0.4 x 22 mm	581280
30G x 7/8 in, 0.3 x 22 mm	581273



Refractive

LASIK Cannula

Intended to irrigate and/or manipulate corneal tissue. 25 Degrees, formed 7 mm from end.




585278

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585278

LASIK Irrigating Cannula (Solomon)

Intended to irrigate and/or manipulate corneal tissue. Thin wall cannula, angled tip features three 0.25 mm ports on 8 mm extension.



585280

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585280

LASIK Irrigating Cannula (Lindstrom)

Intended to irrigate and/or manipulate corneal tissue. Features three 0.25 mm ports located on 7 mm extension with an 11 mm radius, 45 degrees.




585250

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585250

LASIK Flap Spatula/Irrigator (Hersh)


Intended to irrigate and/or manipulate corneal tissue. 45 Degrees, 9 mm from spatulated tip.



585252

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585252



Anesthesia

Retrobulbar (Atkinson)

Intended for injecting anaesthetic agents into the muscle cone.




581638

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 1-1/2 in, 0.6 x 38 mm	581638
● 25G x 1-1/2 in, 0.5 x 38 mm	581637

Retrobulbar/Peribulbar (Atkinson)




585108

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 1-1/4 in, 0.6 x 32 mm	585001
● 25G x 1-1/4 in, 0.5 x 32 mm	585108
● 27G x 1-1/4 in, 0.4 x 32 mm	585109

Retrobulbar (Sharp Point)

Intended for injecting anaesthetic agents into the muscle cone.



581275

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 1-1/2 in, 0.5 x 38 mm	581275

Sub-Tenon's Anesthesia Cannula


Intended for administering anaesthetic agents into the posterior sub-tenon's space.



585176

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 19G x 1 in, 1.1 x 25 mm	585176



ISO Hub Gauge Chart

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

Hydrodissection

Cortical Cleaving Hydrodissector (Flat End)

45 degrees, 11 mm from end. For cleaving cortex from capsular bag.



585158

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G 7/8 in, 0.4 x 22 mm	585158

Micro J Shape Hydrodissector (Pearce)

Hook is 2.25 mm long x 1.75 mm wide.




585029

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585029

Hydrodissection Cannula (Akahoshi Style)

Tapered tip, Thin wall, Angled 1 mm from end.



585198

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G 7/8 in, 0.4 x 22 mm	585198

Nucleus Hydrodissector

Tip is flattened on the horizontal plane.




585155

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G 7/8 in, 0.4 x 22 mm	585155

Micro J Shape Cannula

Hook is 1.5 mm long x 1.75 mm wide. 0.25 mm open end port.




585003

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585003

Nucleus Hydrodissector (Buratto)

Formed, 1 mm bent tip.




585199

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585199

Nucleus Hydrodissector (Chang)

90 Degree bend, 1 mm flattened tip.



585294

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/10 in, 0.4 x 18 mm, Bevel angle 33 degrees	585294
● 27G x 7/10 in, 0.4 x 18 mm, Bevel angle 47 degrees	585295

Nucleus Hydrodissector (Pearce)

Angled 35 degrees, 8 mm from end.



585036

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585037
● 27G x 7/8 in, 0.4 x 22 mm with 8 mm bend	585036

Nucleus Hydrodissector (Sauter)

Tip is flattened on the vertical plane.



585099

Qty: 10/Bx
Brand: Visitec®


Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585099



Capsule Polishers

Capsule Polisher (Kratz)

0.3 mm Top port, 35 degrees, 9 mm from end.




581605

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	581605

Flexible Tip Cannula

Angled 45 degrees, 7 mm from end, 10 mm sleeve.




585017

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585017
● 27G x 7/8 in, 0.4 x 22 mm	585005

Olive Tip Polisher

0.33 mm end port, 30 degrees, 10 mm from end and micro-etched 2 mm from end. Olive shape tip is 2 mm long x 1 mm wide.



585143


Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 1 in, 0.6 x 25 mm	585143

Cystotomes

Irrigating Cystotome

Formed cystotomes are bent to conform with the convexity of the crystalline lens.




581610

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 5/8 in, 0.5 x 16 mm	581610
● 27G x 5/8 in, 0.4 x 16 mm	581618

Irrigating Cystotome

Formed 8.5 mm from bend to tip. Shortened curve designed for deep set eyes or small pupils.



585008

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G 5/8 in, 0.5 x 16 mm	585004
● 27G 5/8 in, 0.4 x 16 mm	585008

Irrigating Cystotome (Blumenthal)

Formed 8.5 mm from bend to tip.



585042

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 5/8 in, 0.4 x 16 mm	585042

Cortex Extractors

Cortex Extractor

0.3 mm end port angled 10 degrees.
Flat end opening has rounded edges.




585106

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 24G x 1/2 in, 0.55 x 12.5 mm	585106

Cortex Extractor

0.4 mm top port.



585193

Qty: 10/Bx
Brand: Visitec®


Description	SKU #
● 23G x 5/8 in, 0.6 x 16 mm	585193



Lacrimal

Lacrimal Cannula

Intended to irrigate and probe the lacrimal system. Curved.




581615

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 21G x 1-3/8 in, 0.8 x 35 mm	581615
● 26G x 1-1/8 in, 0.45 x 28 mm	581276

Lacrimal Cannula

Straight, 13 mm Length.



585068

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 1/2 in, 0.5 x 13 mm	585068



FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT