

# BVI

## A Bright Vision on Illumination

### Light Fibers




Standard  
Light Fibers



Wide View  
Light Fibers



45 Degree Shielded  
Light Fibers



**Vitreco**

Vitreoretinal Surgical Products

[bvimedical.com](http://bvimedical.com)

# Light Fibers

Light Fibers featuring high-end specifications with consistent quality performance for ophthalmic endo-illumination procedures.

## Rigid shaft

- New hardened metal shaft reduces flexibility by 25%
- Reduced bend ability improves surgical control

## Ergonomically designed

- Ergonomic handle design for improved intraocular surgical control and precise surgical maneuvers
- Compatible with most endo-illumination systems due to light source adaptor




## Light fiber tip design

- In Standard, Wide view and 45 Degree Shielded tip designs


## Wide range of sizes

- Available in 23G, 25G and 27G with color coding for identifying gauge size instrument




## Standard Light Fibers

Straight tip for focal endo-illumination.



## Wide View Light Fibers

Wide view tip for global endo-illumination.



## 45 Degree Shielded Light Fibers

45 Degree shielded wide view tip for shielded endo-illumination.

# Light Fibers Options



## Standard Light Fibers

Description	SKU #
● 23G Standard light fibers	LF23.D01
● 25G Standard light fibers	LF25.D01
● 27G Standard light fibers	LF27.D01



## 45 Degree Shielded Light Fibers

Description	SKU #
● 23G 45 Degree shielded light fibers	LF23.D03
● 25G 45 Degree shielded light fibers	LF25.D03
● 27G 45 Degree shielded light fibers	LF27.D03



## Wide View Light Fibers

Description	SKU #
● 23G Wide view light fibers	LF23.D02
● 25G Wide view light fibers	LF25.D02
● 27G Wide view light fibers	LF27.D02

## Gauge Size with Colour Coding



Packaged box/5, Sterile



# Light Source Adaptors

The Light Fibers are compatible with most endo-illumination systems in combination with Light Source Adaptors\*.

Packaged box/1, Sterile

Vitreq® Light Source Adaptor	Description	SKU #
	Constellation® Light Source Adaptor	LF00.R01
	Stellaris® Light Source Adaptor	LF00.R02
	Geuder® Light Source Adaptor	LF00.R03
	Ruck® Light Source Adaptor	LF00.R04
	DORC® Light Source Adaptor	LF00.R05
	Vitreq® Light Fiber connector	LF23.D01

\* All product names, trademarks and registered trademarks are property of their respective owners.

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

Tel: 866-906-8080  
Email: [customersupport@bvimedical.com](mailto:customersupport@bvimedical.com)

