

BVI

Market Leading Portfolio of
Ophthalmic Surgical Devices

Table of Contents

Safety Knives	5
Safety Knives.....	6
Conventional Knives.....	8
Mini-Blades	11
Grindless Beaver® Mini-Blades.....	12
Blades.....	13
Trepines.....	14
Other Specialty Blades	15
Handles and Accessories.....	17
Handles and Accessories Chart	18
 Cannulas	 19
Anesthesia Cannulas.....	20
Anterior Chamber Cannulas.....	21
Anterior Chamber Maintainers.....	22
Hydrodissection Cannulas	23
Cystotome Cannulas.....	26
Irrigation Cannulas	28
Irrigation and Aspiration Cannulas	28
Polishers	30
Glaucoma Revision Picks	31
Lacrimal Cannulas	32
Refractive and FLACS	33
Hubs	34
 Single-Use Instruments	 35
Adjustable Speculums	36
Sprung Speculums	37
Pediatric Speculums.....	38
Fixation Forceps - Toothed.....	38
Fixation Forceps - Notched.....	39
Capsulorhexis Forceps.....	40
General Forceps	41
IOL Forceps	42
Tying Forceps.....	42
Choppers	43
Manipulators.....	44
Dissection Scissors	46
Dressing Scissors.....	46
Sprung Scissors	47
Dressing Forceps.....	48
Needle Holders and Drivers.....	48
Flap Elevators and Separators	49
Epithelium Spatulas, Hoes and Wells.....	51
Markers, Gauges and Calipers	52
Probes and Dilators.....	53
Depressors	54
Clamps and Curettes	54
Fixation Rings	55
Lacrimal Intubation CIS.....	55
IOL Glides	56
 Small Pupil Expansion	 57
I-Ring® Pupil Expander.....	58
 Fluid Control	 59
PVA Eye Spears.....	60
Cellulose Eye Spears.....	61

Sponge Points, Strips and Wicks	62
Corneal Light Shields	62
Full Body Drapes	63
Half Body Drapes	65
Face Drapes	68
Eye Drains	70
Fluid Collection	70
Instrument Wipes	71
Other Fluid Control	72
Non-Ophthalmic	73
Accessories	75
Eye Shields and Bandages	76
Diagnostic Lights	77
General Accessories	78
Custom Packs	80
CustomEyes®	81
Electrosurgery and Electrocautery	83
Cautery Devices	84
Erasers	86
Forceps	88
Dry Eye	89
Permanent Occlusion	90
Mid-Term Occlusion	92
Partial Occlusion	92
Temporary Occlusion	93
Single-Use Punctal Dilators and Probes	93
Other Related Occlusion Options	94
IOLs	95
IPURE®	96
Equipment	97
Endo Optiks®	98
Endo Optiks® Accessories	99
Vitreoretinal	101
Directional Chandelier	101
Micro Forceps and Scissors	103
90° Directional Laser Probes	107
Backflush Instruments	109
Vitreoretinal Consumables	113
Surgical Fluids for Ophthalmic Procedures	119
Surgical Fluids for Ophthalmic Procedures	120
Malosa® Single-Use Packs	121
Phaco Packs	122
LASIK Pack	122
Smile Pack	123
Core Surface Pack	123
Retina Pack	123
Product Index	124



High Quality, Single Handed Safety Activation Knives

Sharpness | Consistency | Safety




Options for Crescent,
MVR, Sideport and Slit

Safety Knives

Safety Slit Knives

Safety Slit Knives are available in a wide range of sizes designed for penetration and accurate incision width.



378224

Single Bevel

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
1.4 mm Slit, Angled, Single bevel	378214
1.5 mm Slit, Angled, Single bevel	378215
1.7 mm Slit, Angled, Single bevel	378217
1.8 mm Slit, Angled, Single bevel	378218
2 mm Slit, Angled, Single bevel	378220
2.2 mm Slit, Angled, Single bevel	378222
2.4 mm Slit, Angled, Single bevel	378224
2.5 mm Slit, Angled, Single bevel	370145
2.6 mm Slit, Angled, Single bevel	370146
2.75 mm Slit, Angled, Single bevel	378227
2.8 mm Slit, Angled, Single bevel	378228
3 mm Slit, Angled, Single bevel	378230
3.2 mm Slit, Angled, Single bevel	378232



378824

Double Bevel

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
2 mm Slit, Angled, Double bevel	378820
2.2 mm Slit, Angled, Double bevel	378822
2.4 mm Slit, Angled, Double bevel	378824
2.5 mm Slit, Angled, Double bevel	378825
2.65 mm Slit, Angled, Double bevel	378826
2.75 mm Slit, Angled, Double bevel	378827
2.8 mm Slit, Angled, Double bevel	378828
2.85 mm Slit, Angled, Double bevel	378829
3 mm Slit, Angled, Double bevel	378830
3.2 mm Slit, Angled, Double bevel	378832



Safety Knives

Safety Crescent Knives

Safety Crescent Knives are angled 55 degrees with the bevel up, designed for optimal tunnel architecture.



378234

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
2.5 mm Crescent, 55 Degrees, Bevel up	378234

Optimum Safety Sideport Knives

Safety Sideport Knives are designed for stab incisions.



378235

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
15 Degree Sideport, Straight	378235
22.5 Degree Sideport, Straight	378236
30 Degree Sideport, Straight	378237
45 Degree Sideport, Straight	378238

Safety Sideport Knives

Safety Sideport Knives are designed for accurate width stab incisions.



378212

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
1 mm Sideport, Angled 45 degrees	378210
1.2 mm Sideport, Angled 45 degrees	378212
0.6 mm 23G Sideport MVR	378233
1.15 mm 20G Sideport, Angled 45 degrees	378231

Kojo Safety Slit Knives

Kojo Safety Slit Knives are slit blades with an arcuate shape.



373025

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
2.5 mm Kojo Slit, 45 Degrees, Arc height is 0.55 mm	373025
3.2 mm Kojo Slit, 45 Degrees, Arc height is 0.71 mm	373032

Conventional Knives

Slit Knives

Slit Knives are available in a wide range of sizes designed for penetration and accurate incision width.



373727

Single Bevel

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
1.8 mm Slit, 45 Degrees, Single bevel up	373718
2.2 mm Slit, 45 Degrees, Single bevel up	373722
2.4 mm Slit, 45 Degrees, Single bevel up	373724
2.5 mm Slit, 45 Degrees, Single bevel up	373725
2.6 mm Slit, 45 Degrees, Single bevel up	373726
2.75 mm Slit, 45 Degrees, Single bevel up	373727
2.8 mm Slit, Straight, Matte	373328
2.8 mm Slit, 45 Degrees, Single bevel up	373728
3 mm Slit, Straight, Matte	373330
3 mm Slit, 45 Degrees, Single bevel up	373730
3.2 mm Slit, Straight, Matte	373332
3.2 mm Slit, 45 Degrees, Single bevel up	373732



370824

Double Bevel

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
2 mm Slit, 45 Degrees, Double bevel up	370820
2.2 mm Slit, 45 Degrees, Double bevel up	370822
2.4 mm Slit, 45 Degrees, Double bevel up	370824
2.5 mm Slit, 45 Degrees, Double bevel up	370825
2.65 mm Slit, 45 Degrees, Double bevel up	370826
2.75 mm Slit, 45 Degrees, Double bevel up	370827
2.8 mm Slit, 45 Degrees, Double bevel up	370828
3 mm Slit, 45 Degrees, Double bevel up	370830

Crescent Knives

Crescent Knives offer the flexibility to create optimal tunnel architecture. Available in straight, bevel up or bevel down.



373809

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
2.5 mm, 55 Degrees, Bevel up	373807
2.5 mm, 55 Degrees, Straight	373808
2.5 mm, 55 Degrees, Bevel down	373809

Sideport Knives

Sideport Knives are designed for accurate width stab incisions.



373712

Qty: 10/Bx
Brand: Beaver®

Description	SKU #
1 mm Sideport Knife, Angled	373710
1.2 mm Sideport Knife, Angled	373712

Conventional Knives

EdgeAhead Sideport and MVR



585232



376620



376630

Qty: Various/Bx
Brand: Beaver®

Description	SKU #
0.6 mm Stiletto, 5/Bx	585232
0.6 mm MVR, Straight, 5/Bx	585231
1.1 mm MVR, Straight, 5/Bx	585240
1.15 mm MVR, 20G, Straight, 10/Bx	376630
1.15 mm Sideport MVR, Angled, 10/Bx	376620

Mini-Blade Knife



Qty: 10/Bx
Brand: Beaver®

Description	SKU #
Surgical blade, 180 Degree cutting surface	374769

Sclerotome Knife



The Sclerotome knife has a smaller cutting edge than the Sclerotome blade (375700) and is designed to have a more precise, controlled dissection.



Qty: 10/Bx
Brand: Beaver®

Description	SKU #
Multi-sided sharp all around	375701



Conventional Knives

Optimum Knives

Optimum Knives are color coded for identification. Offered in a wide range of sizes.



374890

Qty: 10/Bx

Brand: Beaver®

Description	SKU #
15 Degrees, 30 Degree bend left	374890
15 Degrees, Straight, Green	374891
22.5 Degrees, Straight, Slate blue	374892
30 Degrees, Straight, Black	374893
45 Degrees, Straight, Tan	374894

Micro-Unitome Knives

Micro-Unitome Knives are manufactured using a proprietary grindless technology to support a more consistent sharpness. Each knife is color coded for identification.



Qty: 10/Bx

Brand: Beaver®

Description	SKU #
1.5 mm, 15 Degrees, Orange	377512
3 mm, 15 Degrees, Blue	377514
5 mm, 15 Degrees, Green	377516
4 mm, 22.5 Degrees, Light blue	377523
3.5 mm, 30 Degrees, Black	377531
3 mm, 45 Degrees, Tan	377546



Mini-Blades

Available in a variety of designs to meet all small incision and excision needs.

Mini-Blade



Qty: 24/Bx
Brand: Beaver®

Description	SKU #
Orthopedic Mini-Blade, Length 1.255 in (31.88 mm)	374562

Mini-Blade



Qty: 24/Bx
Brand: Beaver®

Description	SKU #
MIS 6400, Narrower version of Beaver6400 for minimally invasive surgery	376401

Mini-Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
Sharp on two sides with pointed tip	376300

Mini-Blade



Qty: 24/Bx
Brand: Beaver®

Description	SKU #
Pointed tip, Smaller alternative to #11 general blade	376500

BeaverGuard Guarded Blades

Preset to ensure uniform depth control.



375435

Qty: 3/Bx
Brand: Beaver®

Description	SKU #
0.35 mm Guarded depth S/Su	375435

Sclerotome Blade

The Sclerotome is a versatile blade suitable for glaucoma and retinal procedures.



Qty: 3/Bx
Brand: Beaver®

Description	SKU #
Multi-sided sharp all around S/Su	375700

Grindless Beaver® Mini-Blades

BVI's Mini-Blades are manufactured using a proprietary grindless technology and are designed to meet all small incision and excision needs. Mini-Blades are preassembled into a unique, patented tip, designed for safe assembly. Each Mini-Blade is color coded for identification.

6400 Mini-Blade



Qty: 20/Bx
Brand: Beaver®

Description	SKU #
Rounded tip, Sharp on one side, Green	BEAVER6400

6700 Mini-Blade



Qty: 20/Bx
Brand: Beaver®

Description	SKU #
Smaller alternative to #15 general blade, Orange	BEAVER6700

6600 Mini-Blade



Qty: 3/Bx
Brand: Beaver®

Description	SKU #
Lamellar, Angled 60 degrees, Smooth undersurface, Sharp all around, Bevel up, Black	BEAVER6600

6900 Mini-Blade



Qty: 20/Bx
Brand: Beaver®

Description	SKU #
Straight, Sharp all around, 180 Degree cutting surface, Double bevel blade, Blue	BEAVER6900

6610 Mini-Blade



Qty: 3/Bx
Brand: Beaver®

Description	SKU #
Lamellar, Straight, Smooth undersurface, Sharp all around, Black	BEAVER6610

6910 Mini-Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
Cleft palate blade, Angled 60 degrees, Sharp all around, 180 Degree cutting surface, Double bevel blade, Blue	BEAVER6910

Blades

Micro-Sharp® Blades

Micro-Sharp Blades are manufactured using a proprietary grindless technology. Each blade is color coded for identification.




Qty: 6/Bx
Brand: Beaver®

Description	SKU #
1.5 mm, 15 Degrees, Orange	377511
3 mm, 15 Degrees, Blue	377513
5 mm, 15 Degrees, Green	377515
4 mm, 22.5 Degrees, Light blue	377522
3.5 mm, 30 Degrees, Black	377530

Optimum Blades

Optimum stab blades are color coded for identification. Offered in a wide range of sizes.



Qty: 5/Bx
Brand: Beaver®

Description	SKU #
15 Degrees, Straight, Green	374881
22.5 Degrees, Straight, Slate blue	374882
30 Degrees, Black	374883



Trephines

Trephine Blades are handcrafted with a sharp, uniform cutting edge, designed for a precise, symmetrical cut of the corneal layers.

Overall length 15.8 mm and provided sterile for single-use.




Qty: 1/Bx
Brand: Beaver®

Description	SKU #
4 mm Trephine Blade	9550026
5 mm Trephine Blade	0009701
5.5 mm Trephine Blade	0009703
6 mm Trephine Blade	0009705
6.25 mm Trephine Blade	0009706
6.5 mm Trephine Blade	0009707
6.75 mm Trephine Blade	0009708
7 mm Trephine Blade	0009709
7.25 mm Trephine Blade	0009719
7.5 mm Trephine Blade	0009711
7.75 mm Trephine Blade	0009722
7.85 mm Trephine Blade	0009710
8 mm Trephine Blade	0009713

Description	SKU #
8.25 mm Trephine Blade	0009727
8.5 mm Trephine Blade	0009715
8.75 mm Trephine Blade	0009723
9 mm Trephine Blade	0009717
9.25 mm Trephine Blade	9550025
9.5 mm Trephine Blade	9550005
10 mm Trephine Blade	0009721
10.5 mm Trephine Blade	9550006
11 mm Trephine Blade	0009725
11.5 mm Trephine Blade	9550013
12 mm Trephine Blade	9550028
13 mm Trephine Blade	9550029
14 mm Trephine Blade	0009728

Other Specialty Blades


Myringotomy Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
Downcutting tip, Offset 45 degrees	377100

Myringotomy Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
Upcutting tip, Offset 45 degrees.	377110

Myringotomy Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
Offset tip, 45 degrees	377120

Myringotomy Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
Narrow shaft, Offset 45 degrees	377121

Myringotomy Knife



Qty: 3/Bx
Brand: Beaver®

Description	SKU #
Spear tip, Offset 45 degrees	377122

Myringotomy Knife



Qty: 3/Bx
Brand: Beaver®

Description	SKU #
Lance tip, Offset 45 degrees	377124

Other Specialty Blades

Beaver® Arthro-Lok® Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
4 mm Rounded tip	376481

Beaver® Arthro-Lok® Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
4 mm Mini-banana, Sharp all around, Curved	376980

Beaver® Arthro-Lok® Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
4 mm Banana, Sharp all around	376984

Beaver® Arthro-Lok® Blade



Qty: 6/Bx
Brand: Beaver®

Description	SKU #
4 mm Mini-meniscus, Flat	379081

Tympanoplasty Blade




Qty: 3/Bx
Brand: Beaver®

Description	SKU #
60 degree angle, Bevel down, Sharp all around, 2.5 mm width	377200

Handles and Accessories

Handles

Stainless steel handles are available in round and hex configurations, in a variety of lengths.



371395



374320



374340



375000

Type: Reusable

Qty: 1/Bx
Brand: Beaver®

Description	SKU #
Round handle, 13 cm Long x 0.5 cm Diameter	371350
Round handle, 13 cm Long x 0.66 cm Diameter	371360
Round handle, 18 cm Long x 4 mm Diameter with wrench	371394
Round handle, 21 cm Long x 4 mm Diameter with wrench	371395
Round handle, 10 cm Long x 0.66 cm Diameter	371410
Flat handle, 17 cm Long	372000
Round handle, 10 cm Long x 0.5 cm Diameter	374310
Round handle, 5 cm Long x 0.5 cm Diameter	374311
Round handle, 7.5 cm Long x 0.5 cm Diameter	374312
Round handle, 15.5 cm Long x 0.5 cm Diameter	374320
Hex handle, 10 cm Long x 0.5 cm Diameter	374330
Hex handle, 15.5 cm Long x 0.5 cm Diameter	374340
Flat handle, 13 cm Long, 1 cm Wide	375000

Wrenches and Collets

All wrenches and collets are stainless.



371301



371305



371415

Qty: 1/Bx
Brand: Beaver®


Description	SKU #
Blade wrench	371301
Collet replacement, 0.5 cm Diameter	371305
Collet replacement, 0.66 cm Diameter	371415

Handles and Accessories Chart

	HANDLES										REPLACEMENT COLLET
	371350	371360	371394	371395	374310	374311	374312	374320	374330	374340	371305
Beaver6400	x	x	x	x	x	x	x	x	x	x	
Beaver6600	x	x	x	x	x	x	x	x	x	x	
Beaver6610	x	x	x	x	x	x	x	x	x	x	
Beaver6700	x	x	x	x	x	x	x	x	x	x	
Beaver6900	x	x	x	x	x	x	x	x	x	x	
Beaver6910	x	x	x	x	x	x	x	x	x	x	
374562	x	x			x	x	x	x	x	x	
375435	x	x			x	x	x	x	x	x	x
375450	x	x			x	x	x	x	x	x	x
375700	x	x			x	x	x	x	x	x	x
376300	x	x			x	x	x	x	x	x	x
376401		x			x	x	x	x	x	x	x
376500	x	x			x	x	x	x	x	x	x
377511	x	x			x	x	x	x	x	x	
377513	x	x			x	x	x	x	x	x	
377515	x	x			x	x	x	x	x	x	
377522	x	x			x	x	x	x	x	x	
377530	x	x			x	x	x	x	x	x	
377120	x	x			x	x	x	x	x	x	x

BLADES





BVI

40 Years of Innovation and Quality Cannulas

Supporting All Areas of Ophthalmic Surgery

- A market leader of ophthalmic cannulas for over 40 years
- Full range of single-use cannulas developed to support all ophthalmic surgical procedures
- A history of innovation to meet your diverse surgical needs



Anesthesia Cannulas

Retrobulbar (Atkinson)

Intended for injecting anaesthetic agents into the muscle cone.




581638

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 1-1/2 in, 0.6 x 38 mm	581638
● 25G x 1-1/2 in, 0.5 x 38 mm	581637

Retrobulbar/Peribulbar (Atkinson)



585108

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 1-1/4 in, 0.6 x 32 mm	585001
● 25G x 1-1/4 in, 0.5 x 32 mm	585108
● 27G x 1-1/4 in, 0.4 x 32 mm	585109

Peribulbar (Atkinson)

Intended for injecting anaesthetic agents around the globe (periocular).




585144

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585144
● 25G x 7/8 in, 0.5 x 22 mm	585027
● 27G x 7/8 in, 0.4 x 22 mm	585145

Retrobulbar (Sharp Atkinson)



581640

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 1-1/2 in, 0.6 x 38 mm	581640
● 25G x 1-1/2 in, 0.5 x 38 mm	581641

Retrobulbar (Sharp Point)

Intended for injecting anaesthetic agents into the muscle cone.




581275

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 1-1/2 in, 0.5 x 38 mm	581275

Retrobulbar, Curved (Atkinson/Straus)

Designed to course around the globe and enter the anterior region of the muscle cone inferotemporal and parallel to the optic nerve. Straight portion 11 mm, curved portion 22 mm.



585019

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 1-3/8 in, 0.5 x 34 mm	585019

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 34.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

Anesthesia Cannulas

Sub-Tenon's Anesthesia Cannula

Intended for administering anaesthetic agents into the posterior sub-tenon's space.



585176

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 19G x 1 in, 1.1 x 25 mm	585176

Anterior Chamber Cannulas

Intended to be used for irrigation/injection of liquid.

Anterior Chamber (Rycroft)

Angled 45 degrees, 4 mm from end.



585025

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 21G x 7/8 in, 0.8 x 22 mm	581281
● 23G x 7/8 in, 0.6 x 22 mm	585025
● 25G x 7/8 in, 0.5 x 22 mm	581639
● 27G x 7/8 in, 0.4 x 22 mm	581280
○* 30G x 7/8 in, 0.3 x 22 mm	581273

Anterior Chamber Cannula

Angled 45 degrees, 6 mm from end.



581284

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	581642
● 27G x 7/8 in, 0.4 x 22 mm	581284
○* 30G x 7/8 in, 0.3 x 22 mm	581285

Anterior Chamber Cannula Straight



585006

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585006

Anterior Chamber Maintainers

Self Retaining ACM (Lewicky)

Intended to aid in maintaining anterior chamber space via irrigation. Flexible tubing 25 cm (10 in). Inner diameter 1 mm (0.04 in); outer diameter 2 mm (0.08 in). Unique annular rings start 0.75 mm from 55 degree bevelled tip for retaining without sutures.



Qty: 3/Bx
Brand: Visitec®

Description	SKU #
20G x 9/64 in, 0.9 x 3.5 mm	585061

ACM Micro Etched (Dahan)

Intended to aid in maintaining anterior chamber space via irrigation. Flexible tubing 15.25 cm (6 in). Inner diameter 0.75 mm (0.03 in), outer diameter 1.6 mm (0.07 in). Special micro etching on tip allows ACM to be self-retaining. Ultra thin wall tip provides high flow rate.



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
20G x 5/32 in, 0.9 x 4 mm, 45 degree beveled tip	585163

Oval ACM (Blumenthal)

Intended to aid in maintaining anterior chamber space via irrigation. 3.5 mm from 40 degree bevelled tip. Flexible tubing 25 cm (10 in). Inner diameter 1 mm (0.04 in), outer diameter 2 mm (0.08 in). Smooth tip allows for straight forward insertion and minimal tissue disturbance.



Qty: 5/Bx
Brand: Visitec®

Description	SKU #
Tapered from 1.25 mm (18G) to 0.9 mm (20G), oval	585149

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 34.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

Hydrodissection Cannulas

Intended to separate the cortex from the capsule via irrigation (hydrodissection).

Nucleus Hydrodissector (Pearce)

Angled 35 degrees, 8 mm from end.




585036

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585037
● 27G x 7/8 in, 0.4 x 22 mm with 8 mm bend	585036
● 27G x 7/8 in, 0.4 x 22 mm	585274

Nucleus Hydrodelineator (P. Koch)

Intended to break up lens contents and facilitate their aspiration via irrigation.



585015

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585015

Cortical Cleaving Hydrodissector (Flat End)

45 degrees, 11 mm from end. For cleaving cortex from capsular bag.



585158

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G 7/8 in, 0.4 x 22 mm	585158

Nucleus Hydrodissector (Sinskey)

35 degrees, 8 mm from end.



585195

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585195

Nucleus Hydrodissector

45 degrees, 7 mm from end.



581592

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
○* 30G x 7/8 in, 0.3 x 22 mm	581592

Nucleus Hydrodissector (Sauter)

Tip is flattened on the vertical plane.



585099

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585099

Hydrodissection Cannulas

Intended to separate the cortex (tissue) from the capsule via irrigation (hydrodissection).

Nucleus Hydrodissector

Tip is flattened on the horizontal plane.



585155

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G 7/8 in, 0.4 x 22 mm	585155

Nucleus Hydrodissector (Chang)

90 Degree bend, 1 mm flattened tip.



585294

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/10 in, 0.4 x 18 mm, Bevel angle 33 degrees	585294
● 27G x 7/10 in, 0.4 x 18 mm, Bevel angle 47 degrees	585295

Tapered Hydrodissector (Blumenthal)

Tip is tapered to 0.3 mm (30G), 45 degrees, 4 mm from end.



585107

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, .5 x 22 mm	585107

Cortical Cleaving Hydrodissector (Round End)

45 degrees, 11 mm from end.



585157

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585157

Nucleus Hydrodissector (Buratto)

Formed, 1 mm bent tip.




585199

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585199

Cortical Cleaving Hydrodissector (Bolger)

Tip is tapered to 0.3 mm (30G), 33 degrees, 1 mm from end.



585159

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585159

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 34.


19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

Hydrodissection Cannulas

Intended to separate the cortex (tissue) from the capsule via irrigation (hydrodissection).

Hydrodissection Cannula (Akahoshi Style)

Tapered tip, thin wall, angled 1 mm from end.



585198

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G 7/8 in, 0.4 x 22 mm	585198

Micro J Shape Cannula

0.25 mm open end port.
Hook is 1.5 mm long x 1.75 mm wide.



585018

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585003
● 27G x 7/8 in, 0.4 x 22 mm	585018

Micro J Shape Cannula (Hessburg)

Hook is 0.85 mm long x 1.7 mm wide.



585148

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 3/4 in, 0.5 x 20 mm	585148

Micro J Shape Hydrodissector (Pearce)

Hook is 2.25 mm long x 1.75 mm wide.



585029

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585029

Cortical Cleaving Hydrodissector (Grabow)

40 degrees, 12 mm from end.



585197

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
○* 30G x 5/8 in, 0.3 x 16 mm	585197

HydroExpression Cannula (Corydon)

1.8 mm from end.



585191

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585091
● 25G x 7/8 in, 0.5 x 22 mm	585090
● 27G x 7/8 in, 0.4 x 22 mm	585191

Cystotome Cannulas

Intended to perform a capsulorhexis. Indicated for use in cataract surgery.

Irrigating Cystotome



581610

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 5/8 in, 0.6 x 16 mm	581599
● 25G x 5/8 in, 0.5 x 16 mm	581610
● 27G x 5/8 in, 0.4 x 16 mm	581618
○* 30G x 5/8 in, 0.3 x 16 mm	581630

Irrigating Cystotome (Blumenthal)



585043

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 5/8 in, 0.4 x 16 mm	585043

Irrigating Cystotome

Formed 8.5 mm from bend to tip.
Shortened curve designed for deep set eyes or small pupils.



585008

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G 5/8 in, 0.5 x 16 mm	585004
● 27G 5/8 in, 0.4 x 16 mm	585008

Irrigating Formed Cystotome

Tip 45 degrees.



585049

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 5/8 in, 0.5 x 16 mm	585049

Irrigating Cystotome



581594

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 5/8 in, 0.6 x 16 mm	581594
● 25G x 5/8 in, 0.5 x 16 mm	581617

Irrigating Side Cutting Cystotome



581602

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 5/8 in, 0.5 x 16 mm	581602

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 34.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

Cystotome Cannulas

Intended to perform a capsulorhexis. Indicated for use in cataract surgery.

Irrigating Cystotome (Blumenthal)

Formed 8.5 mm from bend to tip.



585113

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 5/8 in, 0.5 x 16 mm	585113
● 27G x 5/8 in, 0.4 x 16 mm	585042

Irrigating Cystotome (Simcoe)

45 Degrees and 10.5 mm from tip.



585093

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 1/2 in, 0.4 x 12.5 mm	585093

Blunt Tip Irrigating Cystotome (Kellan)

Tip angled toward bevel. Unique blunt tip creates a wide capsular flap. Useful for pear shape capsulotomies.



585115

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585115

Irrigating Cystotome (Joyce)

Tip 0.5 mm long.



585030

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 5/8 in, 0.5 x 16 mm	585030
● 27G x 5/8 in, 0.4 x 16 mm	585039

Irrigating Cystotome (Reverse)

Tip 90 degrees toward bevel.



585034

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 5/8 in, 0.5 x 16 mm	585034

Irrigation Cannulas

Intended to be used for irrigation/injection of liquid and/or air.

Anterior Chamber Cannula

30 degrees, 6 mm from end.



581277

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
○ 16G x 1 in, 1.6 x 25 mm	581277

Bishop Harmon Cannula

Intended for irrigation of the eye.
40 degrees, 8 mm from end.



585040

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 19G x 1 in, 1.1 x 25 mm	585040

ViscoExpression Cannula (Corydon)



131618

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 5/8 in	131618

Irrigation and Aspiration Cannulas

Intended to inject ophthalmic viscosurgical devices (OVD).

Viscoflow Formed

Formed 11 mm from end, radius is 19 mm.



585007

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585007

Viscoflow Cannula

Angled 45 degrees, 9 mm from end.



585009

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585009
● 25G x 7/8 in, 0.5 x 22 mm	585045
● 27G x 7/8 in, 0.4 x 22 mm	585026
○* 30G x 7/8 in, 0.3 x 22 mm	585046

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 34.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

Irrigation and Aspiration Cannulas

Intended to aspirate the cortex.

Cortex Extractor

0.3 mm top port.



Qty: 10/Bx
Brand: Visitec®

585075

Description	SKU #
● 23G 5/8 in, 0.6 x 16 mm	585075

Cortex Extractor

0.4 mm top port, thin wall cannula,
30 degrees, 2 mm from end, micro-etched
2 mm from end.



Qty: 10/Bx
Brand: Visitec®

585227

Description	SKU #
● 23G x 15/16 in, 0.6 x 23.5 mm	585227

Cortex Extractor

0.3 mm end port angled 10 degrees.
Flat end opening has rounded edges.




Qty: 10/Bx
Brand: Visitec®

585106

Description	SKU #
● 24G x 1/2 in, 0.55 x 12.5 mm	585106

Bimanual Aspiration Cannula

Thin wall cannula, curved 45 degrees,
0.51 mm open end irrigation port.



Qty: 10/Bx
Brand: Visitec®

585222

Description	SKU #
● 23G x 5/8 in, 0.6 x 17 mm	585222

Cortex Extractor

0.4 mm top port.



Qty: 10/Bx
Brand: Visitec®

585193

Description	SKU #
● 23G x 5/8 in, 0.6 x 16 mm	585193

Bimanual Aspiration Cannula

Thin wall cannula, curved 45 degrees,
0.51 mm open end irrigation port.



Qty: 10/Bx
Brand: Visitec®

585223

Description	SKU #
● 21G x 5/8 in, 0.8 x 17 mm	585223

Irrigation and Aspiration Cannulas

Intended for irrigation of the eye and/or aspiration of liquids and body fluids.

I and A Cannula Thinwall (Simcoe)

Tubing: Inner diameter 1 mm (0.04 in), outer diameter 2 mm (0.08 in), 0.3 mm aspiration port.



585047

Qty: 5/Bx
Brand: Visitec®

Description	SKU #
● 23G, 0.6 mm	585047

I and A Cannula (Simcoe)

Tubing: Inner diameter 1 mm (0.04 in), outer diameter 2 mm (0.08 in), 0.4 mm aspiration port.



585110

Qty: 3/Bx
Brand: Visitec®

Description	SKU #
● 23G, 0.6 mm	585110

I and A Cannula (Simcoe)

Tubing: Inner diameter 1 mm (0.04 in), outer diameter 2 mm (0.08 in), 0.3 mm aspiration port.



585086

Qty: 3/Bx
Brand: Visitec®

Description	SKU #
● 23G, 0.6 mm	585086

Polishers

Intended to polish the anterior capsule.

Flexible Tip Cannula

Angled 45 degrees, 7 mm from end, 10 mm sleeve.



585041

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585041
● 25G x 7/8 in, 0.5 x 22 mm	585017
● 27G x 7/8 in, 0.4 x 22 mm	585005

Capsule Polisher

0.5 mm open end port, 30 degrees, 3 mm from end and micro-etched 5 mm from end.



581614

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 21G x 7/8 in, 0.8 x 22 mm	581614

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 34.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

Polishers

Intended to polish the anterior capsule.

Capsule Polisher (Kratz)

0.3 mm Top port, 35 degrees, 9 mm from end.



581604

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 21G x 7/8 in, 0.8 x 22 mm	581604
● 23G x 7/8 in, 0.6 x 22 mm	581605



Olive Tip Polisher

0.33 mm end port, 30 degrees, 10 mm from end and micro-etched 2 mm from end. Olive shape tip is 2 mm long x 1 mm wide.



585143

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 1 in, 0.6 x 25 mm	585143



Posterior Scraper (Simcoe)

Intended to polish the anterior capsule.
0.33 mm open end port, 60 degrees, 7 mm from end.



581645

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	581645



Olive Tip Cannula

0.33 mm end port. Olive shape tip is 2 mm long x 1 mm wide.



585147

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585147



Glaucoma Revision Picks

Intended for bleb needling.

Glaucoma Revision Pick 60 Degree (Lederer)

60 Degrees, 7 mm from end.
Ideal for superior approach.



585133

Qty: 5/Bx
Brand: Visitec®

Description	SKU #
● 27G x 5/8 in, 0.4 x 16 mm	585133



Glaucoma Revision Pick 30 Degree (Lederer)

30 Degrees, 13.5 mm from end.
Ideal for temporal or nasal approach.



585135

Qty: 5/Bx
Brand: Visitec®

Description	SKU #
● 27G x 5/8 in, 0.4 x 16 mm	585135



Lacrimal Cannulas

Lacrimal Cannula

Intended to irrigate and probe the lacrimal system. Curved.



581615

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 21G x 1-3/8 in, 0.8 x 35 mm	581615
● 26G x 1-1/8 in, 0.45 x 28 mm	581276

Lacrimal Cannula (Shahinian)

Intended to irrigate and probe the lacrimal system. 0.3 mm aspiration port.



585033

Qty: 5/Bx
Brand: Visitec®

Description	SKU #
● 23G x 1/2 in, 0.6 x 12 mm	585033

Lacrimal Cannula

Intended to irrigate and probe the lacrimal system. For irrigating punctum or lacrimal duct.



585016

Qty: 5/Bx
Brand: Visitec®

Description	SKU #
19G x 3/4 in (1.1 x 20 mm) shaft attached to 23G x 1/2 in (0.6 x 12 mm) malleable cannula	585016

Lacrimal Cannula

Straight, 13 mm Length.



585068

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 1/2 in, 0.5 x 13 mm	585068

Lacrimal Trephine (Sisler)

Intended for rotary trepanation of an obstructed canaliculus and/or biopsy of the conjunctiva or sclera.



585031

Qty: 3/Bx
Brand: Visitec®

Description	SKU #
● 21G x 1-1/2 in, 0.8 x 38 mm	585031

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 34.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

Refractive and FLACS

LASIK Irrigating Cannula (Lindstrom)

Intended to irrigate and/or manipulate corneal tissue. Features three 0.25 mm ports located on 7 mm extension with an 11 mm radius, 45 degrees.



585250

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585250

LASIK Flap Spatula/Irrigator (Hersh)

Intended to irrigate and/or manipulate corneal tissue. 45 Degrees, 9 mm from spatulated tip.



585252

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585252

LASIK Irrigating Cannula (Solomon)

Intended to irrigate and/or manipulate corneal tissue. Thin wall cannula, angled tip features three 0.25 mm ports on 8 mm extension.



585280

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585280

LASIK Cannula (Daya)

Intended to irrigate and/or manipulate corneal tissue. 45 degrees, 10 mm from end. Flattened 0.3 mm, 1 mm from tip.



585271

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585271

FLACS Cannula (Schweitzer)

Intended to open clear corneal incisions created by the Femtosecond laser, injection of OVD (ophthalmic viscosurgical device) viscoelastic solution into the anterior chamber and to check that capsulorhexis is released for the capsule.



585200

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585200

LASIK Cannula

Intended to irrigate and/or manipulate corneal tissue. 25 Degrees, formed 7 mm from end.



585278

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585278

Hubs


BVI product hubs comply with ISO 6009 Color coding and ISO 594 Conical fittings with a 6% luer taper to assure use with syringes or devices.

COLOR	HUB	GAUGE	ID (IN)	ID (MM)	OD (IN)	OD (MM)
Natural* (Clear)		30G	0.006	0.15	0.012	0.3
Grey		27G	0.008	0.2	0.016	0.4
Brown		26G	0.01	0.25	0.018	0.45
Orange		25G	0.01	0.25	0.02	0.5
Purple		24G	0.012	0.3	0.022	0.55
Blue		23G	0.013	0.33	0.025	0.6
Black		22G	0.016	0.41	0.028	0.7
Green		21G	0.02	0.51	0.032	0.8
Yellow		20G	0.024	0.6	0.036	0.9
Cream		19G	0.028	0.7	0.043	1.1

*Not an ISO color code.

ISO Hub Gauge Chart

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----




BVI

High Quality Single-Use Instruments

Providing Sterile Options for Your Operating Room Procedures

- Single-use instruments offer protection against infection, due to decreased risk of cross-contamination¹
- Surgeons benefit from precision and reliability of brand-new instruments for each patient²



References

¹ Southworth, P.M. Infections and exposures: reported incidents associated with unsuccessful decontamination of reusable surgical instruments. The Journal of Hospital Infection. September 17, 2014

² Scheib, Christian; Single-use instruments deserve consideration for eye surgery: Advantages include cost savings, ease of use and high-quality instruments; December 1, 2016

Adjustable Speculums

Plastic Locking Eye Speculum



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Solid blade, Gear controlled spread, Overall length 83 mm	585297

Mechanical Speculum, Open Blade (PC)



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
5 mm x 14 mm Open blades, 1 mm Wire, Polycarbonate body and screw mechanical speculum	MMSU1290S

Eye Speculum (Lieberman)



Qty: 10/Bx
Brand: Visitec®


Description	SKU #
Open 14 mm wire blades with screw controlled spread	581478US

Mechanical Speculum Solid Blade



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
6 mm x 14.5 mm Solid curved polycarbonate blades, 0.7 mm Thick, 1 mm Wire, Polycarbonate body and screw mechanical speculum	MMSU1488S



Sprung Speculums

Eye Speculum



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Open wire loops, 40 x 30 mm, 14 mm blades, 0.8 mm wire size	581413

Kratz-Barraquer Speculum, Temporal



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
6.5 mm x 14 mm Open blades, 0.9 mm Wire, 30 mm Wide	MMSU1253S

Eye Speculum Open Wire Loops



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Open wire loops. 40 x 30 mm, 14 mm blades, 0.9 mm wire size	581443

Barraquer Speculum, Flat Blade



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
6.5 mm x 14.5 mm Straight closed blades, 1 mm Wire, 30 mm Wide	MMSU1275S

Eye Speculum (Kratz/Barraquer)



Qty: 3/Bx
Brand: Visitec®

Description	SKU #
Open wire loops, 40 x 26.5 mm, Retraction size 13 x 4.5 mm	585094

Kratz-Barraquer Speculum, Curved



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
6.5 mm x 14.5 mm Open curved blades, 0.9 mm Wire, 30 mm Wide	MMSU1287S

Kratz-Barraquer Speculum



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
6.5 mm x 14 mm Open blades, 0.9 mm Wire, 30 mm Wide	MMSU1252S

Barraquer Speculum, Solid Curved (PC)



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
6 mm x 14.5 mm Solid curved polycarbonate blades, 0.7 mm Thick, 0.9 mm Wire, 30 mm	MMSU1474S

Sprung Speculums

Barraquer Speculum, Solid Blade (PC)



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
7 mm x 14.5 mm Solid straight polycarbonate blades, 0.8 mm Thick, 0.9 mm Wire, 30 mm Wide	MMSU1475S

Barraquer Speculum, Solid 18 mm Blade (PC)



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
6 mm x 14.5 mm Solid curved polycarbonate blades, 0.8 mm Thick, 1 mm Wire	MMSU1479S

Pediatric Speculums

Pediatric Eye Speculum



Qty: 3/Bx
Brand: Visitec®

Description	SKU #
Overall size: 34 x 21 mm, Retraction size: 10 x 3.5 mm	585203

Fixation Forceps - Toothed

Fine Toothed Forceps Short (Bonn)



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.3 mm 1 x 2 Teeth, 5 mm Tying platforms, 52 mm Handle	MMSU1123S

Adson Forceps, Toothed



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
1.35 mm x 1.8 mm 1 x 2 Teeth, 122 mm Long	MMSU0238S

Fixation Forceps - Toothed


Toothed Forceps (Bonn)



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.45 mm 1 x 2 Teeth, 5 mm Tying platforms, 80 mm Handle	MMSU1124S

Colibri Toothed Forceps, Short



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.4 mm 1 x 2 Teeth, 5 mm Tying platforms, 52 mm Handle	MMSU1129S

Fixation Forceps - Notched

Fine Notched Forceps, Short (Hoskins)



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.2 mm Notched tips, 5 mm Tying platforms, 52 mm Handle	MMSU1114S

Colibri Notched Forceps, Short



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.2 mm Notched tips, 5 mm Tying platforms, 52 mm Handle	MMSU1209S

Notched Forceps (Hoskins)



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.3 mm Notched tips, 5 mm Tying platforms, 80 mm Handle	MMSU1115S

IOL Removal Forceps



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.1 mm x 0.1 mm Ridges along 4 mm of 15 mm shaft, Angled at 30 degrees, 52 mm Handle	MMSU1460S

Capsulorhexis Forceps

MICS Capsulorhexis Forceps - Curved/Sharp



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.4 mm x 0.3 mm Sharp tips, 12 mm Curved direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1567S

Capsulorhexis Forceps Grasping/Curved Tips - Gen 3

Grasping curved tips, round handle.



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Capsulorhexis Grasping/Curved Tips - Gen 3	581486

MICS Capsulorhexis Forceps - Straight/Sharp




Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.4 mm x 0.3 mm Sharp tips, 12 mm Straight direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1568S

Capsulorhexis Forceps Sharp/Curved Tips - Gen 3

Sharp curved tips, round handle.



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Capsulorhexis Forceps Sharp/Curved Tips - Gen 3	581488

Capsulorhexis Forceps (Utrata) Straight



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Straight, Cystotome tip, 11 mm Long, Flat handle, Overall length 100 mm	581481US

General Forceps

Adson Forceps, Non-Toothed



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.85 mm x 1 mm Tips, 14.5 mm Serrated jaws, 122 mm Long	MMSU0241S

Epilation Forceps



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
5 mm x 7 mm Flat plane tips, 70 mm Handle	MMSU0085S

Jewellers Forceps, Fine Angled



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
7 mm 45 degree angled pointed tips, 80 mm handle	MMSU0283S

Jewellers Forceps, Fine Straight



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
7 mm Pointed tips, 80 mm Handle	MMSU0286S

Moorfields Forceps



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.6 mm tips, 11.5 mm Serrated jaws, 105 mm Long	MMSU0277S

Jewellers Forceps



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
No. 5, Straight forceps, flat handle, Overall length 100 mm	581458

St Martins Forceps



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.9 mm x 0.26 mm 1 x 2 Teeth, 52 mm Handle	MMSU0279S

LASIK Flap-Lifting Forceps



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.3 mm x 11.5 mm Tips, 45 Degree angled shaft, 80 mm Handle	MMSU1169S

General Forceps

Contact Lens Forceps



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.6 mm x 0.4 mm Tips, 9 mm Curved Platforms, 70 mm Handle	MMSU1284S

IOL Forceps

Lens Loading Forceps




Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.6 mm x 8 mm Tips, 45 Degree angled shaft, 80 mm Handle	MMSU1146S

Tying Forceps

Straight Tying Forceps



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.35 mm Tip, 5 mm Tying platforms, 80 mm Handle	MMSU1125S

Angled Tying Forceps (McPherson)




Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.35 mm Tip, 5 mm Tying platforms, 8 mm Shaft angled 45 degree, 80 mm Handle	MMSU1127S

Tying Forceps

Curved Tying Forceps



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.25 mm Tips, 5 mm Curved tying platforms, 80 mm Handle	MMSU1210S

Tying Forceps



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Curved, 6 mm Tying platform, Flat handle, Overall length 96 mm	581411US

Waqar Suture Removal Forceps



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.3 mm x 0.65 mm Pick angled 90 degree, 5 mm Tying platforms, Curved 28 mm shaft, 80 mm Handle	MMSU1583S

Micro Tying Forceps



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Straight, 6 mm Tying platform, Flat handle, Overall length 105 mm	581459US

Choppers

Claw Chopper



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Claw Chopper	MMSU1105S

Greens Chopper, Right Handed



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Greens Chopper, Right Handed	MMSU1119S

Choppers

Phaco Chopper



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Phaco Chopper	MMSU1108S

Nucleus Chopper (P. Koch)



Qty: 3/Bx
Brand: Visitec®

Description	SKU #
45 Degrees, 6 mm from tip, tip angled 90 degrees, 1.5 mm from end	585181

Manipulators

Govan Manipulator



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Govan Manipulator	MMSU1102S

Mushroom Manipulator



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Mushroom Manipulator	MMSU1179S

Sinskey Hook



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Sinskey Hook, 0.2 mm	MMSU1103S

Lester Hook



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Lester Hook	MMSU1207S

Drysdale Manipulator



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Drysdale Manipulator	MMSU1144S

Lebuisson Manipulator




Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Lebuisson Manipulator	MMSU1458S

Manipulators


Phaco Chopper/ Mushroom Manipulator



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Phaco Chopper/Mushroom Manipulator	MMSU1221S


Sinskey Hook - Gen 2



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Hook Sinskey Disposable - Gen 2	581442

Horizontal Chopper/ Kuglen Hook



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Horizontal Chopper/Kuglen Hook	MMSU1231S

IOL Manipulator Blunt Tip



Qty: 3/Bx
Brand: Visitec®

Description	SKU #
Straight, Blunt tip, Fits 0.2 mm diameter holes and larger	585062

ICL Manipulator



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
ICL Manipulator	MMSU1153S

Nucleus Rotator



Qty: 3/Bx
Brand: Visitec®

Description	SKU #
45 Degrees, Blunt forked tip, 9 mm from end, For rotating nucleus during phacoemulsification	585077

Visitec® Micro-Manipulator, Straight



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Micro-Manipulator, Straight	581441

Dissection Scissors

Stevens Tenotomy, Straight, Rounded



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
1.5 mm x 1.5 mm Rounded tips, Straight blades, 110 mm Long	MMSU0214S

Iris Scissors, Straight, Pointed



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.1 mm x 0.5 mm Pointed tips, 32.5 mm Straight blades, 120 mm Long	MMSU0249S

Iris Scissors, Curved, Pointed



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.1 mm x 0.7 mm Pointed tips, 25 mm Curved blades, 115 mm Long	MMSU0244S

Dressing Scissors

Dressing Scissors, Straight, Pointed



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.3 mm Pointed tips, 24.5 mm Straight blades, 102 mm Long, Plastic handle	MMSU0454S

Dressing Scissors, Straight, Rounded



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.8 mm Rounded tips, 24.5 mm Straight blades, 102 mm Long, Plastic handle	MMSU0455S

Sprung Scissors

Vannas Scissors



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Curved, 9.5 mm Blades, with sharp pointed tips, Overall length 90 mm	581422US

Westcott Scissors



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
23 mm Blades with blunt tips, Overall length 105 mm	581437US

Vannas Scissors, Curved



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Pointed tips on 8 mm curved blades, 1.5 mm when closed, 34 mm Handle	MMSU1140S

Westcott Scissors, Curved, Pointed



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.5 mm x 0.6 mm Pointed tips, 17.5 mm Curved blades, 36 mm Long	MMSU1211S

Vannas Scissors, Straight



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Pointed tips on 8 mm straight blades, 1.5 mm when closed, 34 mm Handle	MMSU1141S

Westcott Scissors, Straight, Pointed



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Pointed tips, 17 mm Straight blades, 120 mm Long	MMSU1214S

Westcott Scissors Curved Dolphin Nose





Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.5 mm x 0.6 mm Rounded tips, 22 mm Curved 'Dolphin Nose' blades, 120 mm Long	MMSU1213S

Dressing Forceps


Straight Rampley



Qty: 10/Bx
Brand: Malosa®

Description	SKU #
11 mm x 23 mm Serrated tips, 180 mm Long, Nylon plastic	MMSU0068S

Rampley Short



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
17 mm x 24 mm Serrated tips, 125 mm Long, Nylon plastic	MMSU1578S

Needle Holders and Drivers

Corneal Needle-Holder



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
9 mm Fine curved jaws, 0.45 mm Tip, 106 mm Long	MMSU1110S

Barraquer Needle-Holder



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.42 mm x 0.45 mm Tips, 4 mm Fine curved jaws with 5 mm shaft, 138 mm Long	MMSU3015S

Troutman Needle-Holder



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
4 mm Fine curved jaws, 0.55 mm Tip, 76 mm Flat handle	MMSU1243S

Flap Elevators and Separators

LASIK Flap Separator/Elevator (Seibel)



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.16 mm x 0.28 mm x 1 mm Tip, 20 mm Straight shaft, 0.5 mm x 10.5 mm pointed hook, 52 mm Handle	MMSU1172S

LASIK Flap Separator and Elevator, Binder



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.1 mm x 0.3 mm Rounded tip 0.7 mm long, 10 mm Curved 60 degree angled shaft, 75 mm Handle	MMSU1206S

LASIK Flap Elevator Curved 8 mm



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.2 mm x 0.15 mm Tip, Curved 8 mm tapered shaft at 90 degrees, 75 mm Handle	MMSU1177S

LASIK Flap Separator/Elevator Curved



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.15 mm x 0.3 mm x 8 mm Curved separator at 90 degree to shaft, 0.15 mm x 0.3 mm x 1 mm Pocketing hook, 52 mm Handle	MMSU1294S


Y Shape Corneal Marker



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
2 mm x 8 mm Marking blades, 75 mm Handle	MMSU1134S

Reinstein Lenticule Separator




Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.22 mm x 2 mm Pocketing hook, 1.2 mm x 2 mm Spatula tip, Curved separator at 90 degrees to shaft, 52 mm Handle	MMSU1297S

Flap Elevators and Separators

LASIK Flap Separator/ Elevator Straight




Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.15 mm x 0.3 mm x 8 mm Straight separator at 90 degree to shaft, 0.15 mm x 0.3 mm x 1 mm Pocketing hook, 52 mm Handle	MMSU1295S

LASIK Flap Roller

Used after realignment. Acts as a squeegee to remove excess fluid from the interface; smooths flap and removes folds.



Qty: 5/Bx
Brand: Visitec®

Description	SKU #
14 mm Sleeve	585233



Epithelium Spatulas, Hoes and Wells

Corneal Epithelium Spatula Flat



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
5 mm Spatula with 0.9 mm flat blade, 0.8 mm x 1.7 mm Round bladed tip, 75 mm Handle	MMSU1113S

Epithelium Well



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
3 mm x 8.5 mm ID well, Angled 45 degrees, 75 mm Handle	MMSU1197S
3 mm x 9 mm ID well, Angled 45 degrees, 75 mm Handle	MMSU1196S

Corneal Epithelium Spatula Round



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
2.65 mm x 5.5 mm x 0.35 mm Rounded blade, 75 mm Handle	MMSU1187S

Spatula PRK



Qty: 5/Bx
Brand: Visitec®

Description	SKU #
Spatula tip and top are slightly sharpened to facilitate more aggressive scraping.	585213



Markers, Gauges and Calipers

Mendez Gauge



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
2.5 mm Polycarbonate ring with 13 mm ID, 0-180 Degree marks in increments of 10, 75 mm Handle	MMSU1152S

Pointed IVT Marker 3.5/4, Blue



Qty: 10/Bx
Brand: Malosa®

Description	SKU #
2 x 0.55 mm Tips with 3.5 mm spread, 2 x 0.65 mm Tips with 4 mm spread, Polycarbonate, 108 mm Long	MMSU1277S

Braunstein IVT Marker 3.5/4



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
2 mm x 1.5 mm Tips with 3.5 mm spread, 2 mm x 1.5 mm Tips with 4 mm spread, 90 mm Long	MMSU1229S

Double Ended IVT Guide Pad



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
7 mm Pad with 3.5 guide, 0.8 mm Pad with 4 mm guide, 6 mm 45 Degree angled shafts, 90 mm Long	MMSU1285S

Braunstein IVT Marker 3.5/4, Green



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
2 mm x 1.5 mm Tips with 3.5 mm spread, 2 mm x 1.5 mm Tips with 4 mm spread, Polycarbonate, 108 mm Long	MMSU1276S

Intravitreal Injection Guide



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Polycarbonate lash guard with stainless-steel swage guide tube, Triangular base plate with 3 stabilizing studs, 39 mm, 45 Degree angled wishbone handle	MMSU1577S

Probes and Dilators

Punctal Dilator and Probe



Qty: 3/Bx
Brand: Visitec®

Description	SKU #
25 mm (1 in) medium taper, Blunt, Wilder style dilator mated with a size 0 Bowman lacrimal probe, 0.7 mm x 45 mm (22G x 1 3/4 in)	585124


Golf Club Spud



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.1 mm x 0.8 mm x 2 mm Oval profile, 18 mm Straight shaft, 75 mm Handle	MMSU1417S


Punctum Dilator, Fine/Medium



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.4 mm x 21.5 mm Tip, 0.6 mm x 18 mm Tip, 77 mm Handle	MMSU1181S

Punctum Dilator/ Lacrimal Probe



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.28 mm x 25 mm Dilator, 0.7 mm x 44.5 mm Probe, 70 mm Handle	MMSU1419S

Punctum Dilator/Seeker




Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.45 mm x 21.5 mm Tip, 5 mm Pointed tip, 18 mm Straight shaft, 77 mm Handle	MMSU1193S



Depressors

Schocket Double-Ended Scleral Depressor



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
2.7 mm x 6.5 mm T shaped tip, 5.5 mm x 2.2 mm Rounded tip, 68 mm Handle	MMSU1405S

Ball Tip Scleral Depressor



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
3 mm Ball tip, 21.5 mm Straight shaft, 50 mm Handle	MMSU1408S

Clamps and Curettes


Chalazion Clamp, Oval



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
13 mm x 20 mm Oval profile, 10 mm x 17 mm Aperture, 90 mm Long	MMSU1149S

Chalazion Curette, Double, 1 mm/2 mm



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
1 mm Curette, 1 mm Depth, 2 mm Curette, 1.35 mm Depth, 128 mm Long	MMSU1157S

Chalazion Clamp, Round



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
12 mm Round profile, 9 mm Round aperture, 90 mm Long	MMSU1150S

Chalazion Curette, Double, 2 mm/3 mm



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
2 mm Curette, 1.35 mm Depth, 3 mm Curette, 1.6 mm Depth, 128 mm Long	MMSU1158S

Fixation Rings

Thornton Fixation Ring



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
0.5 mm Thick 6.5 mm radius ring, 2 mm x 0.9 mm Teeth at 5 mm intervals, 75 mm Handle	MMSU1137S

Fixation Ring



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
0.5 mm Thick 13 mm internal diameter ring, 1.5 mm x 0.6 mm Teeth at 3.5 mm intervals, 1.3 mm x 0.4 mm Teeth at 0.9 mm intervals, 75 mm Handle	MMSU1495S

Swivel Fixation Ring



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
13 mm Diameter	372242

Lacrimal Intubation CIS

Guibor Canaliculus Intubation Set (CIS)

Includes gently curved, stainless steel probes with smooth, rounded tips. Unique stepped probe design enables smooth transition in pass-through intubation.




Qty: 5/Bx
Brand: Visitec®

Description	SKU #
Probes - 17.7 cm length, 0.64 mm diameter. Biocompatible silicone tubing - 29 cm length, 0.64 mm diameter	8590350

Canaliculus Intubation Set (CIS)

Includes malleable stainless steel stylette encased in low friction nylon sheath. The stylette handles assist in controlled manipulation during insertion and retraction. The rounded tip, nylon sheath offers a smooth transition to silicone tubing.



Qty: 5/Bx
Brand: Visitec®

Description	SKU #
Probes - 15.5 cm length, 0.51 mm diameter. Nylon sheath - 0.81 mm diameter. Biocompatible silicone tubing - 33 cm length, 0.64 mm diameter	8590450

IOL Glides

Small Pupil IOL Glide (Pearce)

Facilitates placing posterior chamber IOLs into the capsular bag when the pupil is difficult to dilate. Smooth edges and tapered configuration facilitate atraumatic mechanical dilation of the pupil. May also be used in cases with well dilated pupil.



Qty: 50/Bx
Brand: Visitec®

Description	SKU #
33 mm long x 7 mm wide, Tapered to 3 mm, 0.6 mm thick	581001

IOL Glide (Sheets)

Guides IOL into posterior chamber for placement in the bag. Can also be used for Mini-Nuc technique to guide crystalline lens out of anterior chamber.




Qty: 50/Bx
Brand: Visitec®

Description	SKU #
0.6 mm thick	581033

Small Pupil IOL (Bobbitt)

Guides instruments into anterior chamber. Keeps iris flat to avoid possibility of prolapse or trauma. Acts as a channel to facilitate injection of fluid for hydrodissection.




Qty: 50/Bx
Brand: Visitec®

Description	SKU #
20 mm long x 2.75 mm wide, Tapered to 1.75 mm, 0.6 mm thick	581032

IOL Glide (Hessburg)

Phaco glide portion is designed to protect the iris from prolapse or trauma during phacoemulsification and acts as a guide for insertion of instruments. IOL glide portion assists in the implantation of posterior chamber IOLs.



Qty: 50/Bx
Brand: Visitec®

Description	SKU #
18 mm long x 7 mm wide at widest point, 6.25 mm radius, 0.6 mm thick	581034

Surgical Glide (Fichman)

Phaco glide portion is designed to protect the iris from prolapse or trauma during phacoemulsification and acts as a guide for insertion of instruments. IOL glide portion assists in the implantation of posterior chamber IOLs.



Qty: 50/Bx
Brand: Visitec®

Description	SKU #
30 mm x 5 mm, Phaco glide portion: 10 mm long x 2.6 mm wide, tapered to 1.75 mm. IOL glide portion: 20 mm x 5 mm, 2.5 mm radius. 0.6 mm thick	581049



Small Pupil Expansion

Experience the Latest Innovation in Small Pupil Expansion with the I-Ring® Pupil Expander



- Designed to be used without any additional cornea incision
- Created to avoid any iris distortion and damage during surgery and post surgery
- Designed to have 360° expansion on the iris

I-Ring® Pupil Expander

I-Ring® Pupil Expander


The I-Ring® is designed to expand the iris to increase the surgeon's field of view and helps improve access to the lens while protecting the iris. Single-use sterile product.

- Designed to be gentle on iris and other intraocular tissue
- Created to support iris returning to natural shape post surgery
- Straight-forward insertion and removal

Description	SKU #
I-Ring® Pupil Expander, 6/Pk	587000
I-Ring® Pupil Expander, 1/Pk	587001

Positioning Holes

- Designed to ensure Sinskey hook cannot contact iris



Channels

- Designed to not distort or pinch iris when engaged or removed
- Straight-forward engagement and disengagement

Living Hinges

- Created for enhanced flexibility of ring
- Designed for reliable folding during insertion and removal

Sinskey Hook

Used to manipulate nucleus or intraocular lenses.



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Angled 45 degrees, 10 mm from tip, Hook 0.85 mm tip to bend distance, Overall length 130 mm	581442


Micro Iris Hook

Flexible, gold wire hooks designed to retract the iris during surgical procedures.



Qty: 25/Bx
Brand: Visitec®

Description	SKU #
8 mm long x 0.08 mm diameter	585209



BVI

Fluid Control



PVA Eye Spears

Spear overall handle length is 7.16 cm or 2.82 in.

UltraClean™ LASIK PVA Eye Spears

Compressed shape, designed for fast expanding. Made from highly absorbent, virtually lint free PVA material. Designed specifically for rapid recovery of fluids during LASIK and other corneal procedures.



Qty: 100/Bx
Brand: UltraClean™

Description	SKU #
UltraClean™ Eye Spears, 5/Pouch, 20/Pouch per box	40810
UltraClean™ Eye Spears, 10/Pouch, 10/Pouch per box	40815

UltraClean™ LASIK PVA Expanded Eye Spears

Designed for instantaneous wicking, pre-expanded tip. Made from highly absorbent, virtually lint free PVA material. No need to pre-soak and fan.



Qty: Various/Bx
Brand: UltraClean™

Description	SKU #
UltraClean™ PVA Expanded LASIK Eye Spear, 40/Bx	40800
UltraClean™ LASIK Expanded Eye Spear, 100/Bx	40805

Visitec® LASIK PVA Eye Spears

Made from highly absorbent PVA material. Designed for optimal tissue manipulation and fluid absorption to facilitate flap alignment and adhesion.



Qty: 100/Bx
Brand: Visitec®

Description	SKU #
LASIK Eye Spear - Orange	581761

Visitec® LASIK PVA Expanded Eye Spears

Made from highly absorbent PVA material. Pre-expanded for convenience. Designed to provide instantaneous wicking. No need to pre-soak and fan.



Qty: 40/Bx
Brand: Visitec®

Description	SKU #
LASIK Expanded Spear	581709

Ultracell® PVA Eye Spears

Made from highly absorbent PVA material. Designed for rapid recovery of fluids during cataract and other general ophthalmic procedures.



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
PVA Eye Spear with dark blue handle, 100 count, 5 per pouch, 20 pouches per box	40400-8
PVA Eye Spear with dark blue handle, 200 count, 10 per pouch, 20 pouches per box	40405-8

Cellulose Eye Spears

Ultracell® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Designed for fast wicking action, delivers quick absorption of blood and fluid from the surgical area. The flat tip provides precise application at the surgical incision.



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
Ultracell® Cellulose Eye Spear with light blue handle, 144/Bx	40413
Ultracell® Cellulose Eye Spear with light blue handle, 200/Bx	40415

Visitec® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Designed for fast wicking action, delivers quick absorption of fluid from the surgical area. The flat tip provides precise application at the surgical incision. Magenta handle.



Qty: 100/Bx
Brand: Visitec®

Description	SKU #
Sponge Eye, 7 cm (2.75 in) overall length.	581089

Weck-Cel® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Designed for use in delicate, surgical areas and to maintain rigidity during wicking process. Light blue handle.



Qty: Various/Bx
Brand: Weck-Cel®

Description	SKU #
Weck-Cel® Spears, 6 Pack, 144/Bx	0008680
Weck-Cel® Spears, 10 Pack, 180/Bx	0008685

Sponge Points, Strips and Wicks

Weck-Cel® Sponge Points and Strips

Made from highly absorbent, natural, cellulose material.



Qty: Various/Bx
Brand: Weck-Cel®

Description	SKU #
Weck-Cel® Points, 240/Bx	0008661
Weck-Cel® Points and Strips, 24/Bx	0008660

Visitec® Eye Fluid Wick

Wicks absorb and direct fluid from surgical area into fluid collection bag.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Eye fluid wick, 15 cm Long x 0.4 cm Wide	581147

Corneal Light Shields

Ultracell® Corneal Light Shield

Used to protect hinge and posterior side of corneal flap and protect the retina from the operating microscope light. Made from PVA material.



Qty: 20/Bx
Brand: Ultracell®

Description	SKU #
Ultracell® Corneal Light Shield, 7 mm (0.28 in)	40420

Visitec® Corneal Light Shield

PVA absorbent pad. Soft, porous structure designed to protect the retina from operating microscope light. Contains no cellulose. Upon moistening, shield becomes soft and conforms to shape of cornea. Designed to keep cornea moist.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Visitec® Corneal Light Shield, 8 mm (0.3 in) diameter	581062

Full Body Drapes

Full body drapes feature wider widths and a three-sided channel bag for maximum fluid capacity. Full body drapes can be utilized in all ophthalmic procedures.

Full Body Incise Drape

Qty: 10/Bx

Description	SKU #
Non-Woven Incise Drape with 2 Bags	581156
Overall size:	140 cm x 244 cm (55 in x 96 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Fluid collection bag:	25 cm x 15 cm (10 in x 6 in)
Capacity:	2 x 500 mL
Material Type:	Non-Woven Wetlaid



Full Body Aperture Drape

Qty: 10/Bx

Description	SKU #
Full Body Aperture Drape with 2 Bags	581157
Overall size:	140 cm x 244 cm (55 in x 96 in)
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)
Capacity:	2 x 500 mL
Material Type:	Non-Woven Wetlaid



Full Body Incise Drape with Plastic Head and Channel Bag

Qty: 10/Bx

Description	SKU #
Full Body Incise Drape with Plastic Head and Channel Bag	581163
Overall size:	160 cm x 252 cm (63 in x 99 in)
Incise Area:	13 cm x 15 cm (5 in x 6 in)
Channel Bag:	62 cm x 40 cm (24 in x 16 in)
Capacity:	1000 ml
Material Type:	SMS 157 cm (61 in), Plastic 100 cm (39 in)



*Without wick

Full Body Drapes

Non-Woven Full Body Incise Drape

Qty: 10/Bx

Description		SKU #
Non-Woven Full Body Incise Drape		581164
Overall size:	160 cm x 252 cm (60 in x 99 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Channel bag:	62 cm x 40 cm (24 in x 16 in)	
Capacity:	1000 ml	
Material Type:	Non-Woven Wetlaid	

*Without wick



Half Body Drapes

Plastic Incise Drape with Fluid Collection Bag

Qty: 10/Bx

Description	SKU #
Plastic Incise Drape with Fluid Collection Bag	581149
Overall size:	100 cm x 127 cm (39 in x 50 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)
Capacity:	250 ml
Material Type:	Polyethylene Opaque Film



Plastic Aperture Drape with Fluid Collection Bag

Qty: 10/Bx

Description	SKU #
Plastic Aperture Drape with Fluid Collection Bag	581151
Overall size:	100 cm x 127 cm (39 in x 50 in)
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)
Capacity:	250 ml
Material Type:	Polyethylene Opaque Film



Non-Woven Incise Drape with Fluid Collection Bag

Qty: 10/Bx

Description	SKU #
Non-Woven Incise Drape with Fluid Collection Bag	581150
Overall size:	100 cm x 127 cm (39 in x 50 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)
Capacity:	250 ml
Material Type:	Non-Woven Wetlaid




Half Body Drapes

Non-Woven Aperture Drape with Fluid Collection Bag

Qty: 10/Bx


Description	SKU #
Non-Woven Aperture Drape with Fluid Collection Bag	581152
Overall size:	100 cm x 127 cm (39 in x 50 in)
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)
Capacity:	250 ml
Material Type:	Non-Woven Wetlaid



Non-Woven Incise Drape with Channel Bag

Qty: 10/Bx

Description	SKU #
Non-Woven Incise Drape with Channel Bag	581153
Overall size:	140 cm x 127 cm (55 in x 50 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Fluid collection bag with channels:	26 cm x 35 cm (10 in x 14 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid




Half Body Drapes

Non-Woven Aperture Drape with Channel Bag

Qty: 10/Bx

Description		SKU #
Non-Woven Aperture Drape with Channel Bag		581155
Overall size:	140 cm x 127 cm (55 in x 50 in)	
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)	
Fluid collection bag with channels:	26 cm x 35 cm (10 in x 14 in)	
Capacity:	500 ml	
Material Type:	Non-Woven Wetlaid	




Non-Woven Half Body Incise Drape* with Channel Bag

Qty: 10/Bx

Description		SKU #
Non-Woven Half Body Incise Drape* with Channel Bag		581166
Overall size:	160 cm x 180 cm (63 in x 71 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Channel bag:	62 cm x 40 cm (24 in x 16 in)	
Capacity:	1000 ml	
Material Type:	SMS	

*Without wick



Face Drapes

Ideal for refractive or office procedures.


Mini Pre-Split Incise Drape

Adhesive backing.

Qty: 10/Bx

Description	SKU #
Mini Pre-Split Incise Drape	581154
Overall size:	40 cm x 40 cm (16 in x 16 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Pre-split:	5 cm (2 in)
Material Type:	SMS

*Without wick



Mini Plastic Incise Drape

Adhesive backing.

Qty: 10/Bx

Description	SKU #
Mini Plastic Incise Drape	581159
Overall size:	40 cm x 40 cm (16 in x 16 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Material Type:	Polyethylene Opaque Film

*Without wick



Mini Plastic Aperture Drape

Adhesive backing.

Qty: 10/Bx

Description	SKU #
Mini Plastic Aperture Drape	581160
Overall size:	40 cm x 40 cm (16 in x 16 in)
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)
Material Type:	Polyethylene Opaque Film

*Without wick




Face Drapes

Mini Plastic Aperture Drape

Small incise provides complete coverage of lashes, lids and brows. Converts aperture to incise drape.

Qty: 10/Bx

Description		SKU #
Incise Drape		581161
Overall size:	10 cm x 20 cm (4 in x 8 in)	
Adhesive area:	11 cm x 16 cm (4 in x 6 in)	
Material Type:	Polyethylene Opaque Film	



Dual Refractive Drape

Qty: 10/Bx

Description		SKU #
Dual Refractive Drape		581162
Overall size:	40 cm x 26 cm (16 in x 10 in)	
Aperture area:	2.5 cm x 4 cm (1 in x 1.5 in)	
Material Type:	Polyethylene Opaque Film	



*Without wick

Eye Drains

Visitec® LASIK Eye Drain (Chayet)

Designed to reduce stromal bed hydration for a consistent ablation and for nasally hinged flaps.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
LASIK Eye Drain, 10 mm	581098

Fluid Collection

Visitec® Fluid Collection Bags

For procedures where collection and control of fluid is designed to keep surgeon, patient and floor dry. Reduces operating room cleanup and turnaround time. Bags have an adhesive backing to adhere to surgical drape.



581177



581148

Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Fluid Collection Bag: 17 cm x 16 cm (6.75 in x 6 1/4 in), 250 ml capacity	581177
Fluid Collection Bag with Wick Bag: 17 cm x 16 cm (6.75 in x 6.25 in), 250 ml capacity Wick: Includes one, 0.4 cm x 15 cm (0.15625 in x 5.875 in)	581146
Mini Fluid Collection with Wick Bag: 9 cm x 16 cm (3.5 in x 6.25 in), 100 ml capacity Wick: Includes one, 0.4 cm x 15 cm (0.15625 in x 5.875 in)	581148

Instrument Wipes

Visitec® Instrument Wipes



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Instrument Wipe 76.2 mm x 76.2 mm (3 in x 3 in), 2 mm thick	90900

Visitec® Visiwipe Instrument Wipes



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Instrument Wipe 7.3 cm x 7.3 cm (2.875 in x 2.875 in)	581047
Adhesive Wipe 7.6 cm x 4.5 cm (3 in x 1.75 in)	581076

Other Fluid Control

Absorbent Stick

6.6 cm overall length x 5 mm diameter (2.59 in x 0.1875 in).



Qty: 250/Bx
Brand: Visitec®

Description	SKU #
Absorbent Stick 6.3 cm, 10 per pouch, 25 pouches per box.	581060

Visiswab Cotton Swabs

Pointed end may be used as a second instrument to absorb blood and fluid, move sutures and manipulate tissue. Rounded end may be used for scleral depression. Tightly wound ends create less linting than other cotton swabs. 7.6 cm (3 in) overall length. 10 per pouch, 25 pouches per box.



Qty: 250/Bx
Brand: Visitec®

Description	SKU #
Visiswab Cotton 7.6 cm (3 in) length	581069

Non-Ophthalmic

Nasal Packs

Ultracell® brand PVA sponges for postoperative packing following nasal surgery.

- Used following septoplasty, rhinoplasty, and turbinectomy procedures
- May aid in holding septum firmly in place
- Available in a range of sizes: 3 cm pediatric to 10 cm
- Select sizes available in an anatomical shape
- Popular sizes also available in our patented layered design
- Excellent alternative to cumbersome gauze packing
- Latex-free
- Single-use



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
Classic Nasal Pack, 4.5 cm with String, 15 mm Thick, 20/Bx	10100-C
Classic Nasal Pack, 8 cm with String, 20/Bx	10110-C
Classic Nasal Pack, 10 cm with String 15 mm Thick, 20/Bx	10115-C
Anatomical Nasal Pack, 8 cm with String, 15 mm Thick, 10/Bx	10120-C

Epistaxis Products

The Ultracell® brand of epistaxis sponge provides highly absorbent treatment of nosebleeds.

- Fluid absorption causes the packing to expand producing a tamponading effect
- Available in a variety of sizes: 3 cm anterior to 10 cm posterior
- Ideal for use in emergency rooms, clinics, and doctor's offices
- Latex-free
- Single-use



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
Epistaxis Pack, 10 cm with String, 10/Bx	10140-C
Epistaxis Slimpack, 8 cm with String, 20/Bx	10146-C
Epistaxis Pack, 5.5 cm with String, 20/Bx	10150-C
Anterior Epistaxis Pack, 3 cm with String, 20/Bx	10160-C

Non-Ophthalmic

Sinus Packs

A complete line of Ultracell® brand sinus packs specifically designed for use following sinus surgeries.

- Ideal packing following FESS (Functional Endoscopic Sinus Surgery)
- Used for post-operative packing or dressing of sinus cavity
- Available in a range of lengths and thicknesses
- Patented layered design allows for easy manipulation
- Latex-free
- Single-use



Qty: 20/Bx
Brand: Ultracell®

Description	SKU #
Sinus Ultra Slimpack with String 2.5 cm x 0.6 cm x 2 cm	20202-C
Sinus Pack with String 3.5 cm x 1.2 cm x 1.2 cm	20210-C
Sinus Slimpack with String 3.5 cm x 0.9 cm x 1.2 cm	20211-C
Sinus Ultra Slimpack with String 3.5 cm x 0.6 cm x 1.2 cm	20212-C

Ear Wicks

The Ultracell® brand provides both ear wicks and ear packs. Much smaller, compressed form allows straight-forward insertion into swollen ear canal. Pediatric sizes available.

Ear Wicks

- Perfect for use with treatments for otitis externa
- Help maintain antibiotics in the ear canal
- Also available fenestrated
- Latex-free
- Single-use



Qty: Various/Bx
Brand: Ultracell®


Description	SKU #
Ear Wick 9 mm Diameter x 15 mm, 50/Bx	30300-C
Pediatric Ear Wick 7 mm Diameter x 12 mm, 50/Bx	30301-C
Fenestrated Ear Wick 9 mm Diameter x 15 mm, 50/Bx	30305-C
Fenestrated Pediatric Ear Wick 7 mm Diameter x 12 mm, 50/Bx	30306-C
Long Ear Wick 9 mm Diameter x 24 mm, 10/Bx	30310-C
Pediatric Long Ear Wick 7 mm Diameter x 20 mm, 10/Bx	30311-C
Pediatric Long Fenestrated Ear Wick 7 mm Diameter x 20 mm, 10/Bx	30316-C



Quality Accessories to Support You and Your Patients Needs



Eye Shields



Pen Lights

Eye Shields and Bandages

Eye Shield (Left/Right)

White lightweight plastic eye shield with ventilation holes. Flexible eye shield, shaped to fit the right or left eye for protection. Smooth edges do not require a cover. Packaged sterile.



Qty: 50/Bx
Brand: Visitec®

Description	SKU #
Eye Shield, Left/Right 8.5 cm (3.3 in)	581019

Universal Pediatric Eye Shield

Clear lightweight plastic and features smooth edges and ventilation holes. For use on left or right eye. Packaged sterile.



Qty: 25/Bx
Brand: Visitec®

Description	SKU #
Pediatric Eye Shield	581708

Universal Eye Shield

Clear lightweight plastic, for use on left or right eye, that features smooth edges and ventilation holes. Packaged sterile.



Qty: 50/Bx
Brand: Visitec®

Description	SKU #
Universal Eye Shield, 7.5 cm x 6.5 cm (3 in x 2.5 in)	581040

Visitec® Guibor Eye Bandage

Disposable self-adhesive eye dressing designed to protect the eye during drying conditions, such as "dry eye". Straightforward to apply. Packaged sterile.



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Eye Bandage, Guibor bubble bandage	8590300

Diagnostic Lights

Doctor's Six-Pack Diagnostic Penlights

Long-lasting disposable diagnostic penlight features convenient clip-activated switch.



Qty: 6/Bx
Brand: BVI®

Description	SKU #
Diagnostic Penlight, Doctors Six-Pack	8310060

Maxi-Lite® Diagnostic Penlight

Long-lasting disposable diagnostic penlight features twist-on backswitch.



Qty: 1/Ea
Brand: BVI®

Description	SKU #
Diagnostic Penlight, Maxi-Lite®	8310201

Nurse-Mate® Six-Pack Diagnostic Penlights

Long-lasting disposable diagnostic penlight features convenient clip-activated switch and 2 mm to 9 mm pupil gauge. Imprinted with medical caduceus.



Qty: 6/Ea
Brand: Nurse-Mate®

Description	SKU #
Diagnostic Penlight, Nurse-Mate® Six-Pack	8310080

Blu-Spot® Diagnostic Light

Has diffuse cobalt light beam and twist-on backswitch.



Qty: 1/Bx
Brand: Blu-Spot®

Description	SKU #
Diagnostic Penlight	8320150

General Accessories

Superflex Extension Set



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Superflex Extension Set	581007

Irrigation Handle

Standard male luer fitting.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Length: 9.3 cm (3.5 in)	581012

Corneal Marking Pad

Elongated, flat pad allows easy transfer of ink onto markers.



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Gentian violet, Pre-inked, Ready to use	581057

Visitec® Bulk High Flow Filter

Filters out particulates and microbiological contamination of 0.2 micron and larger.



Qty: 50/Bx
Brand: Visitec®

Description	SKU #
0.2 micron Filter	581030

Visitec® Double Male Luer Adapter

Standard male luer fittings on each end for connecting phacoemulsification and I/A machine tubing to instruments.



Qty: 50/Bx
Brand: Visitec®

Description	SKU #
3.78 cm (1-1/2 in) Long	581056



Experience the Evolution

Reimagine Your Operating Room



Custom Packs

CustomEyes®



Knowledgeable Support

- Experienced ophthalmic sales team
- We understand your operating room
- Information at the fingertips



Reliable Quality

- Top-notch quality
- Access to BVI's wide range of single-use instruments



Personalized Packs

- Individualization to each of your surgeon's needs
- Quickly adapt to staff or technique changes




Simplified Surgical Prep

- 3D visualization of pack layers
- Pack assembly to match your operations



Easy Order Process

- Real-time quote and pricing
- Sample in 2 weeks



Custom Packs

CustomEyes® Configurator

The CustomEyes Configurator allows our team of ophthalmic specialists to easily create the right pack for you.



Simple Product Selection

Allows you to choose the product characteristics and select the best product combination for your needs.



Visual Pack Assembly

Optimize the layout of the pack by seeing a virtual 3D layering including order and orientation of each component.



Instant Feedback

Start building a pack together with your rep to get instant feedback and customize quickly.



Actualize

Receive a quote, order a sample and initiate production directly from the app.



CustomEyes®

Custom Procedure Packs

CUSTOM PACKS



Knowledgeable Support



Reliable Quality



Personalized Packs



Simplified Surgical Prep



Easy Order Process



Experience the Evolution



Electrosurgery and Electrocautery



Cautery Devices

Accu-Temp® Cauteries provide a broad selection of tip styles for cauterization in a wide variety of ophthalmic procedures. All products are single-use and provided sterile.

Accu-Temp® High Temperature Cautery

Qty: 10/Bx
Brand: Accu-Temp®



8441500

Description	SKU #
1,150°C (2,102°F), 1/2 in (12.7 mm) Shaft with elongated vascator fine tip	8441500



8442000

Description	SKU #
1,150°C (2,102°F), 1/2 in (12.7 mm) Shaft with fine tip	8442000



8443000

Description	SKU #
1,150°C (2,102°F), 1/2 in (12.7 mm) Shaft with loop tip	8443000



8444000

Description	SKU #
1,150°C (2,102°F), 2 in (50.8 mm) Shaft with loop tip	8444000



8445000

Description	SKU #
1,150°C (2,102°F), 2 in (50.8 mm) Shaft with fine tip	8445000

Accu-Temp® Low Temperature Cautery

Qty: 10/Bx
Brand: Accu-Temp®



8441000

Description	SKU #
612°C (1,134°F), 1/2 in (12.7 mm) Shaft with fine tip	8441000

Cautery Devices

Accu-Temp® Variable Temperature Cautery - Dialed Controlled

Qty: 10/Bx
Brand: Accu-Temp®



Low Temperature 8440000

Description	SKU #
150° - 612°C (302° - 1,134°F), 1/2 in (12.7 mm) Shaft with fine tip	8440000



High Temperature 8446000

Description	SKU #
506° - 1,150°C (943° - 2,102°F), 1/2 in (12.7 mm) Shaft with fine tip	8446000

Accu-Temp® Reusable Handles

The reusable cautery handles are used with the single-use tips.

Qty: 10/Bx
Brand: Accu-Temp®



8461000

Description	SKU #
Reusable Variable Low Temp Handle, 150° - 612°C (302° - 1,134°F)	8461000



8461020

Description	SKU #
Reusable High Temp Handle, 1,150°C (2,102°F)	8461020

Accu-Temp® Single-Use Tips

Each single-use tip is packaged sterile and includes a protective plastic sheath that is designed to protect the reusable handles from contamination during use. Sterile and single-use.

Single-Use Fine Tip



Qty: 25/Bx
Brand: Accu-Temp®

Description	SKU #
Single-use fine tip	8461050

Single-Use Loop Tip



Qty: 25/Bx
Brand: Accu-Temp®

Description	SKU #
Single-use loop tip	8461060

Eraser

Wet-Field® Erasers intended for coagulation in vitreoretinal procedures, glaucoma filtration and cataract incision.
Sterile and single-use.

Wet-Field® Hemostatic Eraser Bipolar Bevel Tip

Qty: 10/Bx
Brand: Wet-Field®



221250

Description	SKU #
18G Bevel tip, Straight	221250

Wet-Field® Hemostatic Eraser Bipolar Blunt Tip

10 per box unless otherwise noted.

Qty: 10/Bx
Brand: Wet-Field®



221251



221260135



221260150



221260

Description	SKU #
18G Blunt tip, Angled 45 degrees	221251
20G Blunt tip, Angled 135 degrees	221260135
20G Blunt tip, Angled 150 degrees	221260150
20G Blunt tip, Straight	221260

Wet-Field® Hemostatic Eraser Bipolar Fine Tip

Qty: 5/Bx
Brand: Wet-Field®



221268



221268135



221268150

Description	SKU #
23G Fine tip, Straight	221268
23G Fine tip, Angled 135 degrees	221268135
23G Fine tip, Angled 150 degrees	221268150
23G Fine tip, Angled 135 degrees	221269135
25G Fine tip, Straight	221267
25G Fine tip, Angled 135 degrees	221267135
25G Fine tip, Angled 150 degrees	221267150

Erasers

Wet-Field® Hemostatic Eraser Bipolar 20G - 23G Tapered Tip

Qty: 5/Bx
Brand: Wet-Field®



Description	SKU #
20 - 23G Tapered fine tip, Straight	221265



Description	SKU #
20 - 23G Tapered blunt tip, Straight	221266

Wet-Field® Osher Eraser

Qty: 5/Bx
Brand: Wet-Field®



Description	SKU #
25G, 115 Degrees	221270

Cables

Qty: 10/Bx
Brand: Wet-Field®



22135010



22155010

Description	SKU #
Wet-Field® Disposable Cable	22135010
Wet-Field® Disposable Cable, Shrouded	22155010

Forceps

Wet-Field® Bipolar Forceps

Bipolar forceps are designed to deliver coagulation with no electrical arcing between the polished instrument tips for minimum coagulum adhesion. The insulated forcep blades are designed to reduce the possibility of electrical shock during operation. Non-sterile and reusable.


Qty: 1/Bx
Brand: Wet-Field®



Description	SKU #
Standard Forceps	221204



Description	SKU #
Straight Iris (Jeweler) Forceps	221206



Description	SKU #
Curved Iris (Jeweler) Forceps	221207



Description	SKU #
Adson Forceps	221208



Description	SKU #
Coaptation Forceps	221209



Description	SKU #
Tenzel Forceps	221211



Providing Relief for Patients with Dry Eye or Lacrimal Disorders


Parasol® | Micro Flow™ Punctal Occlusion



Permanent
Parasol®



Long-Term
EXTEND 180®



Mid-Term
EXTEND®



Partial
Micro Flow™



Short-Term
Collagen

Permanent Occlusion

Dry Eye products that offer predictable retention¹, straight-forward sizing and insertion. All occluder products are made of medical grade silicone and are reversible.

Parasol® Pre-Loaded Punctal Occluders Sterile

Occluders are packaged in a convenient one pair tray, are sterile and pre-loaded on a single-use inserter.



A14-201



Qty: 2/Bx

Parasol® Pre-Loaded Punctal Occluders Sterile	SKU #
X-Small, 0.25 - 0.35 mm	A14-200
Small, 0.35 - 0.65 mm	A14-201
Medium, 0.6 - 0.85 mm	A14-202
Large, 0.9 mm+	A14-203

Parasol® Quickload Punctal Occluders Sterile

Occluders are sterile and pre-mounted on a stand for Quickload convenience. Reusable inserters are sold separately.



A14-201QL

Qty: 16/Bx

Parasol® Quickload Punctal Occluders Sterile	SKU #
Small, 0.35 - 0.65 mm	A14-201QL
Medium, 0.6 - 0.85 mm	A14-202QL
Large, 0.9 mm+	A14-203QL
Combo small and medium	A14-205QL

Parasol® Bulk Punctal Occluders Non-Sterile

Occluders are non-sterile and packaged two occluders per poly bag. Reusable inserters are sold separately.



Qty: 20/Bx

Parasol® Bulk Punctal Occluders Non-Sterile	SKU #
Small, 0.35 - 0.65 mm	A14-201BN
Medium, 0.6 - 0.85 mm	A14-202BN
Large, 0.9 mm+	A14-203BN
Combo small and medium	A14-209BN

¹ C. McCabe, MD, PhD, FACS Plugs Reduce Dry-Eye Symptoms, Improve Vision. November 2009, Review of Ophthalmology.

Permanent Occlusion

Plus Pre-Loaded Punctal Occluders Sterile

The Plus design is the result of a micro-molding manufacturing process using a soft medical grade silicone.

- Sharply angled nose for insertion
- Ultra-smooth surface finish
- Delicate contoured edges
- Designed for consistent retention

Occluders are packaged in a convenient one pair tray. Each occluder is sterile and pre-loaded on a high quality disposable inserter.



A14-204



Qty: 2/Bx

Plus Pre-Loaded Punctal Occluders Sterile	SKU #
Small, 0.4 - 0.55 mm	A14-204
Medium, 0.55 - 0.7 mm	A14-205
Large, 0.7 - 0.85 mm	A14-206
X-Large, 0.85 mm+	A14-207

Plus Punctal Occluders Non-Sterile Bulk

Occluders are supplied non-sterile and have two occluders per poly bag, with 10 bags/pairs per box, unless otherwise noted. Reusable inserters are sold separately.

Qty: 20/Bx



Plus Punctal Occluders Non-Sterile Bulk	SKU #
Small, 0.4 - 0.55 mm	A14-204BN
Medium, 0.55 - 0.7 mm	A14-205BN
Large, 0.7 - 0.85 mm	A14-206BN
X-Large, 0.85 mm+	A14-207BN

Mid-Term Occlusion

Extend® and Extend® 180 Synthetic Punctal Occluders

Provides occlusion for up to 3 months with Extend® or 6 months with Extend® 180, for patients experiencing short-term dry eye symptoms. These are ideal for the transient dry eye associated with LASIK or other ocular surgeries. Can also be used for allergy sufferers, contact lens discomfort, climate change or to sustain drug contact time and effectiveness.

Qty: 20/Bx



A13-204

Extend® Synthetic Punctal Occluders	SKU #
Variety Pack	A13-200
0.2 mm x 2 mm	A13-202
0.3 mm x 2 mm	A13-203
0.4 mm x 2 mm	A13-204
0.5 mm x 2 mm	A13-205

Qty: 20/Bx



A13-2044

Extend® 180 Synthetic Punctal Occluders	SKU #
Variety Pack	A13-2000
0.3 mm x 2 mm	A13-2033
0.3 mm x 2 mm, Single Pair	A13-2033SP
0.4 mm x 2 mm	A13-2044
0.4 mm x 2 mm, Single Pair	A13-2044SP
0.5 mm x 2 mm	A13-2055
0.5 mm x 2 mm, Single Pair	A13-2055SP

Partial Occlusion

Micro Flow™ Punctal Occluders

Designed to provide a straight-forward, reversible means of partial occlusion. Micro Flow™ Punctal Occluders support lacrimal drainage care. Suggest to use when total occlusion is not ideal.

Qty: 2/Bx



A14-401

Micro Flow™ Punctal Occluders Sterile	SKU #
Sterile, Small, 0.4 - 0.55 mm	A14-401
Sterile, Medium, 0.55 - 0.7 mm	A14-402
Sterile, Large, 0.7 - 0.85 mm	A14-403

Qty: 10/Bx

Micro Flow™ Punctal Occluders Non-sterile	SKU #
Non-Sterile, Bulk, Small, 0.4 - 0.55 mm	A20-401BN
Non-Sterile, Bulk, Medium, 0.55 - 0.7 mm	A20-402BN
Non-Sterile, Bulk, Large, 0.7 - 0.85 mm	A20-403BN

Temporary Occlusion

Collagen Synthetic Punctal Occluders

These temporary occluders provide occlusion for 7-10 days. They are ideal for diagnostic use for determining if permanent punctal occlusion is needed for your patient. Six occluders in each foam sleeve, with 10 foam sleeves per box.



A12-103

Qty: 60/Bx

Collagen Synthetic Punctal Occluders	SKU #
0.2 mm x 2 mm	A12-102
0.3 mm x 2 mm	A12-103
0.4 mm x 2 mm	A12-104
Variety Pack	A12-200

Single-Use Punctal Dilators and Probes

Visitec® Punctal Dilator and Probe, Khan

Dual function instrument allows for single handed transition from dilating to probing without release of lateral lid traction.



Qty: 3/Bx

Description	SKU #
25 mm (1 in) medium taper, Blunt, Wilder style dilator mated with a size 0 Bowman lacrimal probe, 0.7 mm x 45 mm (22G x 1.75 in)	585124

Other Related Occlusion Options

Lacrimal Cannula



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 1/2 in, 0.5 x 13 mm	585068

Lacrimal Cannula



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 26G x 1-1/8 in, 0.45 x 28 mm	581276

Punctum Dilator, Fine/Medium

0.4 mm x 21.5 mm Tip,
0.6 mm x 18 mm Tip,
77 mm Handle




Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Punctum Dilator, Fine/Medium	MMSU1181S

Punctum Dilator/Seeker

0.45 mm x 21.5 mm Tip,
5 mm Pointed tip,
18 mm Straight shaft, 77 mm Handle



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Punctum Dilator/Seeker	MMSU1193S

Curved Tying Forceps

0.25 mm Tips, 5 mm Curved tying platforms,
80 mm Handle



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Curved Tying Forceps	MMSU1210S



Cataract Surgery Streamlined

High contrast and high quality¹ vision.
Seamless, predictable, and efficient
preloaded IOL delivery.²



*THIS PRODUCT IS FDA APPROVED AND ONLY SOLD IN THE US.


References

¹US Patent NO: US8647383. ²Data on file, BVI, 2019.

I PURE
ASPHERIC
MONOFOCAL OPTIC

bvimedical.com

Trusted aspheric monofocal IOL in a fully pre-loaded, single-use IOL injector system.



Model Name	B1PC (Clear) B1PY (Yellow)	B3PC (Clear)
Specification	UV blocking	UV blocking
Optic material	Hydrophobic acrylic	Hydrophobic acrylic
Optic design	Aspheric -0.18 microns patented lens design	Aspheric -0.18 microns patented lens design
Manufacturing	Lathe-cut and tumble polished	Lathe-cut and tumble polished
Haptic material	Hydrophobic acrylic with blue PMMA chemically bonded haptic tips	Blue PMMA chemically bonded
Haptic configuration	Modified C-loop, 5 Degree angulation	Modified C-loop, 5 Degree angulation
Dimension (Optic/OAL)	6 mm/12.5 mm	6 mm/12.5 mm
Power	+6 to +30 D (in 0.5 D increments)	+6 to +30 D (in 0.5 D increments)
Nominal A-constant*	118.4	118.4
Optimized constants**	Haigis a0 = -0.542 a1 = 0.161 a2 = 0.204 Hoffer Q pACD = 5.3 Holladay 1 sf = 1.52 SRK/T A = 118.5	Haigis a0 = -0.093 a1 = -0.023 a2 = 0.208 Hoffer Q (pACD) = 5.30 Holladay 1 sf = 1.54 SRK/T A = 118.6
Front injector tip outer diameter	1.82 mm	1.89 mm
Injector type	Single-use full pre-loaded IOL system	Single-use full pre-loaded IOL system

* The A-constant mentioned above is presented as a guideline only for lens power calculations. It is recommended that the A-constant measurement be customized based on the surgeon's experience and measuring equipment.

** <https://iolcon.org/lensesTable.php> (Accessed data Mar 15, 2019)

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

Some of the products and/or specific features as well as the procedures featured in this document may not be approved in your country and thus may not be available there. Design and specifications are subject to change without prior notice as a result of ongoing technical development. Please contact our regional representative regarding individual availability in your respective market. HOYA is a trademark of HOYA Corporation. All other HOYA Surgical Optics names and logos are owned by HOYA Surgical Optics or its affiliates. ©2019 HOYA Medical Singapore Pte. Ltd. All rights reserved.



Most Versatile Glaucoma Treatment in Ophthalmology

Endoscopic CycloPhotocoagulation (ECP)

- Mild to moderate glaucoma
- Combine with phaco and other MIGS
- Advanced glaucoma



Endo Optiks®

Laser Endoscopy System

The E2 and E4 endoscopy systems provide high resolution video imaging and wide field illumination. All endoscopes and probes are for use with these systems.

E2 Ophthalmic Laser Endoscopy System

The E2 system uniquely combines endoscopic imaging and laser treatment capability. The console houses a high resolution video camera, 175 or 300 watt xenon light source and an 810nm diode laser. It offers versatility like no other ophthalmic surgical instrument with targeted, controlled laser delivery for the treatment of glaucoma and retinal disorders and valuable diagnostic imaging.



Description	SKU #
E2 Compact Microprobe Laser and Endoscopy System, 1/Ea	OME2000

E4 Ophthalmic Endoscopy System

The E4 high resolution endoscopy system can be coupled with various external lasers for triple function laser endoscopy.



Description	SKU #
E4 Ophthalmic Laser Endoscopy System, 1/Ea	OME4000

Microendoscopes and Probes



Triple Function Ophthalmic Laser Endoscopes

SKU #	G	Cannula	High Resolution	Laser	Camera	Light	Custom Laser Connection	E2 Compatibility	E4 Compatibility
OME200SMAVA	19	Curved	✓	✓	✓	✓		✓	
OME200SMAVXA	19	Curved	✓	✓	✓	✓	✓	✓	✓
OME200SMAHRA	19	Straight	✓	✓	✓	✓		✓	
OME200SMAHRVA	19	Curved	✓	✓	✓	✓		✓	
OME200SMAHRXA	19	Straight	✓	✓	✓	✓	✓	✓	✓
OME200SMAHRXVA	19	Curved	✓	✓	✓	✓	✓	✓	✓
OME200SMAEA	20	Straight		✓	✓	✓		✓	
OME200SMAXA	20	Straight		✓	✓	✓	✓	✓	✓
OME230SMAHRXA	23	Straight	✓	✓	✓	✓	✓	✓	✓

Endo Optiks®

Dual Function Ophthalmic Endoscopes

SKU #	G	Cannula	High Resolution	Laser	Camera	Light	Custom Laser Connection	E2 Compatibility	E4 Compatibility
OME200HRA	20	Straight	✓		✓	✓		✓	✓
OME230HRA	23	Straight	✓		✓	✓		✓	✓

Endo Optiks® Accessories

Diode Laser Safety Goggles

810nm diode laser safety goggles for use with/without eyeglasses.



Description	SKU #
Diode laser safety goggles, 1/Ea	LG001

Endoscope Tip Protector

Polycarbonate protector tube for 23G endoscope cannulas (not autoclavable).



Description	SKU #
Tip protector for 23G endoscope, 1/Ea	NC001-A

Endoscope Connector Caps

Protective caps for endoscope light, camera and laser connectors.



Description	SKU #
Endoscope connector caps, 50/Ea Video/light and 10/Ea Laser, 1/Ea	EO-CAP-S

Flight Case

Custom foam padded flight cases with wheels and telescoping handles.



FC2720



FC401655

Description	SKU #
Flight case with built-in wheels and telescoping handle for E4 Unit, 1/Ea	FC2720
Flight case with built-in wheels and telescoping handle for E2 Unit, 1/Ea	FC401655

Endo Optiks® Accessories

Footswitches

Non-slip laser energy control for E2 system. Optional illumination control via side paddles compatible with E2 and E4.



Description	SKU #
Footswitch with illumination and laser control, 1/Ea	SW002
Footswitch for laser only, 1/Ea	SW003

Monitors

7 x 11 in Single layer tray with rubber fingered mat for all endoscopes and probes.



Description	SKU #
Operating room system cart with integrated monitor arm	CT003
26 in Medical grade monitor enhanced brightness and contrast (Order Ct002-FPA with this monitor - no stand included), 1/Ea	MON006
17 in Medical grade high resolution LCD flat screen, 1/Ea	MON007

Video Adapters

These couplers offer optimum zoom and manual focus of the endoscopic image. 810nm laser filters installed in all models with 532nm options available.



Description	SKU #
17K Pixel video adapter (for use with 17K high resolution endoscopes), for use with E2 Ophthalmic Laser Endoscopy System, 1/Ea	OME170Z
17K Pixel video adapter (for use with 17K high res. endoscopes), for use with E4 Ophthalmic Laser Endoscopy System, 1/Ea	OME170ZMG
10K Video adapter with 532Nm and 810Nm filter (for use with 10K endoscopes), 1/Ea	OME300ZMG

Sterilization Tray


Medical grade monitors offered specifically to achieve optimum quality of the endoscopic image and compatibility with E2, E4 and mounting arm.

Description	SKU #
Sterilization tray	VP6250




Addressing Control of Light Beam Direction

29G Spotlight Directional Chandelier



Enabling Bimanual Surgery With:

- Directional control of light-beam for enhanced visualization
- 29G valved entry system for straightforward light fiber entrance, fixation and removal
- Stationary and wide field endo-illumination



Patent Pending

Vitreo

Vitreoretinal Surgical Products

bvimedical.com

29G Spotlight Directional Chandelier


29G Spotlight Directional Chandelier

Sets new standards in light beam direction control, it delivers stationary, diffuse wide-field endo-illumination enabling bimanual surgery.



Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
29G Spotlight Directional Chandelier, Sterile	LF29.D02





Improved Intraocular Surgical Control and Precision in Surgical Maneuvers

Micro Forceps and Scissors



Tip Design

- Shortened tip provides stronger grasping force¹

Shaft Design

- Rigid shaft designed to reduce flexibility for better control
- 25% less flexible than competitors while maintaining 32 mm working shaft length¹



Ergonomically Designed Handle

- Synchronized handle squeeze closes the tip to facilitate grasping and cutting
- Round handle allows for tip alignment

¹ Data on file

Vitreco


Vitreoretinal Surgical Products

Micro Forceps

Gauge Size with Color Coding


Single-Use

All Single-Use Sold in 5/Bx



Reusable

All Reusable Sold in 1/Bx



Serrated Gripping

Serrated Gripping Micro Forceps are designed for grasping strong adherent PVR membranes.



Brand: Vitreq®

Description	SKU #
● 23G Single-Use	FC23.D02
● 25G Single-Use	FC25.D02
● 23G Reusable	FC23.R02
● 25G Reusable	FC25.R02

ILM End-gripping

ILM end-gripping forceps have thin tips for ILM and ERM peeling. Designed to improve visualization for grasping thin membranes.



Brand: Vitreq®

Description	SKU #
● 23G Single-Use	FC23.D03
● 25G Single-Use	FC25.D03
● 27G Single-Use	FC27.D03
● 23G Reusable	FC23.R03
● 25G Reusable	FC25.R03
● 27G Reusable	FC27.R03

Asymmetrical Gripping

Asymmetrical tips allow visualization at grasping area and interior end-gripping platforms are designed for assured grasping. The tip design targets ILM, fine epiretinal membranes and diabetic PVR membranes.



Brand: Vitreq®

Description	SKU #
● 23G Single-Use	FC23.D04
● 25G Single-Use	FC25.D04
● 23G Reusable	FC23.R04
● 25G Reusable	FC25.R04

Micro Forceps

Eckardt Endgripping

Eckardt forceps with wider tip are designed to enhance surface grasping and reduce tearing effect. Shorter tip branches designed to enhance grasping and holding power. Designed for ILM, ERM peeling and general peeling of fine membranes which needs end-gripping action.



FC27.D05

FC23.R05

Brand: Vitreq®

Description	SKU #
● 23G Single-Use	FC23.D05
● 25G Single-Use	FC25.D05
● 27G Single-Use	FC27.D05
● 23G Reusable	FC23.R05
● 25G Reusable	FC25.R05
● 27G Reusable	FC27.R05

Shah X-tra Grip*

Fine tipped forceps with serrated jaws platform for maximum grasping of the tissue. Effective for all types of tissue manipulation, ILM, ERM, PVR and diabetic tissue.



FC23.D07

FC23.R07

Brand: Vitreq®

Description	SKU #
● 23G Single-Use	FC23.D07
● 25G Single-Use	FC25.D07
● 27G Single-Use	FC27.D07
● 23G Reusable	FC23.R07
● 25G Reusable	FC25.R07
● 27G Reusable	FC27.R07

Kuhn 45° Power Gripping*

Kuhn 45 Degree Power Gripping forceps are designed for peeling proliferative intravitreal, epi and sub-retinal membranes, lens capsule, iris, small IOFB, fibrinous membranes, as well as ILM and ERM membranes and posterior hyaloids. With a 45 degree crocodile designed inner surface, these instruments are designed to enable selective, strong grasping, and visibility of the target tissue.



FC23.D06

FC25.R06

Brand: Vitreq®

Description	SKU #
● 23G Single-Use	FC23.D06
● 25G Single-Use	FC25.D06
● 23G Reusable	FC23.R06
● 25G Reusable	FC25.R06

Micro Scissors

Vertical

Sharp-tipped vertical cutting scissors for cutting perpendicular to the retina.



SC23.D01


SC23.R01

Brand: Vitreq®

Description	SKU #
● 23G Single-Use	SC23.D01
● 25G Single-Use	SC25.D01
● 27G Single-Use	SC27.D01
● 23G Reusable	SC23.R01

Straight

Sharp-tipped straight cutting scissors for cutting in the iris plane and cutting peripheral membranes.



SC23.D02

SC23.R02

Brand: Vitreq®

Description	SKU #
● 23G Single-Use	SC23.D02
● 25G Single-Use	SC25.D02
● 23G Reusable	SC23.R02
● 25G Reusable	SC25.R02

Horizontal Curved

Sharp-tipped horizontal curved cutting scissors for cutting membranes adherent to the retinal surface. The tip's curvature of 12mm follows the contour of the retina.



SC23.D03

SC23.R03

Brand: Vitreq®

Description	SKU #
● 23G Single-Use	SC23.D03
● 25G Single-Use	SC25.D03
● 27G Single-Use	SC27.D03
● 23G Reusable	SC23.R03
● 25G Reusable	SC25.R03
● 27G Reusable	SC27.R03

Reusable Cleaning Adapter

Unique cleaning/sterilization adapter included with every reusable instrument option.




Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Reusable Cleaning Adapter	FC00.R03



90° Directional Laser Probes



Illuminated
Directional

Standard
Directional

90° Extended Curved Tip

A New Vision on Flexibility

- 90° curved extendable tip for access to peripheral retina
- Hardened metal shaft reduces flexibility by 25%¹
- Retractable mechanism designed to improve control of instrument tip and reduce risk of accidental removal of cannula
- Available in standard (23G, 25G and 27G) and illuminated (23G and 25G) color coded for easy gauge identification
- Compatible with all major equipment

¹ Data on file

Vitreo

Vitreoretinal Surgical Products

bvimedical.com

Vitreoretinal Surgical Products

90° Directional Laser Probes Designed For Improved Surgical Control. The ergonomic handle was designed for improved intraocular surgical control and the sliding knob to control the motion of the instrument tip.


Standard Directional*



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
● 23G Compatible with Alcon®	LP23.A04
● 23G Compatible with Iridex® and B&L	LP23.I04
● 25G Compatible with Alcon®	LP25.A04
● 25G Compatible with Iridex® and B&L	LP25.I04
● 27G Compatible with Alcon®	LP27.A04
● 27G Compatible with Iridex® and B&L	LP27.I04


Illuminated Directional*




Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
● 23G Compatible with Alcon®	LP23.A05
● 23G Compatible with Iridex® and B&L	LP23.I05
● 25G Compatible with Alcon®	LP25.A05
● 25G Compatible with Iridex® and B&L	LP25.I05

* All product names, trademarks and registered trademarks are property of their respective owners.






BVI

Backflush Instruments

Atraumatic Aspiration and Brushing Procedures
Thanks to Improved Intraocular Surgical Control

- Passive and active aspiration control in one instrument
- High volume backflush reservoir for increased capacity during vitreoretinal surgery
- Long shaft (33 mm) makes them suitable for myopic eyes
- Blunt tip and silicone brush needles for atraumatic aspiration and brushing of retinal surface



Vitreo

Vitreoretinal Surgical Products

bvimedical.com

Disposable Backflush Instruments

Disposable backflush instruments are available in three different styles.

Standard

Blunt or brush tip



23G Backflush Instruments Brush
BF23.D02



Blunt tip



Brush tip

Retractable


Silicone brush retracts into the instrument shaft for smooth and valved cannula entry



23G Retractable Brush
BF23.D05

Extendable

8 mm extendable silicone tip for draining fluid from behind retina and other standard backflush procedures



23G Backflush Instruments Extendable Silicone Tip
BF23.D07

Disposable backflush instruments are shipped with cap and tubing for passive or active use during procedures. The direct connection between reservoir and needle allows no contact between hand piece and intraocular fluids.

Backflush Options

All single-use products are provided Sterile and Ready-for-Use.

Standard



BF23.D02

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G Blunt backflush instrument	BF23.D01
● 23G Brush backflush instrument	BF23.D02
● 25G Blunt backflush instrument	BF25.D01
● 25G Brush backflush instrument	BF25.D02
● 27G Blunt backflush instrument	BF27.D01
● 27G Brush backflush instrument	BF27.D02

Retractable



BF23.D05

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G Retractable brush BF-instrument	BF23.D05
● 25G Retractable brush BF-instrument	BF25.D05

Extendable



BF23.D07

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G Extendable silicone tipped BF-instrument	BF23.D07
● 25G Extendable silicone tipped BF-instrument	BF25.D07



Reusable Backflush Instruments

Featuring an ergonomically designed handle, rigid shaft and high volume backflush reservoir for increased backflush capacity, with Passive and Active Aspiration.

Universal Backflush Instrument



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Universal Backflush Instrument	BF00.R01

Backflush Instrument with Luer-lock Connector




Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Universal Backflush Instrument	BF00.R02

Backflush Consumables

Disposable Universal Backflush Needles




BF23.D51

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G BF Blunt needles	BF23.D51
● 23G BF Brush needles	BF23.D52
● 25G BF Blunt needles	BF25.D51
● 25G BF Brush needles	BF25.D52

Blunt Needles with Luer-Lock

Stainless steel cannula with blunt tip for aspiration or injection with syringe or tubing.



CN23.D05

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G Blunt needle with luer-lock connector	CN23.D05
● 25G Blunt needle with luer-lock connector	CN25.D05
● 27G Blunt needle with luer-lock connector	CN27.D05

Backflush Reservoirs



Passive Reservoir
BF00.D51



Active Reservoir
BF00.D52

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
Backflush reservoirs, Passive	BF00.D51
Backflush reservoirs, Active	BF00.D52

Disposable Tubing for Active Aspiration



Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
Disposable tubing for active aspiration	BF00.D53



Vitreoretinal Consumables



Vitreco
Vitreoretinal Surgical Products

bvimedical.com

Cannulas and Needles

New disposable vitreoretinal cannulas and needles include innovative design features to optimize difficult surgical maneuvers. The products are designed to allow the VR surgeon control over some of the most challenging VR procedures. For example, the range includes PFC needles that facilitate injection with minimized pressing force and have a central bore that enables a high volume of heavy liquids to enter the eye at a minimal velocity.

Subretinal Fluid Cannula (Grizzard)

Straight. Sleeve extends 6 mm beyond end.



585021

Qty: 5/Bx
Brand: Visitec®

Description	SKU #
● 20G, 0.9 mm	585021

Vitreoretinal Cannula

Attaches to handpiece or syringe with luer fitting. Can also be used to aspirate or inject silicone oil.



585205

Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 20G x 1-1/8 in, 0.9 x 28 mm	585205

High Viscosity Injector (Hammer)

Thinwall cannula.



585173

Qty: 5/Bx
Brand: Visitec®

Description	SKU #
● 20G (0.9 mm), 4 mm, 10 cm Tube	585173
● 20G (0.9 mm), 4 mm, 16 cm Tube	585208

High Viscosity Injector (Hammer)

Thinwall cannula.



585184

Qty: 5/Bx
Brand: Visitec®

Description	SKU #
● 20G (0.9 mm), 6 mm, 10 cm Tube	585184
● 20G (0.9 mm), 6 mm, 16 cm Tube	585228

Soft Tip Cannula

Soft tip cannulas have flexible, high-flow silicone tips, allowing gentle access to the retina. They are ideal for aspiration of fluid from on or beneath the surface of the retina.



587002


Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 23G x 1-1/3 in, 0.6 x 33 mm, 0.3 mm tip	585226
● 25G x 1-1/3 in, 0.5 x 33 mm, 0.3 mm tip	585229
● 23G x 1-1/3 in, 0.6 x 33 mm, 0.8 mm tip	587002
● 25G x 1-1/3 in, 0.5 x 33 mm, 0.8 mm tip	587003

Cannulas and Needles

VFI Cannula

Silicone oil injection and removal through entry system. High flow rates due to thin walls. Can be directly connected to the silicone oil syringe or to aspiration device.



CN23.D03

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G VFI cannula	CN23.D03
● 25G VFI cannula	CN25.D03
● 27G VFI cannula	CN27.D03

Silicone Tip Cannula

Stainless steel cannula with 2 mm silicone tip for aspiration and manipulation.



CN25.D01

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G Silicone tipped cannula	CN23.D01
● 25G Silicone tipped cannula	CN25.D01

PFC* Injection Needles

New, optimized, PFC injection needle facilitates injection with minimized pressing force. The central bore permits a high volume of heavy liquid to enter the eye at a minimal velocity under controlled conditions.



MD23.D01

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G PFC injection needle	MD23.D01
● 25G PFC injection needle	MD25.D01

* Perfluorocarbon

Pick Needles

Stainless steel cannula with pick tip. Can be connected to a syringe or other handpiece.



CN25.D04

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G Pick needle	CN23.D04
● 25G Pick needle	CN25.D04

Disposable Universal Backflush Needles



BF23.D51

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G BF Blunt needles	BF23.D51
● 23G BF Brush needles	BF23.D52
● 25G BF Blunt needles	BF25.D51
● 25G BF Brush needles	BF25.D52

Blunt Needles with Luer-Lock

Stainless steel cannula with blunt tip for aspiration or injection with syringe or tubing.



CN23.D05

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
● 23G Blunt needle with luer-lock connector	CN23.D05
● 25G Blunt needle with luer-lock connector	CN25.D05
● 27G Blunt needle with luer-lock connector	CN27.D05

Retractable Tissue Manipulators

FEATURES

- Retractable silicone tip with special ribbed structure
- The green colored silicone tip allows for a clear differentiation to the red retinal surface

Tips

Available in 23G and 25G. Disposable and ready to use.



TM23.D01

Qty: 5/Bx
Brand: VitreQ®

Description	SKU #
● 23G Retractable tissue manipulator	TM23.D01
● 25G Retractable tissue manipulator	TM25.D01



Light Source Adaptors

The Light Fibers are compatible with most endo-illumination systems in combination with Light Source Adaptors*.

Constellation® Light Source Adaptor



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Constellation® Light Source Adaptor	LF00.R01

Ruck® Light Source Adaptor



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Ruck® Light Source Adaptor	LF00.R04

Stellaris® Light Source Adaptor



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Stellaris® Light Source Adaptor	LF00.R02

DORC® Light Source Adaptor



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
DORC® Light Source Adaptor	LF00.R05

Geuder® Light Source Adaptor




Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Geuder® Light Source Adaptor	LF00.R03

Accessories

Williamson Scleral Depressor

The Williamson scleral depressor is specially designed to be used for VR-surgery in combination with small gauge entry and BIOM® systems. First scleral indenter specially designed to be used in combination with entry systems with particular attention to sizing and angulation of the indenting surfaces.



Qty: 1/Bx
Brand: Vitreque®



Description	SKU #
3 mm Ball tip, 45 Degree bent tip	VS00.R02

* All product names, trademarks and registered trademarks are property of their respective owners.



Surgical Fluids for Ophthalmic Procedures



Surgical Fluids for Ophthalmic Procedures

Cohesive OVD

NuVisc™ is designed to create and maintain anterior chamber depth and visibility, protect corneal endothelial cells and other intraocular tissues and minimize interaction between tissues during surgical manipulation. NuVisc™ is to be used during ophthalmic surgeries.

NuVisc™ OVD 1.2% Sodium Hyaluronate

Qty: 1/Bx
Brand: NuVisc™

Description	SKU #
NuVisc™ OVD 1.2% sodium hyaluronate, includes 27G Viscoflow® cannula	585305

NuVisc™ Specifications	
Shelf Life	24 Months
Storage Temperature	2-8° C, (36-46° F), do not freeze
Volume	0.8 ml


Chemical Analysis Test Results	
1.6883 in	6.8 - 7.6
Zero shear viscosity	105,000 cps (105 Pa s)
NaHA concentration	12 mg/ml
Molecular weight	Greater than 1 million daltons
Osmolality	320 milliosmoles





Malosa® Single-Use Packs

High Quality Instruments for Your Operating Room



Phaco Packs

Malosa® Core Phaco Pack 1



Qty: 30/Bx

SKU # MMK833/2

Instruments	SKU #
MICS Capsulorhexis - curved/sharp	MMSU1567
Phaco chopper/mushroom manipulator	MMSU1221
Notched forceps (Hoskins)	MMSU1115
Kratz-Barraquer speculum	MMSU1252
Dressing scissors straight rounded	MMSU0455
Clip - Long handled forceps	MMSU1490
Tray - Phaco tray	MMSU0174

Malosa® Core Phaco Pack 2



Qty: 30/Bx

SKU # MMK834/1

Instruments	SKU #
MICS Capsulorhexis - straight/sharp	MMSU1568
Phaco chopper/mushroom manipulator	MMSU1221
Notched forceps (Hoskins)	MMSU1115
Kratz-Barraquer speculum	MMSU1252
Dressing scissors straight rounded	MMSU0455
Clip - Long handled forceps	MMSU1490
Tray - Phaco tray	MMSU0174

LASIK Pack

Malosa® Core LASIK Pack



Qty: 20/Bx

SKU # MMK100/2

Instruments	SKU #
Mechanical speculum, Solid blade (PC)	MMSU1488
Tray - Phaco tray	MMSU0174
LASIK Flap separator and elevator (Binder)	MMSU1206
Y Shape corneal marker	MMSU1134

Smile Pack

Malosa® Core SMILE Pack 2

Qty: 30/Bx

SKU # MMK967



Instruments	SKU #
Reinstein lenticule separator	MMSU1297
Mechanical speculum, Solid blade (PC)	MMSU1488
Notched forceps (Hoskins)	MMSU1115
IOL Removal forceps	MMSU1460
Tray - Phaco tray	MMSU0174

Core Surface Pack

Malosa® Core Surface Treatment Pack

Qty: 20/Bx

SKU # MMK110/3



Instruments	SKU #
Tray - Phaco tray	MMSU0174
Corneal epithelium spatula, Flat	MMSU1113
Epithelium well, 9 mm Trephining	MMSU1196
Mechanical speculum, Solid blade (PC)	MMSU1488

Retina Pack

Malosa® Paediatric Retinal Exam Pack

Qty: 20/Bx

SKU # MMK484



Instruments	SKU #
Neonatal speculum, Small	MMSU1262
Curved scleral depressor	MMSU1409

Product Index

SKU#	Description	PG	SKU#	Description	PG
8660	Weck-Cel® Sponge Points and Strips	62	371395	Round Handle, 21 cm Long x 4 mm Diameter w/ Wrench	17
8661	Weck-Cel® Sponge Points and Strips	62	371410	Round Handle, 10 cm Long x 0.66 cm Diameter	17
8680	Weck-Cel® Cellulose Eye Spears	61	371415	Collet Replacement, 0.66 cm Diameter	17
8685	Weck-Cel® Cellulose Eye Spears	61	372000	Handle Flat, 17 cm Long	17
9701	5 mm Trephine Blade	14	372242	Swivel Fixation Ring	55
9703	5.5 mm Trephine Blade	14	373025	Kojo Safety Slit Knife, 2.5 mm Kojo Slit, 45 Degrees	7
9705	6 mm Trephine Blade	14	373032	Kojo Safety Slit Knife, 3.2 mm Kojo Slit, 45 Degrees	7
9706	6.25 mm Trephine Blade	14	373328	Slit Knife, 2.8 mm Slit, Straight, Matte	8
9707	6.5 mm Trephine Blade	14	373330	Slit Knife, 3 mm Slit, Straight, Matte	8
9708	6.75 mm Trephine Blade	14	373332	Slit Knife, 3.2 mm Slit, Straight, Matte	8
9709	7 mm Trephine Blade	14	373710	Sideport Knife, 1 mm Sideport Knife, Angled	8
9710	7.85 mm Trephine Blade	14	373712	Sideport Knife, 1.2 mm Sideport Knife, Angled	8
9711	7.5 mm Trephine Blade	14	373718	Slit Knife, 1.8 mm Slit, 45 Degrees, Single bevel up	8
9713	8 mm Trephine Blade	14	373722	Slit Knife, 2.2 mm Slit, 45 Degrees, Single bevel up	8
9715	8.5 mm Trephine Blade	14	373724	Slit Knife, 2.4 mm Slit, 45 Degrees, Single bevel up	8
9717	9 mm Trephine Blade	14	373725	Slit Knife, 2.5 mm Slit, 45 Degrees, Single bevel up	8
9719	7.25 mm Trephine Blade	14	373726	Slit Knife, 2.6 mm Slit, 45 Degrees, Single bevel up	8
9721	10 mm Trephine Blade	14	373727	Slit Knife, 2.75 mm Slit, 45 Degrees, Single bevel up	8
9722	7.75 mm Trephine Blade	14	373728	Slit Knife, 2.8 mm Slit, 45 Degrees, Single bevel up	8
9723	8.75 mm Trephine Blade	14	373730	Slit Knife, 3 mm Slit, 45 Degrees, Single bevel up	8
9725	11 mm Trephine Blade	14	373732	Slit Knife, 3.2 mm Slit, 45 Degrees, Single bevel up	8
9727	8.25 mm Trephine Blade	14	373807	Crescent Knife, 2.5 mm, 55 Degrees, Bevel Up	8
9728	14 mm Trephine Blade	14	373808	Crescent Knife, 2.5 mm, 55 Degrees, Straight	8
40413	Ultracell® Cellulose Eye Spears	61	373809	Crescent Knife, 2.5 mm, 55 Degrees, Bevel down	8
40415	Ultracell® Cellulose Eye Spears	61	374310	Round Handle, 10 cm Long x 0.5 cm Diameter	17
40420	Ultracell® Corneal Light Shield	62	374311	Round Handle, 5 cm Long x 0.5 cm Diameter	17
40800	UltraClean™ LASIK PVA Expanded Eye Spears	60	374312	Round Handle, 7.5 cm Long x 0.5 cm Diameter	17
40805	UltraClean™ LASIK PVA Expanded Eye Spears	60	374320	Round Handle, 15.5 cm Long x 0.5 cm Diameter	17
40810	UltraClean™ LASIK PVA Eye Spears	60	374330	Hex Handle, 10 cm Long x 0.5 cm Diameter	17
40815	UltraClean™ LASIK PVA Eye Spears	60	374340	Hex Handle, 15.5 cm Long x 0.5 cm Diameter	17
90900	Visitec® Instrument Wipe	71	374562	Mini-Blade	11
131618	ViscoExpression Cannula (Corydon), 27G x 5/8 in	28	374769	Mini-Blade Knife	9
221204	Standard Forceps	88	374881	Optimum Blade, 15 Degrees, Straight, Green	13
221206	Straight Iris (Jeweler) Forceps	88	374882	Optimum Blade, 22.5 Degrees, Straight, Slate blue	13
221207	Curved Iris (Jeweler) Forceps	88	374883	Optimum Blade, 30 Degrees, Black	13
221208	Adson Forceps	88	374890	Optimum Knife, 15 Degrees, 30 Degree bend left	10
221209	Coaptation Forceps	88	374891	Optimum Knife, 15 Degrees, Straight, Green	10
221211	Tenzel Forceps	88	374892	Optimum Knife, 22.5 Degrees, Straight, Slate blue	10
221250	Wet-Field® Hemostatic Eraser Bipolar Bevel Tip	86	374893	Optimum Knife, 30 Degrees, Straight, Black	10
221251	Wet-Field® Hemostatic Eraser Bipolar Blunt Tip	86	374894	Optimum Knife, 45 Degrees, Straight, Tan	10
221260	Wet-Field® Hemostatic Eraser Bipolar Blunt Tip	86	375000	Handle Flat, 13 cm Long 1 cm Wide	17
221265	Wet-Field® Hemostatic Eraser Bipolar 20G - 23G	87	375435	BeaverGuard Guarded Blade, 0.35 mm Guarded depth	11
221266	Wet-Field® Hemostatic Eraser Bipolar 20G - 23G	87	375700	Sclerotome Blade	11
221267	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	86	375701	Sclerotome Knife	9
221268	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	86	376300	Mini-Blade	11
221270	Wet-Field® Osher Eraser	87	376401	Mini-Blade	11
370145	Safety Slit Knife, 2.5 mm Slit, Angled, Single bevel	6	376481	Beaver® Arthro-Lok® Blade, 4 mm Rounded tip	16
370146	Safety Slit Knife, 2.6 mm Slit, Angled, Single bevel	6	376500	Mini-Blade	11
370820	Slit Knife, 2 mm Slit, 45 Degrees, Double bevel up	8	376620	EdgeAhead Sideport MVR, 1.15 mm Sideport MVR, Angled	9
370822	Slit Knife, 2.2 mm Slit, 45 Degrees, Double bevel up	8	376630	EdgeAhead Sideport MVR, 1.15 mm MVR, Straight	9
370824	Slit Knife, 2.4 mm Slit, 45 Degrees, Double bevel up	8	376980	Beaver® Arthro-Lok® Blade, 4 mm Mini-banana	16
370825	Slit Knife, 2.5 mm Slit, 45 Degrees, Double bevel up	8	376984	Beaver® Arthro-Lok® Blade, 4 mm Banana	16
370826	Slit Knife, 2.65 mm Slit, 45 Degrees, Double bevel up	8	377100	Myringotomy Blade, Downcutting tip	15
370827	Slit Knife, 2.75 mm Slit, 45 Degrees, Double bevel up	8	377110	Myringotomy Blade, Upcutting tip	15
370828	Slit Knife, 2.8 mm Slit, 45 Degrees, Double bevel up	8	377120	Myringotomy Blade, Offset tip	15
370830	Slit Knife, 3 mm Slit, 45 Degrees, Double bevel up	8	377121	Myringotomy Blade, Narrow shaft	15
371301	Blade Wrench	17	377122	Myringotomy Knife, Spear tip	15
371305	Collet Replacement, 0.5 cm Diameter	17	377124	Myringotomy Knife, Spear tip	15
371350	Round Handle, 13 cm Long x 0.5 cm Diameter	17	377200	Tympanoplasty Blade	16
371360	Round Handle, 13 cm Long x 0.66 cm Diameter	17	377511	Micro-Sharp® Blade, 1.5 mm, 15 Degrees, Orange	13
371394	Round Handle, 18 cm Long x 4 mm Diameter w/ Wrench	17	377512	Micro-Unitome Knife, 1.5 mm, 15 Degrees, Orange	10

Product Index

SKU#	Description	PG	SKU#	Description	PG
377513	Micro-Sharp® Blade, 3 mm, 15 Degrees, Blue	13	581148	Visitec® Fluid Collection Bag	70
377514	Micro-Unitome Knife, 3 mm, 15 Degrees, Blue	10	581149	Plastic Incise Drape with Fluid Collection Bag	65
377515	Micro-Sharp® Blade, 5 mm, 15 Degrees, Green	13	581150	Non-Woven Incise Drape with Fluid Collection Bag	65
377516	Micro-Unitome Knife, 5 mm, 15 Degrees, Green	10	581151	Plastic Aperture Drape with Fluid Collection Bag	65
377522	Micro-Sharp® Blade, 4 mm, 22.5 Degrees, Light blue	13	581152	Non-Woven Aperture Drape with Fluid Collection Bag	66
377523	Micro-Unitome Knife, 4 mm, 22.5 Degrees, Light blue	10	581153	Non-Woven Incise Drape with Channel Bag	65
377530	Micro-Sharp® Blade, 3.5 mm, 30 Degrees, Black	13	581154	Mini Pre-Split Incise Drape	68
377531	Micro-Unitome Knife, 3.5 mm, 30 Degrees, Black	10	581155	Non-Woven Aperture Drape with Channel Bag	67
377546	Micro-Unitome Knife, 3 mm, 45 Degrees, Tan	10	581156	Full Body Incise Drape	63
378210	Safety Sideport Knife, 1 mm Sideport MVR,	7	581157	Full Body Aperture Drape	63
378212	Safety Sideport Knife, 1.2 mm Sideport MVR	7	581159	Mini Plastic Incise Drape	68
378214	Safety Slit Knife, 1.4 mm Slit, Angled, Single bevel	6	581160	Mini Plastic Aperture Drape	68
378215	Safety Slit Knife, 1.5 mm Slit, Angled, Single bevel	6	581161	Mini Plastic Aperture Drape	68
378217	Safety Slit Knife, 1.7 mm Slit, Angled, Single bevel	6	581162	Dual Refractive Drape	69
378218	Safety Slit Knife, 1.8 mm Slit, Angled, Single bevel	6	581163	Full Body Incise Drape with Plastic Head and Channel Bag	63
378220	Safety Slit Knife, 2 mm Slit, Angled, Single bevel	6	581164	Non-Woven Full Body Incise Drape	64
378222	Safety Slit Knife, 2.2 mm Slit, Angled, Single bevel	6	581166	Non-Woven Half Body Incise Drape* with Channel Bag	67
378224	Safety Slit Knife, 2.4 mm Slit, Angled, Single bevel	6	581177	Visitec® Fluid Collection Bag	70
378227	Safety Slit Knife, 2.75 mm Slit, Angled, Single bevel	6	581273	Anterior Chamber (Rycroft), 30G x 7/8 in, 0.3 x 22 mm	21
378228	Safety Slit Knife, 2.8 mm Slit, Angled, Single bevel	6	581275	Retrobular (Sharp Point), 25G x 1-1/2 in, 0.5 x 38 mm	20
378230	Safety Slit Knife, 3 mm Slit, Angled, Single bevel	6	581276	Lacrimal Cannula, 26G x 1 1/8 in, 0.45 x 28 mm	32
378231	Safety Sideport Knife, 1.15 mm 20G Sideport	7	581277	Anterior Chamber (Rycroft), 16G x 1 in, 1.6 x 25 mm	28
378232	Safety Slit Knife, 3.2 mm Slit, Angled, Single bevel	6	581280	Anterior Chamber (Rycroft), 27G x 7/8 in, 0.4 x 22 mm	21
378233	Safety Sideport Knife, 0.6 mm 23G Sideport	7	581281	Anterior Chamber (Rycroft), 21G x 7/8 in, 0.8 x 22 mm	21
378234	Safety Crescent Knife, 2.5 mm Crescent, 55 Degrees	7	581284	Anterior Chamber Cannula, 27G x 7/8 in, 0.4 x 22 mm	21
378235	Optimum Safety Sideport Knife, 15 Degree Sideport	7	581285	Anterior Chamber Cannula, 30G x 7/8 in, 0.3 x 22 mm	21
378236	Optimum Safety Sideport Knife, 22.5 Degree Sideport	7	581411US	Tying Forceps	43
378237	Optimum Safety Sideport Knife, 30 Degree Sideport	7	581413	Eye Speculum	37
378238	Optimum Safety Sideport Knife, 45 Degree Sideport	7	581422US	Vannas Scissors	47
378820	Safety Slit Knife, 2 mm Slit, Angled, Double bevel	6	581437US	Westcott Scissors	47
378822	Safety Slit Knife, 2.2 mm Slit, Angled, Double bevel	6	581441	Visitec® Micro-Manipulator	45
378824	Safety Slit Knife, 2.4 mm Slit, Angled, Double bevel	6	581442	Sinskey Hook Gen 2	45, 58
378825	Safety Slit Knife, 2.5 mm Slit, Angled, Double bevel	6	581443	Eye Speculum Open Wire Loops	37
378826	Safety Slit Knife, 2.65 mm Slit, Angled, Double bevel	6	581458	Jewellers Forceps	41
378827	Safety Slit Knife, 2.75 mm Slit, Angled, Double bevel	6	581459US	Micro Tying Forceps	43
378828	Safety Slit Knife, 2.8 mm Slit, Angled, Double bevel	6	581478US	Eye Speculum (Lieberman)	36
378829	Safety Slit Knife, 2.85 Slit, Angled, Double bevel	6	581481US	Capsulorhexis Forceps (Utrata) Straight	40
378830	Safety Slit Knife, 3 mm Slit, Angled, Double bevel	6	581486	Capsulorhexis Forceps Grasping/Curved Tips - Gen 3	40
378832	Safety Slit Knife, 3.2 mm Slit, Angled, Double bevel	6	581488	Capsulorhexis Forceps Sharp/Curved Tips - Gen 3	40
379081	Beaver® Arthro-Lok® Blade, 4 mm Mini-menisus	16	581592	Nucleus Hydrodissector, 30G x 7/8 in, 0.3 x 22 mm	23
581001	Small Pupil IOL Glide (Pearce)	56	581594	Irrigating Cystotome, 23G x 5/8 in, 0.6 x 16 mm	26
581007	Superflex Extension Set	78	581599	Irrigating Cystotome, 23G x 5/8 in, 0.6 x 16 mm	26
581012	Irrigation Handle	78	581602	Irrigating Side Cutting Cystotome	26
581019	Eye Shield (Left/Right)	76	581604	Capsule Polisher (Kratz), 21G x 7/8 in, 0.8 x 22 mm	31
581030	Visitec® Bulk High Flow Filter	78	581605	Capsule Polisher (Kratz), 23G x 7/8 in, 0.6 x 22 mm	31
581032	Small Pupil IOL (Bobbitt)	56	581610	Irrigating Cystotome, 25G x 5/8 in, 0.5 x 16 mm	26
581033	IOL Glide (Sheets)	56	581614	Capsule Polisher, 21G x 7/8 in, 0.8 x 22 mm	30
581034	IOL Glide (Hessburg)	56	581615	Lacrimal Cannula, 21G x 1 3/8 in, 0.8 x 35 mm	32
581040	Universal Eye Shield	76	581617	Irrigating Cystotome, 25G x 5/8 in, 0.5 x 16 mm	26
581047	Visitec® Visiwipe Instrument Wipe	71	581618	Irrigating Cystotome, 27G x 5/8 in, 0.4 x 16 mm	26
581049	Surgical Glide (Fichman)	56	581630	Irrigating Cystotome, 30G x 5/8 in, 0.3 x 16 mm	26
581056	Visitec® Double Male Luer Adapter	78	581637	Retrobular (Atkinson), 25G x 1-1/2 in, 0.5 x 38 mm	20
581057	Corneal Marking Pad	78	581638	Retrobular (Atkinson), 23G x 1-1/2 in, 0.6 x 38 mm	20
581060	Absorbent Stick	72	581639	Anterior Chamber (Rycroft), 25G x 7/8 in, 0.5 x 22 mm	21
581062	Visitec® Corneal Light Shield	62	581640	Retrobular (Sharp Atkinson), 23G x 1-1/2 in, 0.6 x 38 mm	20
581069	Visiswab Cotton Swabs	72	581641	Retrobular (Sharp Atkinson), 25G x 1-1/2 in, 0.5 x 38 mm	20
581076	Visitec® Visiwipe Instrument Wipe	71	581642	Anterior Chamber Cannula, 25G x 7/8 in, 0.5 x 22 mm	21
581089	Visitec® Cellulose Eye Spears	61	581645	Posterior Scraper (Simcoe), 23G x 7/8 in, 0.6 x 22 mm	31
581098	Visitec® LASIK Eye Drain (Chayet)	70	581708	Universal Pediatric Eye Shield	76
581146	Fluid Collection Bag with Wick	70	581709	Visitec® LASIK PVA Expanded Eye Spear	60
581147	Visitec® Eye Fluid Wick	62	581761	Visitec® LASIK PVA Eye Spears	60

Product Index

SKU#	Description	PG	SKU#	Description	PG
585001	Retrobulbar/Peribulbar (Atkinson)	20	585157	Cortical Cleaving Hydrodissector (Round End)	24
585003	Micro J Shape Cannula, 25G x 7/8 in, 0.5 x 22 mm	25	585158	Cortical Cleaving Hydrodissector (Flat End)	23
585004	Irrigating Cystotome, 25G 5/8 in, 0.5 x 16 mm	26	585159	Cortical Cleaving Hydrodissector (Bolger)	24
585005	Flexible Tip Cannula, 27G x 7/8 in, 0.4 x 22 mm	30	585163	ACM Micro Etched (Dahan)	22
585006	Anterior Chamber Cannula Straight, 27G	21	585173	High Viscosity Injector (Hammer)	114
585007	Viscoflow Formed, 27G x 7/8 in, 0.4 x 22 mm	28	585176	Sub-Tenon's Anesthesia Cannula	21
585008	Irrigating Cystotome, 27G 5/8 in, 0.4 x 16 mm	26	585181	Nucleus Chopper (P. Koch)	44
585009	Viscoflow Cannula, 23G x 7/8 in, 0.6 x 22 mm	28	585184	High Viscosity Injector (Hammer)	114
585015	Nucleus Hydrodelineator (P. Koch)	23	585191	HydroExpression Cannula (Corydon)	25
585016	Lacrimal Cannula, 19G x 3/4 in, 1.1 x 20 mm	32	585193	Cortex Extractor, 23G x 5/8 in, 0.6 x 16 mm	29
585017	Flexible Tip Cannula, 25G x 7/8 in, 0.5 x 22 mm	30	585195	Nucleus Hydrodissector (Sinskey)	23
585018	Micro J Shape Cannula, 27G x 7/8 in, 0.5 x 22 mm	25	585197	Cortical Cleaving Hydrodissector (Grabow)	25
585019	Retrobulbar, Curved (Atkinson/Straus), 25G	20	585198	Hydrodissection Cannula (Akaoshi Style)	25
585021	Subretinal Fluid Cannula (Grizzard), 20G, 0.9 mm	114	585199	Nucleus Hydrodissector (Buratto)	24
585025	Anterior Chamber (Rycroft), 23G x 7/8 in, 0.6 x 22 mm	21	585200	FLACS Cannula (Schweitzer)	33
585026	Viscoflow Cannula, 27G x 7/8 in, 0.4 x 22 mm	28	585203	Pediatric Eye Speculum	38
585027	Peribulbar (Atkinson), 25G x 7/8 in, 0.5 x 22 mm	20	585205	Vitreoretinal Cannula, 20G x 1-1/8 in, 0.9 x 28 mm	114
585029	Micro J Shape Hydrodissector (Pearce), 27G	25	585208	High Viscosity Injector (Hammer)	114
585030	Irrigating Cystotome (Joyce), 25G x 5/8 in, 0.5 x 16 mm	25	585209	Micro Iris Hook	58
585031	Lacrimal Trephine (Sisler), 21G x 1-1/2 in, 0.8 x 38 mm	32	585213	Spatula PRK	51
585033	Lacrimal Cannula (Shahinian), 23G x 1/2 in, 0.6 x 12 mm	32	585222	Bimanual Aspiration Cannula, 23G	29
585034	Irrigating Cystotome (Reverse), 25G x 5/8 in, 0.5 x 16 mm	27	585223	Bimanual Aspiration Cannula, 21G	29
585036	Nucleus Hydrodissector (Pearce) with 8 mm bend	23	585226	Soft Tip Cannula, 23G x 1-1/3 in, 0.6 x 33 mm, 0.3 mm tip	114
585037	Nucleus Hydrodissector (Pearce), 25G x 7/8 in, 0.5 x 22 mm	23	585227	Cortex Extractor, 23G x 15/16 in, 0.60 x 23.5 mm	27
585039	Irrigating Cystotome (Joyce), 27G x 5/8 in, 0.4 x 16 mm	27	585228	High Viscosity Injector (Hammer)	114
585040	Bishop Harmon Cannula, 19G x 1 in, 1.1 x 25 mm	28	585229	Soft Tip Cannula, 25G x 1-1/3 in, 0.5 x 33 mm, 0.3 mm tip	114
585041	Flexible Tip Cannula, 23G x 7/8 in, 0.6 x 22 mm	30	585231	EdgeAhead Sideport MVR, 0.6 mm MVR, Straight	9
585042	Irrigating Cystotome (Blumenthal), 27G x 5/8 in, 0.4 x 16 mm	27	585232	EdgeAhead Sideport MVR, 0.6 mm Stiletto	9
585043	Irrigating Cystotome (Blumenthal), 27G x 5/8 in, 0.4 x 16 mm	26	585233	LASIK Flap Roller	50
585045	Viscoflow Cannula, 25G x 7/8 in, 0.5 x 22 mm	28	585240	EdgeAhead Sideport MVR, 1.1 mm MVR, Straight	9
585046	Viscoflow Cannula, 30G x 7/8 in, 0.3 x 22 mm	28	585250	LASIK Irrigating Cannula (Lindstrom)	33
585047	I and A Cannula Thinwall (Simcoe), 23G, 0.6 mm	30	585252	LASIK Flap Spatula/Irrigator (Hersh)	33
585049	Irrigating Formed Cystotome, 25G x 5/8 in, 0.5 x 16 mm	26	585271	LASIK Cannula (Daya), 23G x 7/8 in, 0.6 x 22 mm	33
585061	Self Retaining ACM (Lewicky), 20G x 9/64 in, 0.9 x 3.5 mm	22	585274	Nucleus Hydrodissector (Pearce), 27G x 7/8 in, 0.4 x 22 mm	23
585062	IOL Manipulator Blunt Tip	45	585278	LASIK Cannula, 27G x 7/8 in, 0.4 x 22 mm	33
585068	Lacrimal Cannula, 25G x 1/2 in, 0.5 x 13 mm	32, 94	585280	LASIK Irrigating Cannula (Solomon)	33
585075	Cortex Extractor, 23G 5/8 in, 0.6 x 16 mm	29	585294	Nucleus Hydrodissector (Chang), 27G x 7/10 in, 0.4 x 18 mm	24
585077	Nucleus Rotator	45	585295	Nucleus Hydrodissector (Chang), 27G x 7/10 in, 0.4 x 18 mm	24
585086	I and A Cannula (Simcoe), 23G, 0.6 mm	30	585297	Plastic Locking Eye Speculum	36
585090	HydroExpression Cannula (Corydon), 25G	25	585305	NuVisc™ OVD 1.2% Sodium Hyaluronate	120
585091	HydroExpression Cannula (Corydon), 23G	25	587000	I-Ring® Pupil Expander, 6/Pk	58
585093	Irrigating Cystotome (Simcoe), 27G	27	587001	I-Ring® Pupil Expander, 1/Pk	58
585094	Eye Speculum (Kratz/Barraquer)	37	587002	Soft Tip Cannula, 23G x 1-1/3 in, 0.6 x 33 mm, 0.8 mm tip	114
585099	Nucleus Hydrodissector (Sauter)	23	587003	Soft Tip Cannula, 25G x 1-1/3 in, 0.5 x 33 mm, 0.8 mm tip	114
585106	Cortex Extractor, 24G x 1/2 in, 0.55 x 12.5 mm	29	8310060	Doctor's Six-Pack Diagnostic Penlight	77
585107	Tapered Hydrodelineator (Blumenthal)	24	8310080	Nurse-Mate® Six-Pack Diagnostic Penlight	77
585108	Retrobulbar/Peribulbar (Atkinson), 25G x 1-1/4 in	20	8310201	Maxi-Lite® Diagnostic Penlight	77
585109	Retrobulbar/Peribulbar (Atkinson), 27G x 1-1/4 in	20	8320150	Blu-Spot® Diagnostic Light	77
585110	I and A Cannula (Simcoe), 23G, 0.6 mm	30	8440000	Accu-Temp® Variable Temperature Cautery	85
585113	Irrigating Cystotome (Blumenthal), 25G	27	8441000	Accu-Temp® Low Temperature Cautery	84
585115	Blunt Tip Irrigating Cystotome (Kellan)	27	8441500	Accu-Temp® High Temperature Cautery	84
585124	Punctal Dilator and Probe	53, 93	8442000	Accu-Temp® High Temperature Cautery	84
585133	Glaucoma Revision Pick 60 Degree (Lederer)	31	8443000	Accu-Temp® High Temperature Cautery	84
585135	Glaucoma Revision Pick 30 Degree (Lederer)	31	8444000	Accu-Temp® High Temperature Cautery	84
585143	Olive Tip Polisher, 23G x 1 in, 0.6 x 25 mm	31	8445000	Accu-Temp® High Temperature Cautery	84
585144	Peribulbar (Atkinson), 23G x 7/8 in, 0.6 x 22 mm	20	8446000	Accu-Temp® Variable Temperature Cautery	85
585145	Peribulbar (Atkinson), 27G x 7/8 in, 0.4 x 22 mm	20	8461000	Accu-Temp® Reusable Handle	85
585147	Olive Tip Cannula, 23G x 7/8 in, 0.6 x 22 mm	31	8461020	Accu-Temp® Reusable Handle	85
585148	Micro J Shape Cannula (Hessburg)	25	8461050	Single-use fine tip	85
585149	Oval ACM (Blumenthal)	22	8461060	Single-use loop tip	85
585155	Nucleus Hydrodissector, 27G 7/8 in, 0.4 x 22 mm	24	8590300	Visitec® Guibor Eye Bandage	76

Product Index

SKU#	Description	PG	SKU#	Description	PG
8590350	Guibor Canaliculus Intubation Set (CIS)	55	A14-202QL	Parasol® Quickload Punctal Occluder Sterile	90
8590450	Canaliculus Intubation Set (CIS)	55	A14-203	Parasol® Pre-Loaded Punctal Occluder Sterile	90
9550005	9.5 mm Trephine Blade	14	A14-203BN	Parasol® Bulk Punctal Occluder Non-Sterile	90
9550006	10.5 mm Trephine Blade	14	A14-203QL	Parasol® Quickload Punctal Occluder Sterile	90
9550013	11.5 mm Trephine Blade	14	A14-204	Plus Pre-Loaded Punctal Occluder Sterile	91
9550025	9.25 mm Trephine Blade	14	A14-204BN	Plus Punctal Occluder Non-Sterile, Bulk	91
9550026	4 mm Trephine Blade	14	A14-205	Plus Pre-Loaded Punctal Occluder Sterile	91
9550028	12 mm Trephine Blade	14	A14-205BN	Plus Punctal Occluder Non-Sterile, Bulk	91
9550029	13 mm Trephine Blade	14	A14-205QL	Parasol® Quickload Punctal Occluder Sterile	90
22135010	Wet-Field® Disposable Cable	87	A14-206	Plus Pre-Loaded Punctal Occluder Sterile	91
22155010	Wet-Field® Disposable Cable Shrouded	87	A14-206BN	Plus Punctal Occluder Non-Sterile, Bulk	91
221260135	Wet-Field® Hemostatic Eraser Bipolar Blunt Tip	86	A14-207	Plus Pre-Loaded Punctal Occluder Sterile	91
221260150	Wet-Field® Hemostatic Eraser Bipolar Blunt Tip	86	A14-207BN	Plus Punctal Occluder Non-Sterile, Bulk	91
221267135	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	86	A14-209BN	Parasol® Bulk Punctal Occluder Non-Sterile	90
221267150	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	86	A14-401	Micro Flow™ Punctal Occluder Sterile	92
221268135	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	86	A14-402	Micro Flow™ Punctal Occluder Sterile	92
221268150	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	86	A14-403	Micro Flow™ Punctal Occluder Sterile	92
221269135	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	86	A20-401BN	Micro Flow™ Punctal Occluder Non-sterile	92
10100-C	Classic Nasal Pack, 4.5 cm with String	73	A20-402BN	Micro Flow™ Punctal Occluder Non-sterile	92
10110-C	Classic Nasal Pack, 8 cm with String	73	A20-403BN	Micro Flow™ Punctal Occluder Non-sterile	92
10115-C	Classic Nasal Pack, 10 cm with String	73	BEAVER6400	Grindless Beaver Mini-Blades	12
10120-C	Anatomical Nasal Pack, 8 cm with String	73	BEAVER6600	Grindless Beaver Mini-Blades	12
10140-C	Epistaxis Pack, 10 cm with String	73	BEAVER6610	Grindless Beaver Mini-Blades	12
10146-C	Epistaxis Slimpack, 8 cm with String	73	BEAVER6700	Grindless Beaver Mini-Blades	12
10150-C	Epistaxis Pack, 5.5 cm with String	73	BEAVER6900	Grindless Beaver Mini-Blades	12
10160-C	Anterior Epistaxis Pack, 3 cm with String	73	BEAVER6910	Grindless Beaver Mini-Blades	12
20202-C	Sinus Ultra Slimpack with String	74	BF00.D51	Backflush reservoirs, Passive	112
20210-C	Sinus Pack with String	74	BF00.D52	Backflush reservoirs, Active	112
20211-C	Sinus Slimpack with String	74	BF00.D53	Disposable tubing for active aspiration	112
20212-C	Sinus Ultra Slimpack with String	74	BF00.R01	Universal Backflush Instrument	112
30300-C	Ear Wick	74	BF00.R02	Universal Backflush Instrument with Luer-lock Connector	112
30301-C	Pediatric Ear Wick	74	BF23.D01	Standard 23G Blunt backflush instrument	111
30305-C	Fenestrated Ear Wick	74	BF23.D02	Standard 23G Brush backflush instrument	111
30306-C	Fenestrated Pediatric Ear Wick	74	BF23.D05	Retractable 23G Retractable brush backflush instrument	111
30310-C	Long Ear Wick	74	BF23.D07	Extendable 23G Silicone tipped backflush instrument	111
30311-C	Pediatric Long Ear Wick	74	BF23.D51	23G BF Blunt needles	112, 115
30316-C	Pediatric Long Fenestrated Ear Wick	74	BF23.D52	23G BF Brush needles	112, 115
40400-8	Ultracell® PVA Eye Spears	60	BF25.D01	Standard 25G Blunt backflush instrument	111
40405-8	Ultracell® PVA Eye Spears	60	BF25.D02	Standard 25G Brush backflush instrument	111
A12-102	Collagen Synthetic Punctal Occluder	93	BF25.D05	Retractable 25G Retractable brush backflush instrument	111
A12-103	Collagen Synthetic Punctal Occluder	93	BF25.D07	Extendable 25G Silicone tipped backflush instrument	111
A12-104	Collagen Synthetic Punctal Occluder	93	BF25.D51	25G BF Blunt needles	112, 115
A12-200	Collagen Synthetic Punctal Occluder	93	BF25.D52	25G BF Brush needles	112, 115
A13-200	Extend® Synthetic Punctal Occluder	92	BF27.D01	Standard 27G Blunt backflush instrument	111
A13-2000	Extend® 180 Synthetic Punctal Occluder	92	BF27.D02	Standard 27G Brush backflush instrument	111
A13-202	Extend® Synthetic Punctal Occluder	92	CN23.D01	23G Silicone tipped cannula	115
A13-203	Extend® Synthetic Punctal Occluder	92	CN23.D03	23G VFI cannula	115
A13-2033	Extend® 180 Synthetic Punctal Occluder	92	CN23.D04	23G Pick needle	115
A13-2033SP	Extend® 180 Synthetic Punctal Occluder	92	CN23.D05	23G Blunt needle with luer-lock connector	112, 115
A13-204	Extend® Synthetic Punctal Occluder	92	CN25.D01	25G Silicone tipped cannula	115
A13-2044	Extend® 180 Synthetic Punctal Occluder	92	CN25.D03	25G VFI cannula	115
A13-2044SP	Extend® 180 Synthetic Punctal Occluder	92	CN25.D04	25G Pick needle	115
A13-205	Extend® Synthetic Punctal Occluder	92	CN25.D05	25G Blunt needle with luer-lock connector	112, 115
A13-2055	Extend® 180 Synthetic Punctal Occluder	92	CN27.D03	27G VFI cannula	115
A13-2055SP	Extend® 180 Synthetic Punctal Occluder	92	CN27.D05	27G Blunt needle with luer-lock connector	112, 115
A14-200	Parasol® Pre-Loaded Punctal Occluder Sterile	90	CT003	Monitor	100
A14-201	Parasol® Pre-Loaded Punctal Occluder Sterile	90	EO-CAP-S	Endoscope Connector Caps	99
A14-201BN	Parasol® Bulk Punctal Occluder Non-Sterile	90	FC00.R03	Reusable Cleaning Adapter	106
A14-201QL	Parasol® Quickload Punctal Occluder Sterile	90	FC23.D02	Serrated Gripping 23G Single-Use	104
A14-202	Parasol® Pre-Loaded Punctal Occluder Sterile	90	FC23.D03	ILM End-gripping 23G Single-Use	104
A14-202BN	Parasol® Bulk Punctal Occluder Non-Sterile	90	FC23.D04	Asymmetrical Gripping 23G Single-Use	104

Product Index

SKU#	Description	PG	SKU#	Description	PG
FC23.D05	Eckardt Endgripping 23G Single-Use	105	MMSU0277S	Moorfields Forceps	41
FC23.D06	Kuhn 45° Power Gripping 23G Single-Use	105	MMSU0279S	St Martins Forceps	41
FC23.D07	Shah X-tra Grip 23G Single-Use	105	MMSU0283S	Jewellers Forceps, Fine Angled	41
FC23.R02	Serrated Gripping 23G Reusable	104	MMSU0286S	Jewellers Forceps, Fine Straight	41
FC23.R03	ILM End-gripping 23G Reusable	104	MMSU0454S	Dressing Scissors, Straight, Pointed	46
FC23.R04	Asymmetrical Gripping 23G Reusable	104	MMSU0455S	Dressing Scissors, Straight, Rounded	46
FC23.R05	Eckardt Endgripping 23G Reusable	105	MMSU1102S	Govan Manipulator	44
FC23.R06	Kuhn 45° Power Gripping 23G Reusable	105	MMSU1103S	Sinsky Hook 0.2 mm	44
FC23.R07	Shah X-tra Grip 23G Reusable	105	MMSU1105S	Claw Chopper	43
FC25.D02	Serrated Gripping 25G Single-Use	104	MMSU1108S	Phaco Chopper	44
FC25.D03	ILM End-gripping 25G Single-Use	104	MMSU1110S	Corneal Needle-Holder	48
FC25.D04	Asymmetrical Gripping 25G Single-Use	104	MMSU1113S	Corneal Epithelium Spatula Flat	51
FC25.D05	Eckardt Endgripping 25G Single-Use	105	MMSU1114S	Fine Notched Forceps, Short (Hoskins)	39
FC25.D06	Kuhn 45° Power Gripping 25G Single-Use	105	MMSU1115S	Notched Forceps (Hoskins)	39
FC25.D07	Shah X-tra Grip 25G Single-Use	105	MMSU1119S	Greens Chopper, Right Handed	43
FC25.R02	Serrated Gripping 25G Reusable	104	MMSU1123S	Fine Toothed Forceps Short (Bonn)	38
FC25.R03	ILM End-gripping 25G Reusable	104	MMSU1124S	Toothed Forceps (Bonn)	39
FC25.R04	Asymmetrical Gripping 25G Reusable	104	MMSU1125S	Straight Tying Forceps	42
FC25.R05	Eckardt Endgripping 25G Reusable	105	MMSU1127S	Angled Tying Forceps (McPherson)	42
FC25.R06	Kuhn 45° Power Gripping 25G Reusable	105	MMSU1129S	Colibri Toothed Forceps, Short	39
FC25.R07	Shah X-tra Grip 25G Reusable	105	MMSU1134S	Y Shape Corneal Marker	49
FC27.D03	ILM End-gripping 27G Single-Use	104	MMSU1137S	Thornton Fixation Ring	55
FC27.D05	Eckardt Endgripping 27G Single-Use	105	MMSU1140S	Vannas Scissors, Curved	47
FC27.D07	Shah X-tra Grip 27G Single-Use	105	MMSU1141S	Vannas Scissors, Straight	47
FC27.R03	ILM End-gripping 27G Reusable	104	MMSU1144S	Drysdale Manipulator	44
FC27.R05	Eckardt Endgripping 27G Reusable	105	MMSU1146S	Lens Loading Forceps	42
FC27.R07	Shah X-tra Grip 27G Reusable	105	MMSU1149S	Chalazion Clamp, Oval	54
FC2720	Flight Case	99	MMSU1150S	Chalazion Clamp, Round	54
FC401655	Flight Case	99	MMSU1152S	Mendez Gauge	52
LF00.R01	Constellation® Light Source Adaptor	117	MMSU1153S	ICL Manipulator	45
LF00.R02	Stellaris® Light Source Adaptor	117	MMSU1157S	Chalazion Curette, Double, 1 mm/2 mm	54
LF00.R03	Geuder® Light Source Adaptor	117	MMSU1158S	Chalazion Curette, Double, 2 mm/3 mm	54
LF00.R04	Ruck® Light Source Adaptor	117	MMSU1169S	LASIK Flap-Lifting Forceps	41
LF00.R05	DORC® Light Source Adaptor	117	MMSU1172S	LASIK Flap Separator/Elevator (Seibel)	49
LF29.D02	29G Spotlight Directional Chandelier	102	MMSU1177S	LASIK Flap Elevator Curved 8 mm	49
LG001	Diode Laser Safety Goggles	99	MMSU1179S	Mushroom Manipulator	44
LP23.A04	Laser Probe 23G Compatible with Alcon®	108	MMSU1181S	Punctum Dilator, Fine/Medium	53, 94
LP23.A05	Illuminated Laser Probe 23G Compatible with Alcon®	108	MMSU1187S	Corneal Epithelium Spatula Round	51
LP23.I04	Laser Probe 23G Comp. with Iridex® and B&L	108	MMSU1193S	Punctum Dilator/Seeker	53, 94
LP23.I05	Illuminated Laser Probe 23G Comp. with Iridex® and B&L	108	MMSU1196S	Epithelium Well	51
LP25.A04	Laser Probe 25G Compatible with Alcon®	108	MMSU1197S	Epithelium Well	51
LP25.A05	Illuminated Laser Probe 25G Comp. with Alcon®	108	MMSU1206S	LASIK Flap Separator and Elevator, Binder	49
LP25.I04	Laser Probe 25G Compatible with Iridex® and B&L	108	MMSU1207S	Lester Hook	44
LP25.I05	Illuminated Laser Probe 25G Comp. with Iridex® and B&L	108	MMSU1209S	Colibri Notched Forceps, Short	39
LP27.A04	Laser Probe 27G Comp. with Alcon®	108	MMSU1210S	Curved Tying Forceps	43, 94
LP27.I04	Laser Probe 27G Comp. with Iridex® and B&L	108	MMSU1211S	Westcott Scissors, Curved, Pointed	47
MD23.D01	23G PFC injection needle	115	MMSU1213S	Westcott Scissors Curved Dolphin Nose	47
MD25.D01	25G PFC injection needle	115	MMSU1214S	Westcott Scissors, Straight, Pointed	47
MMK833/2	Malosa® Core Phaco Pack 1	122	MMSU1221S	Phaco Chopper/Mushroom Manipulator	45
MMK834/1	Malosa® Core Phaco Pack 2	122	MMSU1229S	Braunstein IVT Marker 3.5/4	52
MMK100/2	Malosa® Core LASIK Pack	122	MMSU1231S	Horizontal Chopper/Kuglen Hook	45
MMK967	Malosa® Core SMILE Pack 2	123	MMSU1243S	Troutman Needle-Holder	48
MMK110/3	Malosa® Core Surface Treatment Pack	123	MMSU1252S	Kratz-Barraquer Speculum	37
MMK484	Malosa® Paediatric Retinal Exam Pack	123	MMSU1253S	Kratz-Barraquer Speculum, Temporal	37
MMSU0068S	Straight Rambley	48	MMSU1275S	Barraquer Speculum, Flat Blade	37
MMSU0085S	Epilation Forceps	41	MMSU1276S	Braunstein IVT Marker 3.5/4, Green	52
MMSU0214S	Stevens Tenotomy, Straight, Rounded	46	MMSU1277S	Pointed IVT Marker 3.5/4, Blue	52
MMSU0238S	Adson Forceps, Toothed	38	MMSU1284S	Contact Lens Forceps	42
MMSU0241S	Adson Forceps, Non-Toothed	41	MMSU1285S	Double Ended IVT Guide Pad	52
MMSU0244S	Iris Scissors, Curved, Pointed	46	MMSU1287S	Kratz-Barraquer Speculum, Curved	37
MMSU0249S	Iris Scissors, Straight, Pointed	46	MMSU1290S	Mechanical Speculum, Open Blade (PC)	36

Product Index

SKU#	Description	PG
MMSU1294S	LASIK Flap Separator/Elevator Curved	49
MMSU1295S	LASIK Flap Separator/Elevator Straight	50
MMSU1297S	Reinstein Lenticule Separator	49
MMSU1405S	Schocket Double-Ended Scleral Depressor	54
MMSU1408S	Ball Tip Scleral Depressor	54
MMSU1417S	Golf Club Spud	53
MMSU1419S	Punctum Dilator/Lacrimal Probe	53
MMSU1458S	Lebuisson Manipulator	44
MMSU1460S	IOL Removal Forceps	39
MMSU1474S	Barraquer Speculum, Solid Curved (PC)	37
MMSU1475S	Barraquer Speculum, Solid Blade (PC)	38
MMSU1479S	Barraquer Speculum, Solid 18 mm Blade (PC)	38
MMSU1488S	Mechanical Speculum Solid Blade	36
MMSU1495S	Fixation Ring	55
MMSU1567S	MICS Capsulorhexis Forceps - Curved/Sharp	40
MMSU1568S	MICS Capsulorhexis Forceps - Straight/Sharp	40
MMSU1577S	Intravitreal Injection Guide	52
MMSU1578S	Ramley Short	48
MMSU1583S	Waqar Suture Removal Forceps	43
MMSU3015S	Barraquer Needle-Holder	48
MON006	Monitor	100
MON007	Monitor	100
NC001-A	Endoscope Tip Protector	99
OME170Z	Video Adapter	100
OME170ZMG	Video Adapter	100
OME2000	E2 Ophthalmic Laser Endoscopy System	98
OME200HRA	Dual Function Ophthalmic Endoscope	99
OME200SMAEA	Triple Function Ophthalmic Laser Endoscope	98
OME200SMAHRA	Triple Function Ophthalmic Laser Endoscope	98
OME200SMAHRVA	Triple Function Ophthalmic Laser Endoscope	98
OME200SMAHRXA	Triple Function Ophthalmic Laser Endoscope	98
OME200SMAHRXVA	Triple Function Ophthalmic Laser Endoscope	98
OME200SMAVA	Triple Function Ophthalmic Laser Endoscope	98
OME200SMAXA	Triple Function Ophthalmic Laser Endoscope	98
OME200SMAXVA	Triple Function Ophthalmic Laser Endoscope	98
OME230HRA	Dual Function Ophthalmic Endoscope	99
OME230SMAHRXA	Triple Function Ophthalmic Laser Endoscope	98
OME300ZMG	Video Adapter	100
OME4000	E4 Ophthalmic Endoscopy System	98
SC23.D01	Vertical Micro Scissors 23G Single-Use	106
SC23.D02	Straight Micro Scissors 23G Single-Use	106
SC23.D03	Horizontal Curved Micro Scissors 23G Single-Use	106
SC23.R01	Vertical Micro Scissors 23G Reusable	106
SC23.R02	Straight Micro Scissors 23G Reusable	106
SC23.R03	Horizontal Curved Micro Scissors 23G Reusable	106
SC25.D01	Vertical Micro Scissors 25G Single-Use	106
SC25.D02	Straight Micro Scissors 25G Single-Use	106
SC25.D03	Horizontal Curved Micro Scissors 25G Single-Use	106
SC25.R02	Straight Micro Scissors 25G Reusable	106
SC25.R03	Horizontal Curved Micro Scissors 25G Reusable	106
SC27.D01	Vertical Micro Scissors 27G Single-Use	106
SC27.D03	Horizontal Curved Micro Scissors 27G Single-Use	106
SC27.R03	Horizontal Curved Micro Scissors 27G Reusable	106
SW002	Footswitch with illumination	100
SW003	Footswitch for laser only	100
TM23.D01	23G Retractable tissue manipulator	116
TM25.D01	25G Retractable tissue manipulator	116
VP6250	Sterilization Tray	100
VS00.R02	Williamson Scleral Depressor	118

BVI® is refocusing the future of vision.

As one of the fastest-growing, diversified surgical ophthalmic businesses in the world, our purpose-built portfolio spans more than 115 countries.

We've set our sights on touching the lives of millions of patients affected by conditions such as cataracts, refractive error, glaucoma, retinal disease, and dry eye.

Unburdened by legacy or bureaucracy, we have developed our strategy around a simple concept — taking pride in delivering innovative solutions for our physicians and patients, based on their needs.

We trust and empower our associates to make decisions and solve problems because collaboration drives us. Valuing agility, simplicity, and transparency, we stay committed to listening to our customers, delivering for our patients, and keeping the future in focus.



THE FUTURE IN FOCUS



Beaver®
Knives and Blades

Visitec®
Cannulas

Malosa®
Single-Use Instruments

CustomEyes®
Custom Procedure Packs

Vitreo
Vitreoretinal Surgical Products

EndoOptiks®
Microendoscopy

FINEVISION®
Trifocal IOLs

ISOPURE®
Premium Monofocal IOLs

R-EVOLUTION
Ophthalmic Equipment

R-EVOSmart
Ophthalmic Equipment

BVI
Email: customersupport@bvimedical.com

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2021-2023



1594865-01