



BVI

Endo Optiks® Ophthalmic Microendoscopes

See. Diagnose. Treat.

Endo Optiks ophthalmic microendoscopes and probes offer a new perspective for treatment options across many surgical specialties. Ophthalmic microendoscopy is a powerful surgical solution for combined phaco-ECP, vitreoretinal and minimally invasive glaucoma surgery (MIGS).

Microendoscopes and probes can be used with the E2 or E4 ophthalmic laser and microendoscopy system. Microendoscopes and probes are available in a variety of sizes and resolutions to meet all of your surgical needs.

Key Features

- All microendoscopes and probes are autoclavable
- Sizes range from 25 gauge to 18 gauge
- High resolution image display
- Ultra wide field imaging and illumination
- Straight or curved tip configurations



Ophthalmic Microendoscopes and Probes

Our patented ophthalmic laser microendoscopes and probes are available in triple, dual and single function configurations. This delivers a combination of visualized laser application, high resolution imaging and wide field illumination, through one single microprobe.

Triple Function Ophthalmic Laser Microendoscopes									
Item #	Gauge	Cannula	High Resolution	Laser	Camera	Light	Custom Laser Connector	E2 Compatibility	E4 Compatibility
OME230SMAHRXA	23	Straight	√	√	√	√	√	√	√
OME230SMAXA	23	Straight	—	√	√	√	√	√	√
OME200SMAEA	20	Straight	—	√	√	√	—	√	—
OME200SMAXA	20	Straight	—	√	√	√	√	√	√
OME200SMAHRA	19	Straight	√	√	√	√	—	√	—
OME200SMAHRVA	19	Curved	√	√	√	√	—	√	—
OME200SMAHRXA	19	Straight	√	√	√	√	√	√	√
OME200SMAHRXVA	19	Curved	√	√	√	√	√	√	√
OME200SMAVA	19	Curved	—	√	√	√	—	√	—
OME200SMAXVA	19	Curved	—	√	√	√	√	√	√
Dual Function Ophthalmic Endoscopes									
Item #	Gauge	Cannula	High Resolution	Laser	Camera	Light	Custom Laser Connector	E2 Compatibility	E4 Compatibility
OME230LA	23	Straight	—	—	√	√	—	√	√
OME230HRA	23	Straight	√	—	√	√	—	√	√
OME200LA	20	Straight	—	—	√	√	—	√	√
OME200HRA	20	Straight	√	—	√	√	—	√	√
OME200LVA	19	Curved	—	—	√	√	—	√	√
Dual Function Ophthalmic Laser Probes									
Item #	Gauge	Cannula	High Resolution	Laser	Camera	Light	Custom Laser Connector	E2 Compatibility	E4 Compatibility
VPH230LA	23	Straight	—	√	—	√	—	√	√
OME20LA	20	Straight	—	√	—	√	—	√	√
Single Function Ophthalmic Probes									
Item #	Gauge	Cannula	High Resolution	Laser	Camera	Light	Custom Laser Connector	E2 Compatibility	E4 Compatibility
VPH250A	25	Straight	—	√	—	—	—	√	—
VPH230A	23	Straight	—	√	—	—	—	√	—
VP230WA	23	Straight	—	—	—	√	—	√	√
VP200WA	20	Straight	—	—	—	√	—	√	√