

# BVI

## Beaver<sup>®</sup> Safety Knives

World Class Leader in High Quality Surgical Blades

- Single handed activation
- Consistent blade quality
- Safety options include Crescent, MVR, Sideport and Slit



**Beaver<sup>®</sup>**  
Knives and Blades

# Complete Portfolio




## Phaco Slit/Keratome Knives


These knives are designed to make a precise incision size to accommodate the phaco tip.

**Safety Slit** knives are designed for penetration and accurate incision width.



- Available in a wide range of sizes (1.4 – 3.2 mm)
- Single bevel and double bevel blade configurations
- Reference line guides a multi-planar incision



Single Bevel Incision



Double Bevel Incision




Knife Sizes	Reference Line From Tip
1.8 mm and 2 mm	1.75 mm
> 2 mm	2 mm



## Safety

- True no-look, single-handed activation
- Integrated retractable shield prevents sharps injuries during handling and use
- Activated with a slider that provides tactile, visual and audible feedback
- Safety shield does not interfere with grip, offering better control during use



## Paracentesis/Sideport Knives

These knives are designed for accurate width stab incisions.

The angled **Safety Sideport** knives create a consistent incision size.

**Optimum® Safety Sideport** knives enable the surgeon to control the incision width by varying the depth of the incision.



## Crescent Knives

These knives are designed to create optimal scleral tunnel architecture. They can also be used for trabeculectomy and lamellar dissections.

Available in angled, bevel up configuration



## Consistency

- Advanced proprietary manufacturing process, delivers consistency in sharpness
- Safety shield protects blade against inadvertent damage ensuring a consistent incision

## Sharpness

- Enables surgeon to maneuver the blade and optimize incision architecture
- Designed to deliver precision

# Beaver<sup>®</sup> Safety Knives



## Safety Slit Knives

### Single Bevel

Qty: 10/Bx

Description	SKU #
1.4 mm Slit, Angled, Single bevel	378214
1.5 mm Slit, Angled, Single bevel	378215
1.7 mm Slit, Angled, Single bevel	378217
1.8 mm Slit, Angled, Single bevel	378218
2 mm Slit, Angled, Single bevel	378220
2.2 mm Slit, Angled, Single bevel	378222
2.4 mm Slit, Angled, Single bevel	378224
2.5 mm Slit, Angled, Single bevel	370145
2.6 mm Slit, Angled, Single bevel	370146
2.75 mm Slit, Angled, Single bevel	378227
2.8 mm Slit, Angled, Single bevel	378228
3 mm Slit, Angled, Single bevel	378230
3.2 mm Slit, Angled, Single bevel	378232
2.5 mm Kojo Slit, 45 Degrees, Arc height is 0.55 mm	373025
3.2 mm Kojo Slit, 45 Degrees, Arc height is 0.71 mm	373032

### Double Bevel

Qty: 10/Bx

Description	SKU #
2 mm Slit, Angled, Double bevel	378820
2.2 mm Slit, Angled, Double bevel	378822
2.4 mm Slit, Angled, Double bevel	378824
2.5 mm Slit, Angled, Double bevel	378825
2.65 mm Slit, Angled, Double bevel	378826
2.75 mm Slit, Angled, Double bevel	378827
2.8 mm Slit, Angled, Double bevel	378828
2.85 mm Slit, Angled, Double bevel	378829
3 mm Slit, Angled, Double bevel	378830
3.2 mm Slit, Angled, Double bevel	378832



## Optimum Safety Sideport Knives

Qty: 10/Bx

Description	SKU #
15 Degree Sideport, Straight	378235
22.5 Degree Sideport, Straight	378236
30 Degree Sideport, Straight	378237
45 Degree Sideport, Straight	378238



## Safety Sideport Knives

Qty: 10/Bx

Description	SKU #
1 mm Sideport MVR, Angled 45 degrees	378210
1.2 mm Sideport MVR, Angled 45 degrees	378212
0.6 mm 23G Sideport MVR	378233
1.15 mm 20G Sideport MVR, Angled 45 degrees	378231



## Safety Crescent Knives

Qty: 10/Bx

Description	SKU #
2.5 mm Crescent, 55 Degrees, Bevel up	378234

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

**Beaver<sup>®</sup>**  
Knives and Blades

Questions: [bvimedical.com/customer-support](http://bvimedical.com/customer-support)

