

Market Leading Portfolio of Ophthalmic Surgical Devices

Providing Solutions for All Steps of Your Cataract Surgery Pathway



Endoscopy
93



Micro Incisions
9



OVD Injection
118



Pupil Expansion
62



Capsulorhexis
32



Hydrodissection/ Hydrodelineation
29



Anesthesia
26



IOL Implantation
119



Wound Hydration
63



Eye Protection
74



Patient Draping
66

Providing Solutions
For Glaucoma Surgery



Drapes and Covers
66



Knives and Blades
9



Cannulas
25



Single-Use Instruments
41



Fluid Control
63



Laser Endoscopy System
93



OVDs
118

Providing Solutions for All Steps of Your Vitreoretinal Surgery Pathway



Patient Preparation
66




Endoscopy Imaging
93



Pathology Management
99



Endoillumination
97



Coagulation
81



Fluid Control
63



Eye Protection
74

Imagine Better. Together.



An evolving purpose-built portfolio for the next level of patient care.

Table of Contents

Knives and Blades 9

Safety Knives 10
Conventional Knives 12
Mini-Blades® 15
Grindless Mini-Blades® 16
Blades 17
Trephines 18
Other Specialty Blades 19
Handles and Accessories 23
Handles and Accessories Chart 24

Cannulas 25

Anesthesia Cannulas 26
Anterior Chamber Cannulas 27
Anterior Chamber Maintainers 28
Hydrodissection Cannulas 29
Cystotome Cannulas 32
Irrigation Cannulas 34
Irrigation and Aspiration Cannulas 34
Polishers 36
Glaucoma Revision Picks 37
Lacrimal Cannulas 38
Refractive and FLACS 39
Hubs 40

Single-Use Instruments 41

Adjustable Speculums 42
Sprung Speculums 42
Fixation Forceps - Toothed 44
Fixation Forceps - Notched 44
Capsulorhexis Forceps 45
General Forceps 46
IOL Forceps 47
Tying Forceps 47
Choppers 48
Double-Ended Choppers and Manipulators 49
Manipulators 50
Dissection Scissors 52
Dressing Scissors 52
Sprung Scissors 53
Dressing Forceps 54
Needle Holders and Drivers 54
Refractive Surgery Instrumentation 55
Lacrimal Intubation CIS 57
Clamps and Curettes 57
Markers and Gauges 58
Probes and Dilators 59
Fixation Rings 59
Depressors 60
IOL Glides 60

Small Pupil Expansion 61

I-Ring® Pupil Expander 62

Fluid Control 63
PVA Eye Spears 64
Cellulose Eye Spears 65
Sponge Points, Strips and Wicks 65
Corneal Light Shields 65
Full Body Drapes 66
Half Body Drapes 67
Face Drapes 70
Eye Drains 71
Fluid Collection 71
Instrument Wipes 72
Other Fluid Control 72

Accessories 73
Eye Shields and Bandages 74
Diagnostic Lights 75
General Accessories 76

Procedure Packs 77
CustomEyes® 78
CustomEyes® Configurator 80

Electrosurgery and Electrocautery 81
Cautery Devices 82
Erasers 84
Forceps 86

Dry Eye and Lacrimal Disorders 87
Permanent Occlusion 88
Mid-Term Occlusion 90
Partial Occlusion 90
Temporary Occlusion 91
Single-Use Punctal Dilators and Probes 91
Other Related Occlusion Options 92

Equipment 93
Leos® Laser Endoscopy Ophthalmic System 94
VueProbe® Ophthalmic Laser Endoscope 95
Accessories 96

Vitreoretinal 97
Directional Chandelier 97
Micro Forceps and Micro Scissors 99
90° Directional Laser Probes 103
Backflush Instruments 105
Light Fibers 109

Vitreoretinal Consumables 111

Surgical Fluids for Ophthalmic Procedures 117
NuVisc™ OVD 1.2% Sodium Hyaluronate 118

IOLs 119
FINEVISION HP 120

Product Index 121

High-Quality, Single-Handed Activation Safety Knives


Sharpness | Consistency | Safety



Options for Crescent, MVR, Sideport and Slit

Safety Knives

Safety Slit Knives



Safety Slit Knives are available in a wide range of sizes designed for penetration and accurate incision width.



378224



378824

Single Bevel **Qty:** 10/Bx

Description	SKU #
1.4 mm Slit, Angled, Single bevel	378214
1.5 mm Slit, Angled, Single bevel	378215
1.7 mm Slit, Angled, Single bevel	378217
1.8 mm Slit, Angled, Single bevel	378218
2 mm Slit, Angled, Single bevel	378220
2.2 mm Slit, Angled, Single bevel	378222
2.4 mm Slit, Angled, Single bevel	378224
2.5 mm Slit, Angled, Single bevel	370145
2.6 mm Slit, Angled, Single bevel	370146
2.75 mm Slit, Angled, Single bevel	378227
2.8 mm Slit, Angled, Single bevel	378228
3 mm Slit, Angled, Single bevel	378230
3.2 mm Slit, Angled, Single bevel	378232

Double Bevel **Qty:** 10/Bx

Description	SKU #
2.2 mm Slit, Angled, Double bevel	378822
2.4 mm Slit, Angled, Double bevel	378824
2.5 mm Slit, Angled, Double bevel	378825
2.65 mm Slit, Angled, Double bevel	378826
2.75 mm Slit, Angled, Double bevel	378827
2.8 mm Slit, Angled, Double bevel	378828
2.85 mm Slit, Angled, Double bevel	378829
3 mm Slit, Angled, Double bevel	378830
3.2 mm Slit, Angled, Double bevel	378832



Safety Knives

Safety Crescent Knives



Safety Crescent Knives are angled 55° with the bevel up, designed for optimal tunnel architecture.



378234

Qty: 10/Bx

Description	SKU #
2.5 mm Crescent, 55°, Bevel up	378234

Optimum Safety Sideport Knives



Safety Sideport Knives are designed for stab incisions.



378235

Qty: 10/Bx

Description	SKU #
15° Sideport, Straight	378235
22.5° Sideport, Straight	378236
30° Sideport, Straight	378237
45° Sideport, Straight	378238

Safety Sideport Knives



Safety Sideport Knives are designed for accurate width stab incisions.



378212

Qty: 10/Bx

Description	SKU #
1 mm Sideport, Angled 45°	378210
1.2 mm Sideport, Angled 45°	378212
0.6 mm 23G Sideport MVR	378233
1.15 mm 20G Sideport, Angled 45°	378231

Kojo Safety Slit Knives



Kojo Safety Slit Knives are slit blades with an arcuate shape.




373025

Qty: 10/Bx

Description	SKU #
2.5 mm Kojo Slit, 45°, Arc height is 0.55 mm	373025
3.2 mm Kojo Slit, 45°, Arc height is 0.71 mm	373032

Conventional Knives

Slit Knives



Slit Knives are available in a wide range of sizes designed for penetration and accurate incision width.



373727

Single Bevel **Qty:** 10/Bx

Description	SKU #
1.8 mm Slit, 45°, Single bevel up	373718
2.2 mm Slit, 45°, Single bevel up	373722
2.4 mm Slit, 45°, Single bevel up	373724
2.5 mm Slit, 45°, Single bevel up	373725
2.6 mm Slit, 45°, Single bevel up	373726
2.75 mm Slit, 45°, Single bevel up	373727
2.8 mm Slit, Straight, Single bevel up	373328
2.8 mm Slit, 45°, Single bevel up	373728
3 mm Slit, Straight, Single bevel up	373330
3 mm Slit, 45°, Single bevel up	373730
3.2 mm Slit, Straight, Single bevel up	373332
3.2 mm Slit, 45°, Single bevel up	373732



370824

Double Bevel **Qty:** 10/Bx

Description	SKU #
2 mm Slit, 45°, Double bevel	370820
2.2 mm Slit, 45°, Double bevel	370822
2.4 mm Slit, 45°, Double bevel	370824
2.5 mm Slit, 45°, Double bevel	370825
2.65 mm Slit, 45°, Double bevel	370826
2.75 mm Slit, 45°, Double bevel	370827
2.8 mm Slit, 45°, Double bevel	370828
3 mm Slit, 45°, Double bevel	370830

Crescent Knives



Crescent Knives offer the flexibility to create optimal tunnel architecture. Available in straight, bevel up or bevel down.



373809

Qty: 10/Bx

Description	SKU #
2.5 mm, 55°, Bevel up	373807
2.5 mm, Straight	373808
2.5 mm, 55°, Bevel down	373809

Sideport Knives



Sideport Knives are designed for accurate width stab incisions.



373712

Qty: 10/Bx

Description	SKU #
1 mm Sideport Knife, Angled	373710
1.2 mm Sideport Knife, Angled	373712


Conventional Knives

EdgeAhead Sideport and MVR





585232



376620



376630

Qty: Various/Bx

Description	SKU #
0.6 mm Stiletto, 5/Bx	585232
0.6 mm MVR, Straight, 5/Bx	585231
1.1 mm MVR, Straight, 5/Bx	585240
1.15 mm MVR, 20G, Straight, 10/Bx	376630
1.15 mm Sideport MVR, Angled, 10/Bx	376620

Mini-Blade® Knife





Qty: 10/Bx

Description	SKU #
Surgical blade, 180° cutting surface	374769

Sclerotome Knife



The Sclerotome knife has a smaller cutting edge than the Sclerotome blade (375700) and is designed to have a more precise, controlled dissection.



Qty: 10/Bx

Description	SKU #
Multisided sharp all around	375701



Conventional Knives

Optimum Knives



Optimum Knives are color coded for identification. Offered in a wide range of sizes.


374890


374891


374892



374893


374894

Qty: 10/Bx

Description	SKU #
15°, 30° bend left	374890
15°, Straight, Green	374891
22.5°, Straight, Slate blue	374892
30°, Straight, Black	374893
45°, Straight, Tan	374894

Micro-Unitome Knives



Micro-Unitome Knives are manufactured using a proprietary grindless technology to support a more consistent sharpness. Each knife is color coded for identification.


377512


377514


377516


377523


377531


377546

Qty: 10/Bx

Description	SKU #
1.5 mm, 15°, Orange	377512
3 mm, 15°, Blue	377514
5 mm, 15°, Green	377516
4 mm, 22.5°, Light blue	377523
3.5 mm, 30°, Black	377531
3 mm, 45°, Tan	377546



Mini-Blades®

Available in a variety of designs to meet all small incision and excision needs.

Mini-Blade®



Qty: 24/Bx

Description	SKU #
Orthopedic Mini-Blade®, Length 1.255 in (31.88 mm)	374562

Mini-Blade®



Qty: 24/Bx

Description	SKU #
MIS 6400, Narrower version of Beaver6400 for minimally invasive surgery	376401

Mini-Blade®



Qty: 3/Bx

Description	SKU #
Overall length 1.31 inches	376100

Mini-Blade®



Qty: 6/Bx

Description	SKU #
Sharp on two sides with pointed tip	376300

Mini-Blade®



Qty: 24/Bx

Description	SKU #
Pointed tip, Smaller alternative to #11 general blade	376500

Sclerotome Blade

The Sclerotome is a versatile blade suitable for glaucoma and retinal procedures.



Qty: 3/Bx

Description	SKU #
Multisided sharp all around S/Su	375700

BeaverGuard™ Guided Blades

Preset to ensure uniform depth control.



375435

Qty: 3/Bx

Description	SKU #
0.35 mm Guarded depth S/Su	375435
0.5 mm Guarded depth S/Su	375450

Grindless Mini-Blades®

BVI's Mini-Blades® are manufactured using a proprietary grindless technology and are designed to meet all small incision and excision needs. Mini-Blades® are preassembled into a unique, patented tip, designed for safe assembly. Each Mini-Blade® is color coded for identification.

6400 Mini-Blade®



Qty: 20/Bx

Description	SKU #
Rounded tip, Sharp on one side, Green	BEAVER6400

6700 Mini-Blade®



Qty: 20/Bx

Description	SKU #
Smaller alternative to #15 general blade, Orange	BEAVER6700

6600 Mini-Blade®



Qty: 3/Bx

Description	SKU #
Lamellar, Angled 60°, Smooth undersurface, Sharp all around, Bevel up, Black	BEAVER6600

6900 Mini-Blade®



Qty: 20/Bx

Description	SKU #
Straight, Sharp all around, 180° cutting surface, Double bevel blade, Blue	BEAVER6900

6610 Mini-Blade®



Qty: 3/Bx

Description	SKU #
Lamellar, Straight, Smooth undersurface, Sharp all around, Black	BEAVER6610

6910 Mini-Blade®



Qty: 6/Bx

Description	SKU #
Cleft palate blade, Angled 60°, Sharp all around, 180° cutting surface, Double bevel blade, Blue	BEAVER6910

Blades

Micro-Sharp® Blades



Micro-Sharp Blades are manufactured using a proprietary grindless technology. Each blade is color coded for identification.

377511


377513


377515


377522


377530


Qty: 6/Bx

Description	SKU #
1.5 mm, 15°, Orange	377511
3 mm, 15°, Blue	377513
5 mm, 15°, Green	377515
4 mm, 22.5°, Light blue	377522
3.5 mm, 30°, Black	377530

Optimum Blades



Optimum stab blades are color coded for identification. Offered in a wide range of sizes.

374881


374882


374883


Qty: 5/Bx

Description	SKU #
15°, Straight, Green	374881
22.5°, Straight, Slate blue	374882
30°, Black	374883



bvimedical.com

Trephines

Trephine Blades are handcrafted with a sharp, uniform cutting edge, designed for a precise, symmetrical cut of the corneal layers.

Overall length 15.8 mm and provided sterile for single use.

Qty: 1/Bx

Description	SKU #
4 mm Trephine Blade	9550026
5 mm Trephine Blade	0009701
5.5 mm Trephine Blade	0009703
6 mm Trephine Blade	0009705
6.25 mm Trephine Blade	0009706
6.5 mm Trephine Blade	0009707
6.75 mm Trephine Blade	0009708
7 mm Trephine Blade	0009709
7.25 mm Trephine Blade	0009719
7.5 mm Trephine Blade	0009711
7.75 mm Trephine Blade	0009722
7.85 mm Trephine Blade	0009710
8 mm Trephine Blade	0009713

Description	SKU #
8.25 mm Trephine Blade	0009727
8.5 mm Trephine Blade	0009715
8.75 mm Trephine Blade	0009723
9 mm Trephine Blade	0009717
9.25 mm Trephine Blade	9550025
9.5 mm Trephine Blade	9550005
10 mm Trephine Blade	0009721
10.5 mm Trephine Blade	9550006
11 mm Trephine Blade	0009725
11.5 mm Trephine Blade	9550013
12 mm Trephine Blade	9550028
13 mm Trephine Blade	9550029
14 mm Trephine Blade	0009728



Other Specialty Blades


Myringotomy Blade



Qty: 6/Bx

Description	SKU #
Downcutting tip, Offset 45°	377100

Myringotomy Blade



Qty: 6/Bx

Description	SKU #
Upcutting tip, Offset 45°.	377110


Myringotomy Blade



Qty: 6/Bx

Description	SKU #
Offset tip, 45°	377120

Myringotomy Blade



Qty: 6/Bx

Description	SKU #
Narrow shaft, Offset 45°	377121

Myringotomy Knife



Qty: 3/Bx

Description	SKU #
Spear tip, Offset 45°	377122

Myringotomy Knife




Qty: 3/Bx

Description	SKU #
Lance tip, Offset 45°	377124

Other Specialty Blades

Arthro-Lok® Blade



Qty: 6/Bx

Description	SKU #
4 mm Rounded tip	376481

Arthro-Lok® Blade



Qty: 6/Bx

Description	SKU #
4 mm Rounded tip, Curved right, Sharp on one side	376482

Arthro-Lok® Blade



Qty: 6/Bx

Description	SKU #
4 mm Pointed tip, Curved left, Sharp on one side	376780


Arthro-Lok® Blade



Qty: 6/Bx

Description	SKU #
4 mm Mini-banana, Sharp all around, Curved	376980

Arthro-Lok® Blade



Qty: 6/Bx

Description	SKU #
4 mm Mini-banana, Flat	376981

Arthro-Lok® Blade



Qty: 6/Bx

Description	SKU #
4 mm Banana, Sharp all around	376984

Arthro-Lok® Blade



Qty: 6/Bx

Description	SKU #
4 mm Mini-meniscus, Flat	379081

Arthro-Lok® Blade



Qty: 6/Bx

Description	SKU #
4 mm Mini-meniscus, Curved	379082

Other Specialty Blades

DeBakey Blades



374710



374810

Qty: 24/Bx

Description	SKU #
8 mm Cutting width	374710
10 mm Cutting width	374810

Needle Blades

375910



375920



375930



375940



Qty: 3/Bx


Description	SKU #
Needle Blade, 3 mm Cutting edge	375910
Needle Blade, 5 mm Cutting edge	375920
Needle Blade, 7 mm Cutting edge	375930
Needle Blade, 1.5 mm Cutting edge	375940





Other Specialty Blades


Tympanoplasty Blade



Qty: 3/Bx

Description	SKU #
60° angle, Bevel down, Sharp all around, 2.5 mm width	377200

Tympanoplasty Blade



Qty: 3/Bx

Description	SKU #
2.5 mm width, Straight, Sharp all around	377210

Sickle Edge Blade




Qty: 3/Bx

Description	SKU #
Straight, 84 mm length, 5 mm cutting surface	377300

Ring Cutter and Blades



374100




374105 374105 374150 374151

Qty: 1/Bx

Description	SKU #
Ring cutter	374100
Ring cutter hub, Regular speed steel	374105
Ring cutter washer, Regular speed steel	374106
Ring cutter saw blade, Regular speed steel	374150
Ring cutter saw blade, High speed steel	374151

Handles and Accessories

Handles



Stainless steel handles are available in round and hex configurations, in a variety of lengths.



371395



374320



374340



375000

Type: Reusable

Qty: 1/Bx

Description	SKU #
Round handle, 13 cm Long x 0.5 cm Diameter	371350
Round handle, 13 cm Long x 0.66 cm Diameter	371360
Round handle, 18 cm Long x 4 mm Diameter with wrench	371394
Round handle, 21 cm Long x 4 mm Diameter with wrench	371395
Round handle, 10 cm Long x 0.66 cm Diameter	371410
Round handle, 10 cm Long x 0.5 cm Diameter	374310
Round handle, 5 cm Long x 0.5 cm Diameter	374311
Round handle, 7.5 cm Long x 0.5 cm Diameter	374312
Round handle ,15.5 cm Long x 0.5 cm Diameter	374320
Hex handle, 10 cm Long x 0.5 cm Diameter	374330
Hex handle, 15.5 cm Long x 0.5 cm Diameter	374340
Flat handle, 13 cm Long, 1 cm Wide	375000

Wrenches and Collets

All wrenches and collets are stainless.



371301



371305



371415

Qty: 1/Bx

Description	SKU #
Blade wrench	371301
Collet replacement, 0.5 cm Diameter	371305
Collet replacement, 0.66 cm Diameter	371415

Handles and Accessories Chart

BLADES	HANDLES										REPLACEMENT COLLET	
	371350	371360	371394	371395	374310	374311	374312	374320	374330	374340	375000	371305
Beaver6400	x	x	x	x	x	x	x	x	x	x		
Beaver6600	x	x	x	x	x	x	x	x	x	x		
Beaver6610	x	x	x	x	x	x	x	x	x	x		
Beaver6700	x	x	x	x	x	x	x	x	x	x		
Beaver6900	x	x	x	x	x	x	x	x	x	x		
Beaver6910	x	x	x	x	x	x	x	x	x	x		
374562	x	x			x	x	x	x	x	x		x
374710											x	
364810											x	
374580			x	x								
375435	x	x			x	x	x	x	x	x		x
375450	x	x			x	x	x	x	x	x		x
375700	x	x			x	x	x	x	x	x		x
375910	x	x			x	x	x	x	x	x		x
375920	x	x			x	x	x	x	x	x		x
375930	x	x			x	x	x	x	x	x		x
375940	x	x			x	x	x	x	x	x		x
376100	x	x			x	x	x	x	x	x		x
376300	x	x			x	x	x	x	x	x		x
376401		x			x	x	x	x	x	x		x
376481			x	x								
376482			x	x								
376500	x	x			x	x	x	x	x	x		x
376780			x	x								
376800	x	x			x	x	x	x	x	x		x
376980			x	x								
376981			x	x								
376984			x	x								
377100	x	x			x	x	x	x	x	x		x
377110	x	x			x	x	x	x	x	x		x
377120	x	x			x	x	x	x	x	x		x
377121	x	x			x	x	x	x	x	x		x
377511	x	x			x	x	x	x	x	x		
377513	x	x			x	x	x	x	x	x		
377515	x	x			x	x	x	x	x	x		
377522	x	x			x	x	x	x	x	x		
377530	x	x			x	x	x	x	x	x		
377120	x	x			x	x	x	x	x	x		x
377200	x	x			x	x	x	x	x	x		x
377210	x	x			x	x	x	x	x	x		x
377300	x	x			x	x	x	x	x	x		x
379081			x	x								
379082			x	x								

40 Years of Innovation and Quality Cannulas

Supporting All Areas of Ophthalmic Surgery

- A market leader of ophthalmic cannulas for over 40 years
- Full range of single-use cannulas developed to support all ophthalmic surgical procedures
- A history of innovation to meet your diverse surgical needs





Anesthesia Cannulas

Retrobulbar (Atkinson)



Intended for injecting anaesthetic agents into the muscle cone.



581638

Qty: 10/Bx

Description	SKU #
23G x 1-1/2 in, 0.6 x 38 mm	581638
25G x 1-1/2 in, 0.5 x 38 mm	581637

Retrobulbar/Peribulbar (Atkinson)





585108

Qty: 10/Bx

Description	SKU #
23G x 1-1/4 in, 0.6 x 32 mm	585001
25G x 1-1/4 in, 0.5 x 32 mm	585108
27G x 1-1/4 in, 0.4 x 32 mm	585109

Peribulbar (Atkinson)



Intended for injecting anaesthetic agents around the globe (periocular).




585144

Qty: 10/Bx

Description	SKU #
23G x 7/8 in, 0.6 x 22 mm	585144
25G x 7/8 in, 0.5 x 22 mm	585027
27G x 7/8 in, 0.4 x 22 mm	585145

Retrobulbar (Sharp Atkinson)





581640

Qty: 10/Bx

Description	SKU #
23G x 1-1/2 in, 0.6 x 38 mm	581640
25G x 1-1/2 in, 0.5 x 38 mm	581641

Retrobulbar (Sharp Point)



Intended for injecting anaesthetic agents into the muscle cone.



581275

Qty: 10/Bx

Description	SKU #
25G x 1-1/2 in, 0.5 x 38 mm	581275

Retrobulbar, Curved (Atkinson/Straus)



Designed to course around the globe and enter the anterior region of the muscle cone inferotemporal and parallel to the optic nerve. Straight portion 11 mm, curved portion 22 mm.



585019

Qty: 10/Bx

Description	SKU #
25G x 1-3/8 in, 0.5 x 34 mm	585019

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 40.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G

Anesthesia Cannulas

Sub-Tenon's Anesthesia Cannula



Intended for administering anaesthetic agents into the posterior sub-tenon's space.



585176

Qty: 10/Bx

Description	SKU #
19G x 1 in, 1.1 x 25 mm	585176

Anterior Chamber Cannulas

Intended to be used for irrigation/injection of liquid.

Anterior Chamber (Rycroft)



Angled 45°, 4 mm from end.



585025

Qty: 10/Bx

Description	SKU #
21G x 7/8 in, 0.8 x 22 mm	581281
23G x 7/8 in, 0.6 x 22 mm	585025
25G x 7/8 in, 0.5 x 22 mm	581639
27G x 7/8 in, 0.4 x 22 mm	581280
30G x 7/8 in, 0.3 x 22 mm	581273

Anterior Chamber Cannula



Angled 45°, 6 mm from end.



581284

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	581642
27G x 7/8 in, 0.4 x 22 mm	581284
30G x 7/8 in, 0.3 x 22 mm	581285

Anterior Chamber Cannula Straight





585006

Qty: 10/Bx

Description	SKU #
27G x 7/8 in, 0.4 x 22 mm	585006

Anterior Chamber Maintainers

Self-Retaining ACM (Lewicky)



Intended to aid in maintaining anterior chamber space via irrigation. Flexible tubing 25 cm (10 in). Inner diameter 1 mm (0.04 in); outer diameter 2 mm (0.08 in). Unique annular rings start 0.75 mm from 55° bevelled tip for retaining without sutures.




Qty: 3/Bx

Description	SKU #
20G x 9/64 in, 0.9 x 3.5 mm	585061

ACM Micro Etched (Dahan)



Intended to aid in maintaining anterior chamber space via irrigation. Flexible tubing 15.25 cm (6 in). Inner diameter 0.75 mm (0.03 in), outer diameter 1.6 mm (0.07 in). Special micro etching on tip allows ACM to be self-retaining. Ultra thin wall tip provides high flow rate.




Qty: 10/Bx

Description	SKU #
20G x 5/32 in, 0.9 x 4 mm, 45° beveled tip	585163

Oval ACM (Blumenthal)



Intended to aid in maintaining anterior chamber space via irrigation. 3.5 mm from 40° bevelled tip. Flexible tubing 25 cm (10 in). Inner diameter 1 mm (0.04 in), outer diameter 2 mm (0.08 in). Smooth tip allows for straight forward insertion and minimal tissue disturbance.



Qty: 5/Bx

Description	SKU #
Tapered from 1.25 mm (18G) to 0.9 mm (20G), oval	585149

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 40.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G

Hydrodissection Cannulas

Intended to separate the cortex from the capsule via irrigation (hydrodissection).

Nucleus Hydrodissector (Pearce)



Angled 35°, 8 mm from end.



585036

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585037
27G x 7/8 in, 0.4 x 22 mm with 8 mm bend	585036
27G x 7/8 in, 0.4 x 22 mm	585274

Nucleus Hydrodelineator (P. Koch)



Intended to break up lens contents and facilitate their aspiration via irrigation.



585015

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585015

Cortical Cleaving Hydrodissector (Flat End)



45°, 11 mm from end. For cleaving cortex from capsular bag.



585158

Qty: 10/Bx

Description	SKU #
27G 7/8 in, 0.4 x 22 mm	585158

Nucleus Hydrodissector (Sinskey)



35°, 8 mm from end.



585195

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585195

Nucleus Hydrodissector



45°, 7 mm from end.



581592

Qty: 10/Bx

Description	SKU #
30G x 7/8 in, 0.3 x 22 mm	581592

Nucleus Hydrodissector (Sauter)



Tip is flattened on the vertical plane.



585099

Qty: 10/Bx

Description	SKU #
27G x 7/8 in, 0.4 x 22 mm	585099

Hydrodissection Cannulas

Intended to separate the cortex (tissue) from the capsule via irrigation (hydrodissection).

Nucleus Hydrodissector

Tip is flattened on the horizontal plane.





585155 **Qty:** 10/Bx

Description	SKU #
27G 7/8 in, 0.4 x 22 mm	585155

Nucleus Hydrodissector (Chang)

90° bend, 1 mm flattened tip.





585294 **Qty:** 10/Bx

Description	SKU #
27G x 7/10 in, 0.4 x 18 mm, Bevel angle 33°	585294
27G x 7/10 in, 0.4 x 18 mm, Bevel angle 47°	585295

Tapered Hydrodelineator (Blumenthal)

Tip is tapered to 0.3 mm (30G), 45°, 4 mm from end.





585107 **Qty:** 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585107

Cortical Cleaving Hydrodissector (Round End)

45°, 11 mm from end.






585157 **Qty:** 10/Bx

Description	SKU #
27G x 7/8 in, 0.4 x 22 mm	585157

Nucleus Hydrodissector (Buratto)

Formed, 1 mm bent tip.





585199 **Qty:** 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585199

Cortical Cleaving Hydrodissector (Bolger)

Tip is tapered to 0.3 mm (30G), 33°, 1 mm from end.





585159 **Qty:** 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585159

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 40.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G

Hydrodissection Cannulas

Intended to separate the cortex (tissue) from the capsule via irrigation (hydrodissection).

Hydrodissection Cannula (Akahoshi Style)



Tapered tip, thin wall, angled 1 mm from end.



585198

Qty: 10/Bx

Description	SKU #
27G 7/8 in, 0.4 x 22 mm	585198

Micro J Shape Cannula



0.25 mm open end port.
Hook is 1.5 mm long x 1.75 mm wide.



585018

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585003
27G x 7/8 in, 0.4 x 22 mm	585018

Micro J Shape Cannula (Hessburg)



Hook is 0.85 mm long x 1.7 mm wide.



585148

Qty: 10/Bx

Description	SKU #
25G x 3/4 in, 0.5 x 20 mm	585148

Micro J Shape Hydrodissector (Pearce)



Hook is 2.25 mm long x 1.75 mm wide.



585029

Qty: 10/Bx

Description	SKU #
27G x 7/8 in, 0.4 x 22 mm	585029

Cortical Cleaving Hydrodissector (Grabow)



40°, 12 mm from end.



585197

Qty: 10/Bx

Description	SKU #
30G x 5/8 in, 0.3 x 16 mm	585197

HydroExpression Cannula (Corydon)



1.8 mm from end.



585197

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585090
27G x 7/8 in, 0.4 x 22 mm	585191

Cystotome Cannulas

Intended to perform a capsulorhexis. Indicated for use in cataract surgery.

Irrigating Cystotome





581610 **Qty:** 10/Bx

Description	SKU #
23G x 5/8 in, 0.6 x 16 mm	581599
25G x 5/8 in, 0.5 x 16 mm	581610
27G x 5/8 in, 0.4 x 16 mm	581618
30G x 5/8 in, 0.3 x 16 mm	581630

Irrigating Cystotome (Blumenthal)





585043 **Qty:** 10/Bx

Description	SKU #
27G x 5/8 in, 0.4 x 16 mm	585043

Irrigating Cystotome

Formed 8.5 mm from bend to tip. Shortened curve designed for deep set eyes or small pupils.





585008 **Qty:** 10/Bx

Description	SKU #
25G 5/8 in, 0.5 x 16 mm	585004
27G 5/8 in, 0.4 x 16 mm	585008

Irrigating Formed Cystotome

Tip 45°.





585049 **Qty:** 10/Bx

Description	SKU #
25G x 5/8 in, 0.5 x 16 mm	585049

Irrigating Cystotome





581594 **Qty:** 10/Bx

Description	SKU #
23G x 5/8 in, 0.6 x 16 mm	581594
25G x 5/8 in, 0.5 x 16 mm	581617

Irrigating Side Cutting Cystotome





581602 **Qty:** 10/Bx

Description	SKU #
25G x 5/8 in, 0.5 x 16 mm	581602

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 40.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G

Cystotome Cannulas

Intended to perform a capsulorhexis. Indicated for use in cataract surgery.

Irrigating Cystotome (Blumenthal)



Formed 8.5 mm from bend to tip.



585113

Qty: 10/Bx

Description	SKU #
25G x 5/8 in, 0.5 x 16 mm	585113
27G x 5/8 in, 0.4 x 16 mm	585042

Irrigating Cystotome (Simcoe)



45° and 10.5 mm from tip.



585093


Qty: 10/Bx

Description	SKU #
27G x 1/2 in, 0.4 x 12.5 mm	585093

Blunt Tip Irrigating Cystotome (Kellan)



Tip angled toward bevel. Unique blunt tip creates a wide capsular flap. Useful for pear shape capsulotomies.



585115

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585115

Irrigating Cystotome (Joyce)



Tip 0.5 mm long.



585030

Qty: 10/Bx

Description	SKU #
25G x 5/8 in, 0.5 x 16 mm	585030
27G x 5/8 in, 0.4 x 16 mm	585039

Irrigating Cystotome (Reverse)



Tip 90° toward bevel.



585034

Qty: 10/Bx

Description	SKU #
25G x 5/8 in, 0.5 x 16 mm	585034

Irrigation Cannulas

Intended to be used for irrigation/injection of liquid and/or air.

Anterior Chamber Cannula



30°, 6 mm from end.



581277

Qty: 10/Bx

Description	SKU #
◯ 16G x 1 in, 1.6 x 25 mm	581277

Bishop Harmon Cannula



Intended for irrigation of the eye.
40°, 8 mm from end.



585040

Qty: 10/Bx

Description	SKU #
19G x 1 in, 1.1 x 25 mm	585040

Irrigation and Aspiration Cannulas

Intended to inject ophthalmic viscosurgical devices (OVD).

Viscoflow Formed



Formed 11 mm from end, radius is 19 mm.



585007

Qty: 10/Bx

Description	SKU #
27G x 7/8 in, 0.4 x 22 mm	585007

Viscoflow Cannula



Angled 45°, 9 mm from end.



585009

Qty: 10/Bx

Description	SKU #
23G x 7/8 in, 0.6 x 22 mm	585009
25G x 7/8 in, 0.5 x 22 mm	585045
27G x 7/8 in, 0.4 x 22 mm	585026
◯* 30G x 7/8 in, 0.3 x 22 mm	585046

ISO Hub Gauge Chart
*BVI Gauge Color Chart Located on Page 40.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G

Irrigation and Aspiration Cannulas

Intended to aspirate the cortex.

Cortex Extractor

0.3 mm top port.





585075 **Qty:** 10/Bx

Description	SKU #
23G 5/8 in, 0.6 x 16 mm	585075

Cortex Extractor

0.4 mm top port, thin wall cannula, 30°, 2 mm from end, micro-etched 2 mm from end.





585227 **Qty:** 10/Bx

Description	SKU #
23G x 15/16 in, 0.6 x 23.5 mm	585227

Cortex Extractor

0.3 mm end port angled 10°.
Flat end opening has rounded edges.





585106 **Qty:** 10/Bx

Description	SKU #
24G x 1/2 in, 0.55 x 12.5 mm	585106

Bimanual Aspiration Cannula

Thin wall cannula, curved 45°, 0.51 mm open end irrigation port.





585222 **Qty:** 10/Bx

Description	SKU #
23G x 5/8 in, 0.6 x 17 mm	585222

Cortex Extractor

0.4 mm top port.






585193 **Qty:** 10/Bx

Description	SKU #
23G x 5/8 in, 0.6 x 16 mm	585193

Bimanual Aspiration Cannula

Thin wall cannula, curved 45°, 0.51 mm open end irrigation port.





585223 **Qty:** 10/Bx

Description	SKU #
21G x 5/8 in, 0.8 x 17 mm	585223

Irrigation and Aspiration Cannulas

Intended for irrigation of the eye and/or aspiration of liquids and body fluids.

I and A Cannula (Simcoe)



Tubing: Inner diameter 1 mm (0.04 in), outer diameter 2 mm (0.08 in), 0.3 mm aspiration port.



585047

Qty: Varies/Bx

Description	SKU #
23G, 0.6 mm (5/box)	585047
23G, 0.6 mm (3/box)	585086

I and A Cannula (Simcoe)



Tubing: Inner diameter 1 mm (0.04 in), outer diameter 2 mm (0.08 in), 0.4 mm aspiration port.



585110

Qty: 3/Bx

Description	SKU #
23G, 0.6 mm	585110

Polishers

Intended to polish the anterior capsule.

Flexible Tip Cannula



Angled 45°, 7 mm from end, 10 mm sleeve.



585041

Qty: 10/Bx

Description	SKU #
23G x 7/8 in, 0.6 x 22 mm	585041
25G x 7/8 in, 0.5 x 22 mm	585017
27G x 7/8 in, 0.4 x 22 mm	585005

Capsule Polisher (Kratz)



0.3 mm Top port, 35°, 9 mm from end.



581604

Qty: 10/Bx

Description	SKU #
21G x 7/8 in, 0.8 x 22 mm	581604
23G x 7/8 in, 0.6 x 22 mm	581605

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 40.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G

Polishers

Intended to polish the anterior capsule.

Olive Tip Polisher



0.33 mm end port, 30°, 10 mm from end and micro-etched 2 mm from end. Olive shape tip is 2 mm long x 1 mm wide.



585143

Qty: 10/Bx

Description	SKU #
23G x 1 in, 0.6 x 25 mm	585143

Posterior Scraper (Simcoe)



Intended to polish the anterior capsule. 0.33 mm open end port, 60°, 7 mm from end.



581645

Qty: 10/Bx

Description	SKU #
23G x 7/8 in, 0.6 x 22 mm	581645

Olive Tip Cannula



0.33 mm end port. Olive shape tip is 2 mm long x 1 mm wide.



585147

Qty: 10/Bx

Description	SKU #
23G x 7/8 in, 0.6 x 22 mm	585147

Glaucoma Revision Picks

Intended for bleb needling.

Glaucoma Revision Pick 60° (Lederer)



60°, 7 mm from end.
Ideal for superior approach.



585133

Qty: 5/Bx

Description	SKU #
27G x 5/8 in, 0.4 x 16 mm	585133

Glaucoma Revision Pick 30° (Lederer)



30°, 13.5 mm from end.
Ideal for temporal or nasal approach.



585135

Qty: 5/Bx

Description	SKU #
27G x 5/8 in, 0.4 x 16 mm	585135

Lacrimal Cannulas

Lacrimal Cannula

Intended to irrigate and probe the lacrimal system. Curved.



581615 **Qty:** 10/Bx

Description	SKU #
21G x 1-3/8 in, 0.8 x 35 mm	581615
26G x 1-1/8 in, 0.45 x 28 mm	581276

Lacrimal Cannula (Shahinian)

Intended to irrigate and probe the lacrimal system. 0.3 mm aspiration port.



585033 **Qty:** 5/Bx

Description	SKU #
23G x 1/2 in, 0.6 x 12 mm	585033

Lacrimal Cannula

Intended to irrigate and probe the lacrimal system. For irrigating punctum or lacrimal duct.



585016 **Qty:** 5/Bx

Description	SKU #
19G x 3/4 in (1.1 x 20 mm) shaft attached to 23G x 1/2 in (0.6 x 12 mm) malleable cannula	585016

Lacrimal Cannula

Straight, 13 mm Length.



585068 **Qty:** 10/Bx

Description	SKU #
25G x 1/2 in, 0.5 x 13 mm	585068

Lacrimal Trephine (Sisler)

Intended for rotary trepanation of an obstructed canaliculis and/or biopsy of the conjuctiva or sclera.



585031 **Qty:** 3/Bx

Description	SKU #
21G x 1-1/2 in, 0.8 x 38 mm	585031

ISO Hub Gauge Chart

*BVI Gauge Color Chart Located on Page 40.

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G

Refractive and FLACS

LASIK Irrigating Cannula (Lindstrom)



Intended to irrigate and/or manipulate corneal tissue. Features three 0.25 mm ports located on 7 mm extension with an 11 mm radius, 45°.



585250

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585250

LASIK Flap Spatula/Irrigator (Hersh)



Intended to irrigate and/or manipulate corneal tissue. 45°, 9 mm from spatulated tip.



585252

Qty: 10/Bx

Description	SKU #
23G x 7/8 in, 0.6 x 22 mm	585252

LASIK Irrigating Cannula (Solomon)



Intended to irrigate and/or manipulate corneal tissue. Thin wall cannula, angled tip features three 0.25 mm ports on 8 mm extension.



585280

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585280

FLACS Cannula (Schweitzer)



Intended to open clear corneal incisions created by the Femtosecond laser, injection of OVD (ophthalmic viscosurgical device) viscoelastic solution into the anterior chamber and to check that capsulorhexis is released for the capsule.



585200

Qty: 10/Bx

Description	SKU #
25G x 7/8 in, 0.5 x 22 mm	585200

LASIK Cannula



Intended to irrigate and/or manipulate corneal tissue. 25°, formed 7 mm from end.



585278

Qty: 10/Bx

Description	SKU #
27G x 7/8 in, 0.4 x 22 mm	585278

Hubs

BVI product hubs comply with ISO 6009 Color coding and ISO 594 Conical fittings with a 6% luer taper to assure use with syringes or devices.

COLOR	HUB	GAUGE	ID (IN)	ID (MM)	OD (IN)	OD (MM)
Natural* (Clear)		30G	0.006	0.15	0.012	0.3
Grey		27G	0.008	0.2	0.016	0.4
Brown		26G	0.01	0.25	0.018	0.45
Orange		25G	0.01	0.25	0.02	0.5
Purple		24G	0.012	0.3	0.022	0.55
Blue		23G	0.013	0.33	0.025	0.6
Black		22G	0.016	0.41	0.028	0.7
Green		21G	0.02	0.51	0.032	0.8
Yellow		20G	0.024	0.6	0.036	0.9
Cream		19G	0.028	0.7	0.043	1.1

*Not an ISO color code.


ISO Hub Gauge Chart

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G

High-Quality Single-Use Instruments

Providing Sterile Options for Your Operating Room Procedures

- Single-use instruments offer protection against infection, due to decreased risk of cross-contamination[1]
- Surgeons benefit from precision and reliability of brand-new instruments for each patient[2]





Available with contour clips. Contour clips create a round handle.

References

1 Southworth, P.M. Infections and exposures: reported incidents associated with unsuccessful decontamination of reusable surgical instruments. The Journal of Hospital Infection. September 17, 2014

2 Scheib, Christian; Single-use instruments deserve consideration for eye surgery: Advantages include cost savings, ease of use and high-quality instruments; December 1, 2016

Adjustable Speculums


Plastic Locking Eye Speculum



Qty: 10/Bx

Description	SKU #
Solid blade, Gear controlled spread, Overall length 83 mm	585297


Mechanical Speculum Solid Blade



Qty: 20/Bx

Description	SKU #
6 mm x 14.5 mm Solid curved polycarbonate blades, 0.7 mm Thick, 1 mm Wire, Polycarbonate body and screw mechanical speculum	MMSU1488S

Mechanical Speculum, Open Blade (PC)




Qty: 20/Bx

Description	SKU #
5 mm x 14 mm Open blades, 1 mm Wire, Polycarbonate body and screw mechanical speculum	MMSU1290S

Sprung Speculums

Kratz-Barraquer Speculum



Qty: 20/Bx

Description	SKU #
6.5 mm x 14 mm Open blades, 0.9 mm Wire, 30 mm Wide	MMSU1252S

Barraquer Speculum, Flat Blade



Qty: 20/Bx

Description	SKU #
6.5 mm x 14.5 mm Straight closed blades, 1 mm Wire, 30 mm Wide	MMSU1275S

Kratz-Barraquer Speculum, Temporal



Qty: 20/Bx

Description	SKU #
6.5 mm x 14 mm Open blades, 0.9 mm Wire, 30 mm Wide	MMSU1253S

Kratz-Barraquer Speculum, Curved




Qty: 20/Bx

Description	SKU #
6.5 mm x 14.5 mm Open curved blades, 0.9 mm Wire, 30 mm Wide	MMSU1287S

Sprung Speculums

Barraquer Speculum, Solid Curved (PC)



Qty: 20/Bx

Description	SKU #
6 mm x 14.5 mm Solid curved polycarbonate blades, 0.7 mm Thick, 0.9 mm Wire, 30 mm	MMSU1474S

Barraquer Speculum, Solid Blade (PC)



Qty: 20/Bx

Description	SKU #
7 mm x 14.5 mm Solid straight polycarbonate blades, 0.8 mm Thick, 0.9 mm Wire, 30 mm Wide	MMSU1475S

Barraquer Speculum, Solid 18 mm Blade (PC)



Qty: 20/Bx

Description	SKU #
6 mm x 14.5 mm Solid curved polycarbonate blades, 0.8 mm Thick, 1 mm Wire	MMSU1479S

Eye Speculum



Qty: 10/Bx

Description	SKU #
Open wire loops, 40 x 30 mm, 14 mm blades, 0.8 mm wire size	581413

Eye Speculum Open Wire Loops



Qty: 10/Bx

Description	SKU #
Open wire loops. 40 x 30 mm, 14 mm blades, 0.9 mm wire size	581443

Fixation Forceps - Toothed

Fine Toothed Forceps Short (Bonn)





Qty: 20/Bx

Description	SKU #
0.3 mm 1 x 2 Teeth, 5 mm Tying platforms, 52 mm Handle	MMSU1123S

Toothed Forceps (Bonn)





SKU# MMSU1124CS

Qty: 20/Bx

Description	SKU #
0.45 mm 1 x 2 Teeth, 5 mm Tying platforms, 80 mm Handle	MMSU1124S

Colibri Toothed Forceps, Short





Qty: 20/Bx

Description	SKU #
0.4 mm 1 x 2 Teeth, 5 mm Tying platforms, 52 mm Handle	MMSU1129S

Fixation Forceps - Notched

Fine Notched Forceps, Short (Hoskins)





Qty: 20/Bx

Description	SKU #
0.2 mm Notched tips, 5 mm Tying platforms, 52 mm Handle	MMSU1114S

Notched Forceps (Hoskins)





SKU# MMSU1115CS

Qty: 20/Bx

Description	SKU #
0.3 mm Notched tips, 5 mm Tying platforms, 80 mm Handle	MMSU1115S

Colibri Notched Forceps, Short





Qty: 20/Bx

Description	SKU #
0.2 mm Notched tips, 5 mm Tying platforms, 52 mm Handle	MMSU1209S

Capsulorhexis Forceps

MICS Capsulorhexis Forceps - Curved/Sharp





SKU# MMSU1567CS **Qty:** 20/Bx

Description	SKU #
0.4 mm x 0.3 mm Sharp tips, 12 mm Curved direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1567S

MICS Capsulorhexis Forceps - Straight/Sharp





SKU# MMSU1568CS **Qty:** 20/Bx

Description	SKU #
0.4 mm x 0.3 mm Sharp tips, 12 mm Straight direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1568S

MICS Capsulorhexis Forceps - Straight/Grasping





SKU# MMSU1571CS **Qty:** 20/Bx

Description	SKU #
0.4 mm x 0.3 mm Blunt tips, 12 mm Straight direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1571S


MICS Capsulorhexis Forceps - Curved/Grasping






SKU# MMSU1573CS **Qty:** 20/Bx

Description	SKU #
0.4 mm x 0.3 mm Blunt tips, 12 mm Curved direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1573S





General Forceps

Epilation Forceps



Qty: 20/Bx

Description	SKU #
5 mm x 7 mm Flat plane tips, 70 mm Handle	MMSU0085S

Jewellers Forceps, Fine Angled



Qty: 20/Bx

Description	SKU #
7 mm 45° angled pointed tips, 80 mm handle	MMSU0283S

Adson Forceps, Toothed



Qty: 20/Bx

Description	SKU #
1.35 mm x 1.8 mm 1 x 2 Teeth, 122 mm Long	MMSU0238S

Jewellers Forceps, Fine Straight



Qty: 20/Bx

Description	SKU #
7 mm Pointed tips, 80 mm Handle	MMSU0286S

Adson Forceps, Non-Toothed



Qty: 20/Bx

Description	SKU #
0.85 mm x 1 mm Tips, 14.5 mm Serrated jaws, 122 mm Long	MMSU0241S

St Martins Forceps



Qty: 20/Bx

Description	SKU #
0.9 mm x 0.26 mm 1 x 2 Teeth, 52 mm Handle	MMSU0279S

Moorfields Forceps



Qty: 20/Bx

Description	SKU #
0.6 mm tips, 11.5 mm Serrated jaws, 105 mm Long	MMSU0277S

Contact Lens Forceps



Qty: 20/Bx

Description	SKU #
0.6 mm x 0.4 mm Tips, 9 mm Curved Platforms, 70 mm Handle	MMSU1284S

IOL Forceps

Lens Loading Forceps



Qty: 20/Bx

Description	SKU #
0.6 mm x 8 mm Tips, 45° angled shaft, 80 mm Handle	MMSU1146S

IOL Removal Forceps



Qty: 20/Bx

Description	SKU #
0.1 mm x 0.1 mm Ridges along 4 mm of 15 mm shaft, Angled at 30°, 52 mm Handle	MMSU1460S

Tying Forceps

Straight Tying Forceps



SKU# MMSU1125CS

Qty: 20/Bx

Description	SKU #
0.35 mm Tip, 5 mm Tying platforms, 80 mm Handle	MMSU1125S

Angled Tying Forceps (McPherson)



Qty: 20/Bx

Description	SKU #
0.35 mm Tip, 5 mm Tying platforms, 8 mm Shaft angled 45°, 80 mm Handle	MMSU1127S

Curved Tying Forceps



Qty: 20/Bx

Description	SKU #
0.25 mm Tips, 5 mm Curved tying platforms, 80 mm Handle	MMSU1210S

Straight Tying Forceps, Short



Qty: 20/Bx

Description	SKU #
0.26 mm Tips, 5 mm Tying platforms, 52 mm Handle	MMSU1219S

Waqar Suture Removal Forceps



Qty: 20/Bx

Description	SKU #
0.3 mm x 0.65 mm Pick angled 90°, 5 mm Tying platforms, Curved 28 mm shaft, 80 mm Handle	MMSU1583S

Choppers

Claw Chopper



Qty: 20/Bx

Description	SKU #
Claw Chopper	MMSU1105S

Rosen Chopper



Qty: 20/Bx

Description	SKU #
Phaco Chopper	MMSU1106S


Phaco Chopper



Qty: 20/Bx

Description	SKU #
Phaco Chopper	MMSU1108S


Greens Chopper, Right Handed



Qty: 20/Bx

Description	SKU #
Greens Chopper, Right Handed	MMSU1119S


Koch Chopper, 1.4 mm



Qty: 20/Bx

Description	SKU #
Koch Chopper, 1.4 mm	MMSU1228S


Kuglen Hook, Angled



Qty: 20/Bx

Description	SKU #
Kuglen Hook, Angled	MMSU1433S

Butcher Rotator



Qty: 20/Bx

Description	SKU #
Butcher Rotator	MMSU1434S

Koch Chopper, 1.4 mm, Thin Shaft




Qty: 20/Bx

Description	SKU #
Koch Chopper, 1.4 mm, Thin Shaft	MMSU1440S

Choppers


Reverse Dialling Hook



Qty: 20/Bx

Description	SKU #
Reverse Dialling Hook	MMSU1480S

Neuhann Nucleus Divider



Qty: 20/Bx

Description	SKU #
Neuhann Nucleus Divider	MMSU1489S

Nucleus Chopper (P. Koch)



Qty: 3/Bx

Description	SKU #
45°, 6 mm from tip, tip angled 90°, 1.5 mm from end	585181

Double-Ended Choppers and Manipulators


Phaco Chopper/ Mushroom Manipulator



Qty: 20/Bx

Description	SKU #
Phaco Chopper/Mushroom Manipulator	MMSU1221S

Horizontal Chopper/ Kuglen Hook



Qty: 20/Bx

Description	SKU #
Horizontal Chopper/Kuglen Hook	MMSU1231S

Diablo Chopper/Flat Button Manipulator




Qty: 20/Bx

Description	SKU #
Diablo Chopper/Flat Button Manipulator	MMSU1239S

Manipulators

Govan Manipulator



Qty: 20/Bx

Description	SKU #
Govan Manipulator	MMSU1102S


Sinskey Hook



Qty: 20/Bx

Description	SKU #
Sinskey Hook, 0.2 mm	MMSU1103S


Drysdale Manipulator



Qty: 20/Bx

Description	SKU #
Drysdale Manipulator	MMSU1144S


Sinskey Hook, Straight



Qty: 20/Bx

Description	SKU #
Sinskey Hook, Straight	MMSU1147S


ICL Manipulator



Qty: 20/Bx

Description	SKU #
ICL Manipulator	MMSU1153S


Lebuisson Manipulator



Qty: 20/Bx

Description	SKU #
Lebuisson Manipulator	MMSU1458S

Mushroom Manipulator



Qty: 20/Bx

Description	SKU #
Mushroom Manipulator	MMSU1179S

Lester Hook



Qty: 20/Bx

Description	SKU #
Lester Hook	MMSU1207S

Manipulators


Kuglen Hook, Straight



Qty: 20/Bx

Description	SKU #
Kuglen Hook, Straight	MMSU1225S

Diablo Manipulator



Qty: 20/Bx

Description	SKU #
Diablo Manipulator	MMSU1428S

Micro-Manipulator, Straight



Qty: 10/Bx

Description	SKU #
Micro-Manipulator, Straight	581441


Sinskey Hook - Gen 2



Qty: 10/Bx

Description	SKU #
Hook Sinskey Disposable - Gen 2	581442

IOL Manipulator Blunt Tip



Qty: 3/Bx

Description	SKU #
Straight, Blunt tip, Fits 0.2 mm diameter holes and larger	585062

Nucleus Rotator



Qty: 3/Bx

Description	SKU #
45°, Blunt forked tip, 9 mm from end, For rotating nucleus during phacoemulsification	585077

Dissection Scissors

Stevens Tenotomy, Straight, Rounded





Qty: 20/Bx

Description	SKU #
1.5 mm x 1.5 mm Rounded tips, Straight blades, 110 mm Long	MMSU0214S

Iris Scissors, Straight, Pointed





Qty: 20/Bx

Description	SKU #
0.1 mm x 0.5 mm Pointed tips, 32.5 mm Straight blades, 120 mm Long	MMSU0249S

Iris Scissors, Curved, Pointed





Qty: 20/Bx

Description	SKU #
0.1 mm x 0.7 mm Pointed tips, 25 mm Curved blades, 115 mm Long	MMSU0244S

Iris Scissors, Straight, Rounded





Qty: 20/Bx

Description	SKU #
1 mm × 1.5 mm Rounded tips, 30 mm Straight blades, 115 mm Long	MMSU0299S

Dressing Scissors

Dressing Scissors, Straight, Pointed





Qty: 20/Bx

Description	SKU #
0.3 mm Pointed tips, 24.5 mm Straight blades, 102 mm Long, Plastic handle	MMSU0454S

Dressing Scissors, Straight, Rounded





Qty: 20/Bx

Description	SKU #
0.8 mm Rounded tips, 24.5 mm Straight blades, 102 mm Long, Plastic handle	MMSU0455S

Sprung Scissors

Vannas Scissors, Curved



SKU# MMSU1140CS

Qty: 20/Bx

Description	SKU #
Pointed tips on 8 mm curved blades, 1.5 mm when closed, 34 mm Handle	MMSU1140S

Vannas Scissors, Straight



SKU# MMSU1141CS

Qty: 20/Bx

Description	SKU #
Pointed tips on 8 mm straight blades, 1.5 mm when closed, 34 mm Handle	MMSU1141S

Westcott Scissors, Curved, Pointed



Qty: 20/Bx

Description	SKU #
0.5 mm x 0.6 mm Pointed tips, 17.5 mm Curved blades, 36 mm Long	MMSU1211S

Westcott Scissors, Straight Rounded



SKU# MMSU1212CS

Qty: 20/Bx

Description	SKU #
0.5 mm x 0.6 mm Rounded tips, 17 mm Straight blades, 120 mm Long	MMSU1212S

Westcott Scissors Curved Dolphin Nose



Qty: 20/Bx

Description	SKU #
0.5 mm x 0.6 mm Rounded tips, 22 mm Curved 'Dolphin Nose' blades, 120 mm Long	MMSU1213S

Westcott Scissors, Straight, Pointed



Qty: 20/Bx

Description	SKU #
Pointed tips, 17 mm Straight blades, 120 mm Long	MMSU1214S

IOL Cutter (Osher)



SKU# MMSU1459CS

Qty: 20/Bx

Description	SKU #
0.2 mm x 0.5 mm Blunt tips, 9.5 mm Serrated blades, 66 mm Handle	MMSU1459S

Dressing Forceps

Straight Rampley



Qty: 10/Bx

Description	SKU #
11 mm x 23 mm Serrated tips, 180 mm Long, Nylon plastic	MMSU0068S

Rampley Short



Qty: 20/Bx

Description	SKU #
17 mm x 24 mm Serrated tips, 125 mm Long, Nylon plastic	MMSU1578S

Needle Holders and Drivers

Corneal Needle-Holder



SKU# MMSU1110CS

Qty: 20/Bx

Description	SKU #
9 mm Fine curved jaws, 0.45 mm Tip, 106 mm Long	MMSU1110S

Troutman Needle-Holder



Qty: 20/Bx

Description	SKU #
4 mm Fine curved jaws, 0.55 mm Tip, 76 mm Flat handle	MMSU1243S

Castroviejo Needle Holder, Curved



Qty: 20/Bx

Description	SKU #
1.5 mm Tips, 12.5 mm Curved jaws with lock, 142 mm Long	MMSU3013S

Barraquer Needle-Holder



Qty: 20/Bx

Description	SKU #
0.42 mm x 0.45 mm Tips, 4 mm Fine curved jaws with 5 mm shaft, 138 mm Long	MMSU3015S

Refractive Surgery Instrumentation

Corneal Epithelium Spatula Flat





Qty: 20/Bx

Description	SKU #
5 mm Spatula with 0.9 mm flat blade, 0.8 mm x 1.7 mm Round bladed tip, 75 mm Handle	MMSU1113S

Y Shape Corneal Marker





Qty: 20/Bx

Description	SKU #
2 mm x 8 mm Marking blades, 75 mm Handle	MMSU1134S

LASIK Flap Separator/Elevator (Seibel)





Qty: 20/Bx

Description	SKU #
0.16 mm x 0.28 mm x 1 mm Tip, 20 mm Straight shaft, 0.5 mm x 10.5 mm pointed hook, 52 mm Handle	MMSU1172S

LASIK Flap Elevator Curved 8 mm





Qty: 20/Bx

Description	SKU #
0.2 mm × 0.15 mm Tip, Curved 8 mm tapered shaft at 90°, 75 mm Handle	MMSU1177S

LASIK Flap-Lifting Forceps





SKU# MMSU1169CS

Qty: 20/Bx

Description	SKU #
0.3 mm x 11.5 mm Tips, 45° angled shaft, 80 mm Handle	MMSU1169S

LASIK Flap Separator, Flattened Tip





Qty: 20/Bx

Description	SKU #
0.75 mm x 1.1 mm Flattened tip, 23 mm Straight shaft, 75 mm Handle	MMSU1176S

Corneal Epithelium Spatula Round





Qty: 20/Bx

Description	SKU #
2.65 mm x 5.5 mm x 0.35 mm Rounded blade, 75 mm Handle	MMSU1187S

LASIK Flap Separator and Elevator, Binder






Qty: 20/Bx

Description	SKU #
0.1 mm x 0.3 mm Rounded tip 0.7 mm long, 10 mm Curved 60° angled shaft, 75 mm Handle	MMSU1206S

Refractive Surgery Instrumentation


Epithelium Well



Qty: 20/Bx

Description	SKU #
3 mm x 8.5 mm ID well, Angled 45°, 75 mm Handle	MMSU1197S
3 mm x 9 mm ID well, Angled 45°, 75 mm Handle	MMSU1196S


LASIK Flap Separator/ Elevator Straight



Qty: 20/Bx

Description	SKU #
0.15 mm x 0.3 mm x 8 mm Straight separator at 90° to shaft, 0.15 mm x 0.3 mm x 1 mm Pocketing hook, 52 mm Handle	MMSU1295S

Reinstein Lenticule Separator (Hyperopic)



Qty: 20/Bx

Description	SKU #
0.22 mm x 2 mm Pocketing hook, 1.2 mm x 2 mm Spatula tip, Curved separator at 90° to shaft, 52 mm Handle	MMSU1497S


Spatula PRK



Qty: 5/Bx

Description	SKU #
Spatula tip and top are slightly sharpened to facilitate more aggressive scraping.	585213


LASIK Flap Separator/ Elevator Curved



Qty: 20/Bx

Description	SKU #
0.15 mm x 0.3 mm x 8 mm Curved separator at 90° to shaft, 0.15 mm x 0.3 mm x 1 mm Pocketing hook, 52 mm Handle	MMSU1294S

Reinstein Lenticule Separator



Qty: 20/Bx

Description	SKU #
0.22 mm x 2 mm Pocketing hook, 1.2 mm x 2 mm Spatula tip, Curved separator at 90° to shaft, 52 mm Handle	MMSU1297S

LASIK Flap Roller

Used after realignment. Acts as a squeegee to remove excess fluid from the interface; smooths flap and removes folds.



Qty: 5/Bx

Description	SKU #
14 mm Sleeve	585233

Lacrimal Intubation CIS

Guibor Canaliculus Intubation Set (CIS)



Includes gently curved, stainless steel probes with smooth, rounded tips. Unique stepped probe design enables smooth transition in pass-through intubation.



Qty: 5/Bx

Description	SKU #
Probes - 17.7 cm length, 0.64 mm diameter. **Biocompatible silicone tubing** - 29 cm length, 0.64 mm diameter	8590350

Canaliculus Intubation Set (CIS)



Includes malleable stainless steel stylette encased in low friction nylon sheath. The stylette handles assist in controlled manipulation during insertion and retraction. The rounded tip, nylon sheath offers a smooth transition to silicone tubing.



Qty: 5/Bx

Description	SKU #
Probes - 15.5 cm length, 0.51 mm diameter. **Nylon sheath** - 0.81 mm diameter. **Biocompatible silicone tubing** - 33 cm length, 0.64 mm diameter	8590450

Clamps and Curettes

Chalazion Clamp, Oval





Qty: 20/Bx

Description	SKU #
13 mm x 20 mm Oval profile, 10 mm x 17 mm Aperture, 90 mm Long	MMSU1149S

Chalazion Clamp, Round






Qty: 20/Bx

Description	SKU #
12 mm Round profile, 9 mm Round aperture, 90 mm Long	MMSU1150S

Chalazion Curette, Double, 1 mm/2 mm





Qty: 20/Bx

Description	SKU #
1 mm Curette, 1 mm Depth, 2 mm Curette, 1.35 mm Depth, 128 mm Long	MMSU1157S

Chalazion Curette, Double, 2 mm/3 mm






Qty: 20/Bx

Description	SKU #
2 mm Curette, 1.35 mm Depth, 3 mm Curette, 1.6 mm Depth, 128 mm Long	MMSU1158S

Markers and Gauges


Mendez Gauge



Qty: 20/Bx

Description	SKU #
2.5 mm Polycarbonate ring with 13 mm ID, 0-180° marks in increments of 10, 75 mm Handle	MMSU1152S


Pointed IVT Marker 3.5/4, Blue



Qty: 10/Bx

Description	SKU #
2 x 0.55 mm Tips with 3.5 mm spread, 2 x 0.65 mm Tips with 4 mm spread, Polycarbonate, 108 mm Long	MMSU1277S

Braunstein IVT Marker 3.5/4, Green



Qty: 20/Bx

Description	SKU #
2 mm x 1.5 mm Tips with 3.5 mm spread, 2 mm x 1.5 mm Tips with 4 mm spread, Polycarbonate, 108 mm Long	MMSU1276S

Double Ended IVT Guide Pad



Qty: 20/Bx

Description	SKU #
7 mm Pad with 3.5 guide, 0.8 mm Pad with 4 mm guide, 6 mm 45° angled shafts, 90 mm Long	MMSU1285S

Intravitreal Injection Guide




Qty: 20/Bx

Description	SKU #
Polycarbonate lash guard with stainless-steel swage guide tube, Triangular base plate with 3 stabilizing studs, 39 mm, 45° angled wishbone handle	MMSU1577S

Probes and Dilators

Punctum Dilator, Fine/Medium



Qty: 20/Bx

Description	SKU #
0.4 mm x 21.5 mm Tip, 0.6 mm × 18 mm Tip, 77 mm Handle	MMSU1181S

Punctum Dilator/Seeker



Qty: 20/Bx

Description	SKU #
0.45 mm x 21.5 mm Tip, 5 mm Pointed tip, 18 mm Straight shaft, 77 mm Handle	MMSU1193S


Golf Club Spud



Qty: 20/Bx

Description	SKU #
0.1 mm x 0.8 mm x 2 mm Oval profile, 18 mm Straight shaft, 75 mm Handle	MMSU1417S

Punctum Dilator/ Lacrimal Probe



Qty: 20/Bx

Description	SKU #
0.28 mm x 25 mm Dilator, 0.7 mm x 44.5 mm Probe, 70 mm Handle	MMSU1419S

Fixation Rings


Fixation Ring



Qty: 20/Bx

Description	SKU #
0.5 mm Thick 13 mm internal diameter ring, 1.5 mm x 0.6 mm Teeth at 3.5 mm intervals, 1.3 mm x 0.4 mm Teeth at 0.9 mm intervals, 75 mm Handle	MMSU1495S

Swivel Fixation Ring




Qty: 10/Bx

Description	SKU #
13 mm Diameter	372242

Depressors

Schocket Double-Ended Scleral Depressor



Qty: 20/Bx

Description	SKU #
2.7 mm x 6.5 mm T shaped tip, 5.5 mm x 2.2 mm Rounded tip, 68 mm Handle	MMSU1405S

Ball Tip Scleral Depressor



Qty: 20/Bx

Description	SKU #
3 mm Ball tip, 21.5 mm Straight shaft, 50 mm Handle	MMSU1408S

IOL Glides

IOL Glide (Sheets)

Guides IOL into posterior chamber for placement in the bag. Can also be used for Mini-Nuc technique to guide crystalline lens out of anterior chamber.

Qty: 50/Bx

Description	SKU #
0.6 mm thick	581033

Small Pupil IOL (Bobbitt)

Guides instruments into anterior chamber. Keeps iris flat to avoid possibility of prolapse or trauma. Acts as a channel to facilitate injection of fluid for hydrodissection.

Qty: 50/Bx

Description	SKU #
20 mm long x 2.75 mm wide, Tapered to 1.75 mm, 0.6 mm thick	581032

Surgical Glide (Fichman)

Phaco glide portion is designed to protect the iris from prolapse or trauma during phacoemulsification and acts as a guide for insertion of instruments. IOL glide portion assists in the implantation of posterior chamber IOLs.



Qty: 50/Bx

Description	SKU #
30 mm x 5 mm, **Phaco glide portion:** 10 mm long x 2.6 mm wide, tapered to 1.75 mm. **IOL glide portion:** 20 mm x 5 mm, 2.5 mm radius. 0.6 mm thick	581049

Small Pupil Expansion

Experience the Latest Innovation in Small Pupil Expansion with the I-Ring® Pupil Expander

- Designed to be used without any additional cornea incision
- Created to avoid any iris distortion and damage during surgery and post surgery
- Designed to have 360° expansion on the iris



I-Ring® Pupil Expander


I-Ring® Pupil Expander

The I-Ring® is designed to expand the iris to increase the surgeon's field of view and helps improve access to the lens while protecting the iris. Single-use sterile product.

- Designed to be gentle on iris and other intraocular tissue
- Created to support iris returning to natural shape post surgery
- Straight-forward insertion and removal

Description	SKU #
I-Ring® Pupil Expander, 6/Pk	587000
I-Ring® Pupil Expander, 1/Pk	587001

Positioning Holes
- Designed to ensure Sinskey hook cannot contact iris





I-Ring® Pupil Expander

Channels
- Designed to not distort or pinch iris when engaged or removed
- Straight-forward engagement and disengagement

Living Hinges
- Created for enhanced flexibility of ring
- Designed for reliable folding during insertion and removal

Sinskey Hook

Used to manipulate nucleus or intraocular lenses.





Qty: 10/Bx

Description	SKU #
Angled 45°, 10 mm from tip, Hook 0.85 mm tip to bend distance, Overall length 130 mm	581442

Micro Iris Hook

Flexible, gold wire hooks designed to retract the iris during surgical procedures.



Qty: 25/Bx

Description	SKU #
8 mm long x 0.08 mm diameter	585209

BVI

Fluid Control



bvimedical.com

PVA Eye Spears

Spear overall handle length is 7.16 cm or 2.82 in.

UltraClean™ LASIK PVA Eye Spears

Compressed shape, designed for fast expanding. Made from highly absorbent, virtually lint free PVA material. Designed specifically for rapid recovery of fluids during LASIK and other corneal procedures.



Qty: 100/Bx
Brand: UltraClean™

Description	SKU #
UltraClean™ Eye Spears, 5/Pouch, 20/Pouch per box	40810
UltraClean™ Eye Spears, 10/Pouch, 10/Pouch per box	40815

UltraClean™ LASIK PVA Expanded Eye Spears

Designed for instantaneous wicking, pre-expanded tip. Made from highly absorbent, virtually lint free PVA material. No need to pre-soak and fan.



Qty: Various/Bx
Brand: UltraClean™

Description	SKU #
UltraClean™ PVA Expanded LASIK Eye Spear, 40/Bx	40800

Visitec® LASIK PVA Eye Spears

Made from highly absorbent PVA material. Designed for optimal tissue manipulation and fluid absorption to facilitate flap alignment and adhesion.



Qty: 100/Bx
Brand: Visitec®

Description	SKU #
LASIK Eye Spear - Orange 100 count, 10 per pouch, 10 pouches per box	581761

Visitec® LASIK PVA Expanded Eye Spears

Made from highly absorbent PVA material. Pre-expanded for convenience. Designed to provide instantaneous wicking. No need to pre-soak and fan.



Qty: 40/Bx
Brand: Visitec®

Description	SKU #
LASIK Expanded Spear 40 count, 2 per pouch, 20 pouches per box	581709

Ultracell® PVA Eye Spears

Made from highly absorbent PVA material. Designed for rapid recovery of fluids during cataract and other general ophthalmic procedures.



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
PVA Eye Spear with dark blue handle, 100 count, 5 per pouch, 20 pouches per box	40400-8
PVA Eye Spear with dark blue handle, 200 count, 10 per pouch, 20 pouches per box	40405-8

Cellulose Eye Spears

Weck-Cel® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Designed for use in delicate, surgical areas and to maintain rigidity during wicking process. Light blue handle.



Qty: Various/Bx
Brand: Weck-Cel®

Description	SKU #
Weck-Cel® Spears, 6 Pack, 144/Bx	0008680
Weck-Cel® Spears, 10 Pack, 180/Bx	0008685

Sponge Points, Strips and Wicks

Weck-Cel® Sponge Points and Strips

Made from highly absorbent, natural, cellulose material.



Qty: Various/Bx
Brand: Weck-Cel®

Description	SKU #
Weck-Cel® Points, 240/Bx	0008661
Weck-Cel® Points and Strips, 24/Bx	0008660

Visitec® Eye Fluid Wick

Wicks absorb and direct fluid from surgical area into fluid collection bag.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Eye fluid wick, 15 cm Long x 0.4 cm Wide	581147

Corneal Light Shields

Visitec® Corneal Light Shield

PVA absorbent pad. Soft, porous structure designed to protect the retina from operating microscope light. Contains no cellulose. Upon moistening, shield becomes soft and conforms to shape of cornea. Designed to keep cornea moist.



Qty: 20/Bx

Description	SKU #
Visitec® Corneal Light Shield, 8 mm (0.3 in) diameter	581062


Full Body Drapes

Full body drapes feature wider widths and a three-sided channel bag for maximum fluid capacity.
Full body drapes can be utilized in all ophthalmic procedures.

Full Body Incise Drape

Qty: 10/Bx


Description		SKU #
Non-Woven Incise Drape with 2 Bags		581156
Overall size:	140 cm x 244 cm (55 in x 96 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Fluid collection bag:	25 cm x 15 cm (10 in x 6 in)	
Capacity:	2 x 500 mL	
Material Type:	Non-Woven Wetlaid	



Full Body Aperture Drape

Qty: 10/Bx

Description		SKU #
Full Body Aperture Drape with 2 Bags		581157
Overall size:	140 cm x 244 cm (55 in x 96 in)	
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)	
Capacity:	2 x 500 mL	
Material Type:	Non-Woven Wetlaid	




Full Body Incise Drape with Plastic Head and Channel Bag

Qty: 10/Bx

Description		SKU #
Full Body Incise Drape with Plastic Head and Channel Bag		581163
Overall size:	160 cm x 252 cm (63 in x 99 in)	
Incise Area:	13 cm x 15 cm (5 in x 6 in)	
Channel Bag:	62 cm x 40 cm (24 in x 16 in)	
Capacity:	1000 ml	
Material Type:	SMS 157 cm (61 in), Plastic 100 cm (39 in)	

*Without wick



Full Body Drapes


Non-Woven Full Body Incise Drape

Qty: 10/Bx

Description	SKU #
Non-Woven Full Body Incise Drape	581164

Overall size:	160 cm x 252 cm (60 in x 99 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Channel bag:	62 cm x 40 cm (24 in x 16 in)
Capacity:	1000 ml
Material Type:	Non-Woven Wetlaid

*Without wick



Half Body Drapes

Plastic Incise Drape with Fluid Collection Bag

Qty: 10/Bx

Description	SKU #
Plastic Incise Drape with Fluid Collection Bag	581149

Overall size:	100 cm x 127 cm (39 in x 50 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)
Capacity:	250 ml
Material Type:	Polyethylene Opaque Film



Plastic Aperture Drape with Fluid Collection Bag

Qty: 10/Bx

Description	SKU #
Plastic Aperture Drape with Fluid Collection Bag	581151

Overall size:	100 cm x 127 cm (39 in x 50 in)
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)
Capacity:	250 ml
Material Type:	Polyethylene Opaque Film




Half Body Drapes

Non-Woven Incise Drape with Fluid Collection Bag

Qty: 10/Bx


Description		SKU #
Non-Woven Incise Drape with Fluid Collection Bag		581150
Overall size:	100 cm x 127 cm (39 in x 50 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)	
Capacity:	250 ml	
Material Type:	Non-Woven Wetlaid	



Non-Woven Aperture Drape with Fluid Collection Bag

Qty: 10/Bx


Description		SKU #
Non-Woven Aperture Drape with Fluid Collection Bag		581152
Overall size:	100 cm x 127 cm (39 in x 50 in)	
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)	
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)	
Capacity:	250 ml	
Material Type:	Non-Woven Wetlaid	



Non-Woven Incise Drape with Channel Bag

Qty: 10/Bx

Description		SKU #
Non-Woven Incise Drape with Channel Bag		581153
Overall size:	140 cm x 127 cm (55 in x 50 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Fluid collection bag with channels:	26 cm x 35 cm (10 in x 14 in)	
Capacity:	500 ml	
Material Type:	Non-Woven Wetlaid	




Half Body Drapes

Non-Woven Aperture Drape with Channel Bag

Qty: 10/Bx

Description		SKU #
Non-Woven Aperture Drape with Channel Bag		581155
Overall size:	140 cm x 127 cm (55 in x 50 in)	
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)	
Fluid collection bag with channels:	26 cm x 35 cm (10 in x 14 in)	
Capacity:	500 ml	
Material Type:	Non-Woven Wetlaid	




Non-Woven Half Body Incise Drape* with Channel Bag

Qty: 10/Bx

Description		SKU #
Non-Woven Half Body Incise Drape* with Channel Bag		581166
Overall size:	160 cm x 180 cm (63 in x 71 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Channel bag:	62 cm x 40 cm (24 in x 16 in)	
Capacity:	1000 ml	
Material Type:	SMS	

*Without wick



Face Drapes

Ideal for refractive or office procedures.


Mini Pre-Split Incise Drape

Adhesive backing. **Qty:** 10/Bx

Description	SKU #
Mini Pre-Split Incise Drape	581154

Overall size:	40 cm x 40 cm (16 in x 16 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Pre-split:	5 cm (2 in)
Material Type:	SMS

*Without wick



Mini Plastic Incise Drape

Adhesive backing. **Qty:** 10/Bx

Description	SKU #
Mini Plastic Incise Drape	581159

Overall size:	40 cm x 40 cm (16 in x 16 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Material Type:	Polyethylene Opaque Film

*Without wick



Mini Plastic Aperture Drape

Adhesive backing. **Qty:** 10/Bx

Description	SKU #
Mini Plastic Aperture Drape	581160

Overall size:	40 cm x 40 cm (16 in x 16 in)
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)
Material Type:	Polyethylene Opaque Film

*Without wick




Mini Plastic Aperture Drape

Small incise provides complete coverage of lashes, lids and brows. Converts aperture to incise drape.

Qty: 10/Bx

Description	SKU #
Incise Drape	581161

Overall size:	10 cm x 20 cm (4 in x 8 in)
Adhesive area:	11 cm x 16 cm (4 in x 6 in)
Material Type:	Polyethylene Opaque Film



Face Drapes

Dual Refractive Drape

Qty: 10/Bx

Description		SKU #
Dual Refractive Drape		581162
Overall size:	40 cm x 26 cm (16 in x 10 in)	
Aperture area:	2.5 cm x 4 cm (1 in x 1.5 in)	
Material Type:	Polyethylene Opaque Film	

*Without wick



Eye Drains

LASIK Eye Drain (Chayet)

Designed to reduce stromal bed hydration for a consistent ablation and for nasally hinged flaps.


Qty: 20/Bx

Description	SKU #
LASIK Eye Drain, 10 mm	581098


Fluid Collection

Fluid Collection Bags

For procedures where collection and control of fluid is designed to keep surgeon, patient and floor dry. Reduces operating room cleanup and turnaround time. Bags have an adhesive backing to adhere to surgical drape.



581177



581148

Qty: 20/Bx

Description	SKU #
Fluid Collection Bag: 17 cm x 16 cm (6.75 in x 6.25 in), 250 ml capacity	581177
Fluid Collection Bag with Wick Bag: 17 cm x 16 cm (6.75 in x 6.25 in), 250 ml capacity Wick: Includes one, 0.4 cm x 15 cm (0.15625 in x 5.875 in)	581146
Mini Fluid Collection with Wick Bag: 9 cm x 16 cm (3.5 in x 6.25 in), 100 ml capacity Wick: Includes one, 0.4 cm x 15 cm (0.15625 in x 5.875 in)	581148

Instrument Wipes

Visitec® Instrument Wipes



Qty: 20/Bx

Description	SKU #
Instrument Wipe 76.2 mm x 76.2 mm (3 in x 3 in), 2 mm thick	90900

Visitec® Visiwipe Instrument Wipes

Qty: 20/Bx

Description	SKU #
Instrument Wipe 7.3 cm x 7.3 cm (2.875 in x 2.875 in)	581047
Adhesive Wipe 7.6 cm x 4.5 cm (3 in x 1.75 in)	581076

Other Fluid Control

ProOphta Eye Patch

Clear window dressing protects the eye while allowing the patient to see, permits examination without removing the dressing. Blinking is unimpaired. Indicated for use after cataract, corneal and ptosis operations; also for lagophthalmos and dry eye. Packaged non-sterile.



Qty: 50/Bx

Description	SKU #
ProOphta Eye Patch	581705

Absorbent Stick

6.6 cm overall length x 5 mm diameter (2.59 in x 0.1875 in).



Qty: 250/Bx

Description	SKU #
Absorbent Stick 6.3 cm, 10 per pouch, 25 pouches per box.	581060

Niteye® Bubble Bandage

Clear self-adhesive eye bandage for dry eye which protects eyes from drying conditions by sealing in valuable moisture.

Qty: 336/Case

Description	SKU #
Niteye® Bubble Bandage, Non-Sterile	90305

Quality Accessories to Support You and Your Patients' Needs



Eye Shields



Pen Lights




Eye Bandage

Eye Shields and Bandages

Eye Shield (Left/Right)

White lightweight plastic eye shield with ventilation holes. Flexible eye shield, shaped to fit the right or left eye for protection. Smooth edges do not require a cover. Packaged sterile.



Qty: 50/Bx

Description	SKU #
Eye Shield, Left/Right 8.5 cm (3.3 in)	581019

Universal Pediatric Eye Shield

Clear lightweight plastic and features smooth edges and ventilation holes. For use on left or right eye. Packaged sterile.



Qty: 25/Bx

Description	SKU #
Pediatric Eye Shield	581708

Universal Eye Shield

Clear lightweight plastic, for use on left or right eye, that features smooth edges and ventilation holes. Packaged sterile.



Qty: 50/Bx

Description	SKU #
Universal Eye Shield, 7.5 cm x 6.5 cm (3 in x 2.5 in)	581040

Guibor Eye Bandage

Disposable self-adhesive eye dressing designed to protect the eye during drying conditions, such as “dry eye”. Straightforward to apply. Packaged sterile.



Qty: 10/Bx

Description	SKU #
Eye Bandage, Guibor bubble bandage	8590300

Diagnostic Lights

Maxi-Lite® Diagnostic Penlight

Long-lasting disposable diagnostic penlight features twist-on backswitch.



Qty: 1/Ea

Description	SKU #
Diagnostic Penlight, Maxi-Lite®	8310201

Nurse-Mate® Six-Pack Diagnostic Penlights

Long-lasting disposable diagnostic penlight features convenient clip-activated switch and 2 mm to 9 mm pupil gauge. Imprinted with medical caduceus.




Qty: 6/Ea

Description	SKU #
Diagnostic Penlight, Nurse-Mate® Six-Pack	8310080

Blu-Spot® Diagnostic Light

Has diffuse cobalt light beam and twist-on backswitch.



Qty: 1/Bx

Description	SKU #
Diagnostic Penlight	8320150

General Accessories

Superflex Extension Set





Qty: 20/Bx

Description	SKU #
Superflex Extension Set	581007

Bulk High Flow Filter

Filters out particulates and microbiological contamination of 0.2 micron and larger.



Qty: 50/Bx

Description	SKU #
0.2 micron Filter	581030

Irrigation Handle

Standard male luer fitting.





Qty: 20/Bx

Description	SKU #
Length: 9.3 cm (3.5 in)	581012

Double Male Luer Adapter

Standard male luer fittings on each end for connecting phacoemulsification and I/A machine tubing to instruments.




Qty: 50/Bx

Description	SKU #
3.78 cm (1-1/2 in) Long	581056

Corneal Marking Pad

Elongated, flat pad allows easy transfer of ink onto markers.



Qty: 10/Bx

Description	SKU #
Gentian violet, Pre-inked, Ready to use	581057

BVI

Procedure Packs

Experience
the Evolution



Reimagine Your Operating Room



OPTIMIZE EFFICIENCY, ELEVATE SUSTAINABILITY

Efficiency
- Increase Ordering Efficiency
- Increase Delivery Efficiency
- Improve Transition Between Cases
- Increase OR Efficiency

Sustainability
- Reduce OR Waste
- Reduce Carbon Footprint

Optimization
Maximize a Procedure in a Simple Pack That is Customized to Your Specific Needs

bvimedical.com

CustomEyes®

Custom Procedure Packs

Knowledgeable Support

Reliable Quality

Personalized Packs

Simplified Surgical Prep

Easy Order Process



Experience
the Evolution

CustomEyes®

Knowledgeable Support

- Experienced ophthalmic sales team
- We understand your operating room
- Information at the fingertips

Reliable Quality

- Top-notch quality
- Access to BVI's wide range of single-use instruments

Personalized Packs


- Individualization to each of your surgeon's needs
- Quickly adapt to staff or technique changes

Simplified Surgical Prep

- 3D visualization of pack layers
- Pack assembly to match your operations


Easy Order Process


- Real-time quote and pricing
- Sample in 2 weeks




CustomEyes® Configurator

The CustomEyes® Configurator allows our team of ophthalmic specialists to easily create the right pack for you.







Simple Product Selection

Allows you to choose the product characteristics and select the best product combination for your needs.

Visual Pack Assembly

Optimize the layout of the pack by seeing a virtual 3D layering including order and orientation of each component.


Instant Feedback

Start building a pack together with your rep to get instant feedback and customize quickly.

Actualize

Receive a quote, order a sample and initiate production directly from the app.

Electrosurgery and Electrocautery



Cautery Devices

Accu-Temp® Cauteries provide a broad selection of tip styles for cauterization in a wide variety of ophthalmic procedures. All products are single-use and provided sterile.

Accu-Temp® High Temperature Cautery

Qty: 10/Bx
Brand: Accu-Temp®




8441500

Description	SKU #
1,150°C (2,102°F), 1/2 in (12.7 mm) Shaft with elongated vasector fine tip	8441500




8442000

Description	SKU #
1,150°C (2,102°F), 1/2 in (12.7 mm) Shaft with fine tip	8442000




8443000

Description	SKU #
1,150°C (2,102°F), 1/2 in (12.7 mm) Shaft with loop tip	8443000




8444000

Description	SKU #
1,150°C (2,102°F), 2 in (50.8 mm) Shaft with loop tip	8444000




8445000

Description	SKU #
1,150°C (2,102°F), 2 in (50.8 mm) Shaft with fine tip	8445000

Accu-Temp® Low Temperature Cautery

Qty: 10/Bx
Brand: Accu-Temp®




8441000

Description	SKU #
612°C (1,134°F), 1/2 in (12.7 mm) Shaft with fine tip	8441000

Cautery Devices

Accu-Temp® Variable Temperature Cautery - Dialed Controlled

Qty: 10/Bx
Brand: Accu-Temp®





Low Temperature 8440000

Description	SKU #
150° - 612°C (302° - 1,134°F), 1/2 in (12.7 mm) Shaft with fine tip	8440000





High Temperature 8446000

Description	SKU #
506° - 1,150°C (943° - 2,102°F), 1/2 in (12.7 mm) Shaft with fine tip	8446000

Accu-Temp® Reusable Handles

The reusable cautery handles are used with the single-use tips.

Qty: 10/Bx
Brand: Accu-Temp®



8461000

Description	SKU #
Reusable Variable Low Temp Handle, 150° - 612°C (302° - 1,134°F)	8461000



8461020


Description	SKU #
Reusable High Temp Handle, 1,150°C (2,102°F)	8461020

Accu-Temp® Single-Use Tips

Each single-use tip is packaged sterile and includes a protective plastic sheath that is designed to protect the reusable handles from contamination during use. Sterile and single-use.

Single-Use Fine Tip






Qty: 25/Bx
Brand: Accu-Temp®

Description	SKU #
Single-use fine tip	8461050

Single-Use Loop Tip





Qty: 25/Bx
Brand: Accu-Temp®

Description	SKU #
Single-use loop tip	8461060

Erasers

Wet-Field® Erasers intended for coagulation in vitreoretinal procedures, glaucoma filtration and cataract incision. Sterile and single-use.

Wet-Field® Hemostatic Eraser Bipolar Bevel Tip

Qty: 10/Bx
Brand: Wet-Field®





221250


Description	SKU #
18G Bevel tip, Straight	221250

Wet-Field® Hemostatic Eraser Bipolar Blunt Tip

10 per box unless otherwise noted.

Qty: 10/Bx
Brand: Wet-Field®





221251
221260135
221260150
221260

Description	SKU #
18G Blunt tip, Angled 45°	221251
20G Blunt tip, Angled 135°	221260135
20G Blunt tip, Angled 150°	221260150
20G Blunt tip, Straight	221260

Wet-Field® Hemostatic Eraser Bipolar Fine Tip

Qty: 5/Bx
Brand: Wet-Field®



221268
221268135
221268150

Description	SKU #
23G Fine tip, Straight	221268
23G Fine tip, Angled 150°	221268150
23G Fine tip, Angled 135°	221269135
25G Fine tip, Straight	221267
25G Fine tip, Angled 135°	221267135
25G Fine tip, Angled 150°	221267150

Erasers

Wet-Field® Hemostatic Eraser Bipolar 20G - 23G Tapered Tip

Qty: 5/Bx
Brand: Wet-Field®



Description	SKU #
20 - 23G Tapered fine tip, Straight	221265



Description	SKU #
20 - 23G Tapered blunt tip, Straight	221266

Wet-Field® Osher Eraser

Qty: 5/Bx
Brand: Wet-Field®



Description	SKU #
25G, 115°	221270

Cables

Qty: 10/Bx
Brand: Wet-Field®



22135010



22155010

Description	SKU #
Wet-Field® Disposable Cable	22135010
Wet-Field® Disposable Cable, Shrouded	22155010

Forceps

Wet-Field® Bipolar Forceps

Bipolar forceps are designed to deliver coagulation with no electrical arcing between the polished instrument tips for minimum coagulum adhesion. The insulated forcep blades are designed to reduce the possibility of electrical shock during operation. Non-sterile and reusable.

Qty: 1/Bx
Brand: Wet-Field®




Description	SKU #
Standard Forceps	221204




Description	SKU #
Straight Iris (Jeweler) Forceps	221206




Description	SKU #
Curved Iris (Jeweler) Forceps	221207




Description	SKU #
Adson Forceps	221208




Description	SKU #
Coaptation Forceps	221209

Providing Relief for Patients with Dry Eye or Lacrimal Disorders

Parasol® | Micro Flow™ Punctal Occlusion





Permanent Parasol®



Long-Term EXTEND 180®



Mid-Term EXTEND®



Partial Micro Flow™



Short-Term Collagen

bvimedical.com

Permanent Occlusion

Dry Eye products that offer predictable retention[1], straight-forward sizing and insertion. All occluder products are made of medical grade silicone and are reversible.

Parasol® Pre-Loaded Punctal Occluders Sterile

Occluders are packaged in a convenient one pair tray, are sterile and pre-loaded on a single-use inserter.



A14-201



Qty: 2/Bx

Parasol® Pre-Loaded Punctal Occluders Sterile	SKU #
X-Small, 0.25 - 0.35 mm	A14-200
Small, 0.35 - 0.65 mm	A14-201
Medium, 0.6 - 0.85 mm	A14-202
Large, 0.9 mm+	A14-203

Parasol® Quickload Punctal Occluders Sterile

Occluders are sterile and pre-mounted on a stand for Quickload convenience. Reusable inserters are sold separately.



A14-201QL

Qty: 16/Bx

Parasol® Quickload Punctal Occluders Sterile	SKU #
Small, 0.35 - 0.65 mm	A14-201QL
Medium, 0.6 - 0.85 mm	A14-202QL
Large, 0.9 mm+	A14-203QL
Combo small and medium	A14-205QL

Parasol® Punctal Occluders Non-Sterile

Occluders are non-sterile and packaged two occluders per poly bag. Reusable inserters are sold separately.



Qty: 20/Bx

Parasol® Punctal Occluders Non-Sterile	SKU #
Small, 0.35 - 0.65 mm	A14-201BN
Medium, 0.6 - 0.85 mm	A14-202BN
Large, 0.9 mm+	A14-203BN
Combo small and medium	A14-209BN

[1] C. McCabe, MD, PhD, FACSPlugs Reduce Dry-Eye Symptoms, Improve Vision. November 2009, Review of Ophthalmology.

Permanent Occlusion

Plus Pre-Loaded Punctal Occluders Sterile

The Plus design is the result of a micro-molding manufacturing process using a soft medical grade silicone.

- Sharply angled nose for insertion
- Delicate contoured edges
- Ultra-smooth surface finish
- Designed for consistent retention

Occluders are packaged in a convenient one pair tray. Each occluder is sterile and pre-loaded on a high quality disposable inserter.



A14-204



Qty: 2/Bx

Plus Pre-Loaded Punctal Occluders Sterile	SKU #
Small, 0.4 - 0.55 mm	A14-204
Medium, 0.55 - 0.7 mm	A14-205
Large, 0.7 - 0.85 mm	A14-206
X-Large, 0.85 mm+	A14-207

Plus Punctal Occluders Non-Sterile

Occluders are supplied non-sterile and have two occluders per poly bag, with 10 bags/pairs per box, unless otherwise noted. Reusable inserters are sold separately.



Qty: 20/Bx

Plus Punctal Occluders Non-Sterile	SKU #
Small, 0.4 - 0.55 mm	A14-204BN
Medium, 0.55 - 0.7 mm	A14-205BN
Large, 0.7 - 0.85 mm	A14-206BN
X-Large, 0.85 mm+	A14-207BN

Reusable Inserter for Non-Sterile Plugs



Qty: 2/Bx

Description	SKU #
Reusable inserter	A14-300

Mid-Term Occlusion

Extend® and Extend® 180 Synthetic Punctal Occluders

Provides occlusion for up to 3 months with Extend® or 6 months with Extend® 180, for patients experiencing short-term dry eye symptoms. These are ideal for the transient dry eye associated with LASIK or other ocular surgeries. Can also be used for allergy sufferers, contact lens discomfort, climate change or to sustain drug contact time and effectiveness.



A13-204

Qty: 20/Bx

Extend® Synthetic Punctal Occluders	SKU #
Variety Pack	A13-200
0.2 mm x 2 mm	A13-202
0.3 mm x 2 mm	A13-203
0.4 mm x 2 mm	A13-204
0.5 mm x 2 mm	A13-205



A13-2044

Qty: 20/Bx

Extend® 180 Synthetic Punctal Occluders	SKU #
Variety Pack	A13-2000
0.3 mm x 2 mm	A13-2033
0.3 mm x 2 mm, Single Pair	A13-2033SP
0.4 mm x 2 mm	A13-2044
0.4 mm x 2 mm, Single Pair	A13-2044SP
0.5 mm x 2 mm	A13-2055
0.5 mm x 2 mm, Single Pair	A13-2055SP

Partial Occlusion

Micro Flow™ Punctal Occluders

Designed to provide a straight-forward, reversible means of partial occlusion. Micro Flow™ Punctal Occluders support lacrimal drainage care. Suggest to use when total occlusion is not ideal.



A14-401

Qty: 2/Bx

Micro Flow™ Punctal Occluders Sterile	SKU #
Sterile, Small, 0.4 - 0.55 mm	A14-401
Sterile, Medium, 0.55 - 0.7 mm	A14-402
Sterile, Large, 0.7 - 0.85 mm	A14-403

Qty: 10/Bx

Micro Flow™ Punctal Occluders Non-sterile	SKU #
Non-Sterile, Small, 0.4 - 0.55 mm	A20-401BN
Non-Sterile, Medium, 0.55 - 0.7 mm	A20-402BN
Non-Sterile, Large, 0.7 - 0.85 mm	A20-403BN

Temporary Occlusion

Collagen Synthetic Punctal Occluders

These temporary occluders provide occlusion for 7-10 days. They are ideal for diagnostic use for determining if permanent punctal occlusion is needed for your patient. Six occluders in each foam sleeve, with 10 foam sleeves per box.




A12-103

Qty: 60/Bx

Collagen Synthetic Punctal Occluders	SKU #
0.2 mm x 2 mm	A12-102
0.3 mm x 2 mm	A12-103
0.4 mm x 2 mm	A12-104
Variety Pack	A12-200

Single-Use Punctal Dilators and Probes

Punctum Dilator/Lacrimal Probe






Qty: 20/Bx

Description	SKU #
0.28 mm x 25 mm Dilator, 0.7 mm x 44.5 mm Probe, 70 mm Handle	MMSU1419S

Other Related Occlusion Options

Lacrimal Cannula





Qty: 10/Bx

Description	SKU #
25G x 1/2 in, 0.5 x 13 mm	585068

Lacrimal Cannula





Qty: 10/Bx

Description	SKU #
26G x 1-1/8 in, 0.45 x 28 mm	581276

Punctum Dilator, Fine/Medium

0.4 mm x 21.5 mm Tip,
0.6 mm x 18 mm Tip,
77 mm Handle







Qty: 20/Bx

Description	SKU #
Punctum Dilator, Fine/Medium	MMSU1181S

Punctum Dilator/Seeker

0.45 mm x 21.5 mm Tip,
5 mm Pointed tip,
18 mm Straight shaft, 77 mm Handle







Qty: 20/Bx

Description	SKU #
Punctum Dilator/Seeker	MMSU1193S

Curved Tying Forceps

0.25 mm Tips, 5 mm Curved tying platforms,
80 mm Handle





Qty: 20/Bx

Description	SKU #
Curved Tying Forceps	MMSU1210S

Golf Club Spud





Qty: 20/Bx

Description	SKU #
0.1 mm x 0.8 mm x 2 mm Oval profile, 18 mm Straight shaft, 75 mm Handle	MMSU1417S

Leos®

Laser Endoscopy Ophthalmic System



Leos® Laser Endoscopy Ophthalmic System

Leos is a new surgical device that has been carefully designed to significantly enhance the ECP procedure for both the surgeon and operating room staff.

The portable industrial design of the Leos system along with the new, sterile and ready to use, high resolution laser endoscope offers the unique opportunity to efficiently and seamlessly integrate Leos ECP with the cataract procedure.




Description	SKU #
Leos Laser Endoscopy System Includes Footswitch, Literature Pack, Accessory Kit and Dust Cover. Qty: 1/Ea	OME6000U

VueProbe® Ophthalmic Laser Endoscope

VueProbe® Ophthalmic Laser Endoscope

The VueProbe® laser endoscope delivers a never-seen-before clear, digital, endoscopic view of the intraocular structures, providing consistent clarity and repeatability. The probe features new technological advancements in image orientation, automated focus and illumination adjustments.

Description	SKU #
VueProbe 19G Ophthalmic Laser Endoscope, straight, Single use, Sterile, 5 per box.	ELP19SU5
VueProbe 19G Ophthalmic Laser Endoscope, curved, Single use, Sterile, 5 per box.	ELP19CU5



ELP19SU5



ELP19CU5





Accessories

Leos® Footswitch

Wireless footswitch for laser and image orientation control.



Description	SKU #
Leos footswitch, Qty: 1/Ea	FS6000

Leos Dust Cover

Leos System Vinyl Cover

Description	SKU #
Leos dust cover, 1/Ea	DC6000U

Leos Accessory Kit

Kit includes :
- Power cable
- HDMI cable
- Remote Interlock plug
- Laser key
- USB drive

Description	SKU #
Leos accessory kit, Qty: 1/Ea	AK6000U

Leos Literature Pack

Pack includes:
- Operator's Manual
- Quick Reference Guide
- Laminated OR Laser Sign

Description	SKU #
Leos literature pack, Qty: 1/Ea	LIT6000U

Diode Laser Safety Goggles

810nm diode laser safety goggles for use with/without eyeglasses.



Description	SKU #
Diode laser safety goggles, Qty: 1/Ea	LG001

Sterile Drape


Custom drape for Leos control screen to facilitate a sterile procedure.



Description	SKU #
Sterile drape, Qty: 5/Sp	SD6000U


Addressing Control of Light Beam Direction

29G Spotlight Directional Chandelier



Enabling Bimanual Surgery With:

- Directional control of light-beam for enhanced visualization
- 29G valved entry system for straightforward light fiber entrance, fixation and removal
- Stationary and wide field endo-illumination



Patent Pending

Vitreq
Vitreoretinal Surgical Products

bvimedical.com

29G Spotlight Directional Chandelier


29G Spotlight Directional Chandelier

Sets new standards in light beam direction control, it delivers stationary, diffuse wide-field endo-illumination enabling bimanual surgery.



Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
29G Spotlight Directional Chandelier, Sterile	LF29.D02



Patent Pending

Improved Intraocular Surgical Control and Precision in Surgical Maneuvers

Micro Forceps and Scissors




Tip Design

- Shortened tip provides stronger grasping force[1]

Shaft Design

- Rigid shaft designed to reduce flexibility for better control
- 25% less flexible than competitors while maintaining 32 mm working shaft length[1]



Ergonomically Designed Handle

- Synchronized handle squeeze closes the tip to facilitate grasping and cutting
- Round handle allows for tip alignment

[1] Data on file

Vitreo
Vitreoretinal Surgical Products


bvimedical.com

Micro Forceps

Gauge Size with Color Coding


Single-Use

All Single-Use Sold in 5/Bx



Reusable


All Reusable Sold in 1/Bx



Serrated Gripping



Serrated Gripping Micro Forceps are designed for grasping strong adherent PVR membranes.



Brand: Vitreq®

Description	SKU #
23G Single-Use	FC23.D02
25G Single-Use	FC25.D02
23G Reusable	FC23.R02
25G Reusable	FC25.R02

ILM End-gripping



ILM end-gripping forceps have thin tips for ILM and ERM peeling. Designed to improve visualization for grasping thin membranes.




Brand: Vitreq®

Description	SKU #
23G Single-Use	FC23.D03
25G Single-Use	FC25.D03
27G Single-Use	FC27.D03
23G Reusable	FC23.R03
25G Reusable	FC25.R03
27G Reusable	FC27.R03

Asymmetrical Gripping



Asymmetrical tips allow visualization at grasping area and interior end-gripping platforms are designed for assured grasping. The tip design targets ILM, fine epiretinal membranes and diabetic PVR membranes.



Brand: Vitreq®

Description	SKU #
23G Single-Use	FC23.D04
25G Single-Use	FC25.D04
23G Reusable	FC23.R04
25G Reusable	FC25.R04

Micro Forceps

Eckardt Endgripping



Eckardt forceps with wider tip are designed to enhance surface grasping and reduce tearing effect. Shorter tip branches designed to enhance grasping and holding power. Designed for ILM, ERM peeling and general peeling of fine membranes which needs end-gripping action.



Brand: Vitreq®

Description	SKU #
23G Single-Use	FC23.D05
25G Single-Use	FC25.D05
27G Single-Use	FC27.D05
23G Reusable	FC23.R05
25G Reusable	FC25.R05
27G Reusable	FC27.R05

Shah X-tra Grip*



Fine tipped forceps with serrated jaws platform for maximum grasping of the tissue. Effective for all types of tissue manipulation, ILM, ERM, PVR and diabetic tissue.



Brand: Vitreq®

Description	SKU #
23G Single-Use	FC23.D07
25G Single-Use	FC25.D07
27G Single-Use	FC27.D07
23G Reusable	FC23.R07
25G Reusable	FC25.R07
27G Reusable	FC27.R07

Kuhn 45° Power Gripping*



Kuhn 45° Power Gripping forceps are designed for peeling proliferative intravitreal, epi and sub-retinal membranes, lens capsule, iris, small IOFB, fibrinous membranes, as well as ILM and ERM membranes and posterior hyaloids. With a 45° crocodile designed inner surface, these instruments are designed to enable selective, strong grasping, and visibility of the target tissue.



Brand: Vitreq®

Description	SKU #
23G Single-Use	FC23.D06
25G Single-Use	FC25.D06
23G Reusable	FC23.R06
25G Reusable	FC25.R06

*Patent Pending

Micro Scissors

Vertical



Sharp-tipped vertical cutting scissors for cutting perpendicular to the retina.



Brand: Vitreq®

Description	SKU #
23G Single-Use	SC23.D01
25G Single-Use	SC25.D01
27G Single-Use	SC27.D01
23G Reusable	SC23.R01

Straight



Sharp-tipped straight cutting scissors for cutting in the iris plane and cutting peripheral membranes.



Brand: Vitreq®

Description	SKU #
23G Single-Use	SC23.D02
25G Single-Use	SC25.D02
23G Reusable	SC23.R02
25G Reusable	SC25.R02

Horizontal Curved



Sharp-tipped horizontal curved cutting scissors for cutting membranes adherent to the retinal surface. The tip's curvature of 12mm follows the contour of the retina.



Brand: Vitreq®

Description	SKU #
23G Single-Use	SC23.D03
25G Single-Use	SC25.D03
27G Single-Use	SC27.D03
23G Reusable	SC23.R03
25G Reusable	SC25.R03
27G Reusable	SC27.R03

Reusable Cleaning Adapter


Unique cleaning/sterilization adapter included with every reusable instrument option.



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Reusable Cleaning Adapter	FC00.R03

90° Directional Laser Probes



A New Vision on Flexibility

- 90° curved extendable tip for access to peripheral retina
- Hardened metal shaft reduces flexibility by 25%[1]
- Retractable mechanism designed to improve control of instrument tip and reduce risk of accidental removal of cannula
- Available in standard (23G, 25G and 27G) and illuminated (23G and 25G) color coded for easy gauge identification
- Compatible with all major equipment

[1] Data on file


Vitreq
Vitreoretinal Surgical Products

bvimedical.com

Vitreoretinal Surgical Products

90° Directional Laser Probes Designed For Improved Surgical Control. The ergonomic handle was designed for improved intraocular surgical control and the sliding knob to control the motion of the instrument tip.

Standard Directional*





Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
23G Compatible with Alcon®	LP23.A04
23G Compatible with Iridex® and B&L	LP23.I04
25G Compatible with Alcon®	LP25.A04
25G Compatible with Iridex® and B&L	LP25.I04
27G Compatible with Alcon®	LP27.A04
27G Compatible with Iridex® and B&L	LP27.I04

Illuminated Directional*






Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
23G Compatible with Alcon®	LP23.A05
23G Compatible with Iridex® and B&L	LP23.I05
25G Compatible with Alcon®	LP25.A05
25G Compatible with Iridex® and B&L	LP25.I05


* All product names, trademarks and registered trademarks are property of their respective owners.



Backflush Instruments

Atraumatic Aspiration and Brushing Procedures Thanks to Improved Intraocular Surgical Control

- Passive and active aspiration control in one instrument
- High volume backflush reservoir for increased capacity during vitreoretinal surgery
- Long shaft (33 mm) makes them suitable for myopic eyes
- Blunt tip and silicone brush needles for atraumatic aspiration and brushing of retinal surface



Vitreq
Vitreoretinal Surgical Products

bvimedical.com

Disposable Backflush Instruments

Disposable backflush instruments are available in three different styles.

Standard

Blunt or brush tip



23G Backflush Instruments Brush
BF23.D02



Blunt tip



Brush tip

Retractable

Silicone brush retracts into the instrument shaft for smooth and valved cannula entry






23G Retractable Brush
BF23.D05

Extendable

8 mm extendable silicone tip for draining fluid from behind retina and other standard backflush procedures





23G Backflush Instruments Extendable Silicone Tip
BF23.D07

Disposable backflush instruments are shipped with cap and tubing for passive or active use during procedures. The direct connection between reservoir and needle allows no contact between hand piece and intraocular fluids.

Backflush Options

All single-use products are provided Sterile and Ready-for-Use.

Standard



Qty: 5/Bx
BF23.D02 **Brand:** Vitreq®

Description	SKU #
23G Blunt backflush instrument	BF23.D01
23G Brush backflush instrument	BF23.D02
25G Blunt backflush instrument	BF25.D01
25G Brush backflush instrument	BF25.D02
27G Blunt backflush instrument	BF27.D01
27G Brush backflush instrument	BF27.D02

Retractable



Qty: 5/Bx
BF23.D05 **Brand:** Vitreq®

Description	SKU #
23G Retractable brush BF-instrument	BF23.D05
25G Retractable brush BF-instrument	BF25.D05

Extendable



Qty: 5/Bx
BF23.D07 **Brand:** Vitreq®

Description	SKU #
23G Extendable silicone tipped BF-instrument	BF23.D07
25G Extendable silicone tipped BF-instrument	BF25.D07



Reusable Backflush Instruments

Featuring an ergonomically designed handle, rigid shaft and high volume backflush reservoir for increased backflush capacity, with Passive and Active Aspiration.

Universal Backflush Instrument



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Universal Backflush Instrument	BF00.R01

Backflush Instrument with Luer-lock Connector



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Universal Backflush Instrument	BF00.R02

Backflush Consumables

Disposable Universal Backflush Needles



BF25.D51

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
23G BF Blunt needles	BF23.D51
23G BF Brush needles	BF23.D52
25G BF Blunt needles	BF25.D51
25G BF Brush needles	BF25.D52

Blunt Needles with Luer-Lock

Stainless steel cannula with blunt tip for aspiration or injection with syringe or tubing.



CN23.D05

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
23G Blunt needle with luer-lock connector	CN23.D05
25G Blunt needle with luer-lock connector	CN25.D05
27G Blunt needle with luer-lock connector	CN27.D05

Backflush Reservoirs



Passive Reservoir
BF00.D51



Active Reservoir
BF00.D52

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
Backflush reservoirs, Passive	BF00.D51
Backflush reservoirs, Active	BF00.D52

Disposable Tubing for Active Aspiration




Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
Disposable tubing for active aspiration	BF00.D53

Light Fibers

A Bright Vision on Illumination



Standard Light Fibers

45° Shielded Light Fibers

Light Fibers

Light Fibers featuring high-end specifications with consistent quality performance for ophthalmic endo-illumination procedures.

Standard Light Fibers






LF23.D01

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
23G Standard light fibers	LF23.D01
25G Standard light fibers	LF25.D01
27G Standard light fibers	LF27.D01



Straight tip for focal endo-illumination.

45° Shielded Light Fibers






LF27.D03

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
23G 45° shielded light fibers	LF23.D03
25G 45° shielded light fibers	LF25.D03
27G 45° shielded light fibers	LF27.D03



45° shielded wide view tip for shielded endo-illumination.

Vitreoretinal Consumables



Vitreq
Vitreoretinal Surgical Products

bvimedical.com

Cannulas and Needles

New disposable vitreoretinal cannulas and needles include innovative design features to optimize difficult surgical maneuvers. The products are designed to allow the VR surgeon control over some of the most challenging VR procedures. For example, the range includes PFC needles that facilitate injection with minimized pressing force and have a central bore that enables a high volume of heavy liquids to enter the eye at a minimal velocity.

Vitreoretinal Cannula



Attaches to handpiece or syringe with luer fitting. Can also be used to aspirate or inject silicone oil.



585205 **Qty:** 10/Bx

Description	SKU #
20G x 1-1/8 in, 0.9 x 28 mm	585205

High Viscosity Injector (Hammer)



Thinwall cannula.



585184 **Qty:** 5/Bx

Description	SKU #
20G (0.9 mm), 6 mm, 10 cm Tube	585184
20G (0.9 mm), 6 mm, 16 cm Tube	585228

High Viscosity Injector (Hammer)



Thinwall cannula.



585173 **Qty:** 5/Bx

Description	SKU #
20G (0.9 mm), 4 mm, 10 cm Tube	585173
20G (0.9 mm), 4 mm, 16 cm Tube	585208

Soft Tip Cannula



Soft tip cannulas have flexible, high-flow silicone tips, allowing gentle access to the retina. They are ideal for aspiration of fluid from on or beneath the surface of the retina.



587002 **Qty:** 10/Bx

Description	SKU #
23G x 1-1/3 in, 0.6 x 33 mm, 3 mm tip	585226
25G x 1-1/3 in, 0.5 x 33 mm, 3 mm tip	585229
23G x 1-1/3 in, 0.6 x 33 mm, 0.8 mm tip	587002
25G x 1-1/3 in, 0.5 x 33 mm, 0.8 mm tip	587003

Subretinal Fluid Cannula (Grizzard)



Straight. Sleeve extends 6 mm beyond end.



Qty: 5/Bx

Description	SKU #
20G, 0.9 mm	585021

Cannulas and Needles

VFI Cannula



Silicone oil injection and removal through entry system. High flow rates due to thin walls. Can be directly connected to the silicone oil syringe or to aspiration device.



CN27.D03

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
27G VFI cannula	CN27.D03

Silicone Tip Cannula



Stainless steel cannula with 2 mm silicone tip for aspiration and manipulation.



CN25.D01

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
23G Silicone tipped cannula	CN23.D01
25G Silicone tipped cannula	CN25.D01

PFC* Injection Needles



New, optimized, PFC injection needle facilitates injection with minimized pressing force. The central bore permits a high volume of heavy liquid to enter the eye at a minimal velocity under controlled conditions.



MD23.D01

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
23G PFC injection needle	MD23.D01
25G PFC injection needle	MD25.D01

* Perfluorocarbon

Pick Needles



Stainless steel cannula with pick tip. Can be connected to a syringe or other handpiece.




CN25.D04

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
23G Pick needle	CN23.D04
25G Pick needle	CN25.D04

Disposable Universal Backflush Needles





BF25.D51

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
23G BF Blunt needles	BF23.D51
23G BF Brush needles	BF23.D52
25G BF Blunt needles	BF25.D51
25G BF Brush needles	BF25.D52

Blunt Needles with Luer-Lock



Stainless steel cannula with blunt tip for aspiration or injection with syringe or tubing.



CN23.D05

Qty: 5/Bx
Brand: Vitreq®

Description	SKU #
23G Blunt needle with luer-lock connector	CN23.D05
25G Blunt needle with luer-lock connector	CN25.D05
27G Blunt needle with luer-lock connector	CN27.D05

Retractable Tissue Manipulators

Retractable Tissue Manipulators



Available in 23G and 25G.
Disposable and ready to use.

FEATURES

- Retractable silicone tip with special ribbed structure
- The green colored silicone tip allows for a clear differentiation to the red retinal surface



Qty: 5/Bx
Brand: Vitreq®

TM23.D01

Description	SKU #
23G Retractable tissue manipulator	TM23.D01
25G Retractable tissue manipulator	TM25.D01

Squeezers

Unique pencil style injection device that affords surgeons precise control for injecting stains and perfluorocarbon liquids.



The 'Squeezer' gives surgeons perfect control for injecting stains and perfluorocarbon liquids using a single-handed technique. The Squeezer injection grip is the same as the familiar and safe pencil grip that is used in working with forceps, backflush and other surgical instruments, as well as redundant syringe assisted injection with potential uncontrolled injection risks.

Using the Squeezer, surgeons can control the start and speed of injection, as well as the amount of liquid dispensed to avoid overfill and a jet stream effect. The surgeon can easily guide the needle of the Squeezer through the entry system, get close to the retina, and precisely stain only the target tissue with complete control over the injection speed.

The new Squeezer has been developed in collaboration with Prof. Claus Eckardt.

Dye Squeezers



Qty: 5/Bx
Brand: Vitreq®


Patent: P6060106EP

Description	SKU #
2 ml reservoir, Sterile	SQ00.D01

Light Source Adaptors

The Light Fibers are compatible with most endo-illumination systems in combination with Light Source Adaptors*.


Constellation® Light Source Adaptor



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Constellation® Light Source Adaptor	LF00.R01

Ruck® Light Source Adaptor



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Ruck® Light Source Adaptor	LF00.R04


Stellaris® Light Source Adaptor



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Stellaris® Light Source Adaptor	LF00.R02

DORC® Light Source Adaptor



Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
DORC® Light Source Adaptor	LF00.R05

Geuder® Light Source Adaptor




Qty: 1/Bx
Brand: Vitreq®

Description	SKU #
Geuder® Light Source Adaptor	LF00.R03

* All product names, trademarks and registered trademarks are property of their respective owners.

Accessories

Williamson Scleral Depressor



The Williamson scleral depressor is specially designed to be used for VR-surgery in combination with small gauge entry and BIOM® systems. First scleral indenter specially designed to be used in combination with entry systems with particular attention to sizing and angulation of the indenting surfaces.

Qty: 1/Bx
Brand: Vitreq®



Description	SKU #
3 mm Ball tip, 45° bent tip	VS00.R02

Surgical Fluids for Ophthalmic Procedures



bvimedical.com

Surgical Fluids for Ophthalmic Procedures

Cohesive OVD

NuVisc™ is designed to create and maintain anterior chamber depth and visibility, protect corneal endothelial cells and other intraocular tissues and minimize interaction between tissues during surgical manipulation. NuVisc™ is to be used during ophthalmic surgeries.

NuVisc™ OVD 1.2% Sodium Hyaluronate

Qty: 1/Bx
Brand: NuVisc™

Description	SKU #
NuVisc™ OVD 1.2% sodium hyaluronate, includes 27G Viscoflow® cannula	585305

NuVisc™ Specifications	
Shelf Life	24 Months
Storage Temperature	2-8° C, (36-46° F), do not freeze
Volume	0.8 ml

Chemical Analysis Test Results	
1.6883 in	6.8 - 7.6
Zero shear viscosity	105,000 cps (105 Pa s)
NaHA concentration	12 mg/ml
Molecular weight	Greater than 1 million daltons
Osmolality	320 milliosmoles






FINEVISION HP

Hydrophobic Trifocal






FINEVISIONHP.bvimedical.com

Trifocal Hydrophobic

FINEVISION HP

TRIFOCAL OPTIC





Description

Model	POD F GF	
Material	GFY Hydrophobic Acrylic[1]	
Overall diameter	0.45 inches (11.40mm)	
Optic diameter	0.24 inches (6.00mm)	
Optic	Biconvex Aspheric Trifocal	
Haptic design	Double C-loop with Ridgetech® & Posterior 5° Angulated Haptic	
Filtration	UV & Selective Blue Light Filter (~400-450nm)	
Refractive index	1.53	
Abbe number	42	
Additional power (IOL plane)	+1.75D & +3.50D	
Injection system	Medicel Accuject 2.1	
Spherical power	+10D to +30D (0.5D steps)	
Suggested A constant[2]		**Interferometry**
	SRK/T: A	119.5

INJECTION SYSTEM

Medicel Accuject / Viscoject



Optic	Material	Brand	Commercial Name	Accuject 2.0	Accuject 2.1	Accuject 2.2
TRIFOCAL	HYDROPHOBIC	FINEVISION HP TRIFOCAL OPTIC	POD F GF	≤ 24.5D	≤ 35D	≤ 35D

[1] GFY® is patented since 2010.

[2] Values estimated only: surgeons are recommended to personalize their A-constant based on their surgical techniques and equipment, experience with the lens model and postoperative results.

The FINEVISION HP Trifocal IOL (PODF GF) is an FDA-approved device.

Product Index

SKU#	Description	PG
0008660	Weck-Cel® Sponge Points and Strips	65
0008661	Weck-Cel® Sponge Points and Strips	65
0008680	Weck-Cel® Cellulose Eye Spears	65
0008685	Weck-Cel® Cellulose Eye Spears	65
0009701	5 mm Trephine Blade	18
0009703	5.5 mm Trephine Blade	18
0009705	6 mm Trephine Blade	18
0009706	6.25 mm Trephine Blade	18
0009707	6.5 mm Trephine Blade	18
0009708	6.75 mm Trephine Blade	18
0009709	7 mm Trephine Blade	18
0009710	7.85 mm Trephine Blade	18
0009711	7.5 mm Trephine Blade	18
0009713	8 mm Trephine Blade	18
0009715	8.5 mm Trephine Blade	18
0009717	9 mm Trephine Blade	18
0009719	7.25 mm Trephine Blade	18
0009721	10 mm Trephine Blade	18
0009722	7.75 mm Trephine Blade	18
0009723	8.75 mm Trephine Blade	18
0009725	11 mm Trephine Blade	18
0009727	8.25 mm Trephine Blade	18
0009728	14 mm Trephine Blade	18
40400-8	Ultracell® PVA Eye Spears	64
40405-8	Ultracell® PVA Eye Spears	64
40800	UltraClean™ LASIK PVA Expanded Eye Spears	64
40810	UltraClean™ LASIK PVA Eye Spears	64
40815	UltraClean™ LASIK PVA Eye Spears	64
90900	Visitec® Instrument Wipe	72
90305	Niteye® Bubble Bandage	72
221204	Standard Forceps	86
221206	Straight Iris (Jeweler) Forceps	86
221207	Curved Iris (Jeweler) Forceps	86
221208	Adson Forceps	86
221209	Coaptation Forceps	86
221250	Wet-Field® Hemostatic Eraser Bipolar Bevel Tip	84
221251	Wet-Field® Hemostatic Eraser Bipolar Blunt Tip	84
221260	Wet-Field® Hemostatic Eraser Bipolar Blunt Tip	84
221265	Wet-Field® Hemostatic Eraser Bipolar 20G - 23G	85
221266	Wet-Field® Hemostatic Eraser Bipolar 20G - 23G	85
221267	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	84
221268	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	84
221270	Wet-Field® Osher Eraser	85
370145	Safety Slit Knife, 2.5 mm Slit, Angled, Single bevel	10
370146	Safety Slit Knife, 2.6 mm Slit, Angled, Single bevel	10
370820	Slit Knife, 2 mm Slit, 45°, Double bevel up	12
370822	Slit Knife, 2.2 mm Slit, 45°, Double bevel up	12
370824	Slit Knife, 2.4 mm Slit, 45°, Double bevel up	12
370825	Slit Knife, 2.5 mm Slit, 45°, Double bevel up	12
370826	Slit Knife, 2.65 mm Slit, 45°, Double bevel up	12
370827	Slit Knife, 2.75 mm Slit, 45°, Double bevel up	12
370828	Slit Knife, 2.8 mm Slit, 45°, Double bevel up	12
370830	Slit Knife, 3 mm Slit, 45°, Double bevel up	12
371301	Blade Wrench	23
371305	Collet Replacement, 0.5 cm Diameter	23
371350	Round Handle, 13 cm Long x 0.5 cm Diameter	23
371360	Round Handle, 13 cm Long x 0.66 cm Diameter	23
371394	Round Handle, 18 cm Long x 4 mm Diameter w/ Wrench	23
371395	Round Handle, 21 cm Long x 4 mm Diameter w/ Wrench	23
371410	Round Handle, 10 cm Long x 0.66 cm Diameter	23
371415	Collet Replacement, 0.66 cm Diameter	23
372242	Swivel Fixation Ring	59
373025	Kojo Safety Slit Knife, 2.5 mm Kojo Slit, 45°	11
373032	Kojo Safety Slit Knife, 3.2 mm Kojo Slit, 45°	11
373328	Slit Knife, 2.8 mm Slit, Straight, Matte	12
373330	Slit Knife, 3 mm Slit, Straight, Matte	12
373332	Slit Knife, 3.2 mm Slit, Straight, Matte	12
373710	Sideport Knife, 1 mm Sideport Knife, Angled	12
373712	Sideport Knife, 1.2 mm Sideport Knife, Angled	12
373718	Slit Knife, 1.8 mm Slit, 45°, Single bevel up	12
373722	Slit Knife, 2.2 mm Slit, 45°, Single bevel up	12
373724	Slit Knife, 2.4 mm Slit, 45°, Single bevel up	12
373725	Slit Knife, 2.5 mm Slit, 45°, Single bevel up	12
373726	Slit Knife, 2.6 mm Slit, 45°, Single bevel up	12
373727	Slit Knife, 2.75 mm Slit, 45°, Single bevel up	12
373728	Slit Knife, 2.8 mm Slit, 45°, Single bevel up	12
373730	Slit Knife, 3 mm Slit, 45°, Single bevel up	12
373732	Slit Knife, 3.2 mm Slit, 45°, Single bevel up	12
373807	Crescent Knife, 2.5 mm, 55°, Bevel Up	12
373808	Crescent Knife, 2.5 mm, 55°, Straight	12
373809	Crescent Knife, 2.5 mm, 55°, Bevel down	12
374100	Ring Cutter	22
374105	Ring Cutter Hub, Regular speed steel	22
374106	Ring Cutter Washer, Regular speed steel	22
374150	Ring Cutter Saw Blade, Regular speed steel	22
374151	Ring Cutter Saw Blade, High speed steel	22
374310	Round Handle, 10 cm Long x 0.5 cm Diameter	23
374311	Round Handle, 5 cm Long x 0.5 cm Diameter	23
374312	Round Handle, 7.5 cm Long x 0.5 cm Diameter	23
374320	Round Handle, 15.5 cm Long x 0.5 cm Diameter	23
374330	Hex Handle, 10 cm Long x 0.5 cm Diameter	23
374340	Hex Handle, 15.5 cm Long x 0.5 cm Diameter	23
374562	Mini-Blade®	15
374710	DeBakey Blade, 8 mm Cutting width	21
374810	DeBakey Blade, 10 mm Cutting width	21
374769	Mini-Blade® Knife	13
374881	Optimum Blade, 15°, Straight, Green	17
374882	Optimum Blade, 22.5°, Straight, Slate blue	17
374883	Optimum Blade, 30°, Black	17
374890	Optimum Knife, 15°, 30° bend left	14
374891	Optimum Knife, 15°, Straight, Green	14
374892	Optimum Knife, 22.5°, Straight, Slate blue	14
374893	Optimum Knife, 30°, Straight, Black	14
374894	Optimum Knife, 45°, Straight, Tan	14
375000	Handle Flat, 13 cm Long 1 cm Wide	23
375435	BeaverGuard™ Guarded Blade, 0.35 mm Guarded depth	15
375450	BeaverGuard™ Guarded Blade, 0.5 mm Guarded depth	15
375700	Sclerotome Blade	15
375701	Sclerotome Knife	13
357910	Needle Blade, 3 mm Cutting edge	21
357920	Needle Blade, 5 mm Cutting edge	21
375930	Needle Blade, 7 mm Cutting edge	21
375940	Needle Blade, 1.5 mm Cutting edge	21
376100	Mini-Blade®	15
376300	Mini-Blade®	15
376401	Mini-Blade®	15
376481	Arthro-Lok® Blade, 4 mm Rounded tip	20
376482	Arthro-Lok® Blade, 4 mm Rounded tip	20
376500	Mini-Blade®	15
376620	EdgeAhead Sideport MVR, 1.15 mm Sideport MVR, Angled	13
376630	EdgeAhead Sideport MVR, 1.15 mm MVR, Straight	13
376780	Arthro-Lok® Blade, 4 mm Pointed tip	20

Product Index

SKU#	Description	PG
376980	Arthro-Lok® Blade, 4 mm Mini-banana	20
376981	Arthro-Lok® Blade, 4 mm Mini-banana, Flat	20
376984	Arthro-Lok® Blade, 4 mm Banana	20
377100	Myringotomy Blade, Downcutting tip	19
377110	Myringotomy Blade, Upcutting tip	19
377120	Myringotomy Blade, Offset tip	19
377121	Myringotomy Blade, Narrow shaft	19
377122	Myringotomy Knife, Spear tip	19
377124	Myringotomy Knife, Spear tip	19
377200	Tympanoplasty Blade	22
377210	Tympanoplasty Blade	22
377300	Sickle Edge Blade	22
377511	Micro-Sharp® Blade, 1.5 mm, 15°, Orange	17
377512	Micro-Unitome Knife, 1.5 mm, 15°, Orange	14
377513	Micro-Sharp® Blade, 3 mm, 15°, Blue	17
377514	Micro-Unitome Knife, 3 mm, 15°, Blue	14
377515	Micro-Sharp® Blade, 5 mm, 15°, Green	17
377516	Micro-Unitome Knife, 5 mm, 15°, Green	14
377522	Micro-Sharp® Blade, 4 mm, 22.5°, Light blue	17
377523	Micro-Unitome Knife, 4 mm, 22.5°, Light blue	14
377530	Micro-Sharp® Blade, 3.5 mm, 30°, Black	17
377531	Micro-Unitome Knife, 3.5 mm, 30°, Black	14
377546	Micro-Unitome Knife, 3 mm, 45°, Tan	14
378210	Safety Sideport Knife, 1 mm Sideport MVR,	11
378212	Safety Sideport Knife, 1.2 mm Sideport MVR	11
378214	Safety Slit Knife, 1.4 mm Slit, Angled, Single bevel	10
378215	Safety Slit Knife, 1.5 mm Slit, Angled, Single bevel	10
378217	Safety Slit Knife, 1.7 mm Slit, Angled, Single bevel	10
378218	Safety Slit Knife, 1.8 mm Slit, Angled, Single bevel	10
378220	Safety Slit Knife, 2 mm Slit, Angled, Single bevel	10
378222	Safety Slit Knife, 2.2 mm Slit, Angled, Single bevel	10
378224	Safety Slit Knife, 2.4 mm Slit, Angled, Single bevel	10
378227	Safety Slit Knife, 2.75 mm Slit, Angled, Single bevel	10
378228	Safety Slit Knife, 2.8 mm Slit, Angled, Single bevel	10
378230	Safety Slit Knife, 3 mm Slit, Angled, Single bevel	10
378231	Safety Sideport Knife, 1.15 mm 20G Sideport	11
378232	Safety Slit Knife, 3.2 mm Slit, Angled, Single bevel	10
378233	Safety Sideport Knife, 0.6 mm 23G Sideport	11
378234	Safety Crescent Knife, 2.5 mm Crescent, 55°	11
378235	Optimum Safety Sideport Knife, 15° Sideport	11
378236	Optimum Safety Sideport Knife, 22.5° Sideport	11
378237	Optimum Safety Sideport Knife, 30° Sideport	11
378238	Optimum Safety Sideport Knife, 45° Sideport	11
378822	Safety Slit Knife, 2.2 mm Slit, Angled, Double bevel	10
378824	Safety Slit Knife, 2.4 mm Slit, Angled, Double bevel	10
378825	Safety Slit Knife, 2.5 mm Slit, Angled, Double bevel	10
378826	Safety Slit Knife, 2.65 mm Slit, Angled, Double bevel	10
378827	Safety Slit Knife, 2.75 mm Slit, Angled, Double bevel	10
378828	Safety Slit Knife, 2.8 mm Slit, Angled, Double bevel	10
378829	Safety Slit Knife, 2.85 Slit, Angled, Double bevel	10
378830	Safety Slit Knife, 3 mm Slit, Angled, Double bevel	10
378832	Safety Slit Knife, 3.2 mm Slit, Angled, Double bevel	10
379081	Arthro-Lok® Blade, 4 mm Mini-meniscus	20
379082	Arthro-Lok® Blade, 4 mm Mini-meniscus, Curved	20
581007	Superflex Extension Set	76
581012	Irrigation Handle	76
581019	Eye Shield (Left/Right)	74
581030	Bulk High Flow Filter	33
581032	Small Pupil IOL (Bobbitt)	60
581033	IOL Glide (Sheets)	60
581040	Universal Eye Shield	74
581047	Visitec® Visiwipe Instrument Wipe	72
581049	Surgical Glide (Fichman)	60
581056	Double Male Luer Adapter	76
581057	Corneal Marking Pad	76
581060	Absorbent Stick	72
581062	Visitec® Corneal Light Shield	65
581076	Visitec® Visiwipe Instrument Wipe	72
581098	LASIK Eye Drain (Chayet)	71
581146	Fluid Collection Bag with Wick	71
581147	Visitec® Eye Fluid Wick	65
581148	Visitec® Fluid Collection Bag	33
581149	Plastic Incise Drape with Fluid Collection Bag	67
581150	Non-Woven Incise Drape with Fluid Collection Bag	68
581151	Plastic Aperture Drape with Fluid Collection Bag	67
581152	Non-Woven Aperture Drape with Fluid Collection Bag	68
581153	Non-Woven Incise Drape with Channel Bag	68
581154	Mini Pre-Split Incise Drape	70
581155	Non-Woven Aperture Drape with Channel Bag	69
581156	Full Body Incise Drape	66
581157	Full Body Aperture Drape	66
581159	Mini Plastic Incise Drape	70
581160	Mini Plastic Aperture Drape	70
581161	Mini Plastic Aperture Drape	70
581162	Dual Refractive Drape	71
581163	Full Body Incise Drape with Plastic Head and Channel Bag	66
581164	Non-Woven Full Body Incise Drape	67
581166	Non-Woven Half Body Incise Drape* with Channel Bag	69
581177	Fluid Collection Bag	71
581273	Anterior Chamber (Rycroft), 30G x 7/8 in, 0.3 x 22 mm	27
581275	Retrobulbar (Sharp Point), 25G x 1-1/2 in, 0.5 x 38 mm	26
581276	Lacrimal Cannula, 26G x 1 1/8 in, 0.45 x 28 mm	38
581277	Anterior Chamber (Rycroft), 16G x 1 in, 1.6 x 25 mm	34
581280	Anterior Chamber (Rycroft), 27G x 7/8 in, 0.4 x 22 mm	27
581281	Anterior Chamber (Rycroft), 21G x 7/8 in, 0.8 x 22 mm	27
581284	Anterior Chamber Cannula, 27G x 7/8 in, 0.4 x 22 mm	27
581285	Anterior Chamber Cannula, 30G x 7/8 in, 0.3 x 22 mm	27
581413	Eye Speculum	43
581441	Micro-Manipulator	51
581442	Sinskey Hook Gen 2	51, 62
581443	Eye Speculum Open Wire Loops	43
581592	Nucleus Hydrodissector, 30G x 7/8 in, 0.3 x 22 mm	29
581594	Irrigating Cystotome, 23G x 5/8 in, 0.6 x 16 mm	32
581599	Irrigating Cystotome, 23G x 5/8 in, 0.6 x 16 mm	32
581602	Irrigating Side Cutting Cystotome	32
581604	Capsule Polisher (Kratz), 21G x 7/8 in, 0.8 x 22 mm	36
581605	Capsule Polisher (Kratz), 23G x 7/8 in, 0.6 x 22 mm	36
581610	Irrigating Cystotome, 25G x 5/8 in, 0.5 x 16 mm	32
581615	Lacrimal Cannula, 21G x 1 3/8 in, 0.8 x 35 mm	38
581617	Irrigating Cystotome, 25G x 5/8 in, 0.5 x 16 mm	32
581618	Irrigating Cystotome, 27G x 5/8 in, 0.4 x 16 mm	32
581630	Irrigating Cystotome, 30G x 5/8 in, 0.3 x 16 mm	32
581637	Retrobulbar (Atkinson), 25G x 1-1/2 in, 0.5 x 38 mm	26
581638	Retrobulbar (Atkinson), 23G x 1-1/2 in, 0.6 x 38 mm	26
581639	Anterior Chamber (Rycroft), 25G x 7/8 in, 0.5 x 22 mm	27
581640	Retrobulbar (Sharp Atkinson), 23G x 1-1/2 in, 0.6 x 38 mm	26
581641	Retrobulbar (Sharp Atkinson), 25G x 1-1/2 in, 0.5 x 38 mm	26
581642	Anterior Chamber Cannula, 25G x 7/8 in, 0.5 x 22 mm	27
581645	Posterior Scraper (Simcoe), 23G x 7/8 in, 0.6 x 22 mm	37
581705	ProOphta Eye Patch	72
581708	Universal Pediatric Eye Shield	74
581709	Visitec® LASIK PVA Expanded Eye Spear	64

Product Index

SKU#	Description	PG
581761	Visitec® LASIK PVA Eye Spears	64
585001	Retrobulbar/Peribulbar (Atkinson)	26
585003	Micro J Shape Cannula, 25G x 7/8 in, 0.5 x 22 mm	31
585004	Irrigating Cystotome, 25G 5/8 in, 0.5 x 16 mm	32
585005	Flexible Tip Cannula, 27G x 7/8 in, 0.4 x 22 mm	36
585006	Anterior Chamber Cannula Straight, 27G	27
585007	Viscoflow Formed, 27G x 7/8 in, 0.4 x 22 mm	34
585008	Irrigating Cystotome, 27G 5/8 in, 0.4 x 16 mm	32
585009	Viscoflow Cannula, 23G x 7/8 in, 0.6 x 22 mm	34
585015	Nucleus Hydrodelineator (P. Koch)	29
585016	Lacrimal Cannula, 19G x 3/4 in, 1.1 x 20 mm	38
585017	Flexible Tip Cannula, 25G x 7/8 in, 0.5 x 22 mm	36
585018	Micro J Shape Cannula, 27G x 7/8 in, 0.5 x 22 mm	31
585019	Retrobulbar, Curved (Atkinson/Straus), 25G	26
585021	Subretinal Fluid Cannula (Grizzard)	112
585025	Anterior Chamber (Rycroft), 23G x 7/8 in, 0.6 x 22 mm	27
585026	Viscoflow Cannula, 27G x 7/8 in, 0.4 x 22 mm	34
585027	Peribulbar (Atkinson), 25G x 7/8 in, 0.5 x 22 mm	26
585029	Micro J Shape Hydrodissector (Pearce), 27G	31
585030	Irrigating Cystotome (Joyce), 25G x 5/8 in, 0.5 x 16 mm	25
585031	Lacrimal Trephine (Sisler), 21G x 1-1/2 in, 0.8 x 38 mm	38
585033	Lacrimal Cannula (Shahinian), 23G x 1/2 in, 0.6 x 12 mm	38
585034	Irrigating Cystotome (Reverse), 25G x 5/8 in, 0.5 x 16 mm	33
585036	Nucleus Hydrodissector (Pearce) with 8 mm bend	29
585037	Nucleus Hydrodissector (Pearce), 25G x 7/8 in, 0.5 x 22 mm	29
585039	Irrigating Cystotome (Joyce), 27G x 5/8 in, 0.4 x 16 mm	33
585040	Bishop Harmon Cannula, 19G x 1 in, 1.1 x 25 mm	34
585041	Flexible Tip Cannula, 23G x 7/8 in, 0.6 x 22 mm	36
585042	Irrigating Cystotome (Blumenthal), 27G x 5/8 in, 0.4 x 16 mm	33
585043	Irrigating Cystotome (Blumenthal), 27G x 5/8 in, 0.4 x 16 mm	32
585045	Viscoflow Cannula, 25G x 7/8 in, 0.5 x 22 mm	34
585046	Viscoflow Cannula, 30G x 7/8 in, 0.3 x 22 mm	34
585047	I and A Cannula Thinwall (Simcoe), 23G, 0.6 mm	36
585049	Irrigating Formed Cystotome, 25G x 5/8 in, 0.5 x 16 mm	32
585061	Self Retaining ACM (Lewicky), 20G x 9/64 in, 0.9 x 3.5 mm	28
585062	IOL Manipulator Blunt Tip	51
585068	Lacrimal Cannula, 25G x 1/2 in, 0.5 x 13 mm	38, 92
585075	Cortex Extractor, 23G 5/8 in, 0.6 x 16 mm	35
585077	Nucleus Rotator	51
585086	I and A Cannula (Simcoe), 23G, 0.6 mm	36
585090	HydroExpression Cannula (Corydon), 25G	31
585093	Irrigating Cystotome (Simcoe), 27G	33
585099	Nucleus Hydrodissector (Sauter)	29
585106	Cortex Extractor, 24G x 1/2 in, 0.55 x 12.5 mm	35
585107	Tapered Hydrodelineator (Blumenthal)	30
585108	Retrobulbar/Peribulbar (Atkinson), 25G x 1-1/4 in	26
585109	Retrobulbar/Peribulbar (Atkinson), 27G x 1-1/4 in	26
585110	I and A Cannula (Simcoe), 23G, 0.6 mm	36
585113	Irrigating Cystotome (Blumenthal), 25G	33
585115	Blunt Tip Irrigating Cystotome (Kellan)	33
585133	Glaucoma Revision Pick 60° (Lederer)	37
585135	Glaucoma Revision Pick 30° (Lederer)	37
585143	Olive Tip Polisher, 23G x 1 in, 0.6 x 25 mm	37
585144	Peribulbar (Atkinson), 23G x 7/8 in, 0.6 x 22 mm	26
585145	Peribulbar (Atkinson), 27G x 7/8 in, 0.4 x 22 mm	26
585147	Olive Tip Cannula, 23G x 7/8 in, 0.6 x 22 mm	37
585148	Micro J Shape Cannula (Hessburg)	31
585149	Oval ACM (Blumenthal)	28
585155	Nucleus Hydrodissector, 27G 7/8 in, 0.4 x 22 mm	30
585157	Cortical Cleaving Hydrodissector (Round End)	30
585158	Cortical Cleaving Hydrodissector (Flat End)	29
585159	Cortical Cleaving Hydrodissector (Bolger)	30
585163	ACM Micro Etched (Dahan)	28
585173	High Viscosity Injector (Hammer)	112
585176	Sub-Tenon's Anesthesia Cannula	27
585181	Nucleus Chopper (P. Koch)	49
585184	High Viscosity Injector (Hammer)	112
585191	HydroExpression Cannula (Corydon)	31
585193	Cortex Extractor, 23G x 5/8 in, 0.6 x 16 mm	35
585195	Nucleus Hydrodissector (Sinskey)	29
585197	Cortical Cleaving Hydrodissector (Grabow)	31
585198	Hydrodissection Cannula (Akahoshi Style)	31
585199	Nucleus Hydrodissector (Buratto)	30
585200	FLACS Cannula (Schweitzer)	39
585205	Vitreoretinal Cannula, 20G x 1-1/8 in, 0.9 x 28 mm	112
585208	High Viscosity Injector (Hammer)	112
585209	Micro Iris Hook	62
585213	Spatula PRK	56
585222	Bimanual Aspiration Cannula, 23G	35
585223	Bimanual Aspiration Cannula, 21G	35
585226	Soft Tip Cannula, 23G x 1-1/3 in, 0.6 x 33 mm, 0.3 mm tip	112
585227	Cortex Extractor, 23G x 15/16 in, 0.60 x 23.5 mm	27
585228	High Viscosity Injector (Hammer)	112
585229	Soft Tip Cannula, 25G x 1-1/3 in, 0.5 x 33 mm, 0.3 mm tip	112
585231	EdgeAhead Sideport MVR, 0.6 mm MVR, Straight	13
585232	EdgeAhead Sideport MVR, 0.6 mm Stiletto	13
585233	LASIK Flap Roller	56
585240	EdgeAhead Sideport MVR, 1.1 mm MVR, Straight	13
585250	LASIK Irrigating Cannula (Lindstrom)	39
585252	LASIK Flap Spatula/Irrigator (Hersh)	39
585274	Nucleus Hydrodissector (Pearce), 27G x 7/8 in, 0.4 x 22 mm	29
585278	LASIK Cannula, 27G x 7/8 in, 0.4 x 22 mm	39
585280	LASIK Irrigating Cannula (Solomon)	39
585294	Nucleus Hydrodissector (Chang), 27G x 7/10 in, 0.4 x 18 mm	30
585295	Nucleus Hydrodissector (Chang), 27G x 7/10 in, 0.4 x 18 mm	30
585297	Plastic Locking Eye Speculum	42
585305	NuVisc™ OVD 1.2% Sodium Hyaluronate	118
587000	I-Ring® Pupil Expander, 6/Pk	62
587001	I-Ring® Pupil Expander, 1/Pk	62
587002	Soft Tip Cannula, 23G x 1-1/3 in, 0.6 x 33 mm, 0.8 mm tip	112
587003	Soft Tip Cannula, 25G x 1-1/3 in, 0.5 x 33 mm, 0.8 mm tip	112
8310080	Nurse-Mate® Six-Pack Diagnostic Penlight	75
8310201	Maxi-Lite® Diagnostic Penlight	75
8320150	Blu-Spot® Diagnostic Light	75
8440000	Accu-Temp® Variable Temperature Cautery	83
8441000	Accu-Temp® Low Temperature Cautery	82
8441500	Accu-Temp® High Temperature Cautery	82
8442000	Accu-Temp® High Temperature Cautery	82
8443000	Accu-Temp® High Temperature Cautery	82
8444000	Accu-Temp® High Temperature Cautery	82
8445000	Accu-Temp® High Temperature Cautery	82
8446000	Accu-Temp® Variable Temperature Cautery	83
8461000	Accu-Temp® Reusable Handle	83
8461020	Accu-Temp® Reusable Handle	83
8461050	Single-use fine tip	83
8461060	Single-use loop tip	83
8590300	Guibor Eye Bandage	57
8590350	Guibor Canaliculus Intubation Set (CIS)	57
8590450	Canaliculus Intubation Set (CIS)	57
9550005	9.5 mm Trephine Blade	18
9550006	10.5 mm Trephine Blade	18
9550013	11.5 mm Trephine Blade	18

Product Index

SKU#	Description	PG
9550025	9.25 mm Trephine Blade	18
9550026	4 mm Trephine Blade	18
9550028	12 mm Trephine Blade	18
9550029	13 mm Trephine Blade	18
22135010	Wet-Field® Disposable Cable	85
22155010	Wet-Field® Disposable Cable Shrouded	85
221260135	Wet-Field® Hemostatic Eraser Bipolar Blunt Tip	84
221260150	Wet-Field® Hemostatic Eraser Bipolar Blunt Tip	84
221267135	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	84
221267150	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	84
221268150	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	84
221269135	Wet-Field® Hemostatic Eraser Bipolar Fine Tip	84
40400-8	Ultracell® PVA Eye Spears	64
40405-8	Ultracell® PVA Eye Spears	64
A12-102	Collagen Synthetic Punctal Occluder	91
A12-103	Collagen Synthetic Punctal Occluder	91
A12-104	Collagen Synthetic Punctal Occluder	91
A12-200	Collagen Synthetic Punctal Occluder	91
A13-200	Extend® Synthetic Punctal Occluder	90
A13-2000	Extend® 180 Synthetic Punctal Occluder	90
A13-202	Extend® Synthetic Punctal Occluder	90
A13-203	Extend® Synthetic Punctal Occluder	90
A13-2033	Extend® 180 Synthetic Punctal Occluder	90
A13-2033SP	Extend® 180 Synthetic Punctal Occluder	90
A13-204	Extend® Synthetic Punctal Occluder	90
A13-2044	Extend® 180 Synthetic Punctal Occluder	90
A13-2044SP	Extend® 180 Synthetic Punctal Occluder	90
A13-205	Extend® Synthetic Punctal Occluder	90
A13-2055	Extend® 180 Synthetic Punctal Occluder	90
A13-2055SP	Extend® 180 Synthetic Punctal Occluder	90
A14-200	Parasol® Pre-Loaded Punctal Occluder Sterile	88
A14-201	Parasol® Pre-Loaded Punctal Occluder Sterile	88
A14-201BN	Parasol® Punctal Occluder Non-Sterile	88
A14-201QL	Parasol® Quickload Punctal Occluder Sterile	88
A14-202	Parasol® Pre-Loaded Punctal Occluder Sterile	88
A14-202BN	Parasol® Punctal Occluder Non-Sterile	88
A14-202QL	Parasol® Quickload Punctal Occluder Sterile	88
A14-203	Parasol® Pre-Loaded Punctal Occluder Sterile	88
A14-203BN	Parasol® Punctal Occluder Non-Sterile	88
A14-203QL	Parasol® Quickload Punctal Occluder Sterile	88
A14-204	Plus Pre-Loaded Punctal Occluder Sterile	89
A14-204BN	Plus Punctal Occluder Non-Sterile	89
A14-205	Plus Pre-Loaded Punctal Occluder Sterile	89
A14-205BN	Plus Punctal Occluder Non-Sterile	89
A14-205QL	Parasol® Quickload Punctal Occluder Sterile	88
A14-206	Plus Pre-Loaded Punctal Occluder Sterile	89
A14-206BN	Plus Punctal Occluder Non-Sterile	89
A14-207	Plus Pre-Loaded Punctal Occluder Sterile	89
A14-207BN	Plus Punctal Occluder Non-Sterile	89
A14-209BN	Parasol® Punctal Occluder Non-Sterile	88
A14-300	Reusable inserter	89
A14-401	Micro Flow™ Punctal Occluder Sterile	90
A14-402	Micro Flow™ Punctal Occluder Sterile	90
A14-403	Micro Flow™ Punctal Occluder Sterile	90
A20-401BN	Micro Flow™ Punctal Occluder Non-sterile	90
A20-402BN	Micro Flow™ Punctal Occluder Non-sterile	90
A20-403BN	Micro Flow™ Punctal Occluder Non-sterile	90
AK6000U	Leos® Accessory Kit	96
BEAVER6400	Grindless Mini-Blades®	16
BEAVER6600	Grindless Mini-Blades®	16
BEAVER6610	Grindless Mini-Blades®	16
BEAVER6700	Grindless Mini-Blades®	16
BEAVER6900	Grindless Mini-Blades®	16
BEAVER6910	Grindless Mini-Blades®	16
BF00.D51	Backflush reservoirs, Passive	108
BF00.D52	Backflush reservoirs, Active	108
BF00.D53	Disposable tubing for active aspiration	108
BF00.R01	Universal Backflush Instrument	108
BF00.R02	Universal Backflush Instrument with Luer-lock Connector	108
BF23.D01	Standard 23G Blunt backflush instrument	107
BF23.D02	Standard 23G Brush backflush instrument	107
BF23.D05	Retractable 23G Retractable brush backflush instrument	107
BF23.D07	Extendable 23G Silicone tipped backflush instrument	107
BF23.D51	23G BF Blunt needles	108, 113
BF23.D52	23G BF Brush needles	108, 113
BF25.D01	Standard 25G Blunt backflush instrument	107
BF25.D02	Standard 25G Brush backflush instrument	107
BF25.D05	Retractable 25G Retractable brush backflush instrument	107
BF25.D07	Extendable 25G Silicone tipped backflush instrument	107
BF25.D51	25G BF Blunt needles	108, 113
BF25.D52	25G BF Brush needles	108, 113
BF27.D01	Standard 27G Blunt backflush instrument	107
BF27.D02	Standard 27G Brush backflush instrument	107
CN23.D01	23G Silicone tipped cannula	113
CN23.D04	23G Pick needle	113
CN23.D05	23G Blunt needle with luer-lock connector	108, 113
CN25.D01	25G Silicone tipped cannula	113
CN25.D04	25G Pick needle	113
CN25.D05	25G Blunt needle with luer-lock connector	108, 113
CN27.D03	27G VFI cannula	113
CN27.D05	27G Blunt needle with luer-lock connector	108, 113
DC6000U	Leos® Dust Cover	96
ELP19CU5	VueProbe 19G Ophthalmic Laser Endoscope, Curved	95
ELP19SU5	VueProbe 19G Ophthalmic Laser Endoscope, Straight	95
FC00.R03	Reusable Cleaning Adapter	102
FC23.D02	Serrated Gripping 23G Single-Use	100
FC23.D03	ILM End-gripping 23G Single-Use	100
FC23.D04	Asymmetrical Gripping 23G Single-Use	100
FC23.D05	Eckardt Endgripping 23G Single-Use	101
FC23.D06	Kuhn 45° Power Gripping 23G Single-Use	101
FC23.D07	Shah X-tra Grip 23G Single-Use	101
FC23.R02	Serrated Gripping 23G Reusable	100
FC23.R03	ILM End-gripping 23G Reusable	100
FC23.R04	Asymmetrical Gripping 23G Reusable	100
FC23.R05	Eckardt Endgripping 23G Reusable	101
FC23.R06	Kuhn 45° Power Gripping 23G Reusable	101
FC23.R07	Shah X-tra Grip 23G Reusable	101
FC25.D02	Serrated Gripping 25G Single-Use	100
FC25.D03	ILM End-gripping 25G Single-Use	100
FC25.D04	Asymmetrical Gripping 25G Single-Use	100
FC25.D05	Eckardt Endgripping 25G Single-Use	101
FC25.D06	Kuhn 45° Power Gripping 25G Single-Use	101
FC25.D07	Shah X-tra Grip 25G Single-Use	101
FC25.R02	Serrated Gripping 25G Reusable	100
FC25.R03	ILM End-gripping 25G Reusable	100
FC25.R04	Asymmetrical Gripping 25G Reusable	100
FC25.R05	Eckardt Endgripping 25G Reusable	101
FC25.R06	Kuhn 45° Power Gripping 25G Reusable	101
FC25.R07	Shah X-tra Grip 25G Reusable	101
FC27.D03	ILM End-gripping 27G Single-Use	100
FC27.D05	Eckardt Endgripping 27G Single-Use	101
FC27.D07	Shah X-tra Grip 27G Single-Use	101

Product Index

SKU#	Description	PG
FC27.R03	ILM End-gripping 27G Reusable	100
FC27.R05	Eckardt Endgripping 27G Reusable	101
FC27.R07	Shah X-tra Grip 27G Reusable	101
FS6000	Leos® Wireless Footswitch	96
LF00.R01	Constellation® Light Source Adaptor	115
LF00.R02	Stellaris® Light Source Adaptor	115
LF00.R03	Geuder® Light Source Adaptor	115
LF00.R04	Ruck® Light Source Adaptor	115
LF00.R05	DORC® Light Source Adaptor	115
LF23.D01	23G Standard Light Fiber	110
LF23.D03	23G 45° Shielded Light Fiber	110
LF25.D01	25G Standard Light Fiber	110
LF25.D03	25G 45° Shielded Light Fiber	110
LF27.D01	27G Standard Light Fiber	110
LF27.D03	27G 45° Shielded Light Fiber	110
LF29.D02	29G Spotlight Directional Chandelier	98
LG001	Diode Laser Safety Goggles	96
LIT6000U	Leos® literature pack	96
LP23.A04	Laser Probe 23G Compatible with Alcon®	104
LP23.A05	Illuminated Laser Probe 23G Compatible with Alcon®	104
LP23.I04	Laser Probe 23G Comp. with Iridex® and B&L	104
LP23.I05	Illuminated Laser Probe 23G Comp. with Iridex® and B&L	104
LP25.A04	Laser Probe 25G Compatible with Alcon®	104
LP25.A05	Illuminated Laser Probe 25G Comp. with Alcon®	104
LP25.I04	Laser Probe 25G Compatible with Iridex® and B&L	104
LP25.I05	Illuminated Laser Probe 25G Comp. with Iridex® and B&L	104
LP27.A04	Laser Probe 27G Comp. with Alcon®	104
LP27.I04	Laser Probe 27G Comp. with Iridex® and B&L	104
MD23.D01	23G PFC injection needle	113
MD25.D01	25G PFC injection needle	113
MMSU0068S	Straight Rampley	54
MMSU0085S	Epilation Forceps	46
MMSU0214S	Stevens Tenotomy, Straight, Rounded	52
MMSU0238S	Adson Forceps, Toothed	46
MMSU0241S	Adson Forceps, Non-Toothed	46
MMSU0244S	Iris Scissors, Curved, Pointed	52
MMSU0249S	Iris Scissors, Straight, Pointed	52
MMSU0277S	Moorfields Forceps	46
MMSU0279S	St Martins Forceps	46
MMSU0299S	Iris Scissors, Straight, Rounded	52
MMSU0283S	Jewellers Forceps, Fine Angled	46
MMSU0286S	Jewellers Forceps, Fine Straight	46
MMSU0454S	Dressing Scissors, Straight, Pointed	52
MMSU0455S	Dressing Scissors, Straight, Rounded	52
MMSU1102S	Govan Manipulator	50
MMSU1103S	Sinskey Hook 0.2 mm	50
MMSU1105S	Claw Chopper	48
MMSU1106S	Rosen Chopper	48
MMSU1108S	Phaco Chopper	48
MMSU1110CS	Corneal Needle-Holder with contoured clips	54
MMSU1110S	Corneal Needle-Holder	54
MMSU1113S	Corneal Epithelium Spatula Flat	55
MMSU1114S	Fine Notched Forceps, Short (Hoskins)	44
MMSU1115CS	Notched Forceps (Hoskins) with contoured clips	44
MMSU1115S	Notched Forceps (Hoskins)	44
MMSU1119S	Greens Chopper, Right Handed	48
MMSU1123S	Fine Toothed Forceps Short (Bonn)	44
MMSU1124CS	Toothed Forceps (Bonn) with contoured clips	44
MMSU1124S	Toothed Forceps (Bonn)	44
MMSU1125CS	Straight Tying Forceps with contoured clips	47
MMSU1125S	Straight Tying Forceps	47

SKU#	Description	PG
MMSU1127S	Angled Tying Forceps (McPherson)	47
MMSU1129S	Colibri Toothed Forceps, Short	44
MMSU1134S	Y Shape Corneal Marker	55
MMSU1140CS	Vannas Scissors, Curved with contoured clips	53
MMSU1140S	Vannas Scissors, Curved	53
MMSU1141CS	Vannas Scissors, Straight	53
MMSU1141S	Vannas Scissors, Straight	53
MMSU1144S	Drysdale Manipulator	50
MMSU1146S	Lens Loading Forceps	47
MMSU1147S	Sinskey Hook, Straight	50
MMSU1149S	Chalazion Clamp, Oval	57
MMSU1150S	Chalazion Clamp, Round	57
MMSU1152S	Mendez Gauge	58
MMSU1153S	ICL Manipulator	50
MMSU1157S	Chalazion Curette, Double, 1 mm/2 mm	57
MMSU1158S	Chalazion Curette, Double, 2 mm/3 mm	57
MMSU1169CS	LASIK Flap-Lifting Forceps	55
MMSU1169S	LASIK Flap-Lifting Forceps	55
MMSU1172S	LASIK Flap Separator/Elevator (Seibel)	55
MMSU1176S	LASIK Flap Separator, Flattened Tip	55
MMSU1177S	LASIK Flap Elevator Curved 8 mm	55
MMSU1179S	Mushroom Manipulator	50
MMSU1181S	Punctum Dilator, Fine/Medium	59, 92
MMSU1187S	Corneal Epithelium Spatula Round	55
MMSU1193S	Punctum Dilator/Seeker	59, 92
MMSU1196S	Epithelium Well	56
MMSU1197S	Epithelium Well	56
MMSU1206S	LASIK Flap Separator and Elevator, Binder	55
MMSU1207S	Lester Hook	50
MMSU1209S	Colibri Notched Forceps, Short	44
MMSU1210S	Curved Tying Forceps	47, 92
MMSU1211S	Westcott Scissors, Curved, Pointed	53
MMSU1212CS	Westcott Scissors, Straight Rounded with contoured clips	53
MMSU1212S	Westcott Scissors, Straight Rounded	53
MMSU1213S	Westcott Scissors Curved Dolphin Nose	53
MMSU1214S	Westcott Scissors, Straight, Pointed	53
MMSU1219S	Straight Tying Forceps, Short	47
MMSU1221S	Phaco Chopper/Mushroom Manipulator	49
MMSU1225S	Kuglen Hook, Straight	51
MMSU1228S	Koch Chopper, 1.4 mm	48
MMSU1231S	Horizontal Chopper/Kuglen Hook	49
MMSU1239S	Diablo Chopper/Flat Button Manipulator	49
MMSU1243S	Troutman Needle-Holder	54
MMSU1252S	Kratz-Barraquer Speculum	42
MMSU1253S	Kratz-Barraquer Speculum, Temporal	42
MMSU1275S	Barraquer Speculum, Flat Blade	42
MMSU1276S	Braunstein IVT Marker 3.5/4, Green	58
MMSU1277S	Pointed IVT Marker 3.5/4, Blue	58
MMSU1284S	Contact Lens Forceps	46
MMSU1285S	Double Ended IVT Guide Pad	58
MMSU1287S	Kratz-Barraquer Speculum, Curved	42
MMSU1290S	Mechanical Speculum, Open Blade (PC)	42
MMSU1294S	LASIK Flap Separator/Elevator Curved	56
MMSU1295S	LASIK Flap Separator/Elevator Straight	56
MMSU1297S	Reinstein Lenticule Separator	56
MMSU1405S	Schocket Double-Ended Scleral Depressor	60
MMSU1408S	Ball Tip Scleral Depressor	60
MMSU1417S	Golf Club Spud	59, 92
MMSU1419S	Punctum Dilator/Lacrimal Probe	59, 91
MMSU1428S	Diablo Manipulator	51

Product Index

SKU#	Description	PG
MMSU1433S	Kuglen Hook, Angled	48
MMSU1434S	Butcher Rotator	48
MMSU1440S	Koch Chopper, 1.4 mm, Thin Shaft	48
MMSU1458S	Lebuisson Manipulator	50
MMSU1459CS	IOL Cutter (Osher)	53
MMSU1459S	IOL Cutter (Osher)	53
MMSU1460S	IOL Removal Forceps	47
MMSU1474S	Barraquer Speculum, Solid Curved (PC)	43
MMSU1475S	Barraquer Speculum, Solid Blade (PC)	43
MMSU1479S	Barraquer Speculum, Solid 18 mm Blade (PC)	43
MMSU1480S	Reverse Dialling Hook	49
MMSU1488S	Mechanical Speculum Solid Blade	42
MMSU1489S	Neuhann Nucleus Divider	49
MMSU1495S	Fixation Ring	59
MMSU1497S	Reinstein Lenticule Separator (Hyperopic)	56
MMSU1567CS	MICS Capsulorhexis Forceps - Curved/Sharp	45
MMSU1567S	MICS Capsulorhexis Forceps - Curved/Sharp	45
MMSU1568S	MICS Capsulorhexis Forceps - Straight/Sharp	45
MMSU1568S	MICS Capsulorhexis Forceps - Straight/Sharp	45
MMSU1571CS	MICS Capsulorhexis Forceps - Straight/Grasping	45
MMSU1571S	MICS Capsulorhexis Forceps - Straight/Grasping	45
MMSU1573CS	MICS Capsulorhexis Forceps - Curved/Grasping	45
MMSU1573S	MICS Capsulorhexis Forceps - Curved/Grasping	45
MMSU1577S	Intravitreal Injection Guide	58
MMSU1578S	Rampley Short	54
MMSU1583S	Waqar Suture Removal Forceps	47
MMSU3013S	Castroviejo Needle Holder, Curved	54
MMSU3015S	Barraquer Needle-Holder	54
OME6000U	Leos® Laser Endoscopy System	94
SC23.D01	Vertical Micro Scissors 23G Single-Use	102
SC23.D02	Straight Micro Scissors 23G Single-Use	102
SC23.D03	Horizontal Curved Micro Scissors 23G Single-Use	102
SC23.R01	Vertical Micro Scissors 23G Reusable	102
SC23.R02	Straight Micro Scissors 23G Reusable	102
SC23.R03	Horizontal Curved Micro Scissors 23G Reusable	102
SC25.D01	Vertical Micro Scissors 25G Single-Use	102
SC25.D02	Straight Micro Scissors 25G Single-Use	102
SC25.D03	Horizontal Curved Micro Scissors 25G Single-Use	102
SC25.R02	Straight Micro Scissors 25G Reusable	102
SC25.R03	Horizontal Curved Micro Scissors 25G Reusable	102
SC27.D01	Vertical Micro Scissors 27G Single-Use	102
SC27.D03	Horizontal Curved Micro Scissors 27G Single-Use	102
SC27.R03	Horizontal Curved Micro Scissors 27G Reusable	102
SD6000U	Sterile Drape	96
SQ00.D01	Dye Squeezers	114
TM23.D01	23G Retractable tissue manipulator	114
TM25.D01	25G Retractable tissue manipulator	114
VS00.R02	Williamson Scleral Depressor	116
FINEVISION HP	Trifocal Hydrophobic	120

BVI® is refocusing the future of vision.

As one of the fastest-growing, diversified surgical ophthalmic businesses in the world, our purpose-built portfolio spans more than 90 countries.

We've set our sights on touching the lives of millions of patients affected by conditions such as cataracts, refractive error, glaucoma, retinal disease, and dry eye.

Unburdened by legacy or bureaucracy, we have developed our strategy around a simple concept — taking pride in delivering innovative solutions for our physicians and patients, based on their needs.

We trust and empower our associates to make decisions and solve problems because collaboration drives us. Valuing agility, simplicity, and transparency, we stay committed to listening to our customers, delivering for our patients, and keeping the future in focus.

BVI

THE FUTURE IN FOCUS



BVI
Email: customersupport@bvimedical.com

bvimedical.com

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2026

1925979-01 US